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Review: June/2017

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	4
Equipaments and Operational Inputs	4
Transportation and Storage Conditions.....	4
Special Care	5
Samples	5
Procedure	6
A . Preparation of Working Solutions	6
B . DNA/RNA Extraction	6
B.1 Protocol for sample volume of 200µL	6
1. Samples Lysis	6
2. Binding.....	7
3. Washing	7
4. Elution	7
B.2 Protocol for sample volume of 400µL	8
1. Samples Lysis	8
2. Binding	8
3. Whasing	8
4. Elution	9
Process Limitations	9
Bibliographic References	9
Customer Service	9
Universal Symbology	10

FUNCTION

Product developed for Viral DNA /RNA extraction and purification in serum, plasma and other biological fluids.

The nucleic acids obtained can be used for different applications like PCR, RT-PCR, Hibridization, Sequencing, etc.

PRINCIPLE OF ACTION

The **Bio Gene Extração de DNA/RNA Viral** kit is a product developed for extraction and purification of Viral DNA / RNA in biological samples.

The method used is the extraction by silica membrane. The process is carried out in 4 steps: 1) Cell Lysis: cell disruption for the release of nucleic acids; 2) Binding: selective binding of nucleic acid to silica membrane; 3) Washing: to remove residual impurities; 4) Elution: release the nucleic acid from silica membrane. At the end of the process we have concentrated nucleic acid with high purity.

Reagent	Presentation 50 Preps
R1	25 mL
R2	30 mL
R3	12 mL
R4	13 mL
R5	2 x 300 µg
R6	600 µL
R7	50 unid.
R8	150 unid.
R9	100 unid.

REAGENTS

- R1. Lysis Buffer:** Guanidine hydrochloride.
R2. Wash 1: Guanidine hydrochloride, TRIS-HCl.
R3. Wash 2: TRIS-HCl, SDS.
R4. RNase Free Water.
R5. Carrier RNA: Acid polyadenylic.
R6. Proteinase K: Enzyme proteinase K.
R7. Columns: Polypropylene tube with a silica membrane.
R8. Collectors Tubes (2 mL): Polypropylene tube.
R9. Collectors Tubes (1.5 mL): Polypropylene tube.

EQUIPMENT AND OPERATING INPUTS

Materials contained in the kit:

- Reagents described in the table above.
- Usage instructions (manual).

Materials required, but not included in the kit:

- 1- 96-100% ethanol
- 2- PBS 1X may be required for some samples
- 3- Micropipettes and sterile filter pipette tips (0.5-10 μ L, 10-100 μ L, 100-1000 μ L).
- 4- Microcentrifuge
- 5- Vortex
- 6- Heating block or water bath
- 7- Personal protective equipment (gloves, lab coat, glasses)

TRANSPORTATION AND STORAGE CONDITIONS

The kit may be transported between 15 and 30°C. Protect from light and avoid moisture.

The storage temperature is between 15 and 30°C.

After resuspended is recommended to store the Carrier RNA reagent (R5) at -20°C. After first use, is recommended to store the Proteinase K reagent (R6) -20°C.

SPECIAL CARE

1- For professional use only.

2- Strictly follow the methodology proposed to obtain accurate results.

3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose the waste in accordance with appropriate security procedures.

4- It is required skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.

5- It is necessary to have separate areas for the extraction / preparation of reactions and for the amplification / detection of products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.

6- Avoid repeated thawing and freezing of the reagents.

7- Do not use reagents after expiration date.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Kit for Viral DNA / RNA extraction from serum, plasma or other biological fluids. Samples must be collected and stored according to the laboratory's recommendations for molecular testing.

PROCEDURE

Attention:

The Lysis Buffer (**R1**) and Wash 1 (**R2**) reagents contains guanidine salt, and must be handled using suitable PPEs (coats, gloves and goggles). These reagents must not be disposed of on acid solution or on sodium hypochlorite due to the formation of reactive components.

A. Preparation of Working Solutions

Resuspend **R3** and **R5** reagents as shown in the table below:

Reagent	Resuspension Volume
R3 - Wash 2 (concentrated)	Add 48 mL of ethanol 96-100%
R5 - Carrier RNA*	Add 300µL of RNase Free Water (R4)

Before starting the extraction:

- 1- Note if the reagents Wash 2 (**R3**) and Carrier RNA (**R5**) are prepared.
- 2- Set the heating block or water bath at 56°C.
- 3- Preheat reagent RNase free water (**R4**) at 70°C.

B. DNA/RNA extraction

The product is able to purify samples with a volume of 200µL or 400µL. The following protocols are specific for each sample volume.

B.1 Protocol for sample volume of 200µL

1. Samples Lysis

- 1.1 Pipette 200µL of sample (serum, plasma or other biological fluids) at room temperature in Collector Tube 1.5 mL (**R9**).
- 1.2 Add 5µL of Proteinase K (**R6**)* to the Collector Tube 1.5 mL with sample and mix.
- 1.3 Then, pipette 200µL of Lysis Buffer (**R1**) and shake vigorously (10-15 seconds).
- 1.4 Add 5.6 µL of Carrier RNA (**R5**) prepared and mix.
- 1.5 Incubate for 3 minutes at room temperature (18-25°C).

* After use, store at - 20° C.

1.6 Add 200µL of ethanol (96-100%) to each Collector Tube 1.5 mL with sample and mix by vortexing (10-15 seconds).

1.7 Incubate for 5 minutes at room temperature (18-25°C).

2. Binding

2.1 Identify the Column **(R7)** into Collector Tube according to the sample which is being purified.

2.2 Transfer the total of lysed sample volume (~ 610µL) into the respective Column and then centrifuge for 3 minutes at 4,000 x g. If the Column is not completely dried, it must be centrifuged again (approx. 15,000 x g).

2.3 After centrifugation, transfer the Column to new Collector Tube **(R8)** and discard the Collector Tube with the filtrate.

3. Washing

3.1 Add 400µL of Wash 1 **(R2)** to each Column containing sample.

3.2 Centrifuge for 1 minute at 11,000 x g.

3.3 Transfer the Column to new Collector Tube **(R8)** and discard the Collector Tube with the filtrate.

3.4 Add 400µL of Wash 2 **(R3)** prepared.

3.5 Centrifuge for 1 minute at 11,000 x g.

3.6 Transfer the Column to new Collector Tube **(R8)** and discard the Collector Tube with the filtrate.

3.7 Add 200µL of Wash 2 **(R3)** to Column.

3.8 Centrifuge for 5 minutes at $\geq 15,000$ x g, and proceed to step 4. Elution.

4. Elution

4.1 Transfer the Column to a Collector Tube 1.5 mL **(R9)** and discard the Collector Tube with the filtrate.

4.2 Incubate the Column / Tube for 5 minutes at 56°C with the lids opened.

4.3 Add 30~60µL of RNase Free Water **(R4)** preheated to 70°C. Dispense directly in the center of the membrane.

4.4 Incubate at room temperature for 3 minutes.

- 4.5 Centrifuge for 3 minutes at $\geq 15,000 \times g$, and discard the Column.
- 4.6 Storing the DNA / RNA eluted at -20°C if not used immediately. For long time storage, it is recommended to stock DNA/RNA at -70°C

B.2 Protocol for sample volume of 400 μL

1. Samples Lysis

- 1.1 Pipette 400 μL of sample (serum, plasma or other biological fluids) at room temperature in Collector Tube 1.5 mL (**R9**).
- 1.2 Add 10 μL of Proteinase K (**R6**)* to the Collector Tube 1.5 mL with sample and mix.
- 1.3 Then, pipette 400 μL of Lysis Buffer (**R1**) and stir vigorously (10-15 seconds).
- 1.4 Add 11.2 μL of Carrier RNA (**R5**)* prepared and mix.
- 1.5 Incubate for 3 minutes at room temperature ($18-25^{\circ}\text{C}$).
- 1.6 Add 400 μL of ethanol (96-100%) to each well with sample and mix by vortexing (10-15 seconds).
- 1.7 Incubate for 5 minutes at room temperature ($18-25^{\circ}\text{C}$).

2. Binding

- 2.1 Identify the Column (**R7**) into Collector Tube according to the sample which is being purified.
- 2.2 Transfer the total of lysed sample volume ($\sim 610\mu\text{L}$) into the respective Column and then centrifuge for 3 minutes at $4,000 \times g$.
- 2.3 Discard the filter and replace the Collector Tube in the Column.
- 2.4 Transfer the remaining lysate to the Column and centrifuge again for 3 minutes at $4,000 \times g$. If the Column is not completely dried, it must be centrifuged again (approx. $15,000 \times g$ for 1 minute).
- 2.5 After centrifugation, transfer the Column to new Collector Tube (**R8**) and discard the Collector Tube with the filtrate.

3. Washing

- 3.1 Add 400 μL of Wash 1 (**R2**) to each Column containing sample.
- 3.2 Centrifuge for 1 minute at $11,000 \times g$.

* After use, store at -20°C .

3.3 Transfer the Column to new Collector Tube **(R8)** and discard the Collector Tube with the filtrate.

3.4 Add 400µL of Wash 2 **(R3)** prepared.

3.5 Centrifuge for 1 minute at 11,000 x g.

3.6 Transfer the Column to new Collector Tube **(R8)** and discard the Collector Tube with the filtrate.

3.7 Add 200µL of Wash 2 **(R3)** to Column.

3.8 Centrifuge for 5 minutes at $\geq 15,000 \times g$, and proceed to step 4. Elution.

4. Elution

4.1 Transfer the Column to a Collector Tube 1.5 mL **(R9)** and discard the Collector Tube with the filtrate.

4.2 Incubate the Column / Tube for 5 minutes at 56°C with the lid opened.

4.3 Add 30~60µL of RNase Free Water **(R4)** preheated to 70°C. Dispense directly in the center of the membrane.

4.4 Incubate at room temperature for 3 minutes.

4.5 Centrifuge for 3 minutes at $\geq 15,000 \times g$, and discard the Column.

4.6 Storing the DNA / RNA eluted at -20°C if not used immediately. For long time storage, it is recommended to stock DNA/RNA at -70°C

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

BIBLIOGRAPHIC REFERENCES

1- Ogura T, Tsuchiya A, Minas T, Mizuno S. Methods of high integrity RNA extraction from cell/agarose construct. BMC Res Notes. 2015; 8:644.

CUSTOMER SERVICE


Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for kit: 10269360296

SIMBOLOGIA UNIVERSAL

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	CAUTION		DAINGER

BIO GENE

Bioclin ·· QUIBASA

 Rua Teles de Menezes, 92 - Belo Horizonte - MG - Brasil - CEP: 31565-130
Tel +55 31 3439 5454 - Fax +55 31 3439 5455 - www.bioclin.com.br
FARM. RESP. Silvio Wandalsen Arndt - CRF MG 7422
C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira