BIOCONTROL N

REF K073

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FINALIDADE

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USAGE INSTRUCTIONS
LEVEL N - NORMAL

FUNCTION

Reagent used to monitor the accuracy and precision of analytical testing in manual and automated methodologies. Has the concentration of analytes representing the normal parameters. For *in vitro* diagnostics use only.

REAGENTS

Reagent N°1 - Control Serum - Store between 2 and 8°C. Contains: Control serum containing approximately 30 analytes embedded in a lyophilized human matrix and preservative. **Potentially Infectious.**

Warning: The concentration of analytes varies each lot. See table that accompanies the product.

PRESENTATION

Presentation	Volume
1	1 x 3 mL
2	6 x 3 mL
3	1 x 5 mL
4	5 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Automated and semi-automated biochemistry instruments, watch or stopwatch, pipettes, test tubes, kits and water bath at 37°C. Can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- The reagent should be handled cautiously since it is subject to biological contamination.
- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the product is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- The reagent has been tested for anti-HIV, anti- HCV and HBs Antigen using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially Infectious.**
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

RECONSTITUTION OF CONTROL SERUM

Before reconstitute the Control Serum should be noted the volume indicated on the bottle and proceed as described below.

Presentation 1 and 2: Open carefully the bottle and add exactly **3 mL of distilled water**, as volume described in the label of the bottle.

Presentation 3 and 4: Open carefully the bottle and add exactly **5 mL of distilled water**, as volume described in the label of the bottle.

Mix. Leave the bottle resting for 30 minutes. Every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid formation of foam. After mixing, observe the stability according to the table below.

	-20°C	+4°C	+25°C
Bilirubine (protected from the light)	14 days	6 hours	2 hours
CK and CK-MB	30 days	6 hours	2 hours
Other parameters	30 days	5 days	12 hours

Warning:

- The reconstituted Control Serum should not be thawed more than once. Therefore, for minor routines, fractionate the material into aliquots before freezing.
- The Control Serum should be used in an identical way as patient's serum.
- After pipetting, if there is leftover serum control, do not return to the original bottle.

REFERENCE VALUES

The values the enclosed were obtained through repeated measurements using the indicated method. The values provided should be used for guidelines, each laboratory should establish its own limits of precision. The laboratory's average should, however, be within the acceptable average for the given parameters.

See table that accompanies the product.

BIBLIOGRAPHIC REFERENCES

1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. Richtlinien der Bundesärztekammer zur Qualitätssicherung medizinischen Laboratorien. Deutsches Ärzte-blatt 1988;85: B519-B532.
4. Richtlinie der Bundesärztekammer zur Qualitätssicherung quantitativer laboratoriums-medicinischer Untersuchungen. Deutsches Ärzteblatt 2002;98:A 2747-59.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
 Rua Teles de Menezes, 92 - Santa Branca
 CEP 31565-130 - Belo Horizonte - MG - Brasil
 Phone: +55 (31) 3439.5454 - Fax +55 (31) 3439.5455
 E-mail: bioclin@bioclin.com.br
 CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53
 1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
 Phone.: 0800 0315454
 E-mail: sac@bioclin.com.br

ANVISA registration for the Biocontrol N kit: 10269360125

Review: January/2020

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED