FINALIDADE

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BIOCONTROL P

REF K074

USAGE INSTRUCTIONS LEVEL P - PATHOLOGICAL

FUNCTION

Reagent used to monitor the accuracy and precision of analytical testing in manual and automated methodologies. Has the concentration of analytes representing the pathological parameters. For use in *In vitro* Diagnostics only.

REAGENTS

Number 1 - Control Serum - Store between 2 and 8°C. Contains: Control serum containing approximately 30 analytes embedded in a lyophilized human matrix and preservative. **Potentially infectious.**

Warning: The concentration of analytes varies each lot. See table that accompanies the product.

PRESENTATION

Presentation	Volume
1	1 x 3 mL
2	6 x 3 mL
3	1 x 5 mL
4	5 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Automated or semi-automated biochemistry analyzers, watches or stopwatches, pipettes, test tubes, kits and water bath at 37°C. These items are found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

- For professional *in vitro* diagnostic use only.
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- The reagent should be handled cautiously since it is subject to biological contamination.
- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the product is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- The reagent has been tested for anti-HIV and anti-HCV antibodies and HBs Antigen using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. Potentially Infectious.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

RECONSTITUTION OF CONTROL SERUM

Before reconstitute the Control Serum should be noted the volume indicated on the bottle and proceed as described below.

Presentation 1 and 2: Open carefully the bottle and add exactly 3 mL of distilled water, as volume described in the label of the bottle.

Presentation 3 and 4: Open carefully the bottle and add exactly 5 mL of distilled water, as volume described in the label of the bottle.

Mix. Leave the bottle resting for 30 minutes. Every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid formation of foam.

After mixing, observe the stability according to the table below.

	-20°C	+4°C	+25°C
Bilirubin (protected from the light)	14 days	6 hours	2 hours
CK and CK-MB	30 days	6 hours	2 hours
Other parameters	30 days	5 days	12 hours

Warning:

- The reconstituted Control Serum should not be thawed more than once. Therefore, for minor routines, fractionate the material into aliquots before freezing.

- The Control Serum should be used in an identical way as patient's serum.

- After pipetting, if there is leftover serum control, do not return to the original bottle.

REFERENCE VALUES

The values enclosed were obtained through repeated measurements using the method indicated. The values provided should be used for orientation, each laboratory should establish its own limits of precision. The laboratory's average should, however, be within the acceptable average for the given parameters.

See table that accompanies the product.

BIBLIOGRAPHIC REFERENCES

1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsge-sellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien. Deutsches Ärzteblatt 1988;85: B519-B532.
4. Richtlinie der Bundesärztekammer zur Qualitätssicherung quantitativer laboratoriums-medicinischer Untersuchungen. Deutsches Ärzteblatt 2002;98:A 2747-59.
5. QUIBASA: Dados do departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Biocontrol P kit: 10269360125

Review: August/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED