BIOLISA ANTI HBs

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FINALIDADE

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NÚMERO DE TESTES

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1. Frank Fenner and David O. White, Medical Virology, 4th Edition, Academic Press, 1994.
2. Centers for Disease Control. Viral Hepatitis B Fact Sheet.
3. Richman, D., R. Whitley, F. Hayden.Clinical Virology. New York: Churchill Livingstone Inc., 1997.
4. World Health Organization. World Health Organization Hepatitis B Fact Sheet. N°204.Revised October 2000.
5. World Health Organization. Hepatitis B. 2002.
6. Bioclin - Dados de arquivos.

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1. Frank Fenner and David O. White, Medical Virology, 4th Edition, Academic Press, 1994.
2. Centers for Disease Control. Viral Hepatitis B Fact Sheet.
3. Richman, D., R. Whitley, F. Hayden.Clinical Virology. New York: Churchill Livingstone Inc., 1997.
4. World Health Organization. World Health Organization Hepatitis B Fact Sheet. N°204. Revised October 2000.
5. World Health Organization. Hepatitis B. 2002.
6. Bioclin – Datos de archivos

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BIOLISA ANTI HBs

REF K121
USAGE INSTRUCTIONS

FUNCTION
Test for the quantitative detection of antibodies to surface anti-antigen of Hepatitis B virus in serum or human plasma by enzyme immunoassay, microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION
Methodology: Enzyme immunoassay or immunoenzymatic
BIOLISA Anti HBs kit is a quantitative solid phase immunoassay based on the sandwich principle for the detection of Anti-HBs antibodies including IgG, IgM and IgA in serum or plasma human. The microplate is coated with recombinant HBsAg. During the test, sample and conjugate are added to the antigen coated microplate and then incubated. If the sample contains Anti-HBs antibodies, it will bind to the antigen coating the microplate, and simultaneously bind to the conjugate to form Antigen-Anti HBS-Conjugate complexes. After incubation, the microplate is washed to remove unbound materials. After this step, the Substrates A and B are added and then incubated producing a blue color, indicating the amount of Anti-HBs antibodies in the sample. One Stop Solution is added to stop the reaction having a color change to yellow, measured with a microplate reader.

REAGENTS
1- Reference Standard (A - E) - Store between 2 and 8°C. Five (5) bottles (A - E) of Reference Standards containing antibodies anti virus surface antigens Hepatitis B in different concentrations in buffer solution, stabilizer and preservative. **Potentially infectious.**
Reference Standard (A - E) concentration varies with the lot. See bottle label.
2- Conjugate - Store between 2 and 8°C. Anti-human IgG antibodies linked to Peroxidase, stabilizer and preservative.
3- Sensitized Plate - Store between 2 and 8°C.
4- Concentrated Washing - Store between 2 and 8°C. Buffer solution, surfactant and preservative.
5- Substrate A - Store between 2 and 8°C. Solution containing Tetramethylbenzidine (TMB).
6- Substrate B - Store between 2 and 8°C. Solution containing Urea Peroxide.
7- Stop Solution - Store between 2 and 8°C. Sulfuric Acid 1N.
8- Plate Sealers

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Reference Standards	1 Flask (A-E) x 1 mL	2 Flasks (A-E) x 1 mL	4 Flasks (A-E) x 1 mL
2- Conjugate	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
3- Sensitized Plate	1 Unit (96 cavities)	2 Units (96 cavities)	5 Units (96 cavities)
4- Concentrated Washing	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
5- Substrate A	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
6- Substrate B	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
7- Stop Solution	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
8- Plate Sealers	3 Units	5 Units	10 Units

EQUIPMENTS AND OPERATIONAL INPUTS
Materials in the kit:
- Reagents described in the previous item
- Usage Instructions (manual)

Materials needed but not contained in the Kit:
1- Pipette(s) capable of dispensing volumes of 5, 50 and 100 µL with a lower coefficient of variation than 1,5%.
2- Re-pipettor(s) for repetitive dispensing volumes of 100 µL and 300 µL, with a lower coefficient of variation than 1,5% or multichannel pipette (optional).
3- Microplate washer (optional).
4- ELISA reader capable of absorbance at 450/630 nm wavelength.
5- Adjustable volume pipettes (200 µL to 1000 µL) for preparation of the substrate.
6- Test tubes for dilution of the substrate A and B.
7- Paper towel to dry microcavities.
8- Watch or stopwatch.
9- Container to store the Washing Solution after diluted.
10- Distilled or deionized water.
11- Tools of Quality Control.
12- Incubator of 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS
The storage temperature should be from 2 to 8°C. The transport can be in under ambient temperature (up to 30°C) for up to 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE
1- For professional *in vitro* diagnostic use only.
2- Strictly follow the methodology proposed to obtain accurate results.
3- The envelope containing the strips should be opened only after reaching room temperature. Replace strips from unused microcavities in the aluminum bag, seal and store between 2 and 8°C.
4- The water used in material cleansing must be recent and free of contaminants.
5- Deionized saturated columns release alkaline water, several ions and oxidizing agents and reducers, that can significantly alter the results.
6- All the raw material of product is tested and should be non-reactive for Anti-HIV 1 & 2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manual manipulation of any product containing serum is potentially capable of transmitting diseases. Therefore, we must take due care concerning biosafety while handling these products.
7- Always add reagents in the same order to minimize the difference in reaction time between microcavities.
8- As a safety measure, you should cover the plate during the reaction.
9- Ensure that the bottom of the cavity is clean and dry and there are no bubbles in the liquid surface before reading the plate. Do not let the cavities run dry during the test procedure.
10- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
11- Stop Solution contains Sulfuric Acid, which is a strong acid. Handle it with care.
12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
14- Do not use the product in case of damaged packaging.
15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES
Use serum or plasma (EDTA or Heparin). Highly lipemic or hemolyzed samples should not be used. Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for 30 days at -20°C (freezer).

DESCRIPTION OF PROCESS

PREPARATION OF WORKING REAGENT
Washing Solution
Dilute the contents of the Flask N° 4 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. Can be stored at room temperature. If crystallization appears, heat at 37°C until dissolved.

Substrate - Working Solution
Determine the amount of cavities that are being used to prepare an appropriate volume. Prepare the solution by mixing equal parts of Substrate A and Substrate B, about 15 minutes before your use. Keep it protected from light until used
To each microcavity (test), use:
50 µL of Substrate A + 50 µL of Substrate B
For example: Mix 1 mL of Substrate A to 1 mL of Substrate B for two strips of 8 microcavities (16 tests). Leftover reagent occurs.
Use no later than one (1) hour after preparation.

TECHNIQUE
Before starting the assay, bring all Reagents, Samples and Reference Standards to stabilize at room temperature (15 - 30°C) for at least 40 minutes.
Return unused strips to the original sealed packaging.
1- Select the cavities to be used considering: Reference Standards and Samples (which may be tested in duplicate).
2- Select the first cavity for Blank (OPTIONAL).
3- Pipette 50 µL of Reference Standards (A - E) and Samples in the cavities previously determined.
4- Pipette 50 µL of Conjugate into all cavities except in the cavity Blank (if you made the choice to use Blank).
5- Mix gently for ± 30 seconds. Cover cavities with the sealer.
6- Incubate for 60 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
7- Remove the plate sealer.
8- Discard the contents of the cavities by aspiration (Washer) or by decanting (manual).
Use approximately 350 µL of Washing Solution, **previously prepared***, to perform a total of eight (8) wash cycles.
Or
Use approximately 500 µL of Wash Solution, **previously prepared***, for a total of six (6) wash cycles (exclusive to automatic equipment).
Note: Evaluate the aspiration and dispensing capacity of the equipment. To ensure the drying of the plate at the end of washing, beat it for a few seconds on absorbent paper.

9- ATTENTION Follow one of the following procedures:
A) Pipette 100 µL of **previously prepared*** Substrate - Working Solution (A+ B) in all cavities.
***See PREPARATION OF WORKING REAGENTS**
Or
B) Pipette 50 µL of Reagent A in all cavities.
Pipette 50 µL of Reagent B in all cavities.
10- Mix gently for ± 30 seconds. Cover the cavities with plate sealer.
11- Incubate for 30 minutes ± 1 minute in room temperature (15 - 30°C) and protected from light.
12- Pipette 100 µL of Stop Solution in all cavities.
13- Mix gently for ± 30 seconds.
14- Reading: 450 nm (primary filter) / 630 nm (secondary filter) within 30 minutes (maximum).

CALCULATIONS
A calibration curve is used to determine the concentration of Anti-HBs in unknown samples.

Preparation of Calibration Curve
Record the absorbance obtained in the microplate Reader, as outlined in example 1. Calculate averages of duplicates, if carried out duplicates. Plot the absorbance of each Reference Standard versus the corresponding concentration in mUI/mL on graph paper (before you plot them on the graph). Plot the curve.

Note: The data presented in example 1 are for illustration only and may not be used in lieu the calibration curve, which must be constructed in the laboratory.

STANDARD	ABSORBANCE	AVERAGE ABSORBANCE	CONCENTRATION
A	0,026	0,027	0,0 mUI/mL
A	0,029		
B	0,101	0,102	10,0 mUI/mL
B	0,103		
C	0,660	0,665	100,0 mUI/mL
C	0,670		
D	1,700	1,734	300,0 mUI/mL
D	1,768		
E	2,214	2,217	500,0 mUI/mL
E	2,220		

Note: Samples that have absorbance above the Reference Standard must be pre-diluted and retested. The concentration must be multiplied by the dilution factor. Automatic reading and calculation can be performed through the linear regression function in an appropriate computer program.

INTERPRETATION OF RESULTS

Not Reactive: Samples with concentration lower than 10 mUI/mL are considered negative for Anti-HBs.

Reactive: Samples with concentration higher than 10 mIU/mL are considered positive for Anti-HBs.

Note: In case of results between 9 to 11 mIU/mL, the sample must be retested. The samples that obtain result repeatedly between 9 and 11 mUI/ mL should be retested using an alternative method. If results are confirmed, you must collect a new sample in two weeks. The result of the last sample collected should prevail.
The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

Note: The data presented in the examples are for ilustration only and can be used for calculation of the results.

PROCEDURE LIMITATIONS
The interpretation of a diagnostic test, there should not be based on a single assay. Must include other confirmatory testing before a sample is considered positive. One negative result does not exclude the possibility of exposure. Finally, all results must be interpreted in conjunction with other clinical information available.

INTERNAL QUALITY CONTROL
The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PARAMETER OF QUALITY CONTROL

Absorbance of Reference Standard E (500 mUI /mL) ≥ 1,2.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND SPECIFICITY METHODOLOGY
BIOLISA Anti HBs kit was compared with another commercially available ELISA. 07 analyzes were performed and the results were evaluated. The linear equation obtained was Y = 0,866X + 4,516 and correlation coefficient of 0,999. With these results we can conclude that the kit shows good methodological specificity.

Accuracy

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	401	23	104
Standard Deviation	16,4	2,3	4,8
Coefficient of Variation (%)	4,10	10,07	4,58

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	402,67	24,33	101,67
Standard Deviation	3,79	1,15	2,08
Coefficient of Variation (%)	0,94	4,75	2,05

Analytical Sensitivity

The analytical sensitivity is 1 mUI/mL.

Clinical Sensitivity and Specificity

BIOLISA Anti HBs kit examined clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of BIOLISAAnti HBs kit is > 99,9% and clinical specificity of 99,4%.

	TOTAL EXPECTED	BIOLISA ANTI HBs
Negative Sample	513	510
Positive Sample	213	213

Clinical Sensitivity: > 99,9% (213/213)

Clinical Specificity: 99,4% (513/510)

Linearity

The reaction is able to detect concentrations up to the highest point of concentration on the curve. For samples with higher values dilute it, repeat the dosage and multiply the results by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Hepatitis B virus envelope is a spherical, double-stranded DNA, from Hepadnaviridae family. Hepatitis B is transmitted through sexual contact, blood exposure, transmission from mother to child during childbirth or sharing objects cutting perforation. HBV infection has been associated with a variety of liver diseases from mild to chronic, including cirrhosis and hepatocellular carcinoma. In some cases, the virus can persist throughout life. Annually, 1 million people die of chronic hepatitis active, cirrhosis or primary liver cancer. Hepatitis B infects millions of people worldwide and is considered a public health problem. HBsAg is one of the first markers to appear in blood after infection with Hepatitis B virus. The immune response to infection includes the development antibody specific for HBsAg. These

antibodies appear within weeks of being HBsAg detected in blood, when the virus no longer can be passed to other people. The appearance of the Anti-HBs is associated with the recovery and is considered a marker for immunity. These antibodies also appear as a result of successful immunization. Therefore, detection and monitoring antibodies against HBsAg became an important tool in screening and monitoring of infected individuals, as well as an indicator of success of those vaccinated against HBV.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Frank Fenner and David O. White, Medical Virology, 4th Edition, Academic Press, 1994.
2. Centers for Disease Control. Viral Hepatitis B Fact Sheet.
3. Richman, D., R. Whitley, F. Hayden. Clinical Virology. New York: Churchill Livingstone Inc., 1997.
4. World Health Organization. World Health Organization Hepatitis B Fact Sheet. N°204. Revised October 2000.
5. World Health Organization. Hepatitis B. 2002.
6. Bioclin – Dados de arquivos

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA Anti HBs kit: 10269360289

Review: August/2019

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED