BIOLISA CMV IgM

REF **K123**

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FINALIDADE

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- Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
- Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after “Silent” Congenital Cytomegalovirus Infection. N. Engl. J. Med. (1976) 295:468-470.
- Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovirus. N. Engl. J. Med. (1974) 790:291-296.
- Stern, H. Cytomegalovirus Vaccine: Justification and Problems. In: Waterson AP (ed.) Recent Advances in Clinical Virology (1977) 117-134.
- Bioclin - Dados de arquivos

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- Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
- Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. *School Failure and Deafness after “Silent” Congenital Cytomegalovirus Infection*. N. Engl. J. Med.(1976) 295:468-470.
- Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. *Inapparent Congenital Cytomegalovirus*. N. Engl. J. Med. (1974) 790:291-296.
- Stern, H. Cytomegalovirus Vaccine: Justification and Problems. In: Waterson AP (ed.) Recent Advances in Clinical Virology (1977) 117-134.
- Bioclin – Datos de archivos

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REF K123

USAGE INSTRUCTIONS

FUNCTION
Test for the qualitative determination of IgM antibodies to Cytomegalovirus (CMV) in serum or human plasma by enzyme immunoassay microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION
Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA CMV IgM kit is a solid phase enzyme immunoassay based on the principle of indirect qualitative detection of IgM Antibodies to CMV in human serum or plasma. Anti-CMV IgM Antibodies, present in the samples bind to Antigens coated on the microplate forming complexes Antigen-Antibody IgM. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated human Anti-IgM Antibodies are added to the microplate and then incubated. The enzyme-conjugated Anti-Human IgM Antibodies bind to the antigen-bound IgM Antibodies present. A new wash is performed to remove the surplus. After this step, the Substrate is added and incubated, producing a blue color indicating the amount of Anti-CMV IgM Antibodies present in the samples. The Stop Solution is added to stop the reaction with a change in color from blue to yellow, measured on a microplate reader.

REAGENTS
1- Sensitized Plate - Store between 2 and 8°C.
2- Conjugate - Store between 2 and 8°C. Contains: Human Anti-IgM Antibody linked to Peroxidase and preservative.
3- Concentrated Washing Solution - Store between 2 and 8°C. Contains: Buffer Solution, surfactant and preservative.
4- Sample Diluent - Store between 2 and 8°C. Contains: Buffer Solution and preservative.
5- Substrate - Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
6- Stop Solution - Store between 2 and 8°C. Contains: Chloridric Acid 1 M.
7- Cut-Off Calibrator - Store between 2 and 8°C. Contains: Diluted IgM Antibodies Anti-CMV and preservative. **Potentially infectious.**
8- Negative Control - Store between 2 and 8°C. Contains: Non-reactive IgM Antibodies to CMV and preservative. **Potentially infectious.**
9- Positive Control - Store between 2 and 8°C. Contains: IgM Antibodies Anti-CMV and preservative. **Potentially infectious.**
10- Plate Sealers

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1- Sensitized Plate	1 Unit (96 Cavities)	2 Units (192 Cavities)	5 Units (480 Cavities)
2- Conjugate	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3- Concentrated Washing Solution	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
4- Sample Diluent	1 Flask x 42 mL	2 Flasks x 42 mL	5 Flasks x 42 mL
5- Substrate	1 Flask x 11 mL	2 Flasks x 11 mL	5 Flasks x 11 mL
6- Stop Solution	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
7- Cut-Off Calibrator	1 Flask x 100 µL	2 Flasks x 100 µL	5 Flasks x 100 µL
8- Negative Control	1 Flask x 100 µL	2 Flasks x 100 µL	5 Flasks x 100 µL
9- Positive Control	1 Flask x 100 µL	2 Flasks x 100 µL	5 Flasks x 100 µL
10- Plate sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Usage Instructions (manual)

Required materials not contained in the kit:

- Pipette capable of dispensing volumes of 5, 50 and 100 µL with lower coefficient of variation than 1,5%.
- Re-pipettor for repetitive pipetting volumes of 100 µL and 300 µL, with lower coefficient of variation than 1,5% or multicontrol pipette (Optional).
- Microplate washer (optional).
- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- Adjustable volume pipettes (200 µL to 1000 µL) for sample dilution.
- Test tubes for preparation for sample dilution.
- Paper towel to dry cavities
- Stopwatch or watch.
- Flask to store the Washing Solution after diluted.
- Distilled or deionized water.
- Tools of Quality Control.
- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. The transport can be done under ambient temperature (up to 30°C) for up to 72 (seventy two) hours. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain accurate results.
- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused microcavities in the sachet, seal and store between 2 and 8°C.
- The water used in material cleaning must to be recent and free of contaminants.
- Deionized saturated columns release alkaline water several ions and oxidizing and reducing agents that can significantly alter the results.
- Stop Solution contains Hydrochloric Acid, which is a strong acid. Handle it with care.
- All the raw material of product is tested and should be non-reactive for HBsAg, HIV 1 & 2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manual manipulation of any product containing serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.

- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- As a safety measure, you should cover the plate during the reaction.
- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin).
Hemolyzed or highly lipemic samples should not be used.
Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Washing Solution
Dilute the contents of the Flask N° 3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. Can be stored at room temperature. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

TECHNIQUE

- Before starting the assay, bring all reagents, samples, Standard References and Controls to stabilize at room temperature (15 - 30°C) for at least 40 minutes.
- Select the cavities to be used considering: Cut-Off Calibrator, Controls and Samples (it is recommended to test in duplicate). Return the strips of the plate will not be used for the original sealed packaging.
 - Select the first cavity for Blank (OPTIONAL).
 - Prepare a 1:41 dilution of the Samples, Controls and Cut-Off Calibrator in microtubes adding 5 µL more 200 µL of the Sample Diluent. Homogenize.
 - Pipette 100 µL of diluted Samples and Controls and Cut-Off Calibrator into the previously determined wells.
 - Homogenize gently for ± 30 seconds, cover the wells with plate sealer. Change occurs in color from green to blue in the cavities of the samples. Cover wells with plate sealer.
 - Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
 - Remove the plate sealers of cavities.
 - Discard the contents of the cavities by aspiration (Washer) or by decanting (manual). Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate at the end of the wash, beat the plate for a few seconds on absorbent paper.
 - Note:** Poor washing and drying can cause inadequate results.
 - Pipette 100 µL of Conjugate in all cavities except in the Blank cavity (if you made this option).
 - Mix gently for ± 30 seconds. Cover cavities with plate sealer.

- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- Remove the sealer from plate cavities.
- Repeat item 8.
- Pipette 100 µL of Substrate into all wells.
- Mix gently for ± 30 seconds. Cover wells with plate sealer.
- Incubate for 15 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- Remove the plate sealer from cavities.
- Pipette 100 µL of Stop Solution into all wells.
- Mix gently for ± 30 seconds.
- Read the 450 nm within 15 minutes (maximum).

TECHNIQUE VERIFICATION

Verify if the results obtained by the reading of the Blank, Controls and Cut-Off Calibrator are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0,150
Negative Control	< 0,5
Cut-off Calibrator	0,5 to 1,5
Positive Control	> Cut-Off Calibrator

The absorbancies for the Controls and Cut-Off Calibrator were obtained after the reduction of the Blank absorbance. Consider a Blank limit of < 0,150. If the values are out of the expected values, technique must be repeated.

CALCULATIONS

QUALITATIVE

Cut-Off regarded as the mean absorbance obtained with the Cut-Off Calibrator.
Example:

ITEM	ABSORBANCE
Cut-Off Calibrator	0,845
	0,841
Average absorbance from Cut-Off Calibrator	(0,845 + 0,841) / 2 = 0,843

Calculate the index by dividing the sample absorbance by the Cut-Off value. Example:

ITEM	ABSORBANCE
Sample	1,438
Cut-Off value	0,843
Index = Sample / Cut-Off value	1,438 / 0,843 = 1,71

Note: The data presented in the examples are for illustration only and can not be used for calculations of the results.

RESULTS INTERPRETATION

ITEM	ABSORBANCE
	INDEX
Negative	≤ 0,9
Undetermined*	Between 0,9 and 1,5
Positive	≥ 1,5

* Suspicion of Residual IgM for results between 1,0 and 1,5.

Observación: Due to the high sensitivity of the product, detection of residual IgM antibodies may occur (See Procedure Limitations). In the case of results with an index between 1,0 and 1,5, residual IgM antibodies may be suspected, and therefore the IgG avidity test is recommended.

In case of results with index between 0,9 and 1,0, it is recommended to repeat the examination by another alternative method or to wait 2 weeks for a new analysis.

PROCEDURE LIMITATIONS

IgM antibodies are present in the serum for a short period of time, disappearing up to six months after infection, while IgG antibodies remain present for a long period of time. However, the high sensitivity of the new diagnostic methodologies, in some cases, allows to find very low levels of IgM antibodies, called residual IgM, for a longer period of time. Detection of these residual IgM antibodies hinders clinical interpretation over the period of infection. In these cases, to confirm the result, an IgG avidity test is recommended.

The interpretation of a diagnostic test, there should not be based on a single run. This should include confirmation of other tests before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other clinical information available before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	1,762	1,500	1,227
Standard Deviation	0,092	0,078	0,060
Coefficient of Variation (%)	5,247	5,223	4,894

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	1,722	1,462	1,198
Standard Deviation	0,087	0,073	0,062
Coefficient of Variation (%)	5,060	5,009	5,183

Clinical Sensitivity and Specificity

BIOLISA CMV IgM analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the kit BIOLISA CMV IgM is 90,2% and clinical specificity is 93,5%

BIOLISA CMV IgM x EIA REFERENCE

METHOD		EIA REFERENCIA		Total
		Positive	Negative	
BIOLISA CMV IgM	Positive	120	11	131
	Negative	13	159	172
Total Results		133	170	303

Clinical Sensitivity: 90,2% (120/133)

Clinical Specificity: 93,5% (159/170)

Overall Agreement: 92,1% (279/303)

DIAGNOSTIC SIGNIFICANCE

Cytomegalovirus (CMV) is a member of the Herpes Virus family which includes Herpes Simplex virus (HSV) 1 and 2, Varicella Zoster virus (VZV) and Epstein-Barr virus (EBV). It is a human pathogen transmitted through saliva, sexual contact, perinatal, organ transplantation or blood transfusion. In most cases, the infection remains asymptomatic. However, CMV infection can cause severe illness in newborns and immunocompromised individuals such as AIDS patients, cancer or patients who received organ transplants. During immunosuppressive therapy, often there is a reactivation of latent virus or primary infection. CMV infections can be acquired before birth, during delivery and later in life. Can cause severe congenital anomalies such as microcephaly, motor disabilities and mental retardation. Therefore, determination of primary maternal infection and its distinction of latent infection are of great importance. The presence of IgM antibodies indicates the presence of primary infection, while IgG antibodies indicates the immune status of patients

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
2. Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after "Silent" Congenital Cytomegalovirus Infection. N. Engl. J. Med. (1976) 295:468-470.
3. Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovirus. N. Engl. J. Med. (1974) 790:291-296.
4. Stern, H. Cytomegalovirus Vaccine: Justification and Problems. In: Waterson AP (ed.) Recent Advances in Clinical Virology (1977) 117-134.
5. Bioclin – Dados de arquivos

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA CMV IgM kit: 10269360189

Review: November/2017

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED