Bioclin

BIOLISA COVID-19 IgG

REF **K229-1**

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FINALIDADE

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Departamento de Pesquisa e Desenvolvimento.

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR 2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Departamento de Pesquisa e Desenvolvimento.

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USE INSTRUCTIONS

FUNCTION

Test for qualitative determination of IgG antibodies to SARS-CoV-2 (virus that causes COVID-19 disease) in biological samples (serum, plasma or whole blood on filter paper) through enzyme immunoassay test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA COVID-19 IgG kit is an immunoenzymatic assay in phase based on the principle of indirect qualitative detection of antibodies IgG for SARS-CoV-2 in serum, plasma and whole blood on filter paper human. IgG antibodies present in the sample bind to antigens recombinants of SARS-COV-2 coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Anti-IgG antibodies conjugated to Peroxidase is added to the microplate which is then incubated. The enzyme-conjugated antibodies bind to IgG Anti-SARSCoV- 2 present, connected to the antigen-coated plate. New wash is performed to remove surplus. After this stage, the Substrate is added and incubated, producing a blue color that indicates the amount of IgG Anti-SARS-CoV-2 Antibodies present in the sample. The Solution Stop is added to stop the reaction if there is a change in color from blue to yellow, measured in a microplate reader.

REAGENTS

- 1 - Sensitized Plate** - Store between 2 to 8°C. Contains: Sensitized Plate with recombinant SARS-COV-2 antigen and preservative.
- 2 - Conjugate** - Store between 2 to 8°C. Contains: Buffer Solution containing Peroxidase-linked human anti-IgG antibody, surfactant, stabilizers, dye and preservative.
- 3 - Concentrated Washing** - Store between 2 to 8°C. Contains: Buffer Solution, surfactant and preservative.
- 4 - Sample Diluent** - Store between 2 to 8°C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
- 5 - Substrate** - Store between 2 to 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6 - Stop Solution** - Store at 2 to 8°C. Contains: Hydrochloric Acid 1 M.
- 7 - Negative Control** - Store between 2 to 8°C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infective.**
- 8 - Positive Control** - Store between 2 to 8°C. Contains: Buffer Solution, IgG Anti - SARS-CoV-2 antibodies, dye, stabilizers, surfactant and preservative. **Potentially infective.**
- 9 - Plate Sealers**
- 10- Filter Paper Elution Buffer** - Store between 2 to 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 11- Negative Extraction Control** - Store between 2 to 8°C. Contains: Non-reactive sample for IgG anti-SARS-CoV-2 antibodies impregnated in Filter Paper. **Potentially infective.**
- 12- Positive Extraction Control** - Store between 2 to 8°C. Contains: Reactive sample for IgG anti-SARS-CoV-2 antibodies impregnated on Filter Paper. **Potentially infective.**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 – Sensitized Plate	1 Unit (96 cavities)	2 Units (192 cavities)	5 Units (480 cavities)
2 – Conjugate	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3 – Concentrate Washing	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
4 – Sample Diluent	1 Flask x 100 mL	2 Flasks x 100 mL	5 Flasks x 100 mL
5 - Substrate	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
9 – Stop Solution	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
7 –Negative Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
8 –Positive Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
9 – Plate Sealers	3 Units	6 Units	15 Units
10 – Filter Paper Elution Buffer	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
11 – Negative Extraction Control	1 unit	2 units	5 units
12 – Positive Extraction Control	1 unit	2 units	5 units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Operating instructions (manual)

Required materials not contained in the kit:

- 1-** Pipette capable of dispensing volumes of 5 to 500 µL with lower coefficient of variation than 1.5%.
- 2-** Re-pipettor for repetitive pipetting volumes of 500 µL, with lower coefficient of variation than 1.5% or multichannel pipette (optional).
- 3-** Microplate washer (optional).
- 4-** ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5-** Paper towel to dry cavities
- 6-** Stopwatch or watch.
- 7-** Flask to store the washing solution after dilution.
- 8-** Distilled or deionized water.
- 9-** Tools of Quality Control.
- 10-** Incubator 37°C ± 2°C.
- 11-** Paper shredder (3 mm diameter) for filter paper technique.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- 2-** Strictly follow the methodology proposed to obtain accurate results.
- 3-** The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8°C.
- 4-** The water used in material cleaning must to be recent and free of contaminants.
- 5-** Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6-** Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.

7- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.

- 8-** Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 9-** As a safety measure, you should cover the plate during the reaction.
- 10-** You must ensure that the bottom of the cavity is clean and dry, and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- 11-** Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14-** Do not use the product in case of damaged packaging.
- 15-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin).
Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

Whole Blood on Filter Paper

Whole Blood (Puncture or EDTA).
Samples dried on fi lter paper can be stored at room temperature as long as they are protected from direct sunlight and low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8°C. To storage for more than 2 years must be carried out the temperature of -20°C.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENTS

Washing Solution
Dilute the contents of the Reagent N° 3 (Concentrate Washing) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

STABILITY AFTER OPEN

The results of the stability test show that the Biolisa kit COVID-19 IgG is stable after opening for at least 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

TECHNIQUE

Samples of Serum or Plasma

- Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.
- 1-** Select the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used, for the original sealed packaging.
 - 2-** Select the first cavity for Blank (OPTIONAL).

- 3-** Prepare a 1:101 dilution of Samples and Controls in microtubes by adding 5 µL plus 500 µL of sample diluent. Homogenize.
- Observe the color change of the diluent when adding the sample. A Color change indicates that the sample was successfully added to the microtube.
- 4-** Pipette 100 µL of Diluted Sample and Controls into the previously determined wells. **In the Blank cavity, pipette 100µL of Sample Diluent**
- 5-** Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 6-** Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 7-** Remove the sealing of wells.
- 8-** Discard the contents of the wells by aspiration (Washer) or by decanting (Manual).
- Use approximately 300 µL of Washing Solution, previously prepared to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
- Note:** Poor washing/drying can cause inadequate results.
- 9-** Pipette 100 µL of Conjugate into all wells. **Except in the Blank cavity.**
- 10-** Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 11-** Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 12-** Remove the sealer from plate cavities.
- 13-** Repeat item 8.
- 14-** Pipette 100 µL of Substrate into all wells. **Included in the cavity for Blank.**
- 15-** Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 16-** Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 17-** Remove the sealer from plate cavities.
- 18-** Pipette 100 µL of Stop Solution into all wells. **Included in the cavity for Blank.**
- 19-** Mix gently for ± 30 seconds.
- 20-** Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

Whole Blood Sample on Filter Paper

- Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.
- 1-** Select the wells to be used considering Positive Extraction Control, Negative Extraction Control and Whole Blood Samples on Filter Paper (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.
 - 2-** Select the first cavity for Blank (OPTIONAL).
 - 3-** Add a 3 mm disk of Positive Extraction Control, Negative Extraction Control and Whole Blood Sample on Filter Paper, previously perforated, in the determined wells.
 - 4-** Pipette 200 µL of Filter Paper Elution Buffer into each well. **Included in the Blank cavity.**
 - 5-** Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
 - 6-** Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
 - 7-** Remove the sealing of wells.
 - 8-** Discard the contents of the cavities by aspiration (Washing machine) or by decanting (Manual). Remove the paper discs, if necessary, with the aid of a needle or tip.
 - Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
 - Note:** Poor washing/drying can cause inadequate results.
 - 9-** Pipette 100 µL of Conjugate into each well. **Except in the Blank cavity**
 - 10-** Mix gently for ± 30 seconds. Cover cavities with plate sealer.
 - 11-** Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
 - 12-** Remove the sealer from plate cavities.
 - 13-** Repeat item 8.
 - 14-** Pipette 100 µL of Substrate into each well. **Included in the Blank cavity.**

- 15- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 17- Remove the sealer from plate cavities.
- 18- Pipette 100 µL of Stop Solution into all wells. **Included in the Blank cavity**
- 19- Mix gently for ± 30 seconds.
- 20- Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Samples of Serum, Plasma and Filter Paper Whole Blood

Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0.200
Negative Control	< 0.200
Positive Control	> 1.200
Negative Extraction Control	< 0.300
Positive Extraction Control	> 1.200

If the values are out the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Samples of Serum or Plasma

Calculate Cut Off according to the following formula:

Cut Off = (Average Absorbance Positive Control x 0.1) + 0.100

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N° 8)	A1 = 1.898
	A2 = 1.894
Cut Off = (Average Absorbance Positive Control x 0.1) + 0.100	((1.898 + 1.894) / 2 x 0.1) + 0.100 = 0.289

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1.963
Cut Off Value	0.289
Index = Sample / Cut Off Value	1.963 / 0.289 = 6.79

Whole Blood Sample on Filter Paper

Calculate Cut Off according to the following formula:

Cut Off = (Average Absorbance of Positive Extraction Control x 0.1) + 0.100

Example:

ITEM	ABSORBANCE
Positive Extraction Control (Reagent N° 12)	A1 = 1.898
	A2 = 1.894
Cut Off = (Average Absorbance of Positive Extraction Control x 0.1) + 0.100	((1.898 + 1.894) / 2 x 0.1) + 0.100 = 0.289

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	0.950
Cut Off Value	0.289
Index = Sample / Cut Off Value	0.950 / 0.289 = 3.28

INTERPRETATION OF RESULTS

Samples of Serum, Plasma and Filter Paper Whole Blood

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE FOR SERUM SAMPLE, PLASMA AND TOTAL BLOOD ON FILTER PAPER
	INDEX
Negative	< 0.90
Positive	> 1.1
Undetermined	0.9 – 1.1

Note: In the case of an undetermined result, the sample must be reanalyzed. Samples that obtain results repeatedly indeterminate conditions must be retested using an alternative method. If results remain indeterminate, a new sample in two weeks. The result of the last sample must prevail collected. For a more complete assessment of the diagnosis and clinical correlation, it is suggested that each sample be tested for IgG and IgM. The results provided by this kit must be interpreted responsible medical professional, and it is not the only criterion for the determining the diagnosis and /or treatment of the patient.

Note: The data presented in the examples are for illustration only and cannot be used to calculate results.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be established with basis of a single trial. Other confirmatory tests should be included, before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before the definitive diagnosis of the disease.

INTERFERENT

No interference was observed for Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rematoid Factor 980 IU/mL, C Reactive Protein 41.2 mg/dL and Anti Streptolysin O 1023 IU/ml.

CROSS REACTIVITY

A cross-reactivity study was carried out, evaluating 130 samples negative for SARS-CoV-2, but positive for other infections. Among they 10 samples positive for Influenza, 27 positive for Rhinovirus, 13 positive samples for Respiratory Syncytial Virus, 10 positive samples for Zika, 10 positive samples for Dengue, 10 positive samples for Chikungunya, 10 samples positive for Toxoplasmosis, 10 samples positive for Rubella, 10 HIV positive samples, 10 samples positive for HBV and 10 positive samples for HCV. It was not observed cross-reactivity with positive samples for Influenza, Rhinovirus, Respiratory syncytial virus, Zika, Dengue, Chikungunya, Toxoplasmosis, Rubella, HIV, HBV and HCV. Despite the results found, it is not can completely rule out the possibility of cross-reactivity. O final diagnosis should consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE

QUALITY CONTROL

Accuracy

REPETIBILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPETIBILITY	SAMPLE		
	1	2	3
Average	2.501	1.283	0.137
Standard Deviation	0.075	0.078	0.006
Coefficient of Variation (%)	2.998	6.072	4.721

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	2.487	1.282	0.138
Standard Deviation	0.072	0.071	0.009
Coefficient of Variation (%)	2.899	5.573	6.167

Clinical Sensitivity and Specificity

Serum and Plasma Sample

BIOLISA COVID-19 IgG Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA COVID-19 IgG kit is 95.2% and clinical specificity is 95.7%.

	Expected Result	BIOLISA COVID-19 IgG
Positive Sample	62	59
Negative Sample	70	67
Total Tested Sample	132	126

Clinical Sensitivity: 95.2% (59/62)

Clinical Specificity: 95.7% (67/70)

Whole Blood Samples on Filter Paper

BIOLISA COVID-19 IgG Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA COVID-19 IgG kit is 96.0% and clinical specificity is 95.3%.

	Expected Result	BIOLISA COVID-19 IgG
Positive Sample	50	48
Negative Sample	64	61
Total Tested Sample	114	109

Clinical Sensitivity: 96.0% (48/50)

Clinical Specificity: 95.3% (61/64)

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family, and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are: fever, cough, dyspnoea, myalgia or fatigue. It can cause diseases ranging from the cold common to more serious diseases, such as the Eastern Respiratory Syndrome (MERS-CoV) and Severe Acute Respiratory Syndrome SARS-CoV. In December 2019, a series of cases of pneumonia caused disease appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples was found to be Coronavirus infection, which it was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents have already been reported cases of infection.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
4. Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
5. Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
6. TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
7. QUIBASA: Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA COVID-19 IgG kit: 10269360325

Review: June/2020