Bioclin

BIOLISA HBsAg

REF **K120**

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NÚMERO DE TESTES

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- Blumberg, B.S. The Discovery of Australian Antigen and its Relation to Viral Hepatitis. Vitro. 1971;7:223.
- Krugman, S. Glies J.P. Viral Hepatitis, Type B (MS-2-Strain). Further Observations on Natural History and Prevention. New England Journal of Medicine. 288, 755.
- Krugman, S. Overby L.R, et al. Viral Hepatitis Type B Studies On Natural History and Prevention Re- examined. New England Journal of Medicine. 300, 101.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

BIOLISA HBsAg kit is a third generation test for the qualitative detection of the presence of antigen surface of Hepatitis B virus in human serum or plasma (EDTA or Heparin). The test uses monoclonal antibodies to detect various subtypes of HBsAg in serum or plasma. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
BIOLISA HBsAg kit assay is a solid phase immunoassay based on the principle of "sandwich" for the detection of HBsAg in human serum or plasma. The microplate is coated with monoclonal antibodies specific for several subtypes of HBsAg. During the test, sample and conjugated are added to the microplate coated with antibody and then incubated. If the sample contains HBsAg, this will bind to antibodies that coat the microplate and simultaneously be linked to the conjugate to form complex-conjugate-HBsAg antibody. After incubation the microplate is washed to remove the unbound materials. After this step, the Substrates A and B are added and incubated, producing a blue color indicates that the amount of antigen of Hepatitis B in the samples. The Stop Solution is added to stop the reaction, with a color change to yellow, measured in a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8°C. Microplate coated with Anti-Hbs antibodies.
- 2- Conjugated** - Store between 2 and 8°C. Anti-HBs antibodies bound to Peroxidase, stabilizer and preservative.
- 3- Concentrated Washing** - Store between 2 and 8°C. Buffer Solution, surfactant and preservative.
- 4- Substrate A** - Store between 2 and 8°C. Solution containing Tetramethylbenzidine (TMB).
- 5- Substrate B** - Store between 2 and 8°C. Solution containing Urea Peroxide.
- 6- Stop Solution** - Store between 2 and 8°C. Sulfuric Acid 1N.
- 7- Negative Control** - Store between 2 and 8°C. Solution non-reactive for HBsAg, HCV, HIV-1 and HIV-2, stabilizer and preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8°C. Solution containing HbsAg and negative for HCV, HIV-1 and HIV-2, stabilizer and preservative. **Potentially infectious.**
- 9- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Sensitized Plate	1 Unit (96 Cavities)	2 Units (96 Cavities)	5 Units (96 Cavities)
2- Conjugate	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
3- Concentrated Washing	1 Flask x 50 mL	2 Flasks x 50 mL	4 Flasks x 50 mL
4- Substrate A	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
5- Substrate B	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
6- Stop Solution	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
7- Negative Control	1 Flask x 1 mL	2 Flasks x 1 mL	5 Flasks x 1 mL
8- Positive Control	1 Flask x 1 mL	2 Flasks x 1 mL	5 Flasks x 1 mL
9- Plate Sealers	3 Units	5 Units	10 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Usage Instructions (manual)

Required materials not contained in the kit:

- 1- Pipette capable of dispensing volumes of 50 to 500 µL with a lower coefficient of variation than 1,5%.
- 2- Re-pipettor for repetitive pipetting volumes of 100 µL and 300 µL, with lower coefficient of variation than 1,5% or multicontrol pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Adjustable volume pipettes (200 µL to 1000 µL) for preparation of the Substrate.
- 6- Test tubes for preparation of the Substrates A and B.
- 7- Paper towel to dry cavities
- 8- Stopwatch or watch.
- 9- Flask to store the Washing Solution, after diluted.
- 10- Distilled or deionized water.
- 11- Tools of Quality Control.
- 12- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. The transport can be done in under ambient temperature (up to 30°C) for up to 72 (seventy two) hours. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only

- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The envelope containing the strips should be opened only after it reaches room temperature. Place the strip with unused microcavities in the aluminum bag, seal and store at 2 - 8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized saturated columns release alkaline water, several ions and oxidizing and reducing agents, that can significantly alter the results.
- 6- All the raw material of product is tested and should be non-reactive for Anti HIV 1 & 2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manual manipulation of any product containing serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.
- 7- Always add reagents in the same order to minimize the difference in reaction time between the microcavities.
- 8- As a safety measure, you should cover the plate during the reaction.
- 9- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- 10- Do not expose reagents, especially the substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 11- Stop Solution contains Sulfuric Acid, which is a strong acid. Handle it with care.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin).
Hemolyzed or highly lipemic samples should not be used.
Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of the Flask Nº 3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. Can be stored at room temperature. **In case of crystallization, heat it at 37°C until dissolution.**

Substrate - Working Solution

Determine the amount of cavities to be used for preparation of an appropriate volume. Prepare the solution by mixing equal parts of Substrate A and Substrate B, 15 minutes before use. Keep it protected from light until used.

To each microcavity (test), use:

50 µL of substrate A + 50 µL of substrate B

For example, mix 1 mL of Substrate A and 1 mL of Substrate B to two strips with 8 cavities (16 tests). Leftover reagent occurs.

Use no later than one(1) hour after preparation.

STABILITY AFTER OPEN

The results of the stability test show that the Biolisa kit HBsAg is stable after opening for at least 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

TECHNIQUE

Before starting the assay, bring all Reagents, Samples, and Controls to stabilize at room temperature (15 - 30°C) for at least 40 minutes. Return the strips of the microplate that will not be used to the original sealed packaging.

- 1- Select the microcavities to be used considering: Controls and Samples (They can be tested in duplicate).
- 2- Select the first cavity for Blank (OPTIONAL).
- 3- Pipette 50 µL Negative Control, Positive Control and Samples in the cavities previously determined.
- 4- Pipette 50 µL of Conjugate into all cavities included in the Blank cavity.
- 5- Homogenize it gently for ± 30 seconds. Cover cavities with plate sealer.
- 6- Incubate for 80 minutes ± 2 minutes at 37°C ± 2°C.
- 7- Remove the plate sealers from cavities
- 8- Discard the contents of the cavities by aspiration (Washer) or by decanting (manual).
- Use approximately 350 µL of Washing Solution, **previously prepared***, to perform a total of eight (8) wash cycles.

Or

Use approximately 500 µL of Wash Solution, **previously prepared***, for a total of six (6) wash cycles (exclusive to automatic equipment).

Note: Evaluate the aspiration and dispensing capacity of the equipment.
To ensure the drying of the plate, at the end of washing, beat it for a few seconds on absorbent paper.

9- ATTENTION Follow one of the following procedures:

A) Pipette 100 µL of **previously prepared*** Substrate - Working Solution (A + B) in all cavities.

*** See PREPARATION OF REAGENTS**

Or

B) Pipette 50 µL of substrate A in all cavities

Pipette 50 µL of substrate B in all cavities.

10- Homogenize it gently for ± 30 seconds. Cover the cavities with plate sealer.

11- Incubate for 30 minutes ± 1 minute in room temperature (15 - 30°C) and protected from light.

12- Remove the plate sealer from cavities.

13- Pipette 100 µL of Stop Solution in all cavities.

14- Homogenize gently for ± 30 seconds.

15- Read at 450 nm (primary filter) / 630 nm (secondary filter) 30 minutes maximum.

TECHNIQUE VERIFICATION

Make sure to read the results of Reference Standards are consistent with the figures below:

ITEM	VALIDATION REQUISITES
Blank	< 0,050
Negative Control	< 0,100
Positive Control	> 1,000

The absorbances for the above mentioned Controls were obtained after the reduction of the Blank absorbance. For single filter reading (450 nm) to consider the Blank limit < 0,100. If the values are out of the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Calculate the average absorbance from Negative Control.

Example:

ITEM	ABSORBANCE
Negative Control	A1 = 0,023
	A2 = 0,021
Average absorbance from Negative Control	(0,023 + 0,021) / 2 = 0,022

If the results from Controls are valid, calculate the Cut-Off using the following formula:
Example:

ITEM	ABSORBANCE
Cut - Off = Average absorbance from Negative Control + 0,025	0,022 + 0,025 = 0,047

Calculate the index by dividing the sample absorbance by the value of the Cut-Off.
Example:

ITEM	ABSORBANCE
Sample	0,900
Cut-Off Value	0,047
Index: Sample / Cut-Off Value	0,900 / 0,047 = 19,15

RESULTS

RESULTS	QUALITATIVE
	INDEX
Negative	< 0,9
Positive	> 1,1
Undetermined	0,9 - 1,1

NON-REACTIVE

Sample with absorbance lower than the Cut - Off is considered non-reactive to Hepatitis B antigen and can be considerate negative.

REACTIVE

Sample with absorbance higher than the Cut - Off is initially considered reactive to Hepatitis B antigen. The sample must be retested in duplicate before the final interpretation. The sample who is considered reactive in at least one of the retests is supposed to be reactive and must be confirmed through the confirmatory method.

Note: In case of undetermined results, the sample must be retested. Samples that results have been obtained repeatedly undetermined should be retested using an alternative method. If the results remain uncertain, one must collect a new sample in two weeks. The result of the last sample collected should prevail.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

Note: The data presented in the examples are for ilustration only and can be used for calculation of the results.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test, there should not be based on a single run. This should include confirmation of other tests before a sample is considered positive. A negative result does not exclude the possibility of exposure. Finally, all results should be interpreted in conjunction with other clinical information available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PERFORMANCE CHARACTERISTICS
QUALITY CONTROL

Accuracy

ANALYTICAL SENSITIVITY

The analytical sensitivity of the BIOLISA HBsAg kit was determined using the international standard WHO NISBC 80/549 for HBsAg. The analytical sensitivity was 0,05 UI/mL and it is equivalent to 0,2 ng/mL.

CLINICAL SENSITIVITY AND SPECIFICITY

BIOLISA HBsAg kit analyzed clinical samples (seroconversion) compared with other methods EIA. The results show that the clinical sensitivity of the kit BIOLISA HBsAg is 100% and specificty clinic is 99,58%.

	TOTAL EXPECTED	BIOLISA HBsAg
Negative Sample	5714	5690
Positive Sample	400	400

Specificity: 99,58%
Sensitivity:100%

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	0,7541	2,2766	1,3159
Standard Deviation	0,0857	0,1571	0,1168
Coefficient of Variation (%)	11,36	6,90	8,88

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	0,8960	2,5393	1,5161
Standard Deviation	0,1229	0,2358	0,1767
Coefficient of Variation (%)	13,72	9,28	11,65

DIAGNOSTIC SIGNIFICANCE

HBsAg is an early marker that appears in the blood after infection with the Hepatitis B (HBV). This liver infection is transmitted through sexual contact, blood exposure, the transmission from mother to child during birth or sharing objects cutting perforation. The main HBsAg subtypes ad and includes ad and ay, all sharing the common determinant 'a'. HBV infection causes a variety of liver damage, such as acute self - limiting, fulminate hepatitis, chronic hepatitis with progression to cirrhosis and liver failure, and chronic carrier state asymptomatic. The HBV virus in infected people, persists for the rest of their lives and can be transmitted to others. Thus, Hepatitis B has become a public health problem. Infection HBV results in an apparent number of serological markers and one of the first of these markers is the Surface Antigen of Hepatitis B (HBsAg). HBsAg appears 1-10 weeks after exposure and before biochemical evidence of liver disease or jaundice. Three weeks after the onset of acute hepatitis almost half of patients will still be positive for HBsAg. At the chronic carrier state, the HBsAg persists for 6-12 months without serum conversion to the corresponding antibodies. Therefore, screening for HBsAg is highly recommended for all donors, pregnant women and people in high risk groups.

NUMBER OF TESTS

Presentation 1 - 96 tests
Presentation 2 - 192 tests
Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Blumberg, B.S. The Discovery of Australian Antigen and its Relation to Viral Hepatitis. Vitro. 1971;7:223.
2. Krugman, S. Glies J.P. Viral Hepatitis, Type B (MS-2-Strain). Further Observations on Natural History and Prevention. New England Journal of Medicine. 288, 755.
3. Krugman, S. Overby L.R, et al. Viral Hepatitis Type B Studies On Natural History and Prevention Reexamined. New England Journal of Medicine. 300, 101.
4. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA HBsAg kit: 10269360288

Review: September/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED