BIOLISA HCV

REF K128

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FINALIDADE

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NÚMERO DE TESTES

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- Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
- Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
- Van der Poel, C. L., H.T.M. Cuyper, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Four-antigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
- Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
- Bioclin – Dados de arquivos

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1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Four-antigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
5. Bioclin – Datos de archivos

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USAGE INSTRUCTIONS

FUNCTION

Test for the qualitative determination of total antibodies to Hepatitis C Virus (HCV) in serum or human plasma, for immunoassay, in microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
BIOLISA HCV kit is an direct enzyme immunoassay fourth generation in solid-phase based on the qualitative detection of antibodies to Hepatitis C virus in human serum or plasma.
The microplate is coated with HCV recombinant antigens of the Core region, NS5 and NS3.
During testing, samples are added and then incubated in the microplate.
If the samples contains anti-HCV antibodies, these will bind to the antigen coating the microplate to form immobilized antibody-antigen complex.
After the initial incubation, the microplate is washed to remove unbound materials.
The conjugate prepared previously is added to the wells resulting in the formation of the HCV Antigen - anti-HCV Antibody – Conjugate complex.
Then the plate is again incubated.
After the second incubation, the microplate is washed to remove the unbound materials.
The Substrates A and B are then added and the microplate is again incubated to produce a blue color, indicating the amount of HCV antibodies present in the sample.
The Stop Solution is added to microplate to stop the reaction producing a color change from blue to yellow. The color intensity, which is the amount of HCV antibodies in the sample is measured with a microplates reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8°C. Microplates coated with recombinant HCV antigens.
- 2- Conjugate** - Store between 2 and 8°C. Solution HCV antigen linked to Peroxidase in buffer, preservative and stabilizer.
- 3- Concentrated Washing** - Store between 2 and 8°C. Buffer solution, surfactant and preservative.
- 4- Conjugate Diluent** - Store between 2 and 8°C. Buffer solution, surfactant, stabilizer and preservative.
- 5- Substrate A** - Store between 2 and 8°C. Solution containing Tetramethylbenzidine (TMB).
- 6- Substrate B** - Store between 2 and 8°C. Acid solution containing Urea Peroxide.
- 7- Stop Solution** - Store between 2 and 8°C. Sulfuric Acid solution 1N.
- 8- Negative Control** - Store between 2 and 8°C. Plasma non-reactive for anti-HCV antibody and preservative. **Potentially infectious.**
- 9- Positive Control** - Store between 2 and 8°C. Inactivated plasma containing anti-HCV antibodies and preservative. **Potentially infectious.**
- 10- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1- Sensitized Plate	1 Unit (96 cavities)	2 Units (96 cavities)	5 Units (96 cavities)
2- Conjugate	1 Flask x 1,8 mL	2 Flasks x 1,8 mL	5 Flasks x 1,8 mL
3- Concentrated Washing	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
4- Conjugate Diluent	1 Flask x 24 mL	2 Flasks x 24 mL	5 Flasks x 24 mL
5- Substrato A	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
6- Substrato B	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
7- Stop Solution	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
8- Negative Control	1 Flask x 1,5 mL	2 Flasks x 1,5 mL	5 Flasks x 1,5 mL
9- Positive Control	1 Flask x 1,5 mL	2 Flasks x 1,5 mL	5 Flasks x 1,5 mL
10- Plate Sealers	3 Units	5 Units	10 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table.
- Operating instructions (manual).

Required materials not contained in the kit:

- 1- Pipette (s) capable of dispensing 10, 50 and 100 µL volumes with a lower coefficient of variation than 1,5%.
- 2- Re-pipettor (s) for repetitive dispensing of volumes of 100 µL and 300 µL, with a lower coefficient of variation than 1,5% or multichannel pipette (optional).
- 3- Micro plate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Adjustable volume pipettes (200 µL to 1000 µL) for dilution of the substrate.
- 6- Test tubes for dilution of the substrate A and B.
- 7- Paper towel to dry cavities.
- 8- Plastic wrap or micro plate cover for incubation steps.
- 9- Stopwatches or watch.
- 10- Container to store washing solution.
- 11- Distilled or deionized water.
- 12- Tools of Quality Control.
- 13- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. The transport can be in under ambient temperature (up to 30°C) for up to 72 (seventy two) hours. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The envelope containing the strips should be opened only after it reaches room temperature. Place the strip with unused cavities in the aluminum bag, seal and store between 2 and 8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.

- 6- All the raw material of product is tested and should be nonreactive for HBsAg and HIV 1 & 2. However, these tests do not provide total assurance of the absence of infectious agents. The manual manipulation of any product containing serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.
- 7- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 8- As a safety measure, you should cover the plate during the reaction.
- 9- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- 10- Do not expose reagents, especially the substrate, to strong light or hypochlorite fumes during storage or incubation steps.
- 11- Stop Solution contains Sulfuric Acid, which is a strong acid. Handle it with care.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin). Hemolyzed or highly lipemic samples should not be used. Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).

DESCRIPTION OF PROCESS

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of the Flask N° 3 (Concentrated Washing Solution) in 1000 mL of distilled water or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. Can be stored at room temperature. In case of crystallization, heat it at 37°C until dissolution.

Substrate - Working Solution

Determine the amount of cavities to be used for preparation of an appropriate volume. Prepare the solution by mixing equal parts of Substrate A and Substrate B, 15 minutes before use. Keep it protected from light until used. To each cavity (test), use:
50 µL of substrate A + 50 µL of substrate B
For example, mix 1 mL of Substrate A and 1 mL of Substrate B to two strips with 8 cavities (16 tests). Leftover reagent occurs.
Use no later than one (1) hour after preparation.

Conjugate

Determine the amount of cavities that will be used for the preparation of an appropriate amount of the conjugate.
The preparation of the diluted conjugate should occur 15 minutes before use. Use only clean container to prevent contamination.
To prepare the conjugate by diluting Reagent N°2 (Conjugate) in the ratio 1:21 using Reagent N° 4 (Conjugate Diluent), as shown below.
Homogenize the solution, avoiding the formation of foam.
After use, discard the excess diluted Conjugate solution.

Below the dilution table:

Cavities Number	Volume of R4 Required (Conjugate Diluent)	Volume of R2 (Conjugate)
8	1 mL	50 µL
16	2 mL	100 µL
24	3 mL	150 µL
32	4 mL	200 µL
40	5 mL	250 µL
48	6 mL	300 µL
56	7 mL	350 µL
64	8 mL	400 µL
72 - 80	9 mL	450 µL
81 - 96	10 mL	500 µL

TECHNIQUE

Before starting the assay, bring all Reagents, Samples and Controls to stabilize temperature (15 - 30°C) for at least 40 minutes.
Return unused strips to the original sealed packaging.
1- Select the cavities to be used considering: Controls, Samples and Blank (if there is). They can be tested in duplicates.
2- Select the first cavity for Blank (OPTIONAL).
3- Pipette 100 µL of Negative Controls, Positive Controls and Samples in the cavities previously determined.
4- Homogenize gently for ± 30 seconds.
Cover the cavities with sealer.
5- Incubate for 60 minutes ± 2 minutes in an incubator at 37°C ± 1°C.
6- Remove the sealing of the cavities.
7- Discard the contents of the cavities by aspiration (Washer) or by decanting (manual);
Use approximately 350 µL of Washing Solution, **previously prepared***, in advance to make a total of eight (8) washing cycles.
Or
Use approximately 500 µL of Wash Solution, **previously prepared***, for a total of six (6) wash cycles (exclusive to automatic equipment).
Note: Evaluate the aspiration and dispensing capacity of the equipment.
To ensure the drying of the plate at the end of washing, beat it for a few seconds on absorbent paper.
8- Pipette 100 µL of Conjugate, **previously prepared***, into each cavity except in the Blank (If you have made this choice).
9- Homogenize gently for ± 30 seconds. Cover the cavities with the plate sealers.
10- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 1°C.
11- Remove the plate sealer.
12- Repeat item 7.
13- **ATTENTION: follow one of the following procedures:**
A) Pipette 100 µL of the **previously prepared*** Substrate - Working Solution (A + B) in all cavities.
*** See PREPARATION OF WORKING REAGENTS.**
Or
B) Pipette 50 µL of Substrate A and 50 µL Substrate B in all cavities.
14- Homogenize gently for ± 30 seconds. Cover the cavities with sealer.
15- Incubate for 15 minutes ± 2 minute in an incubator at 37°C ± 1°C.
16- Remove the plate sealer.
17- Pipette 100 µL of Stop Solution to each cavity.
18- Homogenize gently for ± 30 seconds.
19- Perform reading at 450nm (primary filter) / 630 nm (secondary filter) up to 15 minutes maximum.

TECHNIQUE VERIFICATION

Make sure to read the results are consistent with the figures below:

ITEM	ABSORBANCE
Blank Cavity	< 0,100
Negative Control	≤ 0,200
Positive Control	≥ 0,600

The absorbance for each control were obtained after the Blank absorbance. For a Single Filter reading (450 nm) to consider the Blank limit of < 0,100. If the values are out of the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

For Cut-Off calculation, calculate the Negative Control average absorbance.

ITEM	ABSORBANCE
Negative Control	A1 = 0,050
	A2 = 0,060
Average Absorbance from Negative Control	(0,050 + 0,060) / 2 = 0,055

If the results obtained from the Control are valid, calculate the Cut-Off according to the formula:
Example:

ITEM	ABSORBANCE
Cut-Off = Average Absorbance of the Negative Control + 0,100	0,055 + 0,100 = 0,155

Calculate the Index dividing the sample absorbance by the Cut-Off value.
Example:

ITEM	ABSORBANCE
Sample	1,100
Cut-Off Value	0,155
Index: Sample / Cut-Off Value	1,100 / 0,155 = 7,096

INTERPRETATION OF RESULTS

RESULTS	QUALITATIVE
	INDEX
Negative	< 0,9
Positive	> 1,1
Undetermined	0,9 - 1,1

Non-Reactive: Sample with absorbance lower or equal to the Cut-Off is considered non-reactive to Hepatitis C antibodies and can be considerate negative.

Reactive: Sample with absorbance higher than the Cut-Off is considered reactive to Hepatitis C antibodies and can be considerate positive, and must be confirmed using confirmation tests.

Note: In case of indeterminate results, the sample must be retested. Samples that obtain repeatedly indeterminate result should be retested using an alternative method. If the results remain uncertain, one must collect a new sample in two weeks. The result of the last sample collected should prevail. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

Note: The date presented in the examples are for ilustration only and can be used for calculation of the results.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test, there should not be based on a single run. This should include confirmation of other tests before a sample is considered positive. A negative result does not exclude the possibility of exposure. Finally, all results should be interpreted in conjunction with other clinical information available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 different samples, obtaining the following results:

REPEATABILITY	SAMPLE 1	SAMPLE 2	SAMPLE 3
Average	4,49	9,01	21,79
Standard Deviation	0,32	0,54	1,07
Coefficient of Variation (%)	7,22	5,94	4,91

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 different samples, obtaining the following results:

REPRODUCIBILITY	SAMPLE 1	SAMPLE 2	SAMPLE 3
Average	4,57	9,40	21,63
Standard Deviation	0,46	1,10	0,20
Coefficient of Variation (%)	10,07	11,73	0,94

Sensitivity and Clinical Specificity

BIOLISA HCV kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA HCV kit is > 99,9% and clinical specificity is 99,74%.

BIOLISA HCV X REFERENCE EIA

	EXPECTED TOTAL	BIOLISA HCV
Sample Positive	421	421
Sample Negative	5169	5156

Sensitivity: > 99,9% (421/421)

Specificity: 99,74% (5156/5169)

DIAGNOSTIC SIGNIFICANCE

The Hepatitis C virus has a small envelope and is the main cause of parenteral transmission of hepatitis non-A, non-B. HCV infection causes a wide variety of chronic liver disease, cirrhosis and liver cancer. The major route of transmission is through blood transfusion and hemoderivatives, organ transplantation, and sharing of contaminated needles and syringes. Antibodies against HCV are found in more than 80% of patients with hepatitis non-A, non-B. Cloning viral genome has made possible the development of serological tests using recombinant antigens. The BIOLISA HCV kit is a third generation immunoassay for the qualitative detection of the presence of anti-HCV IgG in serum or plasma. The test uses recombinant antigens of HCV encoded by the genes for both structural (nucleocapsid) and non-structural proteins to selectively detect HCV antibodies in serum or plasma.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Four-antigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
5. Bioclin – Dados de arquivos.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA HCV kit:10269360304

Review: August/2019

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED