BIOLISA HIV 1/2/O

REF
K118

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Kit BIOLISA HIV 1/2/O is a third generation immunoassay for the qualitative detection of presence of total antibodies (IgG, IgM and IgA) of HIV-1, HIV-2, and/or subtype O in serum or human plasma. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymetric
Kit BIOLISA HIV 1/2/O is a solid phase enzyme immunoassay based on the principle of "sandwich" for the detection of total antibodies (IgG, IgM and IgA) for HIV-1, HIV-2, and/or subtype O in human serum or plasma.
The microplate is coated with HIV recombinant antigen specific for HIV-1 (p24, gp41), HIV-2 (gp36) and subtype O. Antibodies against HIV-1, HIV-2, and/or subtype O in the samples bind to HIV antigens recombinant coated on the microplate into HIV antigens - antibodies complex. After incubation initial, the microplate is washed to remove unbound materials. Enzymes conjugated with antigens recombinant HIV are added to the microplate and then incubated. The enzyme conjugate with antigens recombinant HIV binds to the antigen-antibody complexes immobilized HIV present. After this step, the Substrates A and B are added and incubated, producing a blue color that indicates the amount HIV antibody present in the samples. The Stop Solution is added to stop the reaction, having a color change to yellow, as in a microplate reader.

REAGENTS

- 1- Sensitive Plate** - Store between 2 and 8°C. Microplate coated with recombinant HIV antigens.
- 2- Conjugate** - Store between 2 and 8°C. Recombinant HIV antigens linked to Peroxidase and preservative.
- 3- Concentrated Washing Solution** - Store between 2 and 8°C. Buffer solution, surfactant and preservative.
- 4- Substrate A** - Store between 2 and 8°C. Buffer solution containing Hydrogen Peroxide and preservative.
- 5- Substrate B** - Store between 2 and 8°C. Buffer containing Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8°C. Chloridric Acid 1M.
- 7- HIV Negative Control** - Store between 2 and 8°C. Solution non-reactive for HBsAg, HCV, HIV-1 and HIV-2 and preservative. **Potentially infectious.**
- 8- HIV-1 Positive Control** - Store between 2 and 8°C. Solution containing HIV-1 and negative for HIV-2, HCV, HBsAg and preservative. **Potentially infectious.**
- 9- HIV-2 Positive Control** - Store between 2 and 8°C. Solution containing HIV-2 and negative for HIV-1, HCV, HBsAg and preservative. **Potentially infectious.**
- 10- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1- Sensitive Plate	1 Unit (96 Cavities)	2 Units (96 Cavities)	5 Units (96 Cavities)
2- Conjugate	1 Flask x12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3- Concentrated Washing Solution	1 Flask x 20 mL	2 Flasks x 20 mL	4 Flasks x 20 mL
4- Substrate A	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
5- Substrate B	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
6- Stop Solution	1 Flask x 8 mL	2 Flasks x 8 mL	4 Flasks x 8 mL
7- Negative Control	1 Flask x 1 mL	2 Flasks x 1 mL	4 Flasks x 1 mL
8- HIV-1 Positive Control	1 Flask x 1 mL	2 Flasks x 1 mL	4 Flasks x 1 mL
9- HIV-2 Positive Control	1 Flask x 1 mL	2 Flasks x 1 mL	4 Flasks x 1 mL
10- Plate sealer	3 Units	5 Units	10 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table.
- Operating instructions (manual).

Required materials not contained in the kit:

- 1- Pipette (s) capable of dispensing 50 and 100 µL volumes with a precision of better than 1,5%.
- 2- Re-pipettor (s) for repetitive dispensing of volumes of 100 µL and 300 µL, with precision greater than 1,5% (optional) or multichannel pipette.
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Adjustable volume pipettes (200 µL to 1000 µL) for dilution of the Substrate.
- 6- Test tubes for dilution of the Substrate A and B.
- 7- Paper towel to dry cavities.
- 8- Stopwatches or watch.
- 9- Container to store Washing Solution, after diluted.
- 10- Distilled or deionized water.
- 11- Tools of Quality Control.
- 12- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. The transport at temperatures between 15 and 30°C must not exceed 72 (seventy two) hours. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The envelope containing the strips should be opened only after it reaches room temperature. Place the strip with unused cavities in the aluminum bag, seal and store between 2 and 8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manual manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.
- 7- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 8- As a safety measure, you should cover the plate during the reaction.
- 9- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- 10- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 11- Stop Solution contains Chloridric Acid, which is a strong acid.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin). Hemolyzed or highly lipemic samples should not be used. Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer). For samples tested in duplicate, the required volume is of 0,010 mL of serum/plasma.

DESCRIPTION OF PROCESS

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of the Flask Nº 3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. Store between 2 and 8°C until expiration date printed on the original bottle. Can be stored at room temperature. In case of crystallization, heat it at 37°C until dissolution.

Substrate - Working Solution

Determine the amount of cavities to be used for preparation of an appropriate volume. Prepare the solution by mixing equal parts of Substrate A and Substrate B, 15 minutes before use. Keep it protected from light until used.
To each cavity (test), use:
50 µL of Substrate A + 50 µL of Substrate B
For example: mix 1 mL of Substrate A and 1 mL of Substrate B to two strips with 8 cavities (16 tests). Leftover reagent occurs.
Use no later than one (1) hour after preparation.

TECHNIQUE

Before starting the assay, bring all Reagents, Samples and Controls to stabilize at room temperature (15 - 30°C) for at least 40 minutes. Return unused strips to the original sealed packaging.
1- Select the cavities to be used considering: Controls and Samples. (They can be tested in duplicates).
2- Select the first cavity for Blank (OPTIONAL).
3- Pipette 100 µL of Negative Control, HIV-1 Positive Control, HIV-2 Positive Control and Samples into the cavities previously determined.
4- Homogenize gently for ± 30 seconds. Cover the cavities with sealer.
5- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
6- Remove the plate sealer of the cavities.
7- Discard the contents of the cavities by aspiration (Washer) or by decanting (manual); Use approximately 300 µL of Washing Solution **previously prepared*** to make a total of five (5) washing cycles.
To ensure the drying of the plate at the end of washing, beat it for a few seconds on absorbent paper.
Note: Washing / drying can cause poor results inadequate.
8- Pipette 100 µL of Conjugate into each cavity except in the Blank (If you have made this choice).
9- Mix gently for ± 30 seconds. Cover the cavities with the plate sealer.
10- Incubate for 20 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
11- Remove the plate sealer of the cavities.
12- Repeat item 7.
13- ATTENTION: follow one of the following procedures:
A) Pipette 100 µL of the **prepared Substrate*** - Working Solution (A + B) in all cavities.
*** See PREPARATION OF REAGENTS.**
Or
B) Pipette 50 µL of Substrate A and 50 µL Substrate B in all cavities.
14- Mix gently for ± 30 seconds. Cover the cavities with sealer.
15- Incubate for 10 minutes ± 1 minute in an incubator at 37°C ± 2°C.
16- Remove the plate sealer.
17- Pipette 50 µL of Stop Solution to each cavity.
18- Mix gently for ± 30 seconds.
19- Perform reading at 450 nm (primary filter) / 630 nm (secondary filter) for 30 minutes maximum.

TECHNIQUE VERIFICATION

Verify if the results obtained for the Controls reading are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0,050
Negative Control	< 0,100
Positive Control HIV-1 / HIV-2	> 0,500

The absorbance for each Control were obtained after the reduction of the Blank absorbance. For a single filter reading (450 nm) to consider the Blank limit of < 0,100. If the values are out of the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Calculate the average absorbance from negative control.

ITEM	ABSORBANCE
Negative Control	A1 = 0,020
	A2 = 0,018
Average Absorbance from Negative Control	(0,020 + 0,018) / 2 = 0,019

If the results obtained from the Control are valid, calculate the Cut-Off according to the formula:
Example:

ITEM	ABSORBANCE
Cut-Off = Negative Control Average Absorbance + 0,160	0,019 + 0,160 = 0,179

Calculate the index by dividing the sample by the Cut-Off value.

Example:

ITEM	ABSORBANCE
Sample	1,320
Cut-Off Value	0,179
Index = Sample / Cut-Off Value	1,320 / 0,179 = 7,37

INTERPRETATION OF RESULTS

RESULTS	QUALITATIVE
	INDEX
Negative	< 0,9
Positive	> 1,1
Undetermined	0,9 - 1,1

Non Reactive: Sample with absorbance less than the Cut-Off is considered non-reactive to antibodies HIV-1, HIV-2, and/or subtype O and may be considered negative.

Reactive: Sample with absorbance above the Cut-Off is considered initially reactive for antibodies to HIV-1, HIV-2, and/or subtype O. The sample must be retested in duplicate before the final interpretation. The sample that is reactive in at least one of the reanalysis is presumed to be reactive and must be confirmed by a confirmatory method IFIs and/or Western blot.

Note: In case of indeterminate results, the sample must be retested in duplicate. The sample results repeatedly indeterminate should be retested using a method alternative. If the results remain uncertain, one must collect a new sample.
The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and / or treatment of the patient.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test, there should not be based on a single run. This should include confirmation of other tests before a sample is considered positive. A negative result does not exclude the possibility of exposure. Finally, all results should be interpreted in conjunction with other clinical information available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
REPEATABILITY
The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

REPEATABILITY	SAMPLE 1	SAMPLE 2	SAMPLE 3
Average	0,692	1,540	2,237
Standard Deviation	0,032	0,067	0,073
Coefficient of Variation (%)	4,608	4,338	3,243

REPRODUCIBILITY
The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

REPRODUCIBILITY	SAMPLE 1	SAMPLE 2	SAMPLE 3
Average	0,683	1,558	2,255
Standard Deviation	0,011	0,019	0,016
Coefficient of Variation (%)	1,605	1,191	0,697

Clinical sensitivity and specificity
BIOLISA HIV 1/2/O kit analyzed clinical samples in comparison with other methods of EIA.
The results show that the clinical sensitivity of the BIOLISA HIV 1/2/O kit is > 99,9% and clinical specificity is 99,8%.

BIOLISA HIV 1/2/O x EIA REFERENCE				
METHOD		REFERENCE EIA		Total
BIOLISA HIV 1/2/O	Results	Positive	Negative	
	Positive	74	2	76
	Negative	0	1043	1043
Total		74	1045	1119

Clinical Sensitivity: > 99,9% (95,1 - 100,0%)*
Overall Agreement: 99,8% (99,4 - 100,0%)*
Clinical Specificity: 99,8% (99,3 - 100,0%)*
*Confidence Interval 95%

DIAGNOSTIC SIGNIFICANCE
HIV is the causative agent of Acquired Immune Deficiency Syndrome (AIDS). The main routes of transmission include exposure to blood and blood products, including the sharing of needles and syringes, sexual contact and transmission from mother to child. The virion is surrounded by a lipid envelope that is derived the host cell membrane. Several viral glycoproteins are in the envelope. Each virus contains two copies of genomic RNA. The HIV-1 was isolated from AIDS patients and healthy people with high potential risk for the development of AIDS. The HIV-1 consists of subtype M and subtype O. Strains highly divergent HIV-1 were recognized in 1990 and provisionally grouped as subtype O, because this variation was similar to markers of HIV-1 glycoprotein, but with a slight variation for the marker protein. Although seldom compared to HIV-1 and HIV-2, infections caused by the subtype O have been identified so far in Africa (Cameroon), France and Germany. HIV-2 was isolated AIDS patients and asymptomatic HIV-positive individuals in West Africa. The HIV-1, HIV-2, and Subtype O, induce an immune response. The detection of HIV antibodies in serum, plasma or whole blood is more efficient and a common way of determining whether an individual was exposed to HIV. Despite differences in their biological characteristics, serological activity and genome sequences, HIV-1, HIV-2 and the subtype O showed strong antigenic cross-reactivity. Most HIV-2 positive sera can be identified by serological tests with HIV-1. The test uses recombinant antigens of HIV to selectively detect antibodies to HIV-1, HIV-2 and/or subtype O in serum or plasma.

NUMBER OF TESTS
Presentation 1 - 96 tests
Presentation 2 - 192 tests
Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Chang, SY, Bowman, BH, Weiss, JB, Garcia, RE and White, TJ. The origin of HIV-1 isolate HTLV-III.B. Nature (1993) 3;363:466-9.
2. Arya, SK, Beaver, B, Jagodzinski, L, Ensoli, B, Kanki, PJ,Albert, J, Fenyo, EM, Biberfeld, G, Zagury, JF and Laure, F. New human and simian HIV-related retroviruses possess functional transactivator (tat) gene. Nature (1987) 328:548-550.
3. Caetano JA Immunologic aspects of HIV infection. Acta Med Port (1991) 4 Suppl 1:52S-58S.
4. Janssen, RS, Satten, GA, Stramer, SL, Rawal, BD, O'Brien, TR, Weiblen, BJ, Hecht, FM, Jack, N, Cleghorn, FR, Kahn, JO, Chesney, MA and Busch MP. New testing strategy to detect early HIV-1 infection for use in incidence estimates and for clinical and prevention purposes. JAMA (1998) 280(1):42-48.
5. Travers, K, Mboup, S, Marlink, R, Gueye-Nidaye, A, Siby, T, Thior, I, Traore, I, Dieng-Sarr, A, Sankale, JL and Mullins, C. Natural protection against HIV-1 infection provided by HIV-2. Science (1995) 268:1612-161. 6. Greenberg, AE, Wiktor, SZ, DeCock, KM, Smith, P, Jaffe HW and Dondero, TJ, Jr. HIV-2 and natural protection against HIV-1 infection. Science (1996) 272:1959-1960.
7. Bioclin – Dados de arquivos

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA HIV 1/2/O kit: 10269360200

Review: January/2016

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED