BIOLISA HIV 1/2/0 ANTÍGENO/ANTICORPO K 119

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FINALIDADE

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TÓXICO

BIOLISA HIV 1/2/0 ANTIGEN/ANTIBODY K119
USAGE INSTRUCTIONS

FUNCTION
BIOLISA HIV 1/2/O Antigen/Antibody kit is a fourth generation immunoassay for the qualitative detection of presence of HIV-1 P24 antigen and total antibodies (IgG, IgM and IgA), HIV-1, HIV-2, and/or subtype O in serum or human plasma. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION
Methodology: Enzyme immunoassay or immunoenzymatic
BIOLISA HIV 1/2/O Antigen/Antibody kit is a solid phase enzyme immunoassay based the principle of “sandwich” for the detection of HIV-1 P24 antigen and total antibodies (IgG, IgM and IgA) for HIV-1, HIV-2, and/or subtype O in human serum or plasma. The microplate is coated with monoclonal antibodies and recombinant HIV antigens specific to HIV-1 (p24, gp41) HIV-2 (gp 36) and subtype O. HIV Antigens HIV-1 P24 and/or antibodies of HIV-1, HIV-2, and/or subtype O, present in the sample bind to Antibodies/Antigens coated on the microplate antigen-antibody complexes formed. After initial incubation, the microplate is washed to remove unbound materials. The enzyme-conjugated polyclonal HIV recombinant antigens and is added to the microplate and then incubated. The enzyme-conjugate HIV polyclonal antibody and recombinant antigen complex binds to the immobilized antigen-antibody present. After the second incubation, the microplate is washed to remove materials not connected. After this step, the substrates A and B are added and incubated, producing a blue color that indicates the amount of antigen/HIV antibody present in the samples. The stop solution is added to stop the reaction, with a color change to yellow, so far as a reader microplates.

REAGENTS
1- Sensitive Plate - Store between 2 and 8°C. Microplates coated with monoclonal antibody and recombinant HIV antigens.
2- Conjugate - Store between 2 and 8°C. Policlonal antibody and recombinant HIV antigens linked to Peroxidase and preservative.
3- Concentrated Washing - Store between 2 and 8°C. Buffer solution, surfactant and preservative.
4- Sample Diluent - Store between 2 and 8°C. Buffer Solution and preservative.
5- Substrate A - Store between 2 and 8°C. Buffer solution containing Hydrogen Peroxide and preservative.
6- Substrate B - Store between 2 and 8°C. Buffer solution containing Tetramethylbenzidine (TMB) and preservative.
7- Stop Solution - Store between 2 and 8°C. Chloridric acid 1M.
8- HIV Negative Control - Store between 2 and 8°C. Non-reactive solution for HBsAg, HCV, HIV-1 and HIV-2 and preservative.
9- HIV-1 Positive Control - Store between 2 and 8°C. Solution containing HIV-1 and HIV-2 negative, 2HCV, HBsAg and preservative.
10- HIV-2 Positive Control - Store between 2 and 8°C. Solution containing HIV-2 and negative for HIV-1, HCV, HBsAg and preservative.
11- HIV-1 P24 Positive Control - Store between 2 and 8°C. Solution containing HIV- P24 and negative for HCV, HBsAg and preservative.
12- Plate Sealers

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1- Sensitive Plate	1 Unit (96 Cavities)	2 units (96 Cavities)	5 units (96 Cavities)
2- Conjugated	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3- Concentrated Washing	1 Flask x 20 mL	2 Flasks x 20 mL	4 Flasks x 20 mL
4- Sample Diluent	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
5- Substrate A	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
6- Substrate B	1 Flask x 8 mL	2 Flasks x 8 mL	5 Flasks x 8 mL
7- Stop Solution	1 Flask x 8 mL	2 Flasks x 8 mL	4 Flasks x 8 mL
8- Negative Control	1 Flask x 1 mL	2 Flasks x 1 mL	4 Flasks x 1 mL
9- HIV-1 Positive Control	1 Flask x 1 mL	2 Flasks x 1 mL	4 Flasks x 1 mL
10- HIV-2 Positive Control	1 Flask x 1 mL	2 Flasks x 1 mL	4 Flasks x 1 mL
11- HIV-1 P24 Positive Control	1 Flask x 1 mL	2 Flasks x 1 mL	4 Flasks x 1 mL
12- Plate sealer	3 Units	5 units	10 Units

EQUIPMENTS AND OPERATIONAL INPUTS
Materials in the kit:
- Reagents described in the above table
- Operating instructions (manual)

Required materials not contained in the Kit:
1- Pipette (s) capable of dispensing 25, 50 and 100 µL volumes with a lower coefficient of variation than 1,5%.
2- Re-pipettor (s) for repetitive dispensing of volumes of 100 µL and 300 µL, with lower coefficient of variation than 1,5% or multichannel pipette (optional).
3- Microplate washer (optional).
4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
5- Adjustable volume pipettes (200 µL to 1000 µL) for dilution of the substrate.
6- Test tubes for dilution of the substrate A and B.
7- Paper towel to dry cavities.
8- Stopwatches or watch.
9- Container to store washing solution after diluted.
10- Distilled or deionized water.
11- Tools of Quality Control.
12- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS
The storage temperature should be 2 to 8°C. The transport can be in under ambient temperature (up to 30°C) for up to 72 (seventy two) hours. **Do not freeze.** Keep away from light and avoid moisture.

SPECIAL CARE
1- For *in vitro* diagnostic use only.
2- Strictly follow the methodology proposed to obtain accurate results.
3- The envelope containing the strips should be opened only after it reaches room temperature. Place the strip with unused cavities in the aluminum bag, seal and store between 2 and 8°C.
4- The water used in material cleaning must to be recent and free of contaminants.
5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
6- All the raw material of product is tested and should be nonreactive for HBsAg and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manual manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.
7- Always add reagents in the same order to minimize the difference in reaction time between the cavities.

8- As a safety measure, you should cover the plate during the reaction.
9- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
10- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
11- Stop Solution contains Hydrochloric Acid, which is a strong acid. Handle it with care.
12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Data Sheet Safety of Chemicals) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

SAMPLES
Use serum or plasma (EDTA or Heparin).
Hemolyzed or highly lipemic samples should not be used.
Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).
For samples that will be tested as duplicates, the required volume is of 0,010 mL of serum/plasma.

DESCRIPTION OF PROCESS

PREPARATION OF WORKING REAGENT
Washing Solution
Dilute the contents of the Flask N° 3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. Can be stored at room temperature. In case of crystallization, heat it at 37°C until dissolution.

Substrate - Working Solution
Determine the amount of cavities to be used for preparation of an appropriate volume. Prepare the solution by mixing equal parts of Substrate A and Substrate B, 15 minutes before use. Keep it protected from light until used. To each cavity (test), use:
50 µL of substrate A + 50 µL of substrate B
For example, mix 1 mL of Substrate A and 1 mL of Substrate B to two strips with 8 cavities (16 tests). Leftover reagent occurs.
Use no later than one (1) hour after preparation.

TECHNIQUE
Before starting the assay, bring all Reagents, Samples and Controls to stabilize temperature (15 - 30°C) for at least 40 minutes.
Return unused strips to the original sealed packaging.
1- Select the cavities to be used considering: Controls and Samples. (They can be tested in duplicates).
2- Select the first cavity for Blank (OPTIONAL).
3- Pipette 50 µL of Sample Diluent into all the cavities.
4- Pipette 50 µL of Negative Control, HIV-1 Positive Control, HIV-2 Positive Control, HIV-1 P24 Positive Control and Samples into the cavities previously determined.
5- Homogenize gently for ± 30 seconds. Cover the cavities with sealer.
6- Incubate for 60 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
7- Remove the sealing from all the cavities.
8- Discard the contents of the cavities by aspiration (Washer) or by decanting (manual).
Use approximately 300 µL of washing solution **previously prepared*** to perform a total of five (5) washing cycles.
To ensure the drying of the plate at the end of washing, beat it for a few seconds on absorbent paper.
Note: Poor Washing and drying can cause inadequate results.
9- Pipette 100 µL of Conjugate into each cavity except in the Blank (If you have made this choice).
10- Mix gently for ± 30 seconds. Cover the cavities with the sealer.
11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
12- Remove the plate sealer.

13- Repeat Item 8.
14- ATTENTION: follow one of the following procedures:
A) Pipette 100 µL of the **previously prepared*** substrate - Working solution (A + B) in all cavity.
*** See PREPARATION OF WORKING REAGENTS.**
Or
B) Pipette 50 µL of Substrate A and 50 µL Substrate B in all cavity.
15- Mix gently for ± 30 seconds. Cover the cavities with sealer.
16- Incubate for 30 minutes ± 1 minute in an incubator at 37°C ± 2°C.
17- Remove the plate sealer.
18- Pipette 50 µL of Stop Solution to each cavity.
19- Homogenize gently for ± 30 seconds.
20- Perform reading at 450 nm (primary filter)/630 nm (secondary filter) up to 30 minutes maximum.

TECHNIQUE VERIFICATION
Verify if the results obtained for the controls reading are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0,050
Negative Control	< 0,200
Positive Control HIV-1 / HIV-2 / HIV-1 P24	> 0,500

The absorbance for each control were obtained after the reduction of the blank absorbance. For a single filter reading (450 nm) to consider the blank limit of < 0,100. If the values are out of the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE
Calculate the average absorbance from Negative Control.
Example:

ITEM	ABSORBANCE
Negative Control	A1 = 0,040
	A2 = 0,038
Average Absorbance from Negative Control	(0,040 + 0,038) / 2 = 0,039

If the results obtained from the Control are valid, calculate the Cut-Off according to the formula:
Example:

ITEM	ABSORBANCE
Cut-Off = Negative Control Average Absorbance + 0,160	0,039 + 0,160 = 0,199

Calculate the index dividing the sample absorbance by the Cut-Off value.
Example:

ITEM	ABSORBANCE
Sample	1,900
Cut-Off Value	0,199
Index = Sample / Cut-Off value	1,900 / 0,199 = 9,55

INTERPRETATION OF RESULTS

RESULTS	QUALITATIVE
	INDEX
Negative	< 0,9
Positive	> 1,1
Undetermined	0,9 - 1,1

Non Reactive: Sample with absorbance less than the Cut-off is considered non-reactive for HIV-1 P24 Antigens and HIV-1, HIV-2 antibodies, and/or subtype O can be considered negative.

Reaction: Sample with absorbance equal to or above the Cut-Off is considered initially reactive for HIV-1 P24 antigens and antibodies to HIV-1, HIV-2, and/or subtype O. The sample must be retested in duplicate before the end of interpretation. The sample that is reactive in at least one of the reanalysis, if presumed to be reactive and must be confirmed by a confirmatory method IFIs and/or Western blot.

Note: In case of indeterminate results, the sample must be retested in duplicate. The Sample results repeatedly indeterminate should be retested using a method alternative. If the results remain uncertain, one must collect a new sample.

Note: The date presented in the examples are for ilustration only and can be used for calculation of the results.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test, should not be based on a single run. This should include confirmation of other tests before a sample is considered positive. A negative result does not exclude the possibility of exposure. Therefore, all results should be interpreted in conjunction with other clinical information available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

Precision

REPEATABILITY

20 successive measurements were performed with three samples, using the same lot, obtaining the following results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	0,697	1,561	2,291
Standard Deviation	0,032	0,062	0,055
Coefficient of Variation (%)	4,523	4,002	2,393

REPRODUCIBILITY

20 measurements were performed for 3 consecutive days with three samples, using the same lot, obtaining the following results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	0,697	1,580	2,298
Standard Deviation	0,001	0,022	0,018
Coefficient of Variation (%)	0,143	1,388	0,790

Clinical sensitivity and specificity

BIOLISA HIV 1/2/O Antigen/Antibody Kit analyzed clinical samples in comparison with other methods of EIA.

The results show that the clinical sensitivity of the BIOLISA HIV 1/2/O Antigen/Antibody kit is > 99,9% and clinical specificity is 99,8%.

method		REFERENCE EIA		Total
BIOLISA HIV 1/2/O Antigen/ Antibody		Positive	Negative	
	Positive	146	2	148
	Negative	0	1240	1240
Total		146	1242	1388

Clinical Sensitivity: > 99,9% (97,5 - 100,0%)*

Overall Agreement: 99,9% (99,5 - 100,0%)*

Clinical Specificity: 99,8% (99,4 - 100,0%)*

*Confidence Interval 95%

DIAGNOSTIC SIGNIFICANCE

HIV is the causative agent of Acquired Immune Deficiency Syndrome (AIDS). The main routes of transmission include exposure to blood and blood products, including the sharing of needles and syringes, sexual contact, transmission from mother to child. The virus is surrounded by a lipid envelope that is derived the host cell membrane. Several viral glycoproteins are in the envelope. Each virus contains two copies of genomic RNA positive sense. The HIV-1 was isolated from AIDS patients and people healthy with high potential risk for the development of AIDS. Infection with HIV-1 is identified by the initial antigenemy to the antigen (Ag) HIV-1 is detectable in the blood. In most cases, the antigen levels are often difficult to detect, however, the increase in failure immune system and increased levels of the virus may re-stimulate detectable levels of antigen. The major internal structural proteins of HIV-1 core protein P24, is a component found in the blood during viral antigenemy. Furthermore, HIV-1 subtype consists of M and subtype O. Highly divergent strains of HIV-1 were recognized in 1990 and provisionally grouped as the subtype O, because this variation was similar to markers of HIV-1 glycoprotein, but with a slight variation for the marker protein. Although seldom compared to HIV-1 and HIV-2 infection caused by the subtype so far have been identified in Africa (Cameroon), France and Germany. HIV-2 was isolated from patients with AIDS and asymptomatic HIV-positive individuals in WestAfrica. HIV-1, HIV-2 and subtype O, induce an immune response. The immunological detection of antigens and antibodies HIV in serum, plasma or whole blood is more efficient and a common way of determining whether an individual was exposed to HIV. Despite differences in their biological, serological and genome sequence activities, HIV-1, HIV-2 subtype O shows strong antigenic cross-reactivity. Antibodies against HIV-1 P24 antigen are also included for the detection of pre-serum conversion infection. The majority of HIV-2 positive sera can be identified by serological tests HIV-1.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Chang, SY, Bowman, BH, Weiss, JB, Garcia, RE and White, TJ. The origin of HIV-1 isolate HTLV-IIIb. Nature (1993) 3;363:466-9.
2. Arya, SK, Beaver, B, Jagodzinski, L, Ensoli, B, Kanki, PJ,Albert, J, Fenyo, EM, Biberfeld, G, Zagury, JF and Laure, F. *New human and simian HIV-related retroviruses possess functional transactivator (tat) gene*. Nature (1987) 328:548-550.
3. Caetano JA *Immunologic aspects of HIV infection*. Acta Med Port (1991) 4 Suppl 1:52S-58S.
4. Janssen, RS, Satten, GA, Stramer, SL, Rawal, BD, O'Brien, TR, Weiblen, BJ, Hecht, FM, Jack, N, Cleghorn, FR, Kahn, JO, Chesney, MA and Busch MP. New testing strategy to detect early HIV-1 infection for use in incidence estimates and for clinical and prevention purposes. JAMA (1998) 280(1):42-48.
5. Travers, K, Mboup, S, Marlink, R, Gueye-Nidaye, A, Siby, T, Thior, I, Traore, I, Dieng-Sarr, A, Sankale, JL and Mullins, C. Natural protection against HIV-1 infection provided by HIV-2. Science (1995) 268:1612-161.
6. Greenberg, AE, Wiktor, SZ, DeCock, KM, Smith, P, Jaffe HW and Dondero, TJ, Jr. HIV-2 and natural protection against HIV-1 infection. Science (1996) 272:1959-1960.
7. Bioclin - Dados de arquivos.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA HIV 1/2/O Antigen/Antibody Kit: 10269360199

Review: February/2017

English 2/2

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON