Bioclin

BIOLISA RUBÉOLA IgM

REF K124

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FINALIDADE

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NÚMERO DE TESTES

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- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology, 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE. A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE. A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Test for qualitative determination of IgM antibodies to Rubéola in human samples of serum, plasma and filter paper whole blood by immunoassay microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA Rubella IgM kit is a solid phase immunoenzymatic assay based on the principle of qualitative detection by capture of IgM Antibodies to Rubella samples of serum, plasma and filter paper whole blood. IgM antibodies present in the sample bind to the anti-IgM Antibody coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Rubella antigens conjugated to Peroxidase are added to the microplate which is then incubated. Antigens conjugated to the enzyme bind to the Anti-Rubella IgM Antibodies present, attached to the plate coated with Anti-IgM Antibodies. New washing is performed to remove surplus. After this step, the Substrate is added and incubated, producing a blue color indicating the amount of Anti-Rubella IgM Antibodies present in the sample. The Stop Solution is added to stop the reaction with a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1 - Sensitized Plate** - Store at 2 to 8°C. Contains: Plate sensitized with Anti-IgM antibodies.
- 2- Conjugate** - Store between 2 and 8°C. Contains: Buffer Solution, Rubella Antigen conjugated to Peroxidase surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** - Store at 2 and 8°C. Contains: Buffer Solution, surfactant and preservative.
- 4- Sample Diluent** - Store at 2 and 8°C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
- 5- Substrate** - Store at 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store at 2 and 8°C. Contains: Hydrochloric Acid 1 M.
- 7- Negative Control** - Store at 2 and 8°C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infecting.**
- 8- Positive Control** - Store between 2 and 8°C. Contains: Antibodies IgM Anti-Rubella, Buffer Solution, dye, stabilizers, surfactant and preservative. **Potentially infectious.**
- 9- Plate Sealers**
- 10- Filter Paper Elution Buffer** - Store at 2 and 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 11- Negative Extraction Control** - Store at 2 and 8°C. Contains: Non-reactive sample for IgM anti-rubella antibodies impregnated in Filter Paper. **Potentially infectious.**
- 12- Positive Extraction Control** - Store at 2 to 8°C. Contains: Reactive sample for IgM Anti-Rubella antibodies impregnated in Filter Paper. **Potentially infectious.**

PRESENTATION

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1- Sensitized Plate	1 unit (96 cavities)	2 units (192 cavities)	5 units (480 cavities)
2- Conjugate	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3- Concentrated Washing	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
4- Sample Diluent	1 Flask x 42 mL	2 Flasks x 42 mL	5 Flasks x 42 mL
5- Substrate	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
6- Stop Solution	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
7- Negative Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
8- Positive Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
9- Plate Sealers	3 units	6 units	15 units
10- Elution Buffer in Filter Paper	1 Flask x 20 mL	2 Flask x 20 mL	5 Flask x 20 mL
11- Positive Extraction Control	1 Unit	2 Units	5 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Operating instructions (manual)

Required materials not contained in the kit:

- 1- Pipette capable of dispensing volumes for 5 to 300 µL with lower coefficient of variation than 1.5%.
- 2- Re-pipettor for repetitive pipetting volumes of 300 µL with lower coefficient of variation than 1.5% or multichannel pipette (Optional).
- 3- Microplate Washer (optional).
- 4- ELISA Reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Paper towel to dry cavities
- 6- Watch or stopwatch.
- 7- Flask to store the washing solution after diluted.
- 8- Distilled or deionized water.
- 9- Tools of Quality Control.
- 10- Incubator 37°C ± 2°C.
- 11- Paper shredder (3 mm diameter) for fi lter paper technique.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. The transport in up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store at 2 to 8°C.

- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Saturated deionized columns release alkaline water several ions and oxidizing and reducing agents that can significantly alter the results;
- 6- Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.
- 7- All the raw material of product is tested and should be nonreactive for HBs Antigen, HIV 1 & 2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases.
- Therefore, we must take due care in handling the bio safety of these products.
- 9- As a safety measure, you should cover the plate during the reaction.
- 10- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- 11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin).
Use serum or plasma (EDTA or Heparin). Hemolyzed or highly lipemic samples should not be used. Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

Filter paper whole blood

Whole Blood (Puncture or EDTA).
Samples dried on filter paper can be stored at room temperature provided they are protected from direct sunlight and low humidity. For storage of up to 2 years, the samples should remain refrigerated at 2 to 8°C.
Storage for more than 2 years must be carried out the temperature of -20°C.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of the Flask N°3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. Can be stored at room temperature. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

STABILITY AFTER OPEN

The results of the stability test show that the Biolisa kit Rubella IgM is stable after opening for at least 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

TECHNIQUE

Sample of Serum or Plasma

- Before starting the assay, bring all Reagents, Samples and Controls to stabilize at room temperature (15 - 30°C) for at least 40 minutes.
- 1- Select the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.
 - 2- Select the first cavity for Blank (OPTIONAL).
 - 3- Pipette 100 µL of Sample Diluent. Even in the cavity for Blank.
 - 4- Pipette 5 µL of Sample and Controls into previously determined wells. Observe the color change of the diluent at the time of sample addition. The color change indicates that the sample was added to the well correctly.
 - 5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
 - 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
 - 7- Remove the sealing of wells.
 - 8- Discard the contents of the wells by aspiration (Washer) or by decanting (manual). Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
 - Note:** Poor washing and drying can cause inadequate results.
 - 9- Pipette 100 µL of Conjugate into each well included in the Blank cavity (if you made this option).
 - 10- Mix gently for ± 30 seconds. Cover cavities with plate sealer.
 - 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
 - 12- Remove the sealer from plate cavities.
 - 13- Repeat item 8.
 - 14- Pipette 100 µL of Substrate into all wells, included in the Blank cavity.
 - 15- Shake gently for ± 30 seconds. Cover wells with plate sealer.
 - 16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
 - 17- Remove the sealer from plate cavities.
 - 18- Pipette 100 µL of Stop Solution into all wells, included in the Blank cavity.
 - 19- Mix gently for ± 30 seconds.
 - 20- Read using double filter: 450nm/630nm within up to 15 minutes (maximum).
- Filter Paper Whole Blood Sample**
Before starting the assay, bring all Reagents, Samples and Controls to stabilize at room temperature (15 - 30°C) for at least 40 minutes.
- 1- Select the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.
 - 2- Select the first cavity for Blank (OPTIONAL).
 - 3- Add a disk 3 mm of the Positive Extraction Control, Negative Extraction Control and Whole Blood Samples on Filter Paper in the previously determined wells.
 - 4- Pipette 100 µL of Elution Buffer in filter paper, included in the Blank cavity
 - 5- Pipette 5 µL of Controls into previously determined wells.
 - 6- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
 - 7- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C;
 - 8- Remove the sealing of wells.
 - 9- Discard the contents of the wells by aspiration (Washer) or by decanting (manual). Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
 - Note:** Poor washing and drying can cause inadequate results.
 - 10- Pipette 100 µL of Conjugate into each well, included in the Blank cavity (if you made this option).
 - 11- Mix gently for ± 30 seconds. Cover cavities with plate sealer.
 - 12- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

- 13- Remove the sealer from plate cavities.
- 14- Repeat item 9.
- 15- Pipette 100 µL of Substrate into all wells, included in the Blank cavity.
- 16- Shake gently for ± 30 seconds. Cover wells with plate sealer.
- 17- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 18- Remove the sealer from plate cavities.
- 19- Pipette 100 µL of Stop Solution into all wells, included in the Blank cavity.
- 20- Mix gently for ± 30 seconds.
- 21- Read using double filter: 450nm/630nm within up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Samples of Serum, Plasma and Whole Blood Filter Paper

Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0.150
Negative Control	< 0.150
Positive Control	> 1.000
Negative Extraction Control	< 0.300
Positive Extraction Control	> 1.000

If the values fall out from the expected values, technique must be repeated.

CALCULATIONS

QUALITATIVE

Sample of Serum or Plasma

Calculate Cut Off according to the formula below:

Cut Off = (Positive Control Mean Abs. x 0.1) + 0.200

Example:

ITEM	ABSORBANCE
Positive Control	A1 = 2.134
	A2= 2.190
Cut Off = (Positive Control Mean Abs. x 0.1) + 0.200	Cut Off = ((2.134 + 2.190)/2 x 0.1) + 0.200 Cut Off = 0.416

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.
Example:

ITEM	ABSORBANCE
Sample	1.279
Cut Off Value	0.416
Index: Sample / Cut Off Value	1.279 / 0.416 = 3.074

Filter Paper Whole Blood Sample

Calculate Cut Off according to the formula below:

Cut Off = (Positive Control Mean Abs. x 0.02) + 0.200

Example:

ITEM	ABSORBANCE
Positive Control	2.134
	2.190
Cut Off = (Positive Control Mean Abs. x 0.02) + 0.120	Cut Off=((2.134+ 2.190)/2 x 0.02) + 0.200 Cut Off = 0.243

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.
Example:

ITEM	ABSORBANCE
Sample	1.045
Cut off Value	0.243
Index: Sample / Cut Off Value	1.045 / 0.243 = 4.3

INTERPRETACION OF RESULTS

Samples of serum plasma and filter paper whole blood.

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (No Reagent)	< 0.8
Undetermined	Between 0.8 and 1.2
Positive (Reagent)	> 1.2

The results provided by this kit should be interpreted by the responsible medical professional and not the only criterion for determining the diagnosis and/or treatment of the patient.

PROCEDURE LIMITATIONS

IgM Antibodies are present in the serum for a short period of time, disappearing up to six months after infection, while IgG Antibodies remain present for a long period of time. However, the high sensitivity of the new diagnostic methodologies, in some cases, allows to find very low levels of IgM Antibodies, called residual IgM, for a longer period of time. Detection of these residual IgM Antibodies hinders clinical interpretation over the period of infection. In these cases, to confirm the result, an IgG avidity test is recommended.

The interpretation of a diagnostic test, there should not be based on a single run. This should include confirmation of other tests before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other clinical information available before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	2.523	0.416	1.038
Standard Deviation	0.057	0.030	0.033
Coefficient of Variation (%)	2.269	7.096	3.150

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	2.458	0.371	0.978
Standard Deviation	0.079	0.041	0.052
Coefficient of Variation (%)	3.211	11.008	5.329

Clinical Sensitivity and Specificity

Serum or Plasma Sample of Serum or Plasma

BIOLISA Rubella IgM kit analyzed clinical samples in comparison with other methods of EIA.The results show that the clinical sensitivity of the BIOLISA Rubella IgM kit is > 99.99% and clinical specificity is 98.35%.

	Expected Result	Biolisa CMV IgG
Positive Sample	59	59
Negative Sample	121	119
Total Tested Samples	180	

Clinical Sensitivity: > 99.99% (59/59)

Clinical Specificity: 98.35% (119/121)

Filter Paper Whole Blood Sample

BIOLISA Rubella IgM kit analyzed clinical samples in comparison with other methods of EIA.The results show that the clinical sensitivity of the BIOLISA Rubella IgM kit is > 99.99% and clinical specificity is > 99.99%.

	Expected Result	Biolisa CMV IgG
Positive Sample	53	53
Negative Sample	67	67
Total Tested Samples	96	

Clinical Sensitivity: > 99.99% (53/53)

Clinical Specificity: > 99.99% (67/67)

DIAGNOSTIC SIGNIFICANCE

Rubella is a RNA virus, spherical, small envelope, which belongs to the Togaviridae family. It is commonly known as german or 3 days measles. The Rubella virus infection is transmitted through droplets of saliva, resulting in mild contagious rash in children or young adults. In childhood, infection is a self-limiting, benign, characterized by low fever, headache, lymphadenopathy, arthralgia, and conjunctivitis. However, infection during pregnancy, especially in the first quarter, can lead to miscarriage, intrauterine infection causes fetal death, or congenital anomalies. Congenital Rubella depends on the period in which infection occurs and can result in complications disorders, including deafness, eye problems, including cataracts and glaucoma, congenital heart disease and mental retardation. IgM antibodies against Rubella are produced initially, reaching detectable levels within 2 - 3 days and peak within 14 - 21 days after onset of symptoms that remain detectable during the next 4 - 8 weeks. The diagnosis of active or recent infection can be obtained by the presence of IgM antibodies in the initial sample. After several days, IgG antibodies appear after IgM, with a peak of 14 - 21 days, varying levels persist for a lifetime. The presence of anti-Rubella IgG antibodies is indicative of preview infection and/or immunity.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
2. Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
3. Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
4. Voller, A, Bidwell, DE, A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
5. Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
6. Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176
- 7- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
8. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentatio packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA Rubella IgM Kit: 10269360190

Review: July/2020

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

MANUFACTURED BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON