BIOLISA TOXOPLASMOSE IgM

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FINALIDADE

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1. Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
2. Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
3. McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
4. Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
5. Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
6. Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med.(1978) 298(10):550-553.
7. Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43..
8. Bioclin – Datos de archivos

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USAGE INSTRUCTIONS

FUNCTION

Test for qualitative determination of IgM antibodies to *Toxoplasma gondii* in serum or human plasma by enzyme immunoassay in microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA TOXOPLASMOSIS IgM kit is a solid phase immunoenzymatic assay based on the principle of qualitative detection by capture of IgM Antibodies for *Toxoplasma gondii* in human serum or plasma. IgM antibodies present in the sample bind to the anti-IgM Antibody coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. *Toxoplasma gondii* antigens conjugated to Peroxidase are added to the microplate which is then incubated. Enzyme-conjugated antigens bind to the Anti-*Toxoplasma gondii* IgM Antibodies present, attached to the plate coated with anti-IgM antibodies. New washing is performed to remove surplus. After this step, the Substrate is added and incubated, producing a blue color indicating the amount of Anti-*Toxoplasma gondii* IgM Antibodies present in the sample. The Stop Solution is added to stop the reaction with a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8°C.
- 2- Conjugate** - Store between 2 and 8°C. Contains: *Toxoplasma gondii* antigens conjugate to Peroxidase and preservative.
- 3- Concentrated Washing Solution** - Store between 2 and 8°C. Contains: Buffer Solution, surfactant and preservative.
- 4- Sample Diluent** - Store between 2 and 8°C. Contains: Buffer Solution and preservative.
- 5- Substrate** - Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8°C. Contains: Chloridric Acid 1M.
- 7- Negative Control** - Store between 2 and 8°C. Contains: Non-reactive IgM antibodies to *Toxoplasma gondii* and preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8°C. Contains: IgM Antibodies Anti-*Toxoplasma gondii* and preservative. **Potentially infectious.**
- 9- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1- Sensitized Plate	1 Unit (96 Cavities)	2 Units (192 Cavities)	5 Units (480 Cavities)
2- Conjugate	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3- Concentrated Washing Solution	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
4- Sample Diluent	1 Flask x 42 mL	2 Flasks x 42 mL	5 Flasks x 42 mL
5- Substrate	1 Flask x 11 mL	2 Flasks x 11 mL	5 Flasks x 11 mL
6- Stop Solution	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
7- Negative Control	1 Flask x 100 µL	2 Flasks x 100 µL	5 Flasks x 100 µL
8- Positive Control	1 Flask x 100 µL	2 Flasks x 100 µL	5 Flasks x 100 µL
9- Plate sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Operating instructions (manual)

Required materials not contained in the Kit:

- 1- Pipette capable of dispensing volumes of 5 and 100 µL with lower coefficient of variation than 1,5%.
- 2- Re-pipettor for repetitive pipetting volumes of 100 µL with lower coefficient of variation than 1,5% or multicontrol pipette (Optional).
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Paper towel to dry cavities
- 6- Stopwatch or watch.
- 7- Flask to store the Washing Solution after diluted.
- 8- Distilled or deionized water.
- 9- Tools of Quality Control.
- 10- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. The transport can be done under ambient temperature (up to 30°C) for up to 72 (seventy two) hours. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store at 2-8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- Stop Solution contains Chloridric Acid, which is a strong acid. Handle it with care.
- 7- All the raw material of product is tested and should be nonreactive for HBsAg, HIV 1 & 2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.
- 8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 9- As a safety measure, you should cover the plate during the reaction.

- 10- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- 11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin). Hemolyzed or highly lipemic samples should not be used. Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).

DESCRIPTION OF PROCESS

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of the Flask N°3 (Concentrated Washing Solution) in 1000 mL of distilled water or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. Can be stored at room temperature. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

TECHNIQUE

- Before starting the assay, bring all reagents, Samples and Controls to stabilize at room temperature (15 - 30°C) for at least 40 minutes.
- 1- Select the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the plate will not be used for the original sealed packaging.
- 2- Select the first cavity for Blank (OPTIONAL).
- 3- Pipette 100 µL of Sample Diluent. Even in the cavity for White.
- 4- Pipette 5 µL of Sample and Controls into previously determined wells. Observe the color change of the diluent at the time of sample addition. The color change indicates that the sample was added to the well correctly.
- 5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 7- Remove the plate sealers of cavities.
- 8- Discard the contents of the cavities by aspiration (Washer) or by decanting (Manual).
- Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate at the end of the wash, beat the plate for a few seconds on absorbent paper.
- Note:** Poor washing and drying can cause inadequate results.
- 9- Pipette 100 µL of Conjugate in all cavities except in the Blank cavity (if you made this option).
- 10- Mix gently for ± 30 seconds. Cover cavities with plate sealer.
- 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 12- Remove the sealer from plate cavities.
- 13- Repeat item 8.
- 14- Pipette 100 µL of Substrate into all wells.
- 15- Mix gently for ± 30 seconds. Cover wells with plate sealer.
- 16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 17- Remove the plate sealer from cavities.
- 18- Pipette 100 µL of Stop Solution into all wells.
- 19- Mix gently for ± 30 seconds.
- 20- Read using double filter: 450 nm / 630 nm within up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank Cavity	< 0,150
Negative Control	< 0,150
Positive Control	> 1,000

If the values are out of the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the formula below:
Cut Off = (Positive Control Mean Abs. * 0,1) + 0,100

Example:

ITEM	ABSORBANCE
Positive Control	A1 = 1,780
	A2= 1,700
Cut Off = (Positive Control Mean Abs. * 0,1) + 0,300	Cut Off = ((1,780 + 1,700/2) * 0,1) + 0,100 Cut Off = 0,274

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1,245
Cut Off Value	0,274
Index: Muestra/ Cut Off Value	1,245 / 0,274 = 2,544

INTERPRETATION OF RESULTS

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (No Reagent)	≤ 0,8
Undetermined	Between 0,8 and 1,2
Positive (Reagent)	≥ 1,2

The results provided by this kit should be interpreted by the responsible medical professional and not the only criterion for determining the diagnosis and/or treatment of the patient.

PROCEDURE LIMITATIONS

IgM antibodies are present in the serum for a short period of time, disappearing up to six months after infection, while IgG antibodies remain present for a long period of time. However, the high sensitivity of the new diagnostic methodologies, in some cases, allows to find very low levels of IgM antibodies, called residual IgM, for a longer period of time. Detection of these residual IgM antibodies hinders clinical interpretation over the period of infection.
In these cases, to confirm the result, an IgG avidity test is recommended.

The interpretation of a diagnostic test, should not be based on a single run. This should include confirmation of other tests before that a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other clinical information available before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE
QUALITY CONTROL

Accuracy

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	1,893	1,658	0,076
Standard Deviation	0,096	0,027	0,008
Coefficient of Variation (%)	5,105	1,642	11,058

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	1,827	1,663	0,079
Standard Deviation	0,016	0,020	0,011
Coefficient of Variation (%)	0,893	1,236	13,947

Clinical sensitivity and specificity

BIOLISA Toxoplasmosis IgM kit analyzed clinical samples in comparison with other kit of EIA. The results show that the clinical sensitivity of the BIOLISA Toxoplasmosis IgM kit is 98,18% and clinical specificity is > 99,99%.

	Expected Result	Biolisa Toxoplasmose IgM
Positive Sample	55	54
Negative Sample	57	57
Total Tested Samples	112	

Clinical Sensitivity: 98,18 % (54/55)

Clinical Specificity: > 99,99% (57/57)

DIAGNOSTIC SIGNIFICANCE

Toxoplasma gondii is the causative agent of toxoplasmosis. It is an obligate intracellular protozoan that has been found in many species of birds, reptiles and mammals. The agent can be transmitted through organ transplantation, transfusion of blood and leukocytes, contact with feces of contaminated cats and ingestion of contaminated raw meats. In adults, the infection is usually benign or asymptomatic. However, symptomatic cases, including fatal cases, occur in immunosuppressed patients who have clinical or laboratory evidence of damage to the central nervous system. Following infection, IgM antibodies appear within 5 days and have reduced levels within a few weeks or months. IgG antibodies usually appear 1 - 2 weeks post infection, reaching peak levels in 6-10 weeks persisting for life. In children, the risk of fetal infection varies according to the time of pregnancy when the mother is infected. In maternal infections that occur during the first trimester, the

probability of infection passing to the fetus is lower. However, if transmission occurs, serious outcomes such as miscarriage and hydrocephalus are more likely. Infections acquired later in pregnancy, where fetal transmissions occur more frequently, tend to be less severe but can still generate congenital manifestations, including cerebral calcifications and poor learning.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
2. Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
3. McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
4. Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
5. Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
6. Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
7. Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
8. Bioclin – Dados de arquivos

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA Toxoplasmosis IgM kit: 10269360194

Review: January/2019

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED