BIOLISA TOXOPLASMOSE IgM

REF K126

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Test for qualitative determination of IgM antibodies to *Toxoplasma gondii* in serum or human plasma by enzyme immunoassay in microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic

The BIOLISA Toxoplasmosis IgM kit is a solid phase enzyme immunoassay based on the principle of indirect qualitative detection of IgM Antibodies to *Toxoplasma gondii* in human serum or plasma. Anti-*Toxoplasma gondii* IgM Antibodies, present in the samples bind to Antigens coated on the microplate forming complexes Antigen-Antibody IgM. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated human Anti-IgM Antibodies are added to the microplate and then incubated. The enzyme-conjugated Anti-Human IgM Antibodies bind to the antigen-bound IgM Antibodies present. A new wash is performed to remove the surplus. After this step, the Substrate is added and incubated, producing a blue color indicating the amount of Anti-*Toxoplasma gondii* IgM Antibodies present in the samples. The Stop Solution is added to stop the reaction with a change in color from blue to yellow, measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8°C.
- 2- Conjugate** - Store between 2 and 8°C. Contains: Human Anti-IgM Antibody linked to Peroxidase and preservative.
- 3- Concentrated Washing Solution** - Store between 2 and 8°C. Contains: Buffer Solution, surfactant and preservative.
- 4- Sample Diluent** - Store between 2 and 8°C. Contains: Buffer Solution and preservative.
- 5- Substrate** - Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8°C. Contains: Chloridric Acid 1M.
- 7- Cut-Off Calibrator** - Store between 2 and 8°C. Contains: Diluted IgM Antibodies Anti-*Toxoplasma gondii* and preservative. **Potentially infectious.**
- 8- Negative Control** - Store between 2 and 8°C. Contains: Non-reactive IgM antibodies to *Toxoplasma gondii* and preservative. **Potentially infectious.**
- 9- Positive Control** - Store between 2 and 8°C. Contains: IgM Antibodies Anti-*Toxoplasma gondii* and preservative. **Potentially infectious.**
- 10- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1- Sensitized Plate	1 Unit (96 Cavities)	2 Units (192 Cavities)	5 Units (480 Cavities)
2- Conjugate	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3- Concentrated Washing Solution	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
4- Sample Diluent	1 Flask x 42 mL	2 Flasks x 42 mL	5 Flasks x 42 mL
5- Substrate	1 Flask x 11 mL	2 Flasks x 11 mL	5 Flasks x 11 mL
6- Stop Solution	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
7- Cut-Off Calibrator	1 Flask x 100 µL	2 Flasks x 100 µL	5 Flasks x 100 µL
8- Negative Control	1 Flask x 100 µL	2 Flasks x 100 µL	5 Flasks x 100 µL
9- Positive Control	1 Flask x 100 µL	2 Flasks x 100 µL	5 Flasks x 100 µL
10- Plate sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Operating instructions (manual)

Required materials not contained in the Kit:

- 1- Pipette capable of dispensing volumes of 5, 50 and 100 µL with lower coefficient of variation than 1,5%.
- 2- Re-pipettor for repetitive pipetting volumes of 100µL and 300 µL, with lower coefficient of variation than 1,5% or multichannel pipette (Optional).
- 3- Micro plate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Adjustable volume pipettes (200 µL to 1000 µL) for sample dilution.
- 6- Test tubes for preparation for sample dilution.
- 7- Paper towel to dry cavities
- 8- Watch or stopwatch.
- 9- Flask to store the Washing Solution after diluted.
- 10- Distilled water.
- 11- Tools of Quality Control.
- 12- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. The transport can be done under ambient temperature (up to 30°C) for up to 72 (seventy two) hours. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store at 2-8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- Stop Solution contains Chloridric Acid, which is a strong acid. Handle it with care.

- 7- All the raw material of product is tested and should be nonreactive for HBsAg, HIV 1 & 2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.
- 8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 9- As a safety measure, you should cover the plate during the reaction.
- 10- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- 11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin). Hemolyzed or highly lipemic samples should not be used. Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples can not be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).

DESCRIPTION OF PROCESS

PREPARATION OF WORKING REAGENTS

Washing Solution
Dilute the contents of the Flask N°3 (Concentrated Washing Solution) in 1000 mL of distilled water or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. Can be stored at room temperature. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

TECHNIQUE

- Before starting the assay, bring all Reagents, Samples, Cut-Off Calibrator and Controls to stabilize at room temperature (15-30°C) for at least 40 minutes.
- 1- Select the wells to be used considering: Cut-Off Calibrator, Controls and Samples (it is recommended to test in duplicate). Return the strips of the plate will not be used for the original sealed packaging.
 - 2- Select the first cavity for Blank (OPTIONAL).
 - 3- Prepare a 1:41 dilution of the Samples, Controls and Cut-Off Calibrator in microtubes, adding 5 µL more 200 µL of the Sample Diluent. Homogenize.
 - 4- Pipette 100 µL of diluted Samples and Controls and Cut-Off Calibrator into the previously determined wells.
 - 5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
 - 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
 - 7- Remove the sealing of wells.
 - 8- Discard the contents of the wells by aspiration (Washer) or by decanting (manual).
 - Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles.
 - To ensure drying of the plate at the end of the wash, hit the board for a few seconds on absorbent paper.
 - Note:** Poor washing and drying can cause inadequate results.
 - 9- Pipette 100 µL of Conjugate into each well except in the Blank cavity (if you made this option).
 - 10- Mix gently for ± 30 seconds. Cover cavities with plate sealer.
 - 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
 - 12- Remove the sealer from plate cavities.
 - 13- Repeat item 8.
 - 14- Pipette 100 µL of Substrate into all wells.
 - 15- Shake gently for ± 30 seconds. Cover wells with plate sealer.

- 16- Incubate for 15 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 17- Remove the sealer from plate cavities.
- 18- Pipette 100 µL of Stop Solution into all wells.
- 19- Mix gently for ± 30 seconds.
- 20- Read the 450 nm within 15 minutes (maximum).

TECHNIQUE VERIFICATION

Verify if the results obtained by the reading of the Blank, Controls and Cut-Off Calibrator are compatible to the values below:

ITEM	ABSORBANCE
Blank Cavity	< 0,150
Negative Control	< 0,5
Cut-Off Calibrator	0,5 to 1,5
Positive Control	> Cut-Off Calibrator

The absorbancies for the Controls and Cut-Off Calibrator were obtained after the reduction of the Blank absorbance. Consider a Blank limit of < 0,150. If the values are out of the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Cut-Off regarded as the mean absorbance obtained with the Cut-Off Calibrator. Example:

ITEM	ABSORBANCE
Cut-Off Calibrator	0,952
	0,885
Average Absorbance from Cut-Off Calibrator	(0,952 + 0,885) / 2 = 0,919

Calculate the Index dividing the sample absorbance by the value of the Cut-Off. Example:

ITEM	ABSORBANCE
Sample	1,347
Cut-Off Value	0,919
Index: Sample / Cut-Off Value	1,347 / 0,919 = 1,465

Note: The data presented in the examples are for illustration only and can not be used for calculations of the results.

INTERPRETATION OF RESULTS

RESULTS	QUALITATIVE
	INDEX
Negative	≤ 0,9
Undetermined*	Between 0,9 and 1,5
Positive	≥ 1,5

* Suspicion of Residual IgM for results between 1,0 and 1,5.

Observación: Due to the high sensitivity of the product, detection of residual IgM antibodies may occur (See Procedure Limitations). In the case of results with an index between 1,0 and 1,5, residual IgM antibodies may be suspected, and therefore the IgG avidity test is recommended. In case of results with index between 0,9 and 1,0, it is recommended to repeat the examination by another alternative method or to wait 2 weeks for a new analysis.

PROCEDURE LIMITATIONS

IgM antibodies are present in the serum for a short period of time, disappearing up to six months after infection, while IgG antibodies remain present for a long period of time. However, the high sensitivity of the new diagnostic methodologies, in some cases, allows to find very low levels of IgM antibodies, called residual IgM, for a longer period of time. Detection of these residual IgM antibodies hinders clinical interpretation over the period of infection. In these cases, to confirm the result, an IgG avidity test is recommended.

The interpretation of a diagnostic test, should not be based on a single run. This should include confirmation of other tests before that a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other clinical information available before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE

QUALITY CONTROL

Accuracy

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	0,990	1,862	0,905
Standard Deviation	0,069	0,148	0,043
Coefficient of Variation (%)	6,943	7,965	4,747

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	0,967	1,817	0,882
Standard Deviation	0,068	0,146	0,042
Coefficient of Variation (%)	7,008	8,052	4,807

Clinical sensitivity and specificity

BIOLISA Toxoplasmosis IgM kit analyzed clinical samples in comparison with other kit of EIA. The results show that the clinical sensitivity of the BIOLISA Toxoplasmosis IgM kit is > 99,9% and clinical specificity is 97%.

BIOLISA Toxoplasmosis IgM X Reference EIA				
METHOD		REFERENCE EIA		TOTAL
		Positive	Negative	
BIOLISA TOXOPLASMOSIS IgM	Positive	19	4	23
	Negative	0	130	130
Total Results		19	134	153

Clinical Sensitivity: > 99,9% (19/19)
Clinical Specificity: 97% (130/134)
Overall Agreement: 97,4%

DIAGNOSTIC SIGNIFICANCE

Toxoplasma gondii is the causative agent of toxoplasmosis. It is an obligate intracellular protozoan that has been found in many species of birds, reptiles and mammals. The agent can be transmitted through organ transplants, blood transfusions and leukocyte, contact with cat feces contaminated and ingestion of contaminated raw meats. In adults, the infection is usually benign or asymptomatic. However, symptomatic cases, including fatal cases occur in patients immunosuppressed patients who have clinical or laboratory evidence of central nervous system damage. In children, the risk of fetal infection varies according to the time of pregnancy when the mother is infected. In maternal infections that occur during the first quarter, the probability of infection through to the fetus is smaller. However, if transmission occurs, severe outcomes such as miscarriage and hydrocephalus are more likely. Infections acquired later in pregnancy, where the transmissions fetal occur most frequently tend to be less severe, but even they can generate congenital manifestations, including brain calcification and impaired learning. After infection, IgM antibodies appear 5 days and had reduced levels within a few weeks or months. IgG antibodies usually appear 1 - 2 weeks after infection, reaching levels peak at 6 - 10 weeks persisting for life.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
2. Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
3. McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
4. Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
5. Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
6. Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
7. Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
8. Bioclin – Dados de arquivos

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA Toxoplasmosis IgM kit: 10269360194

Review: June/2018

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED