BIOPROT U/LCR

REF K108

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FINALIDADE

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1. Johnson AM, Rohlfis EM, Silverman LM. Proteínas. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3 ed. Philadelphia: W.B Saunders Company;1999. p. 477-540.
2. Felgenhauer K. Laboratory diagnosis of neurological diseases. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1308-26.
3. Boege F. Urinary proteins. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998.p. 382-400.
4. Orsonneau JL, Douet P, Massoubre C, Lustenberger P, Bernard S. An improved pyrogallol red-molybdate method for determining total urinary protein. Clin Chem 1989;35:2233-6.
5. Watanabe N, Kamei S, Ohkubo A, Yamanaka M, Ohsawa S, Makino K et al. Urinary protein as measured with a pyrogallol red-molybdate complex. Manually and in a Hitachi 726 automated analyzer. Clin Chem 1986; 32:1551-4.
6. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
7. Bioclin – Dados de arquivo.

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Total Protein in urine and liquor samples. Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Colorimetric Test - Pyrogallol Red

The proteins present in urine react to Pyrogallol Red and the Molybdate forming a red color complex, with maximum absorption in 600 nm (578 - 630 nm). The intensity of the color is directly proportional to the protein concentration in the sample.

REAGENTS

Reagent N° 1 - Color Reagent - Store between 2 and 8°C. Contains: Pyrogallol Red < 60 µmol/L, Sodium Molybdate > 40 µmol/L, surfactant, solubilizer, stabilizer and preservative.

Reagent N° 2 - Standard - Store between 2 and 8°C. Contains: Albumin 65 mg/dL and preservative.

Reagent N° 3 - Control - Store between 2 and 8°C. Contains: Proteins* and preservative.

* Value varies with each batch.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2	Reagent N° 3
1	1 x 25 mL	1 x 3,0 mL	1 x 3,0 mL
2	2 x 25 mL	1 x 3,0 mL	1 x 3,0 mL
3	3 x 25 mL	1 x 3,0 mL	1 x 3,0 mL
4	4 x 25 mL	1 x 3,0 mL	1 x 3,0 mL
5	5 x 20 mL	1 x 3,0 mL	1 x 3,0 mL
6	5 x 40 mL	1 x 3,0 mL	1 x 3,0 mL
7	10 x 40 mL	1 x 3,0 mL	1 x 3,0 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the spectrophotometer or colorimeter, watch and stopwatch, pipettes and test tubes. Can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 2 and 8°C. The transport temperature between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- **For professional *in vitro* diagnostic use only.**

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- The presence of detergent in the used material can alter the results significantly.

6- The calibration should be repeated periodically, to verify any alteration in the colorimeter or spectrophotometer response.

7- Reagents N° 2 and 3 contain Sodium Azide, irritating to skin and mucous membrane. Handle with caution.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Urine (aleatory sample and 24 hours) and Cerebrospinal Fluid (Liquor).

Urine sample is stable for 1 day between 15 and 30°C, 7 days between 2 and 8°C, and 1 month at -20°C.⁶

The sample of Liquor is stable for 1 day between 15 and 30°C, 6 days between 2 and 8°C, and 1 year at -20°C.⁶

Preparation of samples:

Urine - There is no necessity of preservatives. Take an aliquot of 20 mL of urine, centrifuge and proceed the technique using the supernatant.

Liquor - Centrifuge the sample and use the supernatant.

In dosages for therapeutic control the samples must be collected always in the same time.

PROCESS DESCRIPTION

REAGENT PREPARATION

The reagents are ready to use.

TECHNIQUE

Mark 3 tests tubes with the letters B (Blank), S (Standard), A (Sample) and proceed as follows:

	B	S	A
Sample	--	--	20 µL
Standard	--	20 µL	--
Distilled Water	20 µL	--	--
Reagent N° 1	1000 µL	1000 µL	1000 µL

Homogenize and incubate for 5 minutes at 37°C and read the absorbance in a 600 nm filter (578 to 630 nm), against a Blank reagent in no more than 30 minutes.

CALCULATIONS

The determination of total protein concentration in urine and liquor is calculated by the reaction:

Protein (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 65$

Example:

In a total protein dosage in liquor, the following absorbance was found:

Sample Absorbance = 0,089
Standard Absorbance = 0,160
Standard Concentration = 65

Protein (mg/dL) = $\frac{0,089}{0,160} \times 65 = 36,16 \text{ mg/dL}$

According to the Lambert-Beer law, you can use the Calibration Factor:

Calibration Factor = $\frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$

Example:

Cf = $\frac{65}{0,160} = 406,25$

Protein (mg/dL) = Sample Absorbance x Cf
Protein (mg/dL) = 0,089 x 406,25 = 36,16
The results are express the mg/dL.

24 hours Urine

Proceed the 24 hours dosage using a aliquot of the total volume. This aliquot must be centrifuge. The concentration in mg/24h is obtained using the dosage result (mg/dL) multiplying this value by the total urine volume (mL) and divide the result by 100.

24 Hour Urine Protein (mg/24h) = $\frac{\text{Aliquot dosage (mg/dL)} \times \text{24 Hours volume (mL)}}{100}$

Example:

In a dosage of total protein in 24 hour urine, the following result was obtained:

Sample dosage (mg/dL) = 5,0
24 Hours Urine volume (mL) = 2000

24 Hours Urine Protein (mg/24h) = $\frac{5,0 \times 2000}{100} = 100 \text{ mg/24 hours}$

INTERFERENCES

Blood presence in urine or liquor produce wrongly elevated results. The presence of detergent in the used material can altered significantly the results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D.

REFERENCE VALUES

The reference values for this method were obtained through the determination of the Total Protein in healthy populations of male and female.

Random Urine 1 to 15 mg/dL

24 Hours Urine 24 to 141 mg/24h

Liquor 15 to 45 mg/dL

To convert the values of mg/dL to g/L (SI), multiply the obtained results by 0,01. These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The Bioprot U/CRL was compared with other commercially available method for dosage of Total Protein in urine. 41 analyzes were performed and the results were evaluated. The linear equation obtained was Y = 1,000X - 0,080, and correlation coefficient 0,99999. With theses results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/24h)	32,84	15,25	27,97
Standard Deviation (mg/24h)	0,46	0,17	0,17
Coefficient of Variation (%)	1,40	1,13	0,59

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/24h)	31,96	15,44	27,89
Standard Deviation (mg/24h)	0,77	0,16	0,10
Coefficient of Variation (%)	2,42	1,06	0,37

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Total Protein in urine. The average found was 1,484 mg/dL, with standard deviation of 0,093 mg/dL. The sensitivity, wich indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 1,764 mg/dL.

Linearity

The reaction is linear until the concentration of 100 mg/dL. For higher values, dilute the sample with NaCl solution at 0,85% and repeat the dosage. Multiply the results by the applied dilution factor.

DIAGNOSTIC SIGNIFICANCE

The proteinuria in high concentration is found in the majority of the kidney problems, such as primary and secondary nephropathies, that can cause the increase of the glomerular permeability or reduction of the kidney tubules absorption. It could also have post kidney causes, such as infections, bleedings or malignant urinary tract diseases. In other acute situations such as fever or physical or psychological stress, high levels in protein can also be found. In liquor, high levels of protein can be found in cases of increase intracranial pressure due to cerebral tumors, intracerebral hemorrhage or traumatisms, in inflammations (mainly due to bacterial meningitis) and in case of multiple sclerosis. Increase the permeability of the barrier of cephalorachiidian liquid-blood is reflected as the increase proportion of liquor/serum of total protein.

BIBLIOGRAPHIC REFERENCES

1. Johnson AM, Rohlfs EM, Silverman LM . Proteinase. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3 ed. Philadelphia: W.B Saunders Company; 1999. p. 477-540.
2. Felgenhauer K. Laboratory diagnosis of neurological diseases. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1308-26.
3. Boege F. Urinary proteins. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 382-400
4. Orsonneau JL, Douet P, Massoubre C, Lustenberger P, Bernard S. An improved pyrogallol red-molybdate method for determining total urinary protein. Clin Chem 1989;35:2233-6.
5. Watanabe N, Kamei S, Ohkubo A, Yamanaka M, Ohsawa S, Makino K et al. Urinary protein as measured with a pyrogallol red-molybdate complex. Manually and in a Hitachi 726 automated analyzer. Clin Chem 1986; 32:1551-4.
6. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
7. Bioclin – Dados de arquivo.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC

REP

OBELIS S.A.
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Bioprot U/CRL kit: 10269360178

Review: February/2019

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

Σ

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

EC

REP

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED