CMV PCR

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Baldanti, F.; Lilleri, D.; Gerna, G. Monitoring human cytomegalovirus infection in transplant recipients. J Clin Virol 2008; 41:237–41.
3. De Vries. J. J; Van der Eijk. A. A; Wolthersc. K. C. Real-time PCR versus viral culture on urine as a gold standard in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2011:2352.
4. Grosse, S. D; Ross, D. S; Dollard, SC. Congenital cytomegalovirus (CMV) infection as a cause of permanent bilateral hearing loss: a quantitative assessment. J Clin Virol 2008;41:57–62.

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Baldanti, F.; Lilleri, D.; Gerna, G. Monitoring human cytomegalovirus infection in transplant recipients. J Clin Virol 2008; 41:237–41.
3. De Vries, J. J.; Van der Eijk, A. A.; Wolthers, K. C. Real-time PCR versus viral culture on urine as a gold standard in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2011:2352.
4. Grosse, S. D.; Ross, D. S.; Dollard, S. C. Congenital cytomegalovirus (CMV) infection as a cause of permanent bilateral hearing loss: a quantitative assessment. J Clin Virol 2008;41:57–62.

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Review: March/2020

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	4
Equipaments and Operational Inputs	4
Transportation and storage Conditions	4
Special Care	4
Samples	5
Process Description	6
A . DNAExtraction	6
B . Preparation of the reagents	6
C . Dilution of Quantitative Standard	7
D . PCR Preparation	7
E . Thermocycler settings for Real-Time PCR	8
F . Result Validation	10
G. Result Interpretation	11
Process Limitations	12
Product Performance / Quality Control	12
Comparison of Methods and Methodology Specificity	12
Repeatability	12
Reproducibility	13
Clinical Sensitivity	13
Analytical Sensitivity	14
Linearity	14
Diagnostic Significance	14
Bibliographic References	14
Customer Service	14
Universal Symbolology	15

FUNCTION

Quantitative detection of Human cytomegalovirus (CMV) DNA through reverse transcription followed by real-time polymerase chain reaction (PCR). For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

The **Bio Gene CMV PCR** test is a *in vitro* assay based on the quantitative detection of CMV DNA by real-time PCR.

The Real-Time PCR method is used to amplify the pathogen DNA.

A thermocycler for Real-Time PCR is used to amplify and detect the fluorescent probe. Software calculates the CMV DNA concentration, expressed in copies/ μL , using the standard curve generated from the quantitative standard, contained in the kit.

Reagent	Presentation	
	50 Tests	150 Tests
R1	1 x 55 μL *	1 x 165 μL *
R2	1 x 525 μL *	3 x 525 μL *
R3	1 x 525 μL	3 x 525 μL
R4	1 x 55 μL *	1 x 165 μL *
R5	1 x 600 μL *	1 x 600 μL *
R6	1 x 50 μL	1 x 150 μL
R7	1 x 500 μL *	1 x 500 μL *
R8	2 x 1,5 mL	2 x 1,5 mL
R9	1 x 1 mL	1 x 2 mL

*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

REAGENTS

R1. PCR Solution: Primer, probe, TRIS-HCl.

R2. Mix Taq: Polymerases, dNTPs, $MgCl_2$, stabilizers.

R3. Mix Buffer: TRIS-HCl.

R4. IC PCR Solution: Primer, probe, TRIS-HCl.

R5. Internal Control: Plasmid, TRIS-HCl.

R6. Negative Control: TRIS-HCl.

R7. Standard A (2×10^5 copies/ μL): Plasmid, TRIS-HCl, EDTA.

R8. Diluent: TRIS-HCl, EDTA.

R9. Water: DNase/RNase free water.

EQUIPMENTS AND OPERATIONAL INPUTS**Materials in the kit:**

- Reagents described in the previous item.
- Usage instructions (manual).

Materials needed but not contained in the kit:

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- 1.5 ml microcentrifuge tubes.
- 4- PCR reaction tubes or plate.
- 5- Powder free latex gloves or similar material.
- 6- Microcentrifuge (3,000 - 12,000 rpm).
- 7- Vortex.
- 8- Micropipettes and sterile filter pipette tips (0.5-10 μL , 10-100 μL , 100-1000 μL).
- 9- Nucleic acid extraction kit.

TRANSPORTATION AND STORAGE CONDITIONS

The required storage temperature is -20°C (-10 to -30°C). **After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension.** Avoid freezing and thawing.

The kit may be transported between 2 and 30°C for 7 days. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.

- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.
- 4- It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.
- 5- It is necessary to have separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplified products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.
- 6- All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.
- 7- Avoid repeated thawing and freezing of the reagents.
- 8- Store the DNA samples at -20°C in case they will not be used immediately.
- 9- Do not use reagents after expiration date.
- 10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 12- Do not use the product in case of damaged packaging.
- 13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

This kit should be used with DNA samples extracted from serum and whole blood. Samples should be collected according to laboratory recommendations for molecular testing. They must be transported and stored as described below¹:

Serum: at -20 °C for 3 days.

Whole Blood: between 2 and 8 ° C for 3 days.

DNA (extracted): Use DNA samples suitable for PCR amplification with adequate purity and concentration.

For long-term storage, it is recommended to freeze samples at -70°C. Avoid repeated freezing and thawing of the samples.

PROCESS DESCRIPTION

A. DNA Extraction

The nucleic acids (DNA) must be extracted in accord to the procedure of the chosen extraction kit.

For control of the extraction process, the **Internal Control (R5)** must be prepared (see item B.3) and added to the samples during the extraction as described below:

- 1- Add 4 μ L of the **Internal Control (R5)** to each tube containing the resuspended samples in extraction/lysis buffer.
- 2- Complete the extraction procedure according to the usage instructions of the extraction kit.

Note: Never add the Internal Control directly to pure biological sample. This may result in degradation of the sample.

B. Preparation of the Reagents

Note: **R5** and **R7** reagents contain RNA/DNA template. They need to be handle (resuspended) in a separated area.

- 1- Before open, spin (*spin* pulse) each of the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)**, **Internal Control (R5)** and **Standard A (R7)**.
- 2- Resuspend each vial of **Mix Taq (R2)*** reagent with 525 μ L of **Mix Buffer (R3)** reagent.

*For the presentation of 150 tests resuspend the **Mix Taq (R2)** vials as need of use.

*The **Mix Taq (R2)** reagent does not contain passive reference Dye (ROX).

- 3- Resuspend the reagents, **PCR Solution (R1)** and **IC PCR Solution (R4)** with **Water (R9)** reagent, in according to the following table:

Reagent	Presentation	
	50 Tests	150 Tests
PCR Solution (R1)	55 µL	165 µL
IC PCR Solution (R4)	55 µL	165 µL

4- Resuspend the reagents, **Internal Control (R5)** and **Standard A (R7)** with **Diluent (R8)** reagent, in according to the following table:

Reagent	Presentantion	
	50 tests	150 tests
Internal Control (R5)	600 µL	600 µL
Standard A (R7)	500 µL	500 µL

C. Dilution of Quantitative Standard

1- Select 3 microtubes (not supplied in the kit) suitable for serial dilution of resuspended **Standard A (R7)**.

2- Pipette 90µL of the **Diluent (R8)** in each microtube and name them as B, C and D respectively.

3- Pipette 10µL of **Standard A (R7)** in microtube B and mix.

4- Replace the tip and pipette 10µL of microtube B in microtube C and mix.

5- Replace the tip and pipette 10µL of microtube C in microtube D and mix.

6- At the end of dilution procedure, it will be available 4 standards with the following concentrations:

Standard A - 2×10^5 copies/µL

Standard B - 2×10^4 copies/µL

Standard C - 2×10^3 copies/µL

Standard D - 2×10^2 copies/µL

D. PCR Preparation

1- Prior beginning the assay, select sufficient number of microtubes/wells to be used on the reaction, in accord to the number of samples, Controls and Quantitative Standards to be used.

2- Prepare the final PCR solution in accord with reactions number required.

Reagents	1 Reaction	25 Reactions	50 Reactions	100 Reactions
Mix Taq (R2)	10 µL	250 µL	500 µL	1 mL
PCR Solution (R1)	1 µL	25 µL	50 µL	100 µL
IC PCR Solution (R4)	1 µL	25 µL	50 µL	100 µL
Water (R9)	3 µL	75 µL	150 µL	300 µL

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

3- Pipette 15µL of the final PCR solution in microtubes or wells previously determined.

4- Add 5µL of extracted DNA sample or 5µL of Quantitative Standards or 5µL of **Negative Control (R6)** in previously determined tubes.

5- Mix well.

6- The total volume of the reaction is 20µL. Each PCR run should include the relevant controls: Negative Control, Internal Control and Quantitative Standards.

7- Transport the tubes/plate to the thermocycler and follow the instructions in the section E (Thermocycler Settings for Real-Time PCR).

E. Thermocycler Settings for Real-Time PCR

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

1- Set the type of experiment:

Quantitative Test with Standard Curve.

2- Set the fluorescent probes:

Target	Detector	Quencher
CMV	FAM	NFQ-MGB
Internal Control	VIC	

Note: The Negative Control and Quantitative Standards do not present the Internal Control, because it is used for extraction and amplification control of the samples.

3- Set the quantitative standards:

Standard A - 2×10^5 copies/ μ L

Standard B - 2×10^4 copies/ μ L

Standard C - 2×10^3 copies/ μ L

Standard D - 2×10^2 copies/ μ L

4- Set the cycling conditions:

Stage	Temperature	Time	Cycles
1	95°C	2 minutes	1
2	95°C	10 seconds	50
	60°C	60 seconds	

Set "Data Collection" as "stage 2, step 2 (60@0:60)".

F. Result Validation

1- Standard Curve

Standard Curve	Permitted Range	Amplification / Detection
Correlation Coefficient (R^2)	$0,99 \leq R^2 \leq 1,00$	Valid

If the R^2 value does not fall within the limits of range, the result must be considered invalid and the test must be repeated.

2- Negative Control

CT Negative Control		Result	Amplification/ Detection
FAM	VIC		
Undetermined	Undetermined	Negative	Valid

3- Samples

CMV		Result	Detection
FAM	VIC		
Determined concentration	$25 \leq CT \leq 31$	Positive	Valid
	$CT < 25$ ou $CT > 31$	Positive	Invalid*
Undetermined concentration	$25 \leq CT \leq 31$	Negative	Valid
	$CT < 25$ ou $CT > 31$	Negative	Invalid*

*The CT values of Internal Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, the concentration of the samples and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

Example: In some cases, the samples with high numbers of copies of DNA / RNA can inhibit the amplification of the Internal Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

G. Result Interpretation

This kit is able to detect from 10 to 1.000.000 copies of DNA per reaction.

The thermocycler software calculates automatically the samples concentration.

Example: If the program shows a concentration of 2.00E+005, this means that the sample concentration is 2.0×10^5 copies/ μ L.

Sample Result in Copies/μL (FAM)	Copies per Reaction
$\geq 1 \times 10^6$	$> 1.000.000$
$2 \leq \text{Quantity} \leq 9 \times 10^5$	Quantity obtained
< 2	< 10

The non-detection of the pathogen DNA does not exclude the presence of infection when the pathogen title is below the detection limit of this kit.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient.

The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

1- Conversion calculation for copies/mL

$$\text{Result in copies/mL} = \frac{\text{Copies/}\mu\text{L} \times \text{Elution volume (}\mu\text{L)}^{**}}{\text{Extracted sample volume (mL)}^{**}}$$

*Result of quantification of the sample provided by the equipment in copies/ μ L.

**According to the protocol of the extraction kit used.

Example: Elution volume: 30 μ L

Sample volume: 200 μ L (0.2 mL)

Quantification of the sample: 5000 copies/ μ L

$$\text{Result in copies/mL} = \frac{(5000 \text{ copies/}\mu\text{L} \times 30\mu\text{L})}{0.2\text{mL}}$$

$$\text{Result in copies/mL} = 750.000 \text{ copies/mL}$$

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

The heparin use in the samples collection for molecular tests is not recommended, due to the PCR inhibition which can generate false negative results.

The samples with presence of PCR inhibitors should be diluted to reduce their effect on the reaction.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The **Bio Gene CMV PCR** test was compared with another protocol for the quantification of the CMV DNA. 50 analyzes were performed and the results were evaluated. Comparing the two tests, the difference between the results obtained was ≤ 0.5 log. This show that the test has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 positive samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Average Concentration (Log)	5.114	4.319
Standard Deviation (Log)	0.070	0.094
Coefficient of variation (%)	1.374	2.198

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample with different concentrations, obtaining the following results:

	Day 1	Day 2	Day 3
Average Concentration (Log)	4.319	4.307	4.395
Standard Deviation (Log)	0.094	0.059	0.061
Coefficient of variation (%)	2.198	1.373	1.403

Clinical Sensitivity

The **Bio Gene CMV PCR** showed clinical sensitivity of 99.9% and a clinical specificity of 99.9%.

Method		Clinical Data		Total
Bio Gene CMV PCR	Result	Positive	Negative	
	Positive	25	0	25
	Negative	0	25	25
TOTAL		25	25	50

Analytical Sensitivity

It was found in the performance study that the analytical sensitivity of the kit is 2.40 copies/ μ L (0.38 log). This means the assay is capable to detect approximately 12 target molecules in 5 μ L DNA extraction product, used on the amplification reaction.

Linearity

The linearity found in the performance studies was 2.03×10^5 copies/ μ L (5.307 log).

DIAGNOSTIC SIGNIFICANCE

Human cytomegalovirus (CMV) belongs to Herpesviridae family and genus Cytomegalovirus and is now called Human betaherpesvirus 5 (HHV5). The CMV is a human pathogen transmitted through saliva, sexual contact, perinatal, organ transplant or blood transfusion. In most cases, the infection is asymptomatic. However, CMV infection can cause serious illness in neonatal and immunocompromised individuals, such as patients with HIV, cancer or patients who have received organ transplants. During immunosuppressive therapy, often there is a reactivation of a latent virus or a primary infection. CMV infections can be acquired before birth, during birth and later in life. It is may cause serious birth defects such as microcephaly, physical disabilities and mental retardation. Because of that, the primary determination of maternal infection and its distinction of latent infection are so important.

Therefore, the quantitative detection of CMV by real time PCR is very important due to its accuracy and agility.

BIBLIOGRAPHIC REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Baldanti, F.; Lilleri, D.; Gerna, G. Monitoring human cytomegalovirus infection in transplant recipients. J Clin Virol 2008; 41:237-41.
3. De Vries. J. J; Van der Eijk. A. A; Wolthersc. K. C. Real-time PCR versus viral culture on urine as a gold standard in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2011:2352.
4. Grosse, S. D; Ross, D. S; Dollard, SC. Congenital cytomegalovirus (CMV) infection as a cause of permanent bilateral hearing loss: a quantitative assessment. J Clin Virol 2008;41:57-62.

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for Bio Gene CMV PCR: 10269360245

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED

Bioclin · QUIBASA

Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130

Tel +55 31 3439 5454 . www.bioclin.com.br

FARM. RESP. Sílvia Wandalsen Arndt - CRF MG 7422

C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira