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Review: April/2020

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	3
Equipaments and Operational Inputs	4
Transportation and Storage Conditions.....	4
Special Care	4
Samples	5
Procedure	5
A . Preparation of Working Solutions	5
B . DNA/RNA Extraction	5
1. Samples Lysis	5
2. Preciitation.....	7
3. Binding.....	7
4. Washing.....	7
5. Elution.....	7
Process Limitations	8
Bibliographic References	8
Customer Service	8
Universal Symbology	9

FUNCTION

Product developed for Viral DNA /RNA extraction and purification in sputum, cerebrospinal fluid, plasma, whole blood, semen, serum, swab, tissue and urine. The nucleic acids obtained can be used for different applications like PCR, RT-PCR, Hibridization, Sequencing, etc.

PRINCIPLE OF ACTION

The **Bio Gene Extração de DNA/RNA Viral** kit is a product developed for extraction and purification of Viral DNA / RNA in biological samples.

The method used is the extraction by silica membrane. The process is carried out in 4 steps: 1) Cell Lysis: cell disruption for the release of nucleic acids; 2) Binding: selective binding of nucleic acid to silica membrane; 3) Washing: to remove residual impurities; 4) Elution: release the nucleic acid from silica membrane. At the end of the process we have concentrated nucleic acid with high purity.

Reagent	Presentation	
	50 Preps	250 Preps
R1	1 x 11 mL	1 x 55 mL
R2	1 x 22 mL	1 x 110 mL
R3	1 x 7 mL	1 x 33 mL
R4	1 x 13 mL	1 x 30 mL
R5	1 x 300 µg	4 x 400 µg
R6	1 x 300 µL	1 x 1.5 mL
R7	1 x 50 unid	1 x 250 unid
R8	3 x 50 unid	3 x 250 unid
R9	2 x 50 unid	2 x 250 unid

REAGENTS

R1. Lysis Buffer: Guanidine Hydrochloride.

R2. Wash 1: Guanidine Hydrochloride, alcohol.

R3. Wash 2: Buffer solution.

R4. RNase Free Water.

R5. Carrier RNA: Acid Polyadenylic.

R6. Proteinase K: Enzyme Proteinase K.

R7. Columns: Polypropylene tube with a silica membrane.

R8. Collectors Tubes (2 mL): Polypropylene tube.

R9. Collectors Tubes (1.5 mL): Polypropylene tube.

EQUIPMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the table above.
- Usage instructions (manual).

Materials required, but not included in the kit:

- 1- 96-100% Ethanol
- 2- PBS 1X may be required for some samples
- 3- Micropipettes and sterile filter pipette tips (0.5-10µL, 10-100µL, 100-1000µL).
- 4- Microcentrifuge
- 5- Vortex
- 6- Heating block.
- 7- Personal protective equipment (gloves, lab coat, glasses)

TRANSPORTATION AND STORAGE CONDITIONS

The kit may be transported between 15 and 30°C. Protect from light and avoid moisture.

The storage temperature is between 15 and 30°C.

After resuspended is recommended to store the Carrier RNA reagent (R5) at -20°C. After first use, is recommended to store the Proteinase K reagent (R6) -20°C.

SPECIAL CARE

1- For professional use only.

- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose the waste in accordance with appropriate security procedures.
- 4- It is required skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.
- 5- It is necessary to have separate areas for the extraction / preparation of reactions and for the amplification / detection of products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.
- 6- Avoid repeated thawing and freezing of the reagents.
- 7- Do not use reagents after expiration date.
- 8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 10- Do not use the product in case of damaged packaging.
- 11- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Kit for Viral DNA / RNA extraction from sputum, cerebrospinal fluid, plasma, whole blood, semen, serum, swab, tissue and urine. Samples must be collected and stored according to the laboratory's recommendations for molecular testing.

PROCEDURE**Attention:**

The Lysis Buffer (**R1**) and Wash 1 (**R2**) reagents contains Guanidine Salt, and must be handled using suitable PPEs (coats, gloves and goggles). These reagents must not be disposed of on Acid Solution or on Sodium Hypochlorite due to the formation of reactive components.

A. Preparation of Working Solutions

Resuspend **R3** and **R5** reagents as shown in the table below:

Reagent	Presentation	
	50 Preps	250 Preps
R3 - Wash 2 (concentrated)	Add 28 mL of Ethanol 96-100%	Add 132 mL of Ethanol 96-100%
R5 - Carrier RNA*	Add 300µL of RNase Free Water (R4)	Add 400µL of RNase Free Water (R4)

Before starting the extraction:

- 1- Note if the reagents Wash 2 (**R3**) and Carrier RNA (**R5**) are prepared.
- 2- Set the heating block at 56°C.
- 3- Preheat the volume of the RNase free water (**R4**) reagent at 70°C to elution.

B. DNA/RNA extraction**1. Samples Lysis****A. Plasma, Suero, Líquido Cefalorraquídeo y Orina**

- 1.1 Pipette 200µL of sample at room temperature in Collector Tube 1.5 mL (**R9**).
- 1.2 Add 5µL of Proteinase K (**R6**)* to the Collector Tube 1.5 mL with sample and mix.
- 1.3 Then, pipette 200µL of Lysis Buffer (**R1**) and shake vigorously (10-15 seconds).
- 1.4 Continue on step **2. Precipitation**.

B. Whole blood

- 1.1 Pipette 200µL of sample at room temperature in Collector Tube 1.5 mL (**R9**).
- 1.2 Add 5µL of Proteinase K (**R6**)* to the Collector Tube 1.5 mL with sample and mix.
- 1.3 Then, pipette 200µL of Lysis Buffer (**R1**) and shake vigorously (10-15 seconds).
- 1.4 Incubate for 5 minutes at 56°C.
- 1.5 Continue on step **2. Precipitation**.

* After use, store at - 20° C.

C. Swab in transport medium / stabilizer

- 1.1 Press the swab with the tip in the transport medium to remove the biological material (it can be centrifuged in a short spin).
- 1.2 Discard the swab and pipette 200 μ L of the transport medium containing the biological material in a 1.5 mL Collector Tube (**R9**).
- 1.3 Add 5 μ L of Proteinase K (**R6**)* to the Collector Tube 1.5 mL with sample and mix.
- 1.4 Then, pipette 200 μ L of Lysis Buffer (**R1**) and shake vigorously (10-15 seconds).
- 1.5 Continue on step **2. Precipitation**.

D. Viscous samples (Semen and Sputum)

- 1.1 Pipette 200 μ L of sample at room temperature in Collector Tube 1.5 mL (**R9**).
- Suggestion:** Very viscous samples should be diluted 1: 2 in PBS 1X.
- 1.2 Add 5 μ L of Proteinase K (**R6**)* to the Collector Tube 1.5 mL with sample and mix.
 - 1.3 Then, pipette 200 μ L of Lysis Buffer (**R1**) and shake vigorously (10-15 seconds).
 - 1.4 Continue on step **2. Precipitation**.

E. Tissue

Materials needed but not contained in the kit:

- Gral and pistil or disruptor / homogenizer of the tissues.
- Liquid nitrogen.
- Needle Coupled Sterile Plastic Syringe (0.7-0.9mm).

Optional: Reducing Agent - β -Mercaptoethanol or 10 - 20 mM Dithiothreitol (DTT)
Extraction of nucleic acids mainly from tissue samples can promote the release of large amounts of endonucleases (DNase and RNase) during the lysis step. The endonucleases cause reduction in the production and quality of nucleic acids, mainly of RNA molecules. For RNase inactivation a reducing agent (DTT or β -Mercaptoethanol) can be used thus preserving the integrity of the nucleic acids. The extraction can be performed without the presence of the reducing agent, observing that the maceration and lysis steps should be performed in an agile manner avoiding the heating and / or thawing of tissue samples previously stored at -20°C or -80°C.

Protocol

- 1.1 To weigh between 20-30 mg of fresh tissue to be extracted.

Note: It is recommended to store the tissues to be extracted properly (-20°C or -80°C) to maintain the integrity of nucleic acids. Avoid freezing and thawing of samples. If commercial stabilizing solution is used, storage should be performed according to the manufacturer's manual.

- 1.2 Carry out the tissue maceration in gral and pistil using liquid nitrogen.

Note: Tissue maceration or rupture can also be performed by equipment developed for this purpose, according to the manufacturer's recommendations.

* After use, store at - 20° C.

- 1.3 Transfer the macerated sample to a 1.5 mL microtube (not supplied).
- 1.4 Then pipette 200µL of Lysis Buffer **(R1)** and mix thoroughly for 1 minute.

Optional: A reducing agent such as β -Mercaptoethanol or 10-20 mM Dithiothreitol (DTT) may be added to the Lysis Buffer **(R1)** at a rate of 1% of the lysis volume, for example, 2µL of reducing agent for each 200µL of Lysis Buffer **(R1)**.

- 1.5 Homogenize the lysed tissue by passing at least 5 times through a needle (0.7 - 0.9 mm) coupled to a sterile plastic syringe.

Note: Homogenization can also be performed by equipment developed for this purpose, according to the manufacturer's recommendations.

- 1.6 Centrifuge the lysed tissue at 15,000 x g for 3 minutes.
- 1.7 Transfer the supernatant to a new 1.5 mL collection tube **(R9)**.
- 1.8 Continue with step 2. **Precipitation.**

2. Precipitation

- 2.1 Add 5.6µL of Carrier RNA **(R5)*** prepared and mix.
- 2.2 Incubate for 3 minutes at room temperature (18-25°C).
- 2.3 Add 200µL of Ethanol (96-100%) to each Collector Tube 1.5 mL with sample and mix by vortexing (10-15 seconds).
- 2.4 Incubate for 5 minutes at room temperature (18-25°C).

3. Binding

- 3.1 Identify the Column **(R7)** into Collector Tube according to the sample which is being purified.
- 3.2 Transfer the total of lysed sample volume (~ 610µL) into the respective Column and then centrifuge for 3 minutes at 4,000 x g. If the Column is not completely dried, it must be centrifuged again (approx. 15,000 x g).
- 3.3 After centrifugation, transfer the Column to new Collector Tube **(R8)** and discard the Collector Tube with the filtrate.

4. Washing

- 4.1 Add 400µL of Wash 1 **(R2)** to each Column containing sample.
- 4.2 Centrifuge for 1 minute at 11,000 x g.
- 4.3 Transfer the Column to new Collector Tube **(R8)** and discard the Collector Tube with the filtrate.
- 4.4 Add 400µL of Wash 2 **(R3)** prepared.
- 4.5 Centrifuge for 1 minute at 11,000 x g.
- 4.6 Transfer the Column to new Collector Tube **(R8)** and discard the Collector Tube with the filtrate.
- 4.7 Add 200µL of Wash 2 **(R3)** to Column.
- 4.8 Centrifuge for 5 minutes at $\geq 15,000 \times g$.

* After use, store at - 20° C.

5. Elution

5.1 Transfer the Column to a Collector Tube 1.5 mL (**R9**) and discard the Collector Tube with the filtrate.

5.2 Incubate the Column / Tube for 5 minutes at 56°C with the lids opened.

5.3 Add 30~60µL of RNase Free Water (**R4**) preheated to 70°C. Dispense directly in the center of the membrane.

5.4 Incubate at room temperature for 3 minutes.

5.5 Centrifuge for 3 minutes at $\geq 15,000 \times g$, and discard the Column.

5.6 Storing the DNA / RNA eluted at -20°C if not used immediately. For long time storage, it is recommended to stock DNA/RNA at -70°C

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

BIBLIOGRAPHIC REFERENCES

1- Ogura T, Tsuchiya A, Minas T, Mizuno S. Methods of high integrity RNA extraction from cell/agarose construct. BMC Res Notes. 2015; 8:644.

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for kit: 10269360296

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	CAUTION		DANGER

BIO **GENE**

Bioclin · **QUIBASA**

 Rua Teles de Menezes, 92, Belo Horizonte, MG, Brasil, CEP: 31565-130
Tel +55 31 3439 5454, www.bioclin.com.br
FARM. RESP. Silvio Wandalsen Arndt - CRF MG 7422
C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira