BIO GENE

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Review: February/2020

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	4
Equipaments and Operational Inputs	4
Transportation and Storage Conditions.....	4
Special Care	5
Samples	6
Procedure	6
A . Preparation of Working Solutions	6
B . DNA Extraction	6
1. Samples Lysis	6
2. Binding.....	7
3. Washing	7
4. Elution	8
Process Limitations	8
Bibliographic References	8
Customer Service	8
Universal Symbology	9

FUNCTION

Product developed for genomic DNA extraction and purification whole blood and cell culture.

The nucleic acids obtained can be used for different applications like PCR, Hibridization, Sequencing, etc.

PRINCIPLE OF ACTION

The **Bio Gene Extração de gDNA** kit is a product developed for extraction and purification of genomic DNA in biological samples.

The method used is the extraction by silica membrane. The process is carried out in 4 steps: 1) Cell Lysis: cell disruption for the release of nucleic acids; 2) Binding: selective binding of nucleic acid to silica membrane; 3) Washing: to remove residual impurities; 4) Elution: release the nucleic acid from silica membrane. At the end of the process we have concentrated nucleic acid with high purity.

Reagent	Presentation	
	50 preps	250 preps
R1	1 x 11 mL	1 x 55 mL
R2	1 x 22 mL	1 x 110 mL
R3	1 x 7 mL	1 x 33 mL
R4	1 x 13 mL	1 x 30 mL
R5	1 x 1.5 mL	4 x 1.75 mL
R6	1 x 50 unid	1 x 250 unid
R7	3 x 50 unid	3 x 250 unid
R8	2 x 50 unid	2 x 250 unid

REAGENTS

- R1. Lysis Buffer:** Buffer solution, Guanidine Hydrochloride and surfactant.
- R2. Wash 1:** Buffer solution, Guanidine Hydrochloride, alcohol and surfactant.
- R3. Wash 2:** Buffer solution and preservative.
- R4. DNase Free Water.**
- R5. Proteinase K:** Enzyme Proteinase K, surfactant and stabilizers.
- R6. Columns:** Polypropylene tube with a silica membrane.
- R7. Collectors Tubes (2 mL):** Polypropylene tube.
- R8. Collectors Tubes (1.5 mL):** Polypropylene tube.

EQUIPMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the table above.
- Usage instructions (manual).

Materials required, but not included in the kit:

- 1- 96-100% Ethanol
- 2- PBS 1X may be required for some samples
- 3- Micropipettes and sterile filter pipette tips (0.5-10 μ L, 10-100 μ L, 100-1000 μ L).
- 4- Microcentrifuge
- 5- Vortex
- 6- Heating block
- 7- Personal protective equipment (gloves, lab coat, glasses)

TRANSPORTATION AND STORAGE CONDITIONS

The kit may be transported between 15 and 30°C. Protect from light and avoid moisture.

The storage temperature is between 15 and 30°C.

After first use, is recommended to store the Proteinase K reagent (R5) -20°C.

SPECIAL CARE

1- For professional use only.

2- Strictly follow the methodology proposed to obtain accurate results.

3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose the waste in accordance with appropriate security procedures.

4- It is required skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.

5- It is necessary to have separate areas for the extraction / preparation of reactions and for the amplification / detection of products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.

6- Avoid repeated thawing and freezing of the reagents.

7- Do not use reagents after expiration date.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Kit for genomic DNA extraction from whole blood or cell culture. Samples must be collected and stored according to the laboratory's recommendations for molecular testing.

PROCEDURE

Attention:

The Lysis Buffer (**R1**) and Wash 1 (**R2**) reagents contains Guanidine Salt, and must be handled using suitable PPEs (coats, gloves and goggles). These reagents must not be disposed of on Acid Solution or on Sodium Hypochlorite due to the formation of reactive components.

A. Preparation of Working Solutions

Resuspend **R3** reagent as shown in the table below:

Reagent	Presentation	
	50 Preps	250 Preps
R3 - Wash 2 (concentrated)	Add 28 mL of Ethanol 96-100%	Add 132 mL of Ethanol 96-100%

Before starting the extraction:

- 1- Note if the reagent Wash 2 (**R3**) is prepared.
- 2- Set the heating block at 70°C.
- 3- Preheat the volume of the DNase free water (**R4**) reagent at 70°C to elution.

B. DNA extraction

1. Samples Lysis

1.1 Pipette 200µL of sample (whole blood or cell culture) at room temperature in Collector Tube 1.5 mL (**R8**).

- 1.2 Add 25µL of Proteinase K (**R5**)* to the Collector Tube 1.5 mL with sample and mix.
- 1.3 Then, pipette 200µL of Lysis Buffer (**R1**) and shake vigorously (10-15 seconds).
- 1.4 Incubate for 15 minutes at temperature 70°C.
- 1.5 Add 200µL of Ethanol (96-100%) to each Collector Tube 1.5 mL with sample and mix by vortexing (10-15 seconds).

2. Binding

- 2.1 Identify the Column (**R6**) into Collector Tube according to the sample which is being purified.
- 2.2 Transfer the total of lysed sample volume (~ 625µL) into the respective Column and then centrifuge for 1 minute at 11,000 x g. If the Column is not completely dried, it must be centrifuged again (approx. 15,000 x g).
- 2.3 After centrifugation, transfer the Column to new Collector Tube (**R7**) and discard the Collector Tube with the filtrate.

3. Washing

- 3.1 Add 400µL of Wash 1 (**R2**) to each Column containing sample.
- 3.2 Centrifuge for 1 minute at 11,000 x g.
- 3.3 Transfer the Column to new Collector Tube (**R7**) and discard the Collector Tube with the filtrate.
- 3.4 Add 600µL of Wash 2 (**R3**) prepared.
- 3.5 Centrifuge for 1 minute at 11,000 x g.
- 3.6 Transfer the Column to new Collector Tube (**R8**) and discard the Collector Tube with the filtrate.
- 3.7 Centrifuge for 1 minute at $\geq 15,000 \times g$.

* After use, store at - 20° C.

4. Elution

4.1 Transfer the Column to a Collector Tube 1.5 mL (**R8**) and discard the Collector Tube with the filtrate.

4.2 Add 60~100µL of DNase Free Water (**R4**) preheated to 70°C. Dispense directly in the center of the membrane.

4.3 Incubate at room temperature for 1 minute.

4.4 Centrifuge for 1 minute at $\geq 15,000 \times g$, and discard the Column.

4.5 Storing the DNA eluted at -20°C if not used immediately. For long time storage, it is recommended to stock DNA at -70°C

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

BIBLIOGRAPHIC REFERENCES

1.Vingataramin, L.; Frost, H. E. A single protocol for extraction of gDNA from bacteria and yeast. Vol 58. University of Sherbrooke, Canada, 2015.

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 031 5454


E-mail: sac@bioclin.com.br

ANVISA registration for kit: 10269360297

SIMBOLOGIA UNIVERSAL

CATALOG NUMBER

MANUFACTURED BY


BATCH CODE

CONTROL


DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

CAUTION

DANGER

BIO GENE

Bioclin ·· QUIBASA

 Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130
Tel +55 31 3439 5454 . Fax +55 31 3439 5455 . www.bioclin.com.br
FARM. RESP. Silvio Wandalsen Arndt - CRF MG 7422
C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira