HPV Alto Risco PCR

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Review: May/2015

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	4
Equipments and Operational Inputs	4
Transportation and Storage Conditions	4
Special Care	4
Samples	5
Process Description	6
A . DNAExtraction	6
B . Preparation of the Reagents	6
C . Dilution of Quantitative Standard	7
D . PCR Preparation	7
E . Thermocycler Settings for Real-Time PCR	8
F . Result Validation	9
G .Result Interpretation	10
Process Limitations	10
Product Performance / Quality Control	11
Comparison of Methods and Methodology Specificity	11
Repeatability	11
Reproducibility	11
Clinical Sensitivity	12
Analytical Sensitivity	12
Diagnostic Significance	12
Bibliographic References	13
Customer Service	13
Universal Symbology	14

FUNCTION

Quantitative detection of *Human Papillomavirus* type 16 (HPV16) DNA using real-time Polymerase Chain Reaction (PCR) technology. Designed for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

The **Bio Gene HPV Alto Risco PCR** test is a *in vitro* assay based on the quantitative detection of HPV16 DNA by real-time PCR.

The Real Time PCR method is used to amplify the pathogen DNA.

A thermocycler for Real Time PCR is used to amplify and detect the fluorescent probe. Software calculates the HPV16 DNA concentration, expressed in copies/ μ L, using the standard curve generated from the quantitative standard, contained in the kit.

Reagent	Presentation
	50 Tests
R1	1 x 55 μ L*
R2	1 x 525 μ L*
R3	1 x 525 μ L
R4	1 x 55 μ L*
R5	1 x 600 μ L*
R6	1 x 50 μ L
R7	1 x 500 μ L*
R8	2 x 1,5 mL
R9	1 x 1 mL

*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

REAGENTS

- R1. PCR Solution:** Primer, probe, TRIS-HCl.
R2. Mix Taq: DNA Polymerase, dNTPs, MgCl₂.
R3. Mix Buffer: TRIS-HCl.
R4. IC PCR Solution: Primer, probe, TRIS-HCl.
R5. Internal Control: Plasmid, TRIS-HCl.
R6. Negative Control: TRIS-HCl.
R7. Standard A (2×10^5 copies/ μ L): Plasmid, TRIS-HCl, EDTA.
R8. Diluent: TRIS-HCl, EDTA.
R9. Water: DNase/RNase free water.

EQUIPMENTS AND OPERATIONAL INPUTS**Materials in the kit:**

- Reagents described in the previous item
- Usage Instructions (manual)

Materials needed but not contained in the kit:

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- PCR reaction tubes or PCR reaction plate.
- 4- Powder free latex gloves or similar material.
- 5- Microcentrifuge (3,000 - 12,000 rpm).
- 6- Vortex.
- 7- Micropipettes and sterile filter pipette tips (0.5-10 μ L, 10-100 μ L, 100-1000 μ L).
- 8- Nucleic acid extraction kit.

TRANSPORTATION AND STORAGE CONDITIONS

The required storage temperature is -20°C (-10 to -30°C).

After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension. Avoid freezing and thawing.

The kit may be transported between 2-30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.

3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.

4- It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.

5- It is necessary to have separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplified products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.

6- All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.

7- Avoid repeated thawing and freezing of the reagents.

8- Store the DNA samples at -20°C in case they will not be used immediately.

9- Do not use reagents after expiration date.

10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

12- Do not use the product in case of damaged packaging.

13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

This kit should be used with DNA samples extracted from cervical or urethral swab. Other types of samples can be used according to medical recommendations or to laboratory.

Samples must be collected according to the laboratory's recommendations for molecular tests and resuspended in transport medium or saline. They must be transported and stored between 2 and 8°C for up to 3 days¹.

Use DNA samples suitable for PCR amplification with adequate purity and concentration. Avoid repeated freezing and thawing of the samples.

Note: Suggested DNA sample concentration **5ng/μL**.

PROCESS DESCRIPTION

A. DNA Extraction

The nucleic acids (DNA) must be extracted in accord to the procedure of the chosen extraction kit.

For control of the extraction process, the **Internal Control (R5)** must be prepared and added to the samples during the extraction as described below:

- 1- Add 4µL of the **Internal Control (R5)** to each tube containing the resuspended samples in extraction/lysis buffer.
- 2- Complete the extraction procedure according to the usage instructions of the extraction kit.

Note: Never add the Internal Control directly to pure biological sample. This may result in degradation of the sample.

B. Preparation of the Reagents

- 1- Before open, spin (*spin* pulse) each of the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)**, **Internal Control (R5)** and **Standard A (R7)**.
- 2- Resuspend the **Mix Taq (R2)** reagent in 525µL of reagent **Mix Buffer (R3)**.
- 3- Resuspend the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)** and **Internal Control (R5)** with **Water (R9)** reagent, in according to the following table:

Reagent	Volume of Water (R9)
PCR Solution (R1)	55 µL
IC PCR Solution (R4)	55 µL
*Internal Control (R5)	600 µL

*R5 and R7 reagents contain DNA template. They need to be handle (resuspended) in a separated area from the area designed to prepare other reagents in order to prevent contamination of other reagents.

C. Dilution of Quantitative Standard

- 1- Resuspend **Standard A (R7)** with 500µL of **Diluent (R8)**.
- 2- Select 3 microtubes (not supplied in the kit) suitable for serial dilution of **Standard A (R7)**.
- 3- Pipette 90µL of the **Diluent (R8)** in each microtube and name them as B, C and D respectively.
- 4- Pipette 10µL of **Standard A (R7)** in microtube B and mix.
- 5- Replace the tip and pipette 10µL of microtube B in microtube C and mix.
- 6- Replace the tip and pipette 10µL of microtube C in microtube D and mix.
- 7- At the end of dilution procedure, it will be available 4 standards with the following concentrations:

Standard A - 2×10^5 copies/µL

Standard B - 2×10^4 copies/µL

Standard C - 2×10^3 copies/µL

Standard D - 2×10^2 copies/µL

D. PCR Preparation

- 1- Prior beginning the assay, select sufficient number of microtubes to be used on the reaction, in accord to the number of samples and quantitative standards to be used.
- 2- Prepare the final PCR solution for each sample / standard / control to be tested.
- 3- Each reaction should contain:
 - 10µL **Mix Taq (R2)** provided in the kit).
 - 1µL **PCR Solution (R1)** provided in the kit).
 - 1µL **IC PCR Solution (R4)** provided in the kit).
 - 3µL **Water (R9)** provided in the kit).

The volumes indicated above are used for 1 reaction. To prepare more reactions, the volume of each reagent must be multiplied by the number of required reactions.

- 4- Add 5µL of extracted DNA sample (the suggested concentration is 5 ng/µL) or 5µL of Quantitative Standards or 5µL of **Negative Control (R6)** in previously determined tubes.
- 5- Mix well.

6- The total volume of the reaction is 20µL. Each PCR run should include the relevant controls: Negative Control, Internal Control and Quantitative Standards.

7- Transport the tubes to the thermocycler and go to the section E (Thermocycler Settings for Real-time PCR).

E. Thermocycler Settings for Real-Time PCR

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

1- Set the type of experiment:

Quantitative Test with Standard Curve.

2- Set the fluorescent probes:

Detector HPV 16	Detector Internal Control
"reporter 1" = FAM	"reporter 2" = VIC

Note: The Negative Control and Quantitative Standards do not present the Internal Control (reporter 2), because it is used for extraction and amplification control of the samples.

3- Set the quantitative standards:

Standard A - 2×10^5 copies/µL

Standard B - 2×10^4 copies/µL

Standard C - 2×10^3 copies/µL

Standard D - 2×10^2 copies/µL

4- Set the cycling conditions:

	Temperature	Time	Cycles
1	95°C	15 minutes	1
2	95°C	10 seconds	50
	60°C	60 seconds	

Set "Data Collection" as "stage 2, step 2 (60@0:60)".

F. Result Validation

1- Standard Curve

Standard Curve	Permitted Range	Amplification / Detection
Correlation Coeficiente (R^2)	$0,99 \leq R^2 \leq 1,00$	Valid

If the R^2 value does not fall within the limits of range, the result must be considered invalid and the test must be repeated.

2- Negative Control

CT Negative Control		Result	Detection
FAM	VIC		
Undetermined	Undetermined	Negative	Valid

3- Samples

HPV		Result	Detection
FAM	VIC		
Determined concentration	$25 \leq CT \leq 31$	Positive	Valid
	$CT < 25$ ou $CT > 31$	Positive	Invalid*
Undetermined concentration	$25 \leq CT \leq 31$	Negative	Valid
	$CT < 25$ ou $CT > 31$	Negative	Invalid*

*The CT values of Internal Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, the concentration of the samples and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

Example: In some cases, the samples with high numbers of copies of DNA/RNA can inhibit the amplification of the Internal Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

G. Result Interpretation

This kit is able to detect from 10 to 1.000.000 copies of DNA per reaction.

The thermocycler software calculates automatically the samples concentration.

Example: If the program shows a concentration of 2.00E+005, this means that the sample concentration is 2.0×10^5 copies/ μ L.

Sample Result in Copies/ μ L (FAM)	Copies per Reaction
$\geq 1 \times 10^6$	$> 1.000.000$
$2 \leq \text{Quantity} \leq 9 \times 10^5$	Quantity obtained
< 2	< 10

The non-detection of the pathogen DNA does not exclude the presence of infection when the pathogen title is below the detection limit of this kit. The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient.

The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The **Bio Gene HPV Alto Risco PCR** test was compared with another test for the quantification of the HPV16 DNA. 25 analyzes were performed and the results were evaluated. Comparing the two tests, the difference between the results obtained was ≤ 0.5 log. This shows that the test has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 positive samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Average Concentration (Log)	2,063	4,292
Standard Deviation (Log)	0,079	0,070
Coefficient of Variation (%)	3,849	1,642

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample with different concentrations, obtaining the following results:

	Day 1	Day 2	Day 3
Average Concentration (Log)	4,292	4,346	4,367
Standard Deviation (Log)	0,070	0,106	0,118
Coefficient of variation (%)	1,642	2,449	2,704

Clinical Sensitivity

The **Bio Gene HPV Alto Risco PCR** showed clinical sensitivity of 99,9% and a clinical specificity of 99,9%.

Method		Clinical Data		Total
Bio Gene HPV Alto Risco PCR	Result	Positive	Negative	
	Positive	25	0	25
	Negative	0	25	25
TOTAL		25	25	50

Analytical Sensitivity

It was found in the performance study that the analytical sensitivity of the kit is 0,23 log. This means the assay is capable to detect approximately 8 target molecules in 5µL DNA extraction product, used on the amplification reaction.

DIAGNOSTIC SIGNIFICANCE

Human Papillomavirus belongs to the Papillomaviridae family, infects lining tissue cells (skin and mucosa), causing lesions due to irregular cell growth. There are more than 150 different subtypes of HPV. Nearly 35 of them infect the female and male anogenital region and may cause from warts to dysplastic lesions from low and high grade. HPV transmission occurs by direct contact between the infected area in the case of genital HPV through sexual contact. HPV types can be classified as low, intermediate and high risk virus, according to the type of injury caused. The HPV genotypes 16, 18, 31, 33, 45, 52, 56, 58 are classified as high risk virus, since they are related to tumors. Because of that, early detection of these viruses is important to prevent the development of malignant lesions HPV-associated. Because of that, the quantitative detection of HPV by real time PCR is very important since it is a sensitive, precise and responsive method.

BIBLIOGRAPHIC REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Schmitz, M; Scheungraber, C; Herrmann, J. et al. Quantitative multiplex PCR assay for the detection of the seven clinically most relevant high-risk HPV types. J Clin Virol. 2009 Apr;44(4):302-7. Epub 2009 Feb 14. PubMed PMID: 19223232.
3. Seaman, W.T; Andrews, E; Couch, M. Detection and quantitation of HPV in genital and oral tissues and fluids by real time PCR. Virol J. 2010 Aug 19;7:194.
4. Seth, R; Rippin, J; Guo, L; Jenkins, D. Detection and quantitation of HPV gene expression using real-time PCR. Methods Mol Med. 2005;119:61-72.

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for Bio Gene HPV Alto Risco PCR: 10269360244

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE


CORROSIVE


IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHTDO NOT USE IF
PACKAGE IS
DAMAGED

BIO **GEÑE**

Bioclin · **QUIBASA**

Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130
Tel +55 31 3439 5454 . Fax +55 31 3439 5455 . www.bioclin.com.br
FARM. RESP. Sílvia Wandalsen Arndt - CRF MG 7422
C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira