Ureaplasma ur PCR

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Review: May/2015

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	4
Equipments and Operational Inputs	4
Transportation and Storage Conditions	4
Special Care	5
Samples	5
Process Description	6
A . DNAExtraction	6
B . Preparation of the reagents	6
C . Dilution of Quantitative Standard	7
D . PCR Preparation	7
E . Thermocycler settings for Real-Time PCR	8
F . Result Validation	9
G . Result Interpretation	10
Process Limitations	11
Product Performance / Quality Control	11
Comparison of Methods and Methodology Specificity	11
Repeatability	11
Reproducibility	12
Clinical Sensitivity	12
Analytical Sensitivity	12
Diagnostic Significance	12
Bibliographic References	13
Customer Service	13
Universal Symbology	14

FUNCTION

Quantitative detection of *Ureaplasma urealyticum* DNA using real-time Polymerase Chain Reaction (PCR) technology. Designed for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

The **Bio Gene UREAPLASMA UR PCR** test is a *in vitro* assay based on the quantitative detection of *U. urealyticum* DNA by real-time PCR.

The Real Time PCR method is used to amplify the pathogen DNA.

A thermocycler for Real Time PCR is use to amplify and detect the fluorescent probe. Software calculates the *U. urealyticum* DNA concentration, expressed in copies/μL, using the standard curve generated from the quantitative standard, contained in the kit.

Reagent	Presentation
	50 Tests
R1	1 x 55μL*
R2	1 x 525μL*
R3	1 x 525μL
R4	1 x 55μL*
R5	1 x 600μL*
R6	1 x 50μL
R7	1 x 500μL*
R8	2 x 1,5 mL
R9	1 x 1 mL

*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

REAGENTS

R1. PCR Solution: Primer, probe, TRIS-HCl.

R2. Mix Taq: DNA Polymerase, dNTPs, MgCl₂.

R3. Mix Buffer: TRIS-HCl.

R4. IC PCR Solution: Primer, probe, TRIS-HCl.

R5. Internal Control: Plasmid, TRIS-HCl.

R6. Negative Control: TRIS-HCl.

R7. Standard A (2×10^5 copies/ μ L): Plasmid, TRIS-HCl, EDTA.

R8. Diluent: TRIS-HCl, EDTA.

R9. Water: DNase/RNase free water.

EQUIPMENTS AND OPERATIONAL INPUTS**Materials in the kit:**

- Reagents described in the previous item
- Usage instructions (manual)

Materials needed but not contained in the kit:

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- PCR reaction tubes or PCR reaction plate.
- 4- Powder free latex gloves or similar material.
- 5- Microcentrifuge (3,000 - 12,000 rpm).
- 6- Vortex.
- 7- Micropipettes and sterile filter pipette tips (0.5-10 μ L, 10-100 μ L, 100-1000 μ L).
- 8- Nucleic acid extraction kit.

TRANSPORTATION AND STORAGE CONDITIONS

The required storage temperature is -20°C (-10 to -30°C).

After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension. Avoid freezing and thawing. The kit may be transported between 2-30°C. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain accurate results.

3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.

4- It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.

5- It is necessary to have separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplified products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.

6- All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.

7- Avoid repeated thawing and freezing of the reagents.

8- Store the DNA samples at -20°C in case they will not be used immediately.

9- Do not use reagents after expiration date.

10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

12- Do not use the product in case of damaged packaging.

13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

This kit should be used with DNA samples extracted from genital or urethral swab. Other types of samples can be used according to medical recommendations or to laboratory.

Samples must be collected according to the laboratory's recommendations for molecular tests and resuspended in transport medium or saline. They must be transported and stored between 2 and 8°C for up to 3 days¹.

Use DNA samples suitable for PCR amplification with adequate purity and concentration. Avoid repeated freezing and thawing of the samples.

Note: Suggested DNA sample concentration **5ng/μL**.

PROCESS DESCRIPTION

A. DNA Extraction

The nucleic acids (DNA) must be extracted in accord to the procedure of the chosen extraction kit.

For control of the extraction process, the **Internal Control (R5)** must be prepared and added to the samples during the extraction as described below:

- 1- Add 4μL of the **Internal Control (R5)** to each tube containing the resuspended samples in extraction/lysis buffer.
- 2- Complete the extraction procedure according to the usage instructions of the extraction kit.

Note: Never add the Internal Control directly to pure biological sample. This may result in degradation of the sample.

B. Preparation of the Reagents

- 1- Before open, spin (*spin* pulse) each of the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)**, **Internal Control (R5)** and **Standard A (R7)**.
- 2- Resuspend the **Mix Taq (R2)** reagent in 525μL of Reagent **Mix Buffer (R3)**.
- 3- Resuspend the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)** and **Internal Control (R5)** with **Water (R9)** reagent, in according to the following table:

Reagent	Volume of Water (R9)
PCR Solution (R1)	55 μL
IC PCR Solution (R4)	55 μL
*Internal Control (R5)	600 μL

*R5 and R7 reagents contain DNA template. They need to be handle (resuspended) in a separated area from the area designed to prepare other reagents in order to prevent contamination of other reagents.

C. Dilution of Quantitative Standard

- 1-** Resuspend **Standard A (R7)** with 500 µL of **Diluent (R8)**.
- 2-** Select 3 microtubes (not supplied in the kit) suitable for serial dilution of **Standard A (R7)**.
- 3-** Pipette 90µL of the **Diluent (R8)** in each microtube and name them as B, C and D respectively.
- 4-** Pipette 10µL of **Standard A (R7)** in microtube B and mix.
- 5-** Replace the tip and pipette 10µL of microtube B in microtube C and mix.
- 6-** Replace the tip and pipette 10µL of microtube C in microtube D and mix.
- 7-** At the end of dilution procedure, it will be available 4 standards with the following concentrations:

Standard A - 2×10^5 copies/µL

Standard B - 2×10^4 copies/µL

Standard C - 2×10^3 copies/µL

Standard D - 2×10^2 copies/µL

D. PCR Preparation

- 1-** Prior beginning the assay, select sufficient number of microtubes to be used on the reaction, in accord to the number of samples and quantitative standards to be used.
- 2-** Prepare the final PCR solution for each sample / standard / control to be tested.
- 3-** Each reaction should contain:
 - 10µL **Mix Taq (R2)** provided in the kit).
 - 1µL **PCR Solution (R1)** provided in the kit).
 - 1µL **IC PCR Solution (R4)** provided in the kit).
 - 3µL **Water (R9)** provided in the kit).

The volumes indicated above are used for 1 reaction. To prepare more reactions, the volume of each reagent must be multiplied by the number of required reactions.

4- Add 5µL of extracted DNA sample (the suggested concentration is 5 ng/µL) or 5µL of Quantitative Standards or 5µL of **Negative Control (R6)** in previously determined tubes.

5- Mix well.

6- The total volume of the reaction is 20µL. Each PCR run should include the relevant controls: Negative Control, Internal Control and Quantitative Standards.

7- Transport the tubes to the thermocycler and go to the section E (Thermocycler Settings for Real-Time PCR).

E. Thermocycler Settings for Real-Time PCR

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

1- Set the type of experiment:

Quantitative Test with Standard Curve.

2- Set the fluorescent probes:

Detector Ureaplasma ur	Detector Internal Control
“reporter 1” = FAM	“reporter 2” = VIC

Note: The Negative Control and Quantitative Standards do not present the Internal Control (reporter 2), because it is used for extraction and amplification control of the samples.

3- Set the quantitative standards:

Standard A - 2×10^5 copies/µL

Standard B - 2×10^4 copies/µL

Standard C - 2×10^3 copies/µL

Standard D - 2×10^2 copies/µL

4- Set the cycling conditions:

	Temperature	Time	Cycles
1	95°C	15 minutes	1
2	95°C	10 seconds	50
	60°C	60 seconds	

Set “Data Collection” as “stage 2, step 2 (60@0:60)”.

F. Result Validation**1- Standard Curve**

Standard Curve	Permitted Range	Amplification / Detection
Correlation Coeficiente (R^2)	$0,99 \leq R^2 \leq 1,00$	Valid

If the R^2 value does not fall within the limits of range, the result must be considered invalid and the test must be repeated.

2- Negative Control

CT Negative Control		Result	Detection
FAM	VIC		
Undetermined	Undetermined	Negative	Valid

3- Samples

Ureaplasma ur		Result	Detection
FAM	VIC		
Determined concentration	$25 \leq CT \leq 31$	Positive	Valid
	$CT < 25$ ou $CT > 31$	Positive	Invalid*
Undetermined concentration	$25 \leq CT \leq 31$	Negative	Valid
	$CT < 25$ ou $CT > 31$	Negative	Invalid*

*The CT values of Internal Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, the concentration of the samples and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

Example: In some cases, the samples with high numbers of copies of DNA / RNA can inhibit the amplification of the Internal Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

G. Result Interpretation

This kit is able to detect from 10 to 1.000.000 copies of DNA per reaction. The thermocycler software calculates automatically the samples concentration.

Example: If the program shows a concentration of 2.00E+005. This means that the sample concentration is 2.0×10^5 copies/ μ L.

Sample Result in Copies/ μ L (FAM)	Copies per Reaction
$\geq 1 \times 10^6$	$> 1.000.000$
$2 \leq \text{Quantity} \leq 9 \times 10^5$	Quantity obtained
< 2	< 10

The non-detection of the pathogen DNA does not exclude the presence of infection when the pathogen title is below the detection limit of this kit. The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient.

The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The **Bio Gene Ureaplasma ur PCR** test was compared with another test for the quantification of the *U. urealyticum* DNA. 25 analyzes were performed and the results were evaluated. Comparing the two tests, the difference between the results obtained was ≤ 0.5 log. This shows that the test has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 positive samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Average Concentration (Log)	4,337	2,402
Standard Deviation (Log)	0,093	0,082
Coefficient of Variation (%)	2,157	3,415

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample with different concentrations, obtaining the following results:

	Day 1	Day 2	Day 3
Average Concentration (Log)	2,402	2,401	2,395
Standard Deviation (Log)	0,082	0,084	0,090
Coefficient of variation (%)	3,415	3,514	3,781

Clinical Sensitivity

The **Bio Gene Ureaplasma ur PCR** showed clinical sensitivity of 99.9% and a clinical specificity of 99.9%.

Method		Clinical Data		Total
Bio Gene Ureaplasma ur PCR	Result	Positive	Negative	
	Positive	25	0	25
	Negative	0	25	25
TOTAL		25	25	50

Analytical Sensitivity

It was found in the performance study that the analytical sensitivity of the kit is 0,15 log. This means the assay is capable to detect approximately 7 target molecules in 5µL DNA extraction product, used on the amplification reaction.

DIAGNOSTIC SIGNIFICANCE

Ureaplasma urealyticum belongs to Mycoplasmataceae family and is the agent responsible for 20-30% of cases of non-gonococcal urethritis. In men it mainly causes urethritis and in women it can cause genitourinary tract infection, infertility and pregnancy complications such as premature birth. In neonates it can result in meningitis and pneumonia. The transmission occurs

through sexual contact or perinatally. Early diagnosis of *U. urealyticum* allow proper direction of the treatment, avoiding progression of the disease with consequences for the host. Because of that, the quantitative detection of *U. urealyticum* by real time PCR is very important since it is a sensitive, precise and responsive method.

BIBLIOGRAPHIC REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Mallard, K; Schopfer, K; Bodmer, T. Development of real-time PCR for the differential detection and quantification of *Ureaplasma urealyticum* and *Ureaplasma parvum*. J Microbiol Methods. 2005 Jan;60(1):13-9.
3. Xiao, L; Glass, J.I; Paralanov, V. et al. Detection and characterization of human *Ureaplasma* species and serovars by real-time PCR. J Clin Microbiol. 2010 Aug;48(8):2715-23.
4. Yi, J; Yoon, B.H; Kim, E.C. Detection and biovar discrimination of *Ureaplasma urealyticum* by real-time PCR. Mol Cell Probes. 2005 Aug;19(4):255-60.

CUSTOMER SERVICE




Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for Bio Gene *Ureaplasma ur* PCR: 10269360252

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED

BIO **GEÑE**

Bioclin · **QUIBASA**

Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130
Tel +55 31 3439 5454 . Fax +55 31 3439 5455 . www.bioclin.com.br
FARM. RESP. Sílvia Wandalsen Arndt - CRF MG 7422
C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira