Zika Vírus PCR

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Review: November/2020

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	4
Equipaments and Operational Inputs	4
Transportation and storage Conditions	4
Special Care	4
Samples	5
Process Description	6
A . RNAExtraction	6
B . Preparation of the reagents	6
C . Dilution of Quantitative Standard	7
D . PCR Preparation	7
E . Thermocycler settings for Real-Time PCR	8
F . Result Validation	10
G. Result Interpretation	11
Process Limitations	12
Product Performance / Quality Control	12
Comparison of Methods and Methodology Specificity	12
Repeatability	12
Reproducibility	13
Clinical Sensitivity	13
Analytical Sensitivity	14
Linearity	14
Diagnostic Significance	14
Bibliographic References	14
Customer Service	14
Universal Symbolology	15

FUNCTION

Quantitative detection of Zika virus RNA through reverse transcription followed by real-time polymerase chain reaction (RT-PCR). For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

The **Bio Gene Zika Virus PCR** test is a *in vitro* assay based on the quantitative detection of Zika virus RNA by real-time RT-PCR.

The Real-Time RT-PCR method is used to amplify the pathogen RNA. A thermocycler for Real-Time PCR is used to amplify and detect the fluorescent probe. Software calculates the Zika virus RNA concentration, expressed in copies/ μ L, using the standard curve generated from the quantitative standard, contained in the kit.

Reagent	Presentation	
	50 Tests	150 Tests
R1	1 x 55 μ L *	1 x 165 μ L *
R2	1 x 525 μ L *	3 x 525 μ L *
R3	1 x 525 μ L	3 x 525 μ L
R4	1 x 55 μ L *	1 x 165 μ L *
R5	1 x 600 μ L *	1 x 600 μ L *
R6	1 x 50 μ L	1 x 150 μ L
R7	1 x 500 μ L *	1 x 500 μ L *
R8	2 x 1,5 mL	2 x 1,5 mL
R9	1 x 1 mL	1 x 2 mL

*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

REAGENTS

R1. PCR Solution: Primer, probe, TRIS-HCl.

R2. Mix Taq: Polymerases, dNTPs, $MgCl_2$, stabilizers.

R3. Mix Buffer: TRIS-HCl.

R4. IC PCR Solution: Primer, probe, TRIS-HCl.

R5. Internal Control: Plasmid, TRIS-HCl.

R6. Negative Control: TRIS-HCl.

R7. Standard A (2×10^5 copies/ μ L): Plasmid, TRIS-HCl, EDTA.

R8. Diluent: TRIS-HCl, EDTA.

R9. Water: DNase/RNase free water.

EQUIPMENTS AND OPERATIONAL INPUTS**Materials in the kit:**

- Reagents described in the previous item.
- Usage instructions (manual).

Materials needed but not contained in the kit:

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- 1.5 ml microcentrifuge tubes.
- 4- PCR reaction tubes or plate.
- 5- Powder free latex gloves or similar material.
- 6- Microcentrifuge (3,000 - 12,000 rpm).
- 7- Vortex.
- 8- Micropipettes and sterile filter pipette tips (0.5-10 μ L, 10-100 μ L, 100-1000 μ L).
- 9- Nucleic acid extraction kit.

TRANSPORTATION AND STORAGE CONDITIONS

The required storage temperature is -20°C (-10 to -30°C). **After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension.** Avoid freezing and thawing.

The kit may be transported between 2 and 30°C for 7 days. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.

- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.
- 4- It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.
- 5- It is necessary to have separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplified products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.
- 6- All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.
- 7- Avoid repeated thawing and freezing of the reagents.
- 8- Store the RNA samples at -20°C in case they will not be used immediately.
- 9- Do not use reagents after expiration date.
- 10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 12- Do not use the product in case of damaged packaging.
- 13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

This kit should be used with RNA samples extracted from serum, plasma, urine, cerebrospinal fluid and amniotic fluid. Samples should be collected according to laboratory recommendations for molecular testing. They must be transported and stored as described below¹:

Serum and plasma: at -20°C for up to 3 days.

Urine: between 2 to 8°C for up to 3 days.

Cerebrospinal fluid and amniotic fluid: between 2 to 8°C for up to 1 day.

RNA (extracted): Use RNA samples suitable for RT-PCR amplification with adequate purity and concentration.

For long-term storage, it is recommended to freeze samples at -70°C. Avoid repeated freezing and thawing of the samples.

PROCESS DESCRIPTION

A. RNA Extraction

The nucleic acids (RNA) must be extracted in accord to the procedure of the chosen extraction kit.

For control of the extraction process, the **Internal Control (R5)** must be prepared (see item B.4) and added to the samples during the extraction as described below:

- 1- Add 4µL of the **Internal Control (R5)** to each tube containing the resuspended samples in extraction/lysis buffer.
- 2- Complete the extraction procedure according to the usage instructions of the extraction kit.

Note: Never add the Internal Control directly to pure biological sample. This may result in degradation of the sample.

B. Preparation of the Reagents

Note: **R5** and **R7** reagents contain RNA/DNA template. They need to be handle (resuspended) in a separated area.

- 1- Before open, spin (*spin* pulse) each of the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)**, **Internal Control (R5)** and **Standard A (R7)**.
- 2- Resuspend each vial of **Mix Taq (R2)*** reagent with 525µL of **Mix Buffer (R3)** reagent.

*For the presentation of 150 tests resuspend the **Mix Taq (R2)** vials as need of use.

*The **Mix Taq (R2)** reagent does not contain passive reference Dye (ROX).

- 3- Resuspend the reagents, **PCR Solution (R1)** and **IC PCR Solution (R4)** with **Water (R9)** reagent, in according to the following table:

Reagent	Presentation	
	50 Tests	150 Tests
PCR Solution (R1)	55 µL	165 µL
IC PCR Solution (R4)	55 µL	165 µL

4- Resuspend the reagents, **Internal Control (R5)** and **Standard A (R7)** with **Diluent (R8)** reagent, in according to the following table:

Reagent	Presentantion	
	50 tests	150 tests
Internal Control (R5)	600 µL	600 µL
Standard A (R7)	500 µL	500 µL

C. Dilution of Quantitative Standard

1- Select 3 microtubes (not supplied in the kit) suitable for serial dilution of resuspended **Standard A (R7)**.

2- Pipette 90µL of the **Diluent (R8)** in each microtube and name them as B, C and D respectively.

3- Pipette 10µL of **Standard A (R7)** in microtube B and mix.

4- Replace the tip and pipette 10µL of microtube B in microtube C and mix.

5- Replace the tip and pipette 10µL of microtube C in microtube D and mix.

6- At the end of dilution procedure, it will be available 4 standards with the following concentrations:

Standard A - 2×10^5 copies/µL

Standard B - 2×10^4 copies/µL

Standard C - 2×10^3 copies/µL

Standard D - 2×10^2 copies/µL

D. PCR Preparation

1- Prior beginning the assay, select sufficient number of microtubes/wells to be used on the reaction, in accord to the number of samples, Controls and Quantitative Standards to be used.

2- Prepare the final PCR solution in accord with reactions number required.

Reagents	1 Reaction	25 Reactions	50 Reactions	100 Reactions
Mix Taq (R2)	10 µL	250 µL	500 µL	1 mL
PCR Solution (R1)	1 µL	25 µL	50 µL	100 µL
CI PCR Solution (R4)	1 µL	25 µL	50 µL	100 µL
Water (R9)	3 µL	75 µL	150 µL	300 µL

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

3- Pipette 15µL of the final PCR solution in microtubes or wells previously determined.

4- Add 5µL of extracted RNA sample or 5µL of Quantitative Standards or 5µL of **Negative Control (R6)** in previously determined tubes.

5- Mix well.

6- The total volume of the reaction is 20µL. Each PCR run should include the relevant controls: Negative Control, Internal Control and Quantitative Standards.

7- Transport the tubes/plate to the thermocycler and follow the instructions in the section E (Thermocycler Settings for Real-Time PCR).

E. Thermocycler Settings for Real-Time PCR

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

1- Set the type of experiment:

Quantitative Test with Standard Curve.

2- Set the fluorescent probes:

Target	Detector	Quencher
Zika	FAM	NFQ-MGB
Internal Control	VIC	

Note: The Negative Control and Quantitative Standards do not present the Internal Control, because it is used for extraction and amplification control of the samples.

3- Set the quantitative standards:

Standard A - 2×10^5 copies/ μ L

Standard B - 2×10^4 copies/ μ L

Standard C - 2×10^3 copies/ μ L

Standard D - 2×10^2 copies/ μ L

4- Set the cycling conditions:

Stage	Temperature	Time	Cycles
1	55°C	10 minutes	1
2	95°C	2 minutes	1
3	95°C	10 seconds	50
	60°C	60 seconds	

Set "Data Collection" as "stage 3, step 2 (60@0:60)".

F. Result Validation

1- Standard Curve

Standard Curve	Permitted Range	Amplification / Detection
Correlation Coefficient (R^2)	$0,99 \leq R^2 \leq 1,00$	Valid

If the R^2 value does not fall within the limits of range, the result must be considered invalid and the test must be repeated.

2- Negative Control

CT Negative Control		Result	Amplification/ Detection
FAM	VIC		
Undetermined	Undetermined	Negative	Valid

3- Samples

Zika		Result	Detection
FAM	VIC		
Determined concentration	$25 \leq CT \leq 35$	Positive	Valid
	$CT < 25$ ou $CT > 35$	Positive	Invalid*
Undetermined concentration	$25 \leq CT \leq 35$	Negative	Valid
	$CT < 25$ ou $CT > 35$	Negative	Invalid*

*The CT values of Internal Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, the concentration of the samples and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

Example: In some cases, the samples with high numbers of copies of DNA / RNA can inhibit the amplification of the Internal Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

G. Result Interpretation

This kit is able to detect from 10 to 1.000.000 copies of RNA per reaction.

The thermocycler software calculates automatically the samples concentration.

Example: If the program shows a concentration of 2.00E+005, this means that the sample concentration is 2.0×10^5 copies/ μ L.

Sample Result in Copies/ μ L (FAM)	Copies per Reaction
$\geq 1 \times 10^6$	> 1.000.000
$2 \leq \text{Quantity} \leq 9 \times 10^5$	Quantity obtained
< 2	< 10

The non-detection of the pathogen RNA does not exclude the presence of infection when the pathogen title is below the detection limit of this kit.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient.

The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

1- Conversion calculation for copies/mL

$$\text{Result in copies/mL} = \frac{\text{Copies/}\mu\text{L} \times \text{Elution volume } (\mu\text{L})^{**}}{\text{Extracted sample volume (mL)}^{**}}$$

*Result of quantification of the sample provided by the equipment in copies/ μL .

**According to the protocol of the extraction kit used.

Example: Elution volume: 30 μL

Sample volume: 200 μL (0.2 mL)

Quantification of the sample: 5000 copies/ μL

$$\text{Result in copies/mL} = \frac{(5000 \text{ copies/}\mu\text{L} \times 30\mu\text{L})}{0.2\text{mL}}$$

$$\text{Result in copies/mL} = 750.000 \text{ copies/mL}$$

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

The heparin use in the samples collection for molecular tests is not recommended, due to the PCR inhibition which can generate false negative results.

The samples with presence of PCR inhibitors should be diluted to reduce their effect on the reaction.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The **Bio Gene Zika Virus PCR** test was compared with another protocol for the quantification of the Zika virus RNA. 50 analyzes were performed and the results were evaluated. Comparing the two tests, the difference between the results obtained was ≤ 0.5 log. This shows that the test has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 positive samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Average Concentration (Log)	4.298	2.483
Standard Deviation (Log)	0.025	0.032
Coefficient of variation (%)	0.595	1.321

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample with different concentrations, obtaining the following results:

	Day 1	Day 2	Day 3
Average Concentration (Log)	4.298	4.240	4.249
Standard Deviation (Log)	0.025	0.063	0.029
Coefficient of variation (%)	0.595	1.507	0.687

Clinical Sensitivity

The **Bio Gene Zika Virus PCR** showed clinical sensitivity of 99.9% and a clinical specificity of 99.9%.

Method		Clinical Data		Total
Bio Gene Zika Virus PCR	Result	Positive	Negative	
	Positive	25	0	25
	Negative	0	25	25
TOTAL		25	25	50

Analytical Sensitivity

It was found in the performance study that the analytical sensitivity of the kit is 0.0937 copies/ μ L (- 1.028 log). This means the assay is capable to detect approximately 0.468 target molecules in 5 μ L RNA extraction product, used on the amplification reaction. The sensitivity was 28 copies/mL. *

* The sensitivity may vary according to the used extraction kit.

Linearity

The linearity found in the performance studies was 1.56×10^5 copies/ μ L (5.193 log).

DIAGNOSTIC SIGNIFICANCE

The Zika virus belongs to the *Flaviviridae* family. The infection is transmitted through the bite of the *Aedes aegypti* mosquito infected with the virus. Symptoms appear after a short incubation period, usually, it is self-limiting disease, benign, characterized by low-grade fever, headache, conjunctivitis, malaise and rash. However, there is a correlation that Zika virus infection can cause Guillain Barré syndrome and during pregnancy, especially in the first trimester, can lead to serious complications such as microcephaly, which results in severe sequelae in the newborn.

Therefore, the quantitative detection of Zika virus by Real-Time PCR, mainly for diagnosis in pregnant women, is very important due to its accuracy and agility.

BIBLIOGRAPHIC REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Faye O, Faye O, Dupressoir A, Weidmann M, Ndiaye M, Sall AA. One-step RT-PCR for detection of Zika virus. J Clin Virol. 2008;43:96–101. 10.1016/j.jcv.2008.05.005.
3. World Health Organization (WHO). Weekly epidemiological record. 2015 Nov; vol 90; 45 (609-616).

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for Bio Gene Zika Virus PCR: 10269360300

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED

Bioclin · QUIBASA

Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130

Tel +55 31 3439 5454 . www.bioclin.com.br

FARM. RESP. Sílvia Wandalsen Arndt - CRF MG 7422

C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira