BIOLÁTEX ASO

REF K004

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for the qualitative and semiquantitative determination of Antistreptolysin O by the agglutination of latex particles without dilution of the sample. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Latex

This method is based on a particle agglutination reaction of latex coated with Streptolysin O, especially treated to prevent non specific agglutination. Agglutination is visible in samples with concentrations equal to or greater than 200 IU/mL, according to International Standards established by the references of the WHO.

REAGENTS

- 1 - Number 1 - ASO Latex** - Store between 2 and 8°C. **Do not freeze.** Contains: Latex Particles in suspension sensitized with Streptolysin O. Homogenize before using it.
- 2 - Number 2 - Positive Control** - Store between 2 and 8°C. Contains: Latex Agglutinating Solution.
- 3 - Number 3 - Negative Control** - Store between 2 and 8°C. Contains: 0,9% Saline Solution and preservative.

PRESENTATION

Reagent	Presentation	
	1 (Economic Packing)	2 (Normal Packing)
Reagent Nº 1	2 mL	2 mL
Reagent Nº 2	--	1 mL
Reagent Nº 3	--	1 mL

The normal package is accompanied by test cards.

EQUIPMENTS AND OPERATIONAL INPUTS

Cards, slides or dark bottom plates, spatulas, watch or stopwatch. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- 2-** Strictly follow the methodology proposed to obtain exact results.
- 3-** Water used in material cleaning must be be recent and free of contaminants.

- 4-** Do not freeze reagents.
- 5-** Always use reagents from same lot.
- 6-** Do not use lipemic serum. Do not use plasma.
- 7-** For cards reuse, wash immediately after use in distilled or deionized water until all residue is removed. If the washing can not be immediate, wash the card with neutral detergent and rinse thoroughly in distilled or deionized water until all the residue is removed.
- 8-** Handle with care Reagents Nº 3 that contain Sodium Azide.
- 9-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 10-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 11-** Do not use the product in case of damaged packaging.
- 12-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum, with no previous dilution. The analyte is stable for 2 days if kept between 2 and 8°C and 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing.

PROCESS DESCRIPTION
TECHNIQUE

QUALITATIVE TEST

In each circle of the card place the following:

	Circle Nº 1	Circle Nº 2	Circle Nº 3
Negative Control	20 µL	--	--
Positive Control	--	20 µL	--
Sample	--	--	20 µL
Reagent Nº 1*	20 µL	20 µL	20 µL

*Previously homogenized.

Homogenize with the aid of a spatula using the entire length of each circle of the card. Soon after, agitate the card with circular motions for 2 minutes. Make immediately reading with artificial light.

A clear agglutination indicates the presence of Antistreptolysin O at a concentration equal to or greater than 200 IU/mL. In this case, perform the semiquantitative test.

SEMIQUANTITATIVE TEST

- 1- Perform sample dilutions with saline, starting from initial sample (1:2, 1:4, 1:8, 1:16, 1:32, etc).
- 2- Follow the process as described in qualitative test for each of the dilutions.
- It will be considered as title the biggest serum dilution that presents agglutination.

RESULTS

Positive: Clear presence of agglutination.

Negative: Absence of agglutination (homogenized suspension).

CALCULATIONS

Sample	Concentration (UI/mL)
No dilution	200
1/2	400
1/4	800
1/8	1600
1/16	3200

Results can be expressed in titles or in IU/mL.

IU/mL = 200 x title from last dilution (dilution number)

Negative Tests: Express the results with negative or minor than 200 IU/mL.

PROCEDURE LIMITATIONS

Do not use plasma, hemolysate serum or limpemic, as they may produce nonspecific agglutination.

By correlating methods for determining Anti-streptolysin O, to verify the sensitivity of the reagents. The obtained results should only be compared when expressed in IU/mL.

The detergent residue on the card can interfere with the test and generate a false result.

INTERFERENT

No interference was observed for Bilirubin up to 20 mg/dL, Hemoglobin up to 10 g/L and Lipids up to 10 g/L. Rheumatoid Factor up to 300 IU/mL does not interfere with the results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's traceability was determined through the reference material from WHO.

REFERENCE VALUES

Up to 200 IU/mL

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Biolatex ASO kit was compared with another method to measure the presence of Antistreptolysin O. According to results from 40 clinical samples, all methods showed excellent correlation. With these results we can conclude the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Samples	N° of Repetitions	Expected Results	Results Obtained
01	20	Negative	Negative
02	20	Positive (1:8)	Positive (1:8)
03	20	Positive (1:64)	Positive (1:64)

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Samples	N° of Repetitions	Expected Results	Results Obtained
01	20	Negative	Negative
02	20	Positive (1:8)	Positive (1:8)
03	20	Positive (1:64)	Positive (1:64)

Analytical Sensitivity

The analytical sensitivity study of the Kit Biolatex ASO was performed by dilution of a positive sample with a known concentration of 1000 IU/mL for Anti-Streptolysin O. The analytical sensitivity was 200 IU/mL.

High Dose Hook Effect

Hook effect wasn't observed with a high concentration of Antistreptolysin O up to 1500 IU/mL.

DIAGNOSTIC SIGNIFICANCE

The Streptolysin O is a toxic immunogenic exoenzymes produced by many β -hemolytic streptococci group A. Taking into account the large amount of exoenzymes released *in vivo*, the determination the antibody

response against the presence of Streptolysin O has become a routine procedure for diagnosis and treatment of rheumatic fever, acute glomerulonephritis, scarlet fever, tonsillitis, erysipelas, puerperal sepsis and other infections streptococcus group A.

NUMBER OF TESTS

Presentation 1 (Economic Packing) 100 tests/20 μ L of Sample / 20 μ L of Reagent

Presentation 2 (Normal Packing) 100 tests / 20 μ L of Sample / 20 μ L of Reagent

BIBLIOGRAPHIC REFERENCES

- 1 - PLOTZ and SINGER, A. J. Med. 22, 1979.
- 2 - ADAMS, L. E.; HESS, E., J. Amer. Technol, 48, 1978.
- 3 - NORMAUSELL, D., Immunochemistry, 9, 1972.
- 4 - DITO, W., Am. Soc. Clin. Pat., 69, 1976.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Biolatex ASO kit: 10269360137

Review: January/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED