BIOLÁTEX FR

REF K043

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination qualitative and semiquantitative of Rheumatoid Factors (RF), by agglutination of latex particles, without previous dilution of the sample. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Latex

This method is based on the particle agglutination reaction of latex particles sensitized with IgG, especially treated to prevent nonspecific agglutination. Agglutination is visible in samples with RF concentrations equal to or greater than 8 IU/mL, according to the references established by the WHO International Standards.

REAGENTS

Number 1 - RF Latex - Store between 2 and 8°C. **Do not freeze.** Contains: Latex particles in suspension sensitized with human Gamma Globulin. Homogenize before using it.

Number 2 - Positive Control - Store between 2 and 8°C. Contains: Latex Agglutinating Solution, 0,9% Saline Solution, stabilizer and preservative. **Potentially infectious.**

Number 3 - Negative Control - Store between 2 and 8°C. Contains: 0,9% Saline Solution and preservative.

PRESENTATION

Reagent	Presentation	
	1 (Economic Packing)	2 (Normal Packing)
Nº 1	2 mL	2 mL
Nº 2	--	1 mL
Nº 3	--	1 mL

The normal package is accompanied by test cards.

EQUIPMENTS AND OPERATIONAL INPUTS

Cards, slides or dark bottom plates, spatulas, watches or stopwatches. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Do not freeze reagents.

5- Always use reagents from same lot.

6- Do not use lipemic serum. Do not use plasma.

7- For cards reuse, wash immediately after use in distilled or deionized water until all residue is removed. If the washing can not be immediate, wash the card with neutral detergent and rinse thoroughly in distilled or deionized water until all the residue is removed.

8- Handle with care. Reagents N° 2 and 3, they contain Sodium Azide.

9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

10- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

11- Do not use the product in case of damaged packaging.
12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum, with no previous dilution. The analyte is stable for 02 days if kept between 2 and 8°C and 3 months at -20°C⁵.

Samples with presence of fibrin should be centrifuged before testing.

PROCESS DESCRIPTION TECHNIQUE

QUALITATIVE TEST

In each circle of the card place the following:

	Circle Nº1	Circle Nº2	Circle Nº3
Negative Control	20 µL	--	--
Positive Control	--	20 µL	--
Sample	--	--	20 µL
Reagent Nº 1*	20 µL	20 µL	20 µL

*Previously homogenized.

Homogenize with the aid of a spatula using the entire length of each circle of the card. Soon after, agitate the card with circular motions for 2 minutes. Make immediately reading with artificial light.

A clear agglutination indicates the presence of Rheumatoid Factor at a concentration equal to or greater than 8 IU/mL. In this case, perform the semiquantitative test.

SEMIQUANTITATIVE TEST

1- Perform sample dilutions with saline, starting from initial sample (1:2, 1:4, 1:8, 1:16, 1:32, etc).

2- Follow the process as described in qualitative test for each of the dilutions.

It will be considered as title the biggest serum dilution that presents agglutination.

RESULTS

Positive: Clear presence of agglutination

Negative: Absence of agglutination (homogenized suspension)

CALCULATIONS

Results are expressed in IU/mL.

Sample	Concentration (UI/mL)
No dilution	8
1 / 2	16
1 / 4	32
1 / 8	64
1 / 16	128

Results can be expressed in titles or in IU/mL.

IU/mL = 8 x title from last dilution (N° of dilution)

Negative Tests: express the results as negative or minor than 8 IU/mL.

PROCEDURE LIMITATIONS

Do not use plasma, hemolysed or lipemic sera, as they may produce nonspecific agglutination.

By correlating methods for determining Rheumatoid Factors, to verify the sensitivity of the reagents. The obtained results should only be compared when expressed in UI/mL.

The detergent residue on the card can interfere with the test and generate a false result.

INTERFERENT

No interference was observed for Bilirubin up to 20 mg/dL, Hemoglobin up to 10 g/L and Lipids up to 10 g/L.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's traceability was determined through the reference material NIBSC 64/002 (Rheumatoid Arthritis Serum).

REFERENCE VALUES

Up to 8 IU/mL

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE**QUALITY CONTROL****Accuracy****COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY**

The Biolatex RF kit was compared with another method for confirmation of inflammatory process. According to results from 104 clinical samples, all methods showed excellent correlation. With these results we can conclude the kit has good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Results Obtained
01	40	Positive 1/2	Positive 1/2
02	40	Positive 1/32	Positive 1/32
03	40	Negative	Negative

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Results Obtained
01	40	Positive 1/2	Positive 1/2
02	40	Positive 1/32	Positive 1/32
03	40	Negative	Negative

Analytical Sensitivity

The analytical sensitivity study of the Kit Biolatex RF was performed by dilution of a positive sample with a known concentration of 40 IU/mL for Anti-Streptolysin O.

The analytical sensitivity was 8 IU/mL.

High Dose Hook Effect

Hook effect wasn't observed with a high concentration of Rheumatoid Factor up to 1500 IU/mL.

DIAGNOSTIC SPECIFICITY

Rheumatoid Factor refers itself to a group of macro globulins (anti globulins) that react with the Fc fragment from human gamma globulin. In rheumatoid arthritis the findings of Rheumatoid Factor represent the most important sera data, present in 75% patients. Its

plasmatic level may be significantly heighten in old age, conjunctive tissue diseases, chronic liver disease, syphilis, tuberculosis, leprosy, bacterial endocarditis, mononucleosis, sarcoidosis, leishmaniasis, rubella, neoplasms, parasitic infections, blood transfusions, renal transplantation, Sjogren's syndrome.

NUMBER OF TESTS

Presentation 1 (Economic Packing) / 100 tests/20 µL of Sample / 20 µL of Reagent

Presentation 2 (Normal Packing) / 100 tests / 20 µL of Sample / 20 µL of Reagent

BIBLIOGRAPHIC REFERENCES

- 1 - PLOTZ and SINGER, A. J. Med. 22, 1979.
- 2 - ADAMS, L. E.; HESS, E., J. Amer. Technol, 48, 1978.
- 3 - NORMAUSELL, D., Immunochemistry, 9, 1972.
- 4 - DITO, W., Am. Soc. Clin. Pat., 69, 1976.
- 5 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium


CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Biolatex RF kit: 10269360138

Review: January/2020

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED