CÁLCULO RENAL

REF K008

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of the usual components of Kidney Stone. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Bioclin

The physical-chemical analysis of Kidney Stone is of great clinical importance in preventive orientation from new formations. The Bioclin Kidney Stone kit allows the fast, simple and secure identification of Carbonate, Oxalate, Ammonium, Phosphate, Calcium, Magnesium, Urate and Cystine.

Carbonate: Formation of carbonic gas in acidic medium.

Oxalate: Selective precipitation in slightly acidic medium.

Phosphate: Formation of Molybdenum Blue in the presence of Ammonium Molybdate.

Calcium: Selective precipitation in neutral medium.

Magnesium: Formation of violaceous, by complexation with specific ligand.

Ammonium: Formation yellowish with Iodine-Iodide Potassium and Mercury.

Urate: Formation of Blue Tungsten through reduction of Sodium Phosphotungstate.

Cystine: Formation of red color in the presence of Sodium Cianetonitroprussiato.

REAGENTS

All the reagents are ready to be used. Included the set 01 bottle of Reagent Nº 15 - Standard, containing all usual components from Kidney Stones (Oxalate, Calcium, Carbonate, Magnesium, Ammonium, Cystine and Uric Acid), and 01 bottle of Reagent Nº 8 - Negative Control, free of analytes. For comparison purposes, these reagents should be treated as a sample. **All Standard (R15) tests are expected to be positive and all Negative Control (R8) tests result in negative.**

PRESENTATION

Kit composed with 13 Reagents, 1 Standard and 1 Negative Control, enough to perform 10 complete analysis. The kit includes an accessory (dropper) to be used in the technique for measuring reagent Nº 1 and distilled or deionized water.

Reagent	Volume
Nº 1 - Acid Reagent	10,0 mL
Nº 2 - Oxalate	5,0 mL
Nº 3 - Phosphorus I	3,0 mL
Nº 4 - Phosphorus II	3,0 mL
Nº 5 - Alkaline Reagent	2 x 8,0 mL
Nº 6 - Calcium	7,0 mL
Nº 7 - Magnesium	10,0 mL
Nº 8 - Negative Control	580 mg
Nº 9 - Ammonium	5,0 mL
Nº 10 - Uric Acid I	5,0 mL
Nº 11 - Uric Acid II	5,0 mL
Nº 12 - Cystine I	2,0 mL
Nº 13 - Cystine II	2,0 mL
Nº 14 - Cystine III	2,0 mL
Nº 15 - Standard	580 mg

*Includes accessory (dropper)

EQUIPMENTS AND OPERATIONAL INPUTS

Porcelain crucible, Bunsen burner, centrifuge, water bath, watch or stopwatch, pipettes, test tubes. Can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperatures should be between 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- The reagents Nº 1 and Nº 4 contain Hydrochloric Acid and Sulfuric Acid respectively, which are strong acids. Handle them with care.
- 5- The reagents Nº 5 contain Sodium Hydroxide and reagents Nº 9 and Nº 10 contain Potassium Hydroxide, which are strong bases. Handle them with care.
- 6- The reagent Nº 7 contains Sodium Azide, irritant to skin and mucous membranes. Handle with care.
- 7- The reagent Nº 13 contains Cyanide which is toxic. Handle with care.
- 8- It is convenient to proceed with the complete analysis of Standard that follows the kit, with the purpose of familiarizing yourself with the test.

9- If the stones has sent less than 40 mg or the bottom of the Khan tube is not filled in, proceed with the analysis reducing it to half the volume of all reagents utilized in step B.1.

10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

12- Do not use the product in case of damaged packaging.

13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Kidney Stones

PROCESS DESCRIPTION

TECHNIQUE

A) PHYSICAL ANALYSIS

Describe the dimension, weight (if possible), form, color, surface (smooth or irregular) of Kidney Stone. During its pulverization, observe the consistence (petrified, soft).

B) CHEMICAL ANALYSIS

1- Preparing the Sample to be analyzed

ATTENTION: The sample should only be prepared on the day of the tests.

Pulverize the stones. Weigh* approximately 40-50 mg of the homogeneous powder.

Transfer the weighed amount to a test tube, and add 10 drops of Reagent Nº 5 and 10 drops of distilled or deionized water.

Heat in water bath at 56°C for 5 minutes, shaking the tube 2 to 3 times during this period. Centrifuge at 3000 rpm for 3 minutes. Transfer the en tire supernatant to another tube and mark it with the letter S (for use in analysis of the supernatant) and mark the tube containing the precipitate with the letter P (for solubilized use). If the material completely solubilizes, separate 1,0 mL as supernatant and the remaindr as precipitate.

* If the laboratory can not weigh the amount of pulverized the stone required to perform the tests, use a sufficient amount to fill the bottom of the Khan tube (12 x 75 mm).

2- Solubilization of the Precipitate and analysis of Carbonate

CARBONATE: To the tube marked with the letter P, add 10 drops of Reagent Nº 1, observing simultaneously if there is effervescence. If so, the test is positive for Carbonate. Add then 10 drops of distilled water and mix. Heat to the first

sign of boiling (In flame directly, hot plate or any other heating device). Allow to cool and, in case of high turbidity, centrifuge again at 3000 rpm por 3 minutes (to avoid influence on results of Calcium and Oxalate). Use this solution for analysis of Oxalate, Calcium and Magnesium.

3- Analysis of the Solution of the Precipitate

OXALATE: Take 0,1 mL from tube P and add 5 drops of Reagent Nº 2. Wait for 30 seconds. The formation of turbidity or white precipitate indicates the presence of Oxalate.

CALCIUM: Take 0,1 mL from tube P and add 7 drops of Reagent Nº 6. Wait for 30 seconds. The formation of turbidity or white precipitate indicates the presence of Calcium.

MAGNESIUM: Transfer 0,05 mL from tube P to a Erlenmeyer and add 20 mL of distilled or deionized water. Add to this solution 1 drop of Reagent Nº 5. Mix. This will be the diluted sample. Use a test tube of and add 10 drops of Reagent Nº 7. Add 0,05 mL of the diluted sample and mix. The appearance of violaceous or intense pink color indicates the presence of Magnesium.

4- Analysis of Supernatant

URATE: Transfer 0,1 mL from tube S to a test tube. Add 5 drops of Reagent Nº 10 and 5 drops of Reagent Nº 11. Wait 5 minutes. If there is an intense blue color and it persists during the 5 minutes, Urate is present. The appearance of weak blue color not persistent during the 5 minutes should be interpreted as negative.

CYSTINE: Transfer 0,1 mL from tube S to a test tube. Add 1 drop of Reagent Nº 12 and 1 drop of reagent Nº 13. Wait 5 minutes. Add 2 drops of Reagent Nº 14. The appearance of intense red color indicates the presence of Cystine.

Note: The color disappears quickly formed. Caution to manipulate Reagent Nº 13. Contains Cyanide.

AMMONIUM: Transfer 0,1 mL from tube S to a test tube. Add 10 drops of distilled or deionized water. Mix. Add 5 drops of Reagent Nº 9. The formation of precipitate orange - yellow indicates the presence of Ammonium.

Note: Precipitation of the Reagent N° 9 may occur, but it doesn't interfere on the test results. The color formed can acquire a dark color if there is a high concentration of Ammonium in the sample. The formation of a white precipitate should be interpreted as a negative result.

PHOSPHATE: Transfer 0,1 mL from tube S to a test tube. Add 1 mL of deionized or distilled water and 01 drop of Reagent N° 1. Mix. Heat the tube to the first boiling point (direct flame, hot plate or any other heating device) - *Be careful to do not project the liquid.* Optionally, the tube may be heated in boiling water for 5 minutes. Next, add 2 drops of Reagent N° 3. Mix. Add 2 drops of Reagent N° 4. Leave resting in 2 minutes and add 2 drops of Reagent N° 5. Wait for 5 minutes to read the color verification. The appearance of blue color indicates the presence of Phosphate.

Note: Reagent N° 3 may cause a yellowish coloration. With this, the reaction in the presence of phosphate will no longer be blue and will become greenish. However, this change does not affect the quality and reliability of the result.

RESULTS

A model result to be provided is described as follows:

Material: Kidney Stone

1. Physical Analysis

- 1.1. Weight: 50 mg
- 1.2. Dimensions: 5 x 3 mm
- 1.3. Form: Oval
- 1.4. Color: Dark gray
- 1.5. Surface: Irregular
- 1.6. Consistency: Stony

2. Chemical Analysis

- 2.1. Carbonate: Negative
- 2.2. Oxalate: Positive
- 2.3. Phosphate: Positive
- 2.4. Calcium: Positive
- 2.5. Magnesium: Negative
- 2.6. Ammonium: Negative
- 2.7. Urata: Negative
- 2.8. Cystine: Negative

3. Conclusion

The examined material contain Oxalate and Calcium Phosphate.

LABORATORY INTERPRETATION

It is important for the Technical Manager to understand that:

- 1- The calculi can be pure or mixed. Mixed calculi can be formed by Calcium Oxalate, Hydroxyapatite [Ca₁₀(PO₄)₆(OH)₂] and struvite (NH₄MgO₄·6H₂O). The presence of Uric Acid and / or Cystine in mixed stones is less common but may occur.
- 2- In tests where an anion is detected, there will necessarily be the result of a corresponding cation, and vice versa.

- 3- Positive results for Oxalate will be followed by positive results for Calcium, but not always positive results for Calcium are followed by Oxalate. In this case, observe the results for Phosphate and Carbonate.
 - 4- Ammonium and Magnesium usually have coincident results (or are both negative or positive). If this does not occur, it is necessary to repeat the analysis. When positive, they will be followed by Phosphate.
 - 5- Positive results for Phosphate, Ammonium, Magnesium and Calcium and negative for Oxalate are due to the concomitant presence of struvite and apatite, and absence of Calcium Oxalate.
 - 6- Carbonate is generally formed in the calculi during the formation of Phosphate through the absorption of carbon dioxide, so the presence of both anions and Calcium cation is very common.
 - 7- Cystine is usually isolated from the others. In the case of a positive result, the others probably will not be, unless the patient has other conditions for this: diet rich in purines (Uric Acid), excess or malabsorption of Ascorbic Acid (Oxalate) or other pathologies.
- The BIOCLIN KIDNEY STONES KIT ALLOWS THE SECURE IDENTIFICATION OF THE CATIONS AND ANIONS PRESENT IN THE CALCULI. HOWEVER, THE FINAL PATIENT EVALUATION AND ITS DIAGNOSIS MUST BE GIVEN BY THE SPECIALIST PHYSICIAN, BY ANAMNESIS AND OTHER AUXILIARY EXAMS.**

INTERFERENT

The Phosphate Test is sensitive to the presence of orthophosphates and pyrophosphates. Sulfides, sulfites, hexacyanoferrates (II) and tartrates can seriously affect the reaction, generating false-negative. The presence of arsenates may cause a false positive. For the Calcium and Oxalate tests, false-positive results can be obtained in the presence of silver salts. However, none of these interferons are involved in the constitution of common renal stone samples. The presence of cystine stones can cause a slight change in staining in the uric acid test, which disappears within minutes. The wait time to perform the recommended reading by the technique is sufficient for the interfering color to disappear completely.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
METHOD COMPARISON AND METHODOLOGICAL SPECIFICITY
The Kidney Stone kit was compared with other methods for Kidney Stone dosing kit. According to the results of 10 samples, the correlation was 100%. With these results we can conclude that the kit shows good methodological specificity.

Precision
REPEATABILITY
The repeatability was calculated from 10 successive determinations using 3 different samples, and the results were 100% correlated.

REPRODUCIBILITY
The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 different samples, and the results had a correlation of 100%.

DIAGNOSTIC SIGNIFICANCE
Stones are found on the renal pelvis, ureter, and bladder. Small stones can come out with urine. The study of the structure and composition of urinary stones points to important data that has allowed to clarify many questions regarding its origin and formation. It is evident that its origin is not unique, and that there are several factors capable of producing the formation of stones; In each given case, the formation of the calculi may be due to one of these factors alone or the simultaneous presence of several of them. Fluctuations in urinary pH play a key role in the genesis of renal calculi. Normally, the urine pH ranges from 6.0 to 6.5. At pH below this value, around 5.3, uric acid becomes poorly soluble and can crystallize. Calcium oxalate also precipitates at acidic or neutral pH, and struvite at alkaline pH, usually associated with infection by microorganisms such as Proteus. The formation of mixed calculi is favored by the obstruction caused by calcium crystals followed by infection. This infection can lead to accumulation and deposition of ammonium crystals.

NUMBER OF TESTS
Biocal Kidney Stone kit contains reagents for 10 complete analysis.

BIBLIOGRAPHIC REFERENCES
1- HENRY, Richard J .; CANNON, Donald C .; WINKELMAN, James W .. Clinical Chemistry: Principles And Technics. 2. ed. New York: Harper & Row, Publishers, Inc, 1974. 1629 p.
2- HENRY, John Bernard. Clinical Diagnoses: And Treatments By Laboratory Methods. 20. ed. Barueri, Sp: Manole, 2008. 1726 p. Translation: Ida Cristina Gubert.
3- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Kidney Stone kit:
10269360091

Review: January/2020

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED