## CALIBRADOR HbA1c

### K092

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### FINALIDADE

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1 - The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.

2 - LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985- 92.

3 - JEPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.

4 - HOEZEL W, WEYKAMP C et al, Clin. Chem., (2004), 50:1:166-74.

5 - Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

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## HbA1c CALIBRATOR K092 USAGE INSTRUCTIONS

### FUNCTION

Reagent used to monitor the accuracy and precision for quantitative determination of hemoglobin A1c (HbA1c) in automated methodologies in photometrics systems, only for *in vitro* diagnostic.

### REAGENTS

**Number 1 - Calibrator Level 1** - Store between 2 and 8°C. Contains: Lyophilized HbA1c. **Potentially infectious.**

**Number 2 - Calibrator Level 2** - Store between 2 and 8°C. Contains: Lyophilized HbA1c. **Potentially infectious.**

**Number 3 - Calibrator Level 3** - Store between 2 and 8°C. Contains: Lyophilized HbA1c. **Potentially infectious.**

**Number 4 - Calibrator Level 4** - Store between 2 and 8°C. Contains: Lyophilized HbA1c. **Potentially infectious.**

**Number 5 - Calibrator Level 5** - Store between 2 and 8°C. Contains: Lyophilized HbA1c. **Potentially infectious.**

**Warning:** The concentration varies with the lot. See bottle label.

### Presentation

Presentation	Calibrator N°1	Calibrator N°2	Calibrator N°3	Calibrator N°4	Calibrator N°5
1	1 x 0,25 mL	1 x 0,25 mL	1 x 0,25 mL	1 x 0,25 mL	1 x 0,25 mL
2	2 x 0,25 mL	2 x 0,25 mL	2 x 0,25 mL	2 x 0,25 mL	2 x 0,25 mL
3	4 x 0,25 mL	4 x 0,25 mL	4 x 0,25 mL	4 x 0,25 mL	4 x 0,25 mL
4	1 x 1,0 mL	1 x 1,0 mL	1 x 1,0 mL	1 x 1,0 mL	1 x 1,0 mL
5	2 x 1,0 mL	2 x 1,0 mL	2 x 1,0 mL	2 x 1,0 mL	2 x 1,0 mL
6	4 x 1,0 mL	4 x 1,0 mL	4 x 1,0 mL	4 x 1,0 mL	4 x 1,0 mL

### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport, at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture.

### SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Reagent should be handled with caution, for it is subject to biological contamination.
- 5- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the product is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- 6- Reagent was tested for anti-HIV, anti-HCV antibodies and HB antigen. Risk of infection can not be excluded and the reagent must be handled with the same caution observed for the patient's serum. Potentially infecting.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

### PROCESS DESCRIPTION

#### Calibrator Reconstitution

**Before reconstituting the Calibrator, the volume indicated on the bottle label must be observed and proceed as described below:**

Remove the plastic caps from the Calibrator 1 to 5 along with the metal seal, and then carefully remove the rubber caps.

**Presentations 1 to 3 (0.25 mL):** Add exactly **0.25 mL** of distilled water into each bottle.

**Presentations 4 to 6 (1.0 mL):** Add exactly **1.0 mL** of distilled water to each bottle.

Replace the rubber caps and allow to stand for 10 minutes. Homogenize carefully by inversion to prevent material from adhering to the bottle wall. Avoid foaming.

**WARNING:** After reconstitution the calibrators will be ready for use and should NOT be submitted to the recommended hemolysis procedure for the samples.

#### Storage and Stability

Store away from light and heat.

After reconstitution the calibrators are stable for up to 10 days at 2 to 8°C or 30 days at -20°C.

Calibrators should be kept out of their ideal storage temperature only for the time needed for their use in the equipment. The reconstituted calibrator should not be frozen and thawed more than once. Therefore, for smaller routines, fractionate the material into aliquots before freezing.

### TECHNICAL

The use of the Calibrator should be according to the programming of the HbA1c Birreagente - Bioclin kit.

### BIBLIOGRAPHIC REFERENCES

- 1 - The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- 2 - LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL *et al.* Clin. Chem., (2001), 47:1985-92.
- 3 - JEPSSON JO, KOBOLD U, BARR J, FINKE A *et al.* Clin. Chem. Lab. Med., (2002), 40:78-79.
- 4 - HOEZEL W, WEYKAMP C *et al.* Clin. Chem., (2004), 50:1:166-74.
- 5 - Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: [bioclin@bioclin.com.br](mailto:bioclin@bioclin.com.br)  
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium


**CUSTOMER SERVICE**

Customer Advisory Service  
Phone: 0800 0315454.  
E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

ANVISA registration for HbA1c Calibration Kit: 10269360149

**Review:** November/2019

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <math>\leq N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED