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- The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985-92.
- JEPPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.
- HOEZEL W, WEYKAMP C et al, Clin. Chem., (2004), 50:1:166-74.
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1 - The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.

2 - LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985- 92.

3 - JEPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.

4 - HOEZEL W, WEYKAMP C et al, Cin. Chem., (2004), 50:1:166-74.

5 - Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

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USAGE INSTRUCTIONS

FUNCTION

Reagent used to monitor the accuracy and precision for quantitative determination of hemoglobin A1c (HbA1c) in automated methodologies in photometrics systems, only for *in vitro* diagnostic.

CALIBRATOR

Number 1 - Calibrator level 1 - Store between 2 and 8°C.
Contains: Lyophilized HbA1c.

Number 2 - Calibrator level 2 - Store between 2 and 8°C.
Contains: Lyophilized HbA1c.

Number 3 - Calibrator level 3 - Store between 2 and 8°C.
Contains: Lyophilized HbA1c.

Number 4 - Calibrator level 4 - Store between 2 and 8°C.
Contains: Lyophilized HbA1c.

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TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For *in vitro* diagnostic only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Reagent should be handled with caution, for it is subject to biological contamination.
- 5- Reagent was tested for anti-HIV, anti-HCV antibodies and HB antigen. Risk of infection can not be excluded and the reagent must be handled with the same caution observed for the patient's serum. Potentially infecting.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or requested by the SAC (Advisory Service Customer) of Quibasa.

PROCESS DESCRIPTION

Calibrator utilization must obey the specific techniques for equipments used for calibrations.

Preparation: To use, the calibrators should be treated as sample, requiring the completion of hemolysis as instructed use of Glycosylated Hemoglobin Kit K091.

Storage and Stability

Store away from light and heat, 2 to 8°C.
After preparation, the calibrator is stable for up to 30 days if kept between 2 and 8°C.

REFERENCE VALUES

The values presented in these tables were obtained through repeated dosages using the indicated method.

These values should be used as guidance;
Each laboratory should establish its own limits of precision.
The laboratory average, shall therefore be inside the acceptable average for the standards supplied.

See Reference Values table enclosed.

PRODUCT PERFORMANCE

CONTROL QUALITY

Precision

REPEATABILITY

20 successive measurements were performed with four standards of different levels, obtaining the following results:

	Calibra- tor 1	Calibra- tor 2	Calibra- tor 3	Calibra- tor 4
Average Concentration (%)	2,84	7,14	10,07	13,73
Standard Deviation (%)	0,03	0,04	0,05	0,09
Coefficient of Variation (%)	0,93	0,56	0,49	0,65

REPRODUCIBILITY

20 dosages were performed during 3 consecutive days with 3 samples, obtaining the following results:

	Calibra- tor 1	Calibra- tor 2	Calibra- tor 3	Calibra- tor 4
Average Concentration (%)	2,84	7,14	10,08	13,71
Standard Deviation (%)	0,01	0,01	0,01	0,04
Coefficient of Variation (%)	0,20	0,08	0,11	0,26

BIBLIOGRAPHIC REFERENCES

- 1 - The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- 2 - LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL *et al.* Clin. Chem., (2001), 47:1985-92.
- 3 - JEPSSON JO, KOBOLD U, BARR J, FINKE A *et al.* Clin. Chem. Lab. Med., (2002), 40:78-79.
- 4 - HOEZEL W, WEYKAMP C *et al.* Clin. Chem., (2004), 50:1:166-74.
- 5 - Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

MANUFACTURER'S DATA

QUIBASA QUÍMICA BÁSICA Ltda
 Rua Teles de Menezes, 92 - Santa Branca
 CEP 31565-130 - Belo Horizonte - MG - Brasil
 Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
 E-mail: bioclin@bioclin.com.br
 CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
 Phone: 0800 0315454.
 E-mail: sac@bioclin.com.br

ANVISA registration for HbA1c Calibration Kit: 10269360149

Review: August/12

SIMBOLOGIA UNIVERSAL

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON