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REFERENCIAS BIBLIOGRÁFICAS
1. FAIRBANKS, V.F.; KLEE, G.G. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
2. WICH, M.; PINGERRA, W.; LEHMANN, P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
3. GUDER, W.G.; ZAWTA, B. et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46. a.
4. DATI, F.; SHUMANN, G.; THOMAS, L.; AGUZZI, F.; BAUDNER, S.; BIENVENU, J. et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34:517-20.
5. THOMAS L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:273-5.
6- BAKKER AJ, MÜCKE M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9):1240-1243.
7. Bioclin – Dados de arquivo

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USAGE INSTRUCTIONS


FUNCTION

Method for quantitative determination of Latent Binding Capacity of the Iron in serum or plasma samples. Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric test

A known concentration of Iron ions is incubated with serum and will bind specifically with the Transferrin in insaturated Iron-binding sites. The exceeding Iron ions are measured through the reaction with Ferene. The difference between the exceeding Iron and the total quantity of iron initially added to the serum is equivalent to the quantity bounded to the Transferrin, which will be the Latent Binding Capacity of the Iron – CLLF.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Buffer (pH 8,7), Iron Sulphate II and Tyurea.

Number 2 – Color Reagent - Store between 2 and 8°C. Contains: Ascorbic Acid, Ferene and Tyurea.

Reagent N° 3 - Calibrator - Store between 2 and 8°C. Contains: Serum Proteins and 0.9% Sodium Azide.

Potentially Infectious.

Note: The concentration of Serum Proteins varies with the batch.

Reagent N° 4 - Control - Store between 2 and 8°C. Contains: Serum Proteins and 0.9% Sodium Azide.

Potentially Infectious.

Note: The concentration of Serum Proteins varies with the batch.

PRESENTATION

Presentation	Reag. 1	Reag. 2	Reag. 3	Reag. 4
1	1 x 20 mL	1 x 5 mL	1 mL	1 mL
2	2 x 20 mL	2 x 5 mL	1 mL	1 mL
3	3 x 20 mL	3 x 5 mL	1 mL	1 mL
4	4 x 20 mL	4 x 5 mL	1 mL	1 mL
5	5 x 20 mL	5 x 5 mL	1 mL	1 mL
6	1 x 40 mL	1 x 10 mL	1 mL	1 mL
7	2 x 40 mL	2 x 10 mL	1 mL	1 mL
8	4 x 40 mL	4 x 10 mL	1 mL	1 mL
9	2 x 40 mL	1 x 20 mL	1 mL	1 mL
10	4 x 40 mL	2 x 20 mL	1 mL	1 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry auto-chemistry analyzers, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 2 and 8°C. The transport temperature between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Reagent N° 2 must be kept protected from light.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Calibration should be repeated periodically to verify any changes in the response of the colorimeter or spectrophotometer.
- 6- Hemolysis, even discrete, interferes in the dosage.
- 7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 9- Do not use the product in case of damaged packaging.
- 10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum or heparinized plasma free of hemolysis.

Serum samples and heparinized Plasma are stable for 7 days at 2 to 8°C, or 1 month at -20°C.

The sample for therapeutic control should always be collected at the same time.

The sample should not be frozen more than once.

PROCESS DESCRIPTION

REAGENT PREPARATION

The reagents are ready to use.

TECHNIQUE

For the use of the kit, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol P and Biocontrol N Bioclin kits.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

CALCULATIONS

The obtained value corresponds to the CLLF.

CLLF - Latent Iron-Binding Capacity

CTLF - Total Iron-Binding Capacity

IST - Transferrin Saturation Index

CTLF = CLLF + Iron Serum Concentration

$$\text{IST}(\%) = \frac{\text{Serum Iron} \times 100}{\text{CTLF}}$$

INTERFERENCES

There was no sign of interference of Ascorbic Acid up to 30 mg/dL, conjugated and free Bilirubin up to 60 mg/dL, lipemia up to 2000 mg/dL of Triglycerides, RF up to 350 IU/mL, Copper up to 15 mg/dL, Zinc up to 15 mg/dL and Hemoglobin < 200 mg/dL.

In rare cases, samples from patients with gammopathies may show false results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The Iron-Binding Capacity values of the kit's internal calibrator and the Biocal were defined by iron and transferrin dosages. Therefore, they are traceable to the primary iron standards (NIST SRM 937) and transferrin (ERM DA470k / IFCC).

REFERENCE VALUES

The reference values for this method were obtained through the determination of the Iron-Binding Capacity in healthy populations of male and female.

Latent Capacity (CLLF)	140 - 280 $\mu\text{g/dL}$
Total Binding Capacity (CTLF)	250 - 410 $\mu\text{g/dL}$
Transferrin Saturation (IST)	20 - 50%

To convert the values of $\mu\text{g/dL}$ to $\mu\text{mol/L}$ (SI), multiply the obtained results by 0,179.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS

Iron-Binding Capacity Automated kit was compared with other methods for determination of Iron-Binding Capacity dosage commercially available measurement. 07 analyzes were performed and the results were evaluated. The equation linear obtained was $Y = 1,0128X - 2,969$ and the correlation coefficient 0,9987. With these results, we can conclude that the kit shows good methodological specificity.

Precision
REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration ($\mu\text{g/dL}$)	179,05	384,60	487,70
Standard Deviation ($\mu\text{g/dL}$)	1,50	4,03	3,60
Coefficient of Variation (%)	0,84	1,05	0,74

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration ($\mu\text{g/dL}$)	155,65	334,03	422,75
Standard Deviation ($\mu\text{g/dL}$)	1,45	3,07	3,80
Coefficient of Variation (%)	0,93	0,92	0,90

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Iron-Binding Capacity. The average found was $5,81 \mu\text{g/dL}$, with deviation standard of $0,07 \mu\text{g/dL}$. The sensitivity, which indicates the method limit detection, corresponds the average plus 3 times the standard deviation, and is equal to $6,025 \mu\text{g/dL}$.

Linearity

The reaction is linear until the concentration of $750 \mu\text{g/dL}$ or $135 \mu\text{mol/L}$.

For higher values, dilute the sample with NaCl solution at 0,85% and repeat the dosage. Multiply the results by the dilution factor applied.

DIAGNOSTIC SIGNIFICANCE

Several diseases involving Iron in the organism may be evaluated from the determination of the Iron Binding Capacity combined with the determination of Iron Serum. The sum of those two values represents the maximum iron concentration that the proteins are capable of binding. In a iron deficiency, the iron binding capacity is increased, and in severe or malignant inflammatory diseases this binding capacity is reduced.

BIBLIOGRAPHIC REFERENCES

1. FAIRBANKS, V.F.; KLEE, G.G. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
2. WICH, M.; PINGERRA, W.; LEHMANN, P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
3. GUDER, W.G.; ZAWTA, B. et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.
4. DATI, F.; SHUMANN, G.; THOMAS, L.; AGUZZI, F.; BAUDNER, S.; BIENVENU, J. et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34:517-20.
5. THOMAS L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:273-5.
- 6- BAKKER AJ, MÜCKE M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9):1240-1243.
7. Bioclin – Dados de arquivo

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC

REP

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

English 2/2

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Iron-Binding Capacity Automation kit: 10269360180

Review: January/2019

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

EC

REP

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

CONTROL

CONTROL+

POSITIVE CONTROL

CONTROL-

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED