CAPACIDADE LIGADORA DE FERRO

REF K009

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Iron-Binding Capacity. Colorimetric test, only for *in vitro* diagnostic.

PRINCIPLE OF ACTION

Methodology: Colorimetric

The iron is transported in blood by a specific Betaglobulin denominated Trasferrin. Each Transferrin is capable of transporting 2 atoms of ion Fe⁺³ and in physiologic conditions the transferring is found approximately 30% saturated. To determinate the Iron-Binding Capacity (IBC) the transferrin molecule, incubate the testing serum with the Iron Standard. This iron will saturate the available sites of Transferring and the excess of Iron (not bind) will be dosed through a Ferrozine colored complex, presenting the IBC.

REAGENTS

Number 1 - Buffer - Store between 15 and 30°C. Contains: Tris Buffer 2,5 mmol/L, Hydroxilamine 300 mmol/L and surfactant.

Number 2 - Color Reagent - Store between 15 and 30°C. Contains: Ferrozine 30 mmol/L and stabilizer.

Number 3 - Standard - Store between 15 and 30°C. Contains: Iron 500 mcg/dL.

PRESENTATION

Reagent	Volume
Nº 1	50 mL
Nº 2	3,2 mL
Nº 3	20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, water-bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized in laboratories of clinical analysis.

TRANSPORTAION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must to be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- The water level in water-bath must by higher than the reagent level in test tubes.

6- Hemolysis, even discrete, interferes in the dosage.

7- Frequently, the absorbance initial values (in sample with translucid supernatant) are equal to zero. In this case, use only the value of A₂, in calculations.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Data Sheet Safety of Chemicals) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum obtained free of hemolysis. The analytical is stable for 7 days between 2 and 8°C.

In dosages for therapeutic control the samples must be collected always in the same time.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

The glassware must be meticulously clean and absent of the Iron presence.

Mark 3 tess tubes with the letters B (Blank), P (Standard), A (Sample) and proceed as follows:

	Blank	Standard	Sample
Reagent Nº 1	1,0 mL	--	1,5 mL
Distilled Water	1,5 mL	2,0 mL	--
Sample	--	--	500 µL
Reagent Nº 3	--	500 µL	500 µL

Smoothly homogenize and incubate in water-bath at 37°C for 10 minutes. Determinate the absorbance of the Sample in 560 nm (540 - 580 nm), hitting the zero with the Blank. This will be the absorbance A₁.

And then add:

	Blank	Standard	Sample
Reagent Nº 2	1 drop	1 drop	1 drop

Smoothly homogenize and incubate in water-bath at 37°C for 10 minutes. Determinate the absorbance of the sample in 560 nm (540 - 580 nm), hitting the zero with the Blank. This will be the absorbance A₂.

CALCULATIONS

The calculations below are based, also, in the methodology described in Iron Serum Bioclin kit.

CLLF - Latent Iron-Binding Capacity

CTLF - Total Iron-Binding Capacity

IST - Transferrin Saturation Index

CLLF (mcg/dL) = 500 - $\left(\frac{(A_2 - A_1)}{\text{Standard Absorbance}} \times 500 \right)$

CTLF (mcg/dL) = CLLF + Iron Serum

IST (%) = $\frac{\text{Iron Serum}}{\text{Total Capacity (mcg/dL)}} \times 100$

The results are expressed in mcg/dL.

PROCEDURE LIMITATIONS

An insufficient shaking of the reagent can result in wrongly elevated values. The use of iron contaminated water will increase the intensity of Blank's color, leading to wrongly low results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

The reference values in mcg/dL, for this method were obtained through the determination of the Iron-Binding Capacity in healthy populations of male and female. Latent Capacity (CLLF).....140 - 280 mcg/dL Total Binding Capacity (CTLF).....250 - 410 mcg/dL Transferrin Saturation (IST).....20 - 50%

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY
The Iron Binding Capacity kit was compared with other methods for measurement of Iron-Binding Capacity commercially available. 07 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1,003X - 0,671$, and the correlation coefficient 0,999. With these results, we can conclude that the kit shows good methodological specificity.

Precision
REPEATABILITY
The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mcg/dL)	144,35	202,55	297,20
Standard Deviation (mcg/dL)	0,75	0,76	0,77
Coefficient of Variation (%)	0,52	0,37	0,26

REPRODUCIBILITY
The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mcg/dL)	143,67	202,73	297,12
Standard Deviation (mcg/dL)	0,83	0,20	0,08
Coefficient of Variation (%)	0,57	0,10	0,03

Linearity
The reaction is linear until the residual values of the Binding Capacity are equal to 450 mcg/dL. For higher values, dilute the serum with distilled or deionized water and multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE
The iron transport in plasma is made by a specific Transferrin (Betaglobulin), that can be dosed indirectly through the Iron quantity that it can bind. This quantity is denominated the Total Iron-Binding Capacity. (CTLF). For the overload diagnosis or Iron deficiency, the exams must insert a complete hematological profile, dosing the iron, of CTLF and calculating the Transferrin Saturation Index (IST).
The most frequent causes for the CTLF increase are: iron deficiency anemia, pregnancy, oral anti-contraceptive use and hypoxia state.

The capacity is reduced in the chronic inflammatory diseases major, malignant diseases, nephrosis, iron intoxication and thalassemia major.
The IST is increasing the intoxication with iron, nephrosis, thalassemia, pyridoxine, hemolytic anemia and sideroblastic, and reducing the iron deficiency, malignant diseases, chronic infestations, pregnancy and menstrual cycle.

NUMBER OF TESTS
20 Tests / 500 µL of Sample
40 Tests / 250 µL of Sample

BIBLIOGRAPHIC REFERENCES
1 - PETERS T., et coll.: J. Lab and Clin Med. 48:280,1956.
2 - HENRY, R. J., Cannon. D. C. and Winkleman, J. W., eds.: Clinical Chemistry: Principles and Technics, Haper and Row Publishers, 1974.
3 - ZAK, B.: Clin Chim. Acta 3: 328, 1958.
4 - BRITTENHAM, G.: Clin Chim. Acta. 91: 203, 1979.
5 - TONKS, D. B.: Chim. Acta 2: 393, 1983.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Iron Binding Capacity kit: 10269360096

Review: April/2013

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED