CLORETOS COLORIMÉTRICO

REF K050

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FINALIDADE

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
USAGE INSTRUCTIONS

FUNCTION

Method for determination of Chloride. Colorimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

In the presence of Chloride ions, the Mercury Thiocyanate, in an acid environment, forms the Mercuric Chloride and Thiocyanate ions. The Thiocyanate ions react with the Iron ions forming a Iron Thiocyanate, of yellow-orange color, who is proportional to the Chloride concentration in the sample.

REAGENTS

- Number 1 - Color Reagent** - Store between 15 and 30°C. Contains: Mercury Thiocyanate 2 mmol/L, Iron Nitrate 17 mmol/L and Nitric Acid solubilizer 0,2%. Handle carefully - **Toxic reagent**.
- Number 2 - Activator** - Store between 15 and 30°C. Contains: Perchloric Acid 1,50 mol/L.
- Number 3 - Standard** - Store between 15 and 30°C. Contains: Chloride 100 mEq/L and Sodium Azide 7,7 mmol/L.

PRESENTATION

Reagents	Volume
Nº 1	350 mL
Nº 2	10 mL
Nº 3	3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, watch and stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. Can be found at market specialized on Clinical Analysis Laboratories.

TRANSPORTAION AND STORAGE CONDITIONS

The storage and transport temperature shall be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must to be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- The serum or plasma shall be prepared in maximum one hour after collection. A superior time produce CO₂ loss who interferes with the chloride distribution between cells and plasma.
- Urine and liquor must be centrifuge.
- Use heparinized plasma, because the regular anticoagulants interfere in the determination.

- Handle reagent Nº 1 carefully - **toxic reagent**.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Sample or serum (heparinized), urine, sweat, liquor. The analytical is stable for 7 days between 2 and 8°C and 1 year at -20°C. ^o See **special care**.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

TECHNIQUE I – Serum, plasma or liquor:
Mark 3 test tubes with letters B (Blank), A (Sample), P (Standard) and proceed as follows:

	Blank	Sample	Standard
Reagent Nº 1	3,5 mL	3,5 mL	3,5 mL
Reagent Nº 2	100 µL	100 µL	100 µL
Standard	---	10 µL	---
Reagent Nº 3	---	---	10 µL

Homogenize and read the Sample and Standard absorbance at 500 nm, matching the zero with the Blank. The color is stable for 2 hours.

TECHNIQUE II – Urine and sweat:
Dilute the urine 1:4 with distilled or deionized water (0,1 mL of urine + 0,3 mL of distilled or deionized water). Follow the scheme below, eliminating the reagent Nº 2 (Activator). Multiply the obtained result by 4.

	Blank	Sample	Standard
Reagent Nº 1	3,5 mL	3,5 mL	3,5 mL
Diluted urine or seat	---	10 µL	---
Reagent Nº 3	---	---	10 µL

Homogenize. Read the Sample and Standard absorbance at 500 nm, matching the zero with the Blank.
Attention: To dose Chloride in sweat please follow the methodology above (eliminating the Activator), however without diluting the sweat.

CALCULATIONS

Carry out distinct calculations for each technical procedure.

$$\text{Chloride (mEq/L)} = \frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 100$$

As the reaction follows the Lambert Beer Law, it can use the Calibration Factor:

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (100 mEq/L)}}{\text{Standard Absorbance}}$$

$$\text{mEq/L} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

$$\text{Urine (mEq/24 hours)} = \text{mEq/L} \times \text{Volume (liters)}$$

The results are expressed in mEq/L or mmol/L, for the samples (serum, plasma, liquor and sweat) and in mEq/24 hours or mmol/24 hours for urine.

PROCEDURE LIMITATIONS

The use of automatic pipette appropriately calibrated for pipetting the Standard and Sample, in order to minimize problems with inaccurate volume.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 919A.

REFERENCE VALUES

The reference values in mEq/L, for this method, were obtained through the determination of Chloride in healthy populations of male and female.
Serum or plasma (Of all ages): 96 to 107 mEq/L
Urine: 150 to 250 mEq/24 hours
Liquor: 118 to 132 mEq/L
Sweat: 5 to 35 mEq/L

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

RECOVERY

The Chloride Colorimetric kit recovery analysis was made with 5 sample determinations. The accurate was calculated and was found in a good agreement with the reference values, obtaining a recovery between 95% and 103%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

Chloride Colorimetric kit was compared with other method commercially available for dosage of Chloride. 07 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1,022X - 1,0165$ with correlation coefficient equal to 0,9729. With theses results, we can conclude that the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mEq/L)	106,21	91,00	78,30
Standard Deviation (mEq/L)	1,14	2,15	0,98
Coefficient of Variation (%)	1,07	2,36	1,25

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mEq/L)	106,05	90,73	78,37
Standard Deviation (mEq/L)	0,17	0,38	0,12
Coefficient of Variation (%)	0,16	0,42	0,15

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Chlorides. The average found was 0,694 mEq/L, with standard deviation of 0,089 mEq/L. The sensitivity, which indicates the indicative of methods detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,960 mEq/L.

Linearity

The reaction is linear between 20 and 120 mEq/L. For larger values, dilute the serum with distilled or deionized water, repeat the determination and multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Chloride is the major of the extra cellular liquid (LEC) anion, being important regulating the water maintenance and distribution, in an anionic and cationic balance of LEC and in osmotic pressure. Its concentration in blood is regulated by the kidneys, suprarenal gland, lungs, skin, intestine gastric tract and blood pH.

In several clinical situations, can occur reduction (hypochloremia) or increase (hyperchloremia) of the plasmatic concentration. Hypochloremia - In nephritis associated to pyelonephritis, and patients in Addisonian crisis, the values of Chloride are reduced. Extended vomits or persistent gastric secretion can lead to chlorine loss and low plasmatic values. It can also occur hypochloremia in metabolic alkalosis, bromate intoxication, cranial injuries, kidney absorption malfunction.

Hyperchloremia - Pathological causes can be quoted: acute kidney insufficiency, metabolic acidosis, nephrotic syndrome. In insipidus diabetes, in adrenocortical hyperfunction and hyperparathyroidism. In cystic fibrosis, systemic diseases that lead to intestinal obstruction, bile cirrhosis and lung abnormalities, high values of Chloride in sweat are important indicatives of diagnosis.

NUMBER OF TESTS

100 Tests / 10 μ L of sample / 3,5 mL of Reagent

BIBLIOGRAPHIC REFERENCES

- 1 - SCHALES O.; SCHALES S., J. Biol. Chem., 1941, 140, 879.
- 2 - ZALL, D. M.; FISHER, D.; GARNER, M. Q., Anal. Chem., 1.956, 28,1665.
- 3 - BENNINGTON, L. James, Dict & Encycl. of Lab. Med. and Tech., 1.984.
- 4 - BURTIS, A. Carl; Ashwood, R. Edward, Clin., Chem., Tietz Text book of; 2nd. ed., 1.986, 1366.
- 5 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

**OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Chloride Colorimetric kit: 10269360118

Review: April/2018

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED