Bioclin

COLINESTERASE

REF K094

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FINALIDADE

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- 1 - Recommendations of the German Society for Clinical Chemistry. Standardization of methods for the estimation of enzyme activities in biological fluids: Standard method for the determination of Cholinesterase activity. J Clin Chem Clin Biochem 1992; 30: 163 - 170.
- 2 - Thomas L, Clinical laboratory diagnostics. 1st ed Frankfurt: THBooks Verlagsgesellschaft; 1998. p-65-71.
- 3 - Hallbach J, Klinische Chemie für den Einstieg. 1st ed Stuttgart: Thieme; 2001. p. 143-144.
- 4 - BAKKER AJ, MÜCKE M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9):1240-1243.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1 - Recommendations of the German Society for Clinical Chemistry. Standardization of methods for the estimation of enzyme activities in biological fluids: Standard method for the determination of Cholinesterase activity. J Clin Chem Clin Biochem 1992; 30: 163-170.
- 2 - Thomas L, Clinical laboratory diagnostics. 1st ed Frankfurt: THBooks Verlagsgesellschaft; 1998. p-65-71.
- 3 -Hallbach J, Klinische Chemie für den Einstieg. 1st ed Stuttgart: Thieme; 2001. p. 143-144.
- 4 - BAKKER AJ, MÜCKE M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9):1240-1243.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of activity of Cholinesterase. in human samples of serum or plasma. Kinetic enzymatic test, only for *in vitro* diagnostic use.


PRINCIPLE OF ACTION

Methodology: Kinetic Enzymatic (DGKC).

On the catalytic action of Cholinesterase, the Butyrylthiocholine is hydrolyzed to Thiocholine and Butyrate:

The Thiocholine reduces the Hexacyanoferrate (III), yellow, Hexacyanoferrate (II), colorless. The decrease in absorbance is measured at 405 nm.

REAGENTS

Number 1 - Buffer Reagent - Store between 2 and 8°C. Contains: Pyrophosphate < 200 mmol/L, Potassium Hexacyanoferrate (III) < 10 mmol/L, stabilizer and preservative.

Number 2 - Substrate - Store between 2 and 8°C. Contains: Butyrylthiocholine < 100 mmol/L, buffer < 30 mmol/L, stabilizers and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 20 mL	1 x 5 mL
2	2 x 20 mL	2 x 5 mL
3	1 x 40 mL	1 x 10 mL
4	2 x 40 mL	2 x 10 mL
5	2 x 40 mL	1 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated Spectrophotometer, pipettes, tips, test tubes, watch or stopwatch, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

7- Do not use the product in case of damaged packaging.
8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or heparinized plasma or EDTA can be used. Serum and plasma are for up to 1 week if kept at room temperature, 2 weeks if kept between 2 and 8°C and 6 months at -20°C.

PROCESS DESCRIPTION

The calibration stability of the Cholinesterase kit installed on refrigerated equipment is at least 21 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

REACTION CONDITIONS

It is indispensable condition the use of thermostated cuvette at 37°C with 1cm optical path and 405 nm reading.

TECHNIQUE

Bioclin recommends using the kit Biocal as calibrator and as control serum, Biocontrol N and P Bioclin Kits

	Blank	Sample
Sample	-	10 µL
Distilled Water	10 µL	-
Reagent 1	400 µL	400 µL
Homogenize and incubate for approximately 3 minutes at 37°C.		
Reagent 2	100 µL	100 µL

Homogenize again, incubate for exactly 2 minutes in a thermostated cell at 37°C and read the absorbance (405nm). Start the timer and repeat the readings at 1, 2 and 3 minutes.

As long as the limitations of each equipment are observed, the volumes of samples and reagents can be modified proportionally without impairing the test performance.

CALCULATIONS

Considering the optical path of 1 cm and reaction temperature 37°C the calculation will be:

$$\Delta A/\text{min} = [\Delta A/\text{min Sample}] - [\Delta A/\text{min Blank}]$$

$$\text{Activity (U/L)} = \Delta A/\text{min} \times 62000$$

INTERFERENT

Does not produce interference significant: Ascorbic Acid up to 30 mg/dL, Bilirubin up to 45 mg/dL, Hemoglobin up to 1000 mg/dL and lipemia up to 1400 mg/dL Triglyceride equivalents.

In rare cases, samples from patients with gammopathies may show false results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

Calibration of the kit should be done using the BIOCAL calibrator, which is traceable to the reference method DGKC (Deutsche Gesellschaft für Klinische Chemie).

REFERENCE VALUES

The reference values in U/L, for this method, were obtained through the determination of Cholinesterase in healthy populations of male and female.

Men	4620 - 11500 U/L
Women	3930 - 10800 U/L

The results provided by this kit should be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Cholinesterase kit was compared with other method commercially available to measure Cholinesterase. 40 analyzes were conducted and the results were evaluated. The linear equation obtained was $Y = 1.0028X - 170.85$ and correlation coefficient 0.9964. With these results we can conclude the kit shows good methodological specificity.

Precision**REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	2464.27	8602.94	3561.00
Standard Deviation (U/L)	13.88	114.28	32.90
Coefficient of Variation (%)	0.56	1.33	0.92

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	2419.81	8400.17	3506.55
Standard Deviation (U/L)	43.80	224.10	66.36
Coefficient of Variation (%)	1.81	2.67	1.89

Sensitivity

Minimum limit of detection is 50 U/L.

Linearity

Linear test at 20000 U/L.

DIAGNOSTIC SIGNIFICANCE

(ChE) Serum Cholinesterase (pseudocholinesterase II or cholinesterase, EC3.1.1.8) is found primarily in the liver but also in white matter of the brain, pancreas, heart and serum. ChE serum levels may be determined as an indicator of hepatic poisoning by inhalation or skin contact by some organophosphorus compounds (including some insecticides neurological or gases) in liver diseases, or before of anesthesia with succinylcholine thus eliminating a deficiency congenita these enzymes that could lead to an apnea prolonged due to slow degradation of miorelaxant. A decrease of 15 - 25% is observed in mild poisoning, a decrease of 25 - 35% is observed in moderate poisoning, a decrease of 30 - 50% is observed in severe poisoning, but is also observed in patients with severe or chronic hepatitis. A decrease of 50 - 70% occurs in patients with an advanced cirrhosis and carcinoma liver metastasis.

BIBLIOGRAPHIC REFERENCES

- 1 - Recommendations of the German Society for Clinical Chemistry. Standardization of methods for the estimation of enzyme activities in biological fluids: Standart method for the determination of Cholinesterase activity. J Clin Chem Clin Biochem 1992; 30: 163-170.
- 2 - Thomas L, Clinical laboratory diagnostics. 1st ed Frankfurt: THBooks Verlagsgesellschaft; 1998. p-65-71.
- 3 - Hallbach J, Klinische Chemie fur den Einstieg. 1st ed Stuttgart: Thieme ; 2001. p. 143-144.
- 4 - BAKKER AJ, MÜCKE M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9):1240-1243.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for Cholinesterase kit: 10269360159

Review: August/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED