CONTROLE HbA1c

K110

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FINALIDADE

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- 1 - The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- 2 - LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985-92.
- 3 - JEPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.
- 4 - HOEDEL W, WEYKAMP C et al, Clin. Chem., (2004), 50:1:166-74.
- 5 - Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1 - The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- 2 - LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985-92.
- 3 - JEPPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.
- 4 - HOEZEL W, WEYKAMP C et al, Cin. Chem., (2004), 50:1:166-74.
- 5 - Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Control used to monitor accuracy and precision of quantitative determinations of A1c Hemoglobin (HbA1c) in automated methodologies in photometric systems. Only for *in vitro* diagnostic use.

REAGENTS

Number 1 - Control Level 1 - Store between 2 and 8°C.

Contains: Lyophilized HbA1c. **Potentially infectious.**

Number 2 - Control Level 2 - Store between 2 and 8°C.

Contains: Lyophilized HbA1c. **Potentially infectious.**

Warning: The concentration varies with the lot. See bottle label.

PRESENTATION

Control	Presentation 1	Presentation 2
Nº 1	1 x 0,25 mL	1 x 1,0 mL
Nº 2	1 x 0,25 mL	1 x 1,0 mL

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Reagent should be handled with caution, for it is subject to biological contamination.

5- Reagent was tested for anti-HIV antibodies, anti-HCV and HB antigen. Risk of infection can not be excluded and the reagent must be handled with the same caution observed for the patient's serum. Potentially infecting.

6- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the product is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

Reconstitution of Control

Before reconstituting the Control, the volume indicated on the bottle label must be observed and proceed as described below:

Remove the plastic caps from Controls 1 to 2 along with the metal seal, and then carefully remove the rubber caps.

Presentation 1 (0.25 mL): Add exactly **0.25 mL of distilled water** into each bottle.

Presentation 2 (1.0 mL): Add exactly **1.0 mL of distilled water** into each bottle.

Replace the rubber caps and allow to stand for 10 minutes.

Homogenize carefully by inversion to prevent material from adhering to the bottle wall. Avoid foaming.

ATTENTION: After reconstitution, the controls will be ready for use and should NOT be submitted to the recommended hemolysis procedure for the samples.

Storage and Stability

Store away from light and heat.

After reconstitution the controls are stable for up to 10 days at 2 to 8°C or 30 days at -20°C.

The controls should be kept out of their ideal storage temperature only for the time needed for their use in the equipment. The reconstituted control should not be frozen and thawed more than once. Therefore, for smaller routines, fractionate the material into aliquots before freezing.

TECHNICAL

The use of the Control should be according to the programming of the HbA1c Birreagente - Bioclin kit.

BIBLIOGRAPHIC REFERENCES

- 1 - The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- 2 - LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985-92.
- 3 - JEPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.
- 4 - HOESEL W, WEYKAMP C et al, Clin. Chem., (2004), 50:1:166-74.
- 5 - Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
 CEP 31565-130 - Belo Horizonte - MG - Brasil
 Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
 E-mail: bioclin@bioclin.com.br
 CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
 1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
 Phone: 0800 0315454.
 E-mail: sac@bioclin.com.br

ANVISA registration for HbA1c Control Kit: 10269360177

Review: November/2019

SIMBOLOGIA UNIVERSAL

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON