CONTROLE HbA1c

K110

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Control used to monitor accuracy and precision of quantitative determinations of A1c Hemoglobin (HbA1c) in automated methodologies in photometric systems. Only for *in vitro* diagnostic use.

CALIBRATOR

Number 1 - Control Level 1 - Store between 2 and 8°C. Contains liquid stable with hemoglobin base. For concentration see label.

Number 2 - Control Level 2 - Store between 2 and 8°C. Contains liquid stable with hemoglobin base. For concentration see label.

PRESENTATION

Control	Volume
Nº 1	0,25 mL
Nº 2	0,25 mL

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For *in vitro* diagnostic only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Reagent should be handled with caution, for it is subject to biological contamination.
- 5- Reagent was tested for anti-HIV antibodies, anti-HCV and HB antigen. Risk of infection can not be excluded and the reagent must be handled with the same caution observed for the patient's serum. Potentially infecting.
- 6- We recommend applying the local, state and federal regulations for environmental protection, so the disposal of reagents and biological material can be made in accordance with the current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available in the website www.bioclin.com.br or upon request to Bioclin's Customer Services.

PROCESS DESCRIPTION

Control utilization must obey the specific techniques for the used equipment.

Preparation: Controls should be treated as a sample, according to the instructions for using the Bioclin kit K091.

REFERENCE VALUES

The values presented in these tables were obtained through repeated dosages using the indicated method. These values should be used as guidance. Each laboratory should establish its range of reference values, according to the population served.

The average of the laboratory, however, should be within the range acceptable for the given parameters.

See Reference Values table enclosed

PRODUCT PERFORMANCE

QUALITY CONTROL

Precision

REPEATABILITY

20 successive measurements were performed with two controls, obtaining the following results:

	Control 1	Control 2
Average Concentration (%)	3,23	7,14
Standard Deviation (%)	0,03	0,04
Coefficient of Variation (%)	0,93	0,56

REPRODUCIBILITY

20 dosages were performed during 3 consecutive days with 2 controls, obtaining the following results:

	Control 1	Control 2
Average Concentration (%)	2,84	7,14
Standard Deviation (%)	0,01	0,01
Coefficient of Variation (%)	0,20	0,08

BIBLIOGRAPHIC REFERENCES

- 1 - The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- 2 - LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985-92.
- 3 - JEPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.
- 4 - HOEDEL W, WEYKAMP C et al, Clin. Chem., (2004), 50:1:166-74.
- 5 - Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

MANUFACTURER'S DATA

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454.

E-mail: sac@bioclin.com.br

ANVISA registration for HbA1c Control Kit: 10269360177

Review: July/12

SIMBOLOGIA UNIVERSAL

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON