COVID-19 IgG/IgM BIO

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FINALIDADE

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- 1- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- 2- Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- 3- Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
- 4- ANVISA. RDC Nº107, DE 5 DE SETEMBRO DE 2016. DOU nº172, 06 de setembro de 2016. Ministério da Saúde.
- 5- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
4- ANVISA. RDC Nº107, DE 5 DE SETEMBRO DE 2016. DOU nº172, 06 de setembro de 2016. Ministério da Saúde.
5- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Rapid immunochromatographic test for qualitative determination of IgM and IgG antibodies to the COVID-19 virus in samples of whole blood (EDTA, Heparin or Citrate), serum or plasma (EDTA or Heparin). For in vitro diagnostics only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
The COVID-19 IgG/IgM BIO kit is an immunochromatographic test that separately detects antibodies of the IgM and IgG classes to the COVID-19 virus. The kit contains goat anti-human IgG antibodies in the IgG test line, and goat anti-human IgM antibodies in the IgM test line. After adding the sample, antigens from the COVID-19 virus conjugated to colloidal gold particles, present in the kit, bind to the sample antibodies and migrate by capillarity through the membrane. In the case of IgG reagent samples, a red line will be formed in the IgG test region. IgM reagent samples will form a red line in the IgM test region. The sample-antigen mixture continues running through the membrane to the control line, where a red line is formed, for both reactive and non-reactive samples. The formation of the red line in the region of the control line shows that the kit is suitable for use.

REAGENTS

Reagent N° 1: Cassette - Store at 2 to 30°C. **Do not freeze.** Contains:

- A)** 01 plastic cassette (device) with 01 test strip containing:
- Control Region (C): Rabbit IgG antibody goat.
 - Test Region M (1): Goat anti-human IgM antibody.
 - Test Region G (2): Goat anti-human IgG antibody.
 - Conjugate: COVID-19 Virus Antigen conjugated with colloidal gold particles.
- B)** 01 silica sachet

Reagent N° 2: Diluent - Store at 2 to 30°C. **Do not freeze.** Contains: Buffer and preservative.

Collection Accessories:

- Disposable plastic pipette
- Disposable lancet (Sterile)

PRESENTATION

Reagent	Number Cassette per Packaging	Diluent
Cassette + Diluent (Includ pipette)*	10	1.50 mL
	25	3.50 mL
Cassette + Diluent (Includ lancet and pipette)*	10	1.50 mL
	25	3.50 mL

* The number of accessories is proportional to the number of cassettes.

EQUIPMENTS AND OPERATIONAL INPUTS

Pipettes, tips, watches and stopwatches, alcohol 70% (p/p), cotton and waste for disposal. They can be found at markets specialized in laboratories of clinical analysis.

CONDITIONS OF STORAGE AND TRANSPORTATION

Transport can be carried out at temperatures up to 45°C by maximum, 7 days. The storage temperature must be 2 at 30°C. Avoid moisture. **Do not freeze.**

WARNINGS

- 1- For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- The used device has potential infectious risk. The safety routine must be followed strictly and the use of disposable gloves is indispensable.
- 4- Use clean and dry recipients to settle the samples.
- 5- Open the aluminum envelop, to take the cassette off, only for its immediate use.
- 6- The kit must not be frozen or exposed to temperatures above 45°C this expose may cause its deterioration.
- 7- The reagents, as well as the sample, must reach environment temperature before the test begins.
- 8- Do not mix reagents between kits with the same batch or with different batch numbers.
- 9- Do not use kit after the expiration date. The expiry date is indicated on the product label.
- 10- The cassette (test device), the pipette and the lancet are disposable and cannot be reused.
- 11- The kit materials that had contact with the patient's biological material must have a special destination. The lancets and pipettes must be disposed of with sharps. Cassettes (test devices) and vials of diluent must be discarded in biological waste.
- 12- Boxes, instructions for use, aluminum envelope and silica sachet can be disposed of in regular trash.
- 13- Reagent N° 2 contains sodium azide, irritating to skin and mucous membranes. Handle the reagent with care.
- 14- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.
- 15- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by Quibasa's Customer Service Department.
- 16- Do not use the product in case of damage to the packaging.
- 17- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.


SAMPLES

Whole Blood: Can be collected by venipuncture using the suitable anticoagulant (EDTA, Heparin or Citrate) or by digital puncture with the aid of a lancet. Whole blood collected by finger prick should be used for immediate testing. Whole blood collected by venipuncture with anticoagulant can be stored between 2 and 8°C and tested up to 24 hours after collection. **Do not freeze.**


Plasma and Serum: Must be separated from the blood interface immediately after centrifugation, which should occur after the gathering. Plasma or serum can be stored for up to 3 days between 2 and 8°C and 6 months at -20°C. To obtain the plasma, perform collection using EDTA or Heparin as an anticoagulant. Do not use hemolyzed or lipemic samples.

DESCRIPTION OF THE TECHNICAL PROCESS - VENOPUNCTION

- 1- The sample must be at a temperature between 15 and 30°C before start the test.
- 2- Remove the cassette (test device) from the protective packaging place it on a clean, level surface and identify it properly.
- 3- For Whole Blood: Transfer 20 µL of whole blood into the sample window. Hold the bottle vertically and apply 2 drops (60-80 µL) of Diluent (Reagent N° 2) to the sample well.

- 4- For Serum or Plasma: Transfer 10 µL of serum or plasma into the sample window. Hold the bottle vertically and apply 2 drops (60-80 µL) of Diluent (Reagent N° 2) in the sample well.


- 5- Wait for the lines to form. Interpret the results between 10 and 15 minutes. Do not interpret after 15 minutes.

TECHNIQUE - PUNCTURE DIGITAL / LANCET

- 1- Remove the cassette (test device) from the protective packaging, place it on a clean and level surface and identify it properly.
 - 2- Wash and dry the patient's hands. Pass 70% alcohol (p/p) on fingertip for antiseptis of the area used.
- Note:** Due to the absence of 70% alcohol (p/p) in the market, due to the COVID-19 pandemic, it is recommended the use of topical antiseptics, indicating the solution chlorhexidine gluconate alcoholic strength 0.5%.
- The product must be applied in a sufficient quantity to moisten the entire area to be treated with the aid of gauze sterile. After application, completely dry the if necessary, repeat the application and wait for drying complete product.

- 3- Press the tip of the finger, thus promoting a accumulation of blood in this region, so that it will be punctured with the lancet.
- 4- Remove the protective cover from the lancet.
- 5- Position and press the lancet firmly over the area to be punctured. Then, blood will flow out through the perforated area.
- 6- Discard the first drop of blood with the aid of a dry cotton. The first drop of blood may contain fluids interstitial and intracellular cells that can interfere with the test.

- 7- Collect blood using a disposable plastic pipette that comes with the kit. Tighten the plastic pipette above the dash marked and touch the open cavity of the plastic pipette to the drop of blood. Relieve the pressure in the plastic pipette so that the blood is aspirated. Aspirate the blood to the mark marked on the plastic pipette disposable (equivalent to 10 µL of blood).
- 8- In the sample well, dispense the blood by pressing pipette again.
- 9- Repeat the procedures of items 7 and 8 again, to that the equivalent of 20 µL of blood is dispensed in total.
- 10- Hold the bottle vertically and apply 2 drops (60-80 µL) of Diluent (Reagent No. 2) in the sample well.
- 11- Wait for the lines to form. Interpret the results between 10 and 15 minutes. **Do not interpret after 15 minutes.**

INTERPRETATION OF THE RESULT

IgM Reagent Test: Formation of a red line on the test region M (1) and another line in the control region (C) in first 10 to 15 minutes. Do not interpret after 15 minutes.

IgG Reagent Test: Formation of a red line on the test region G (2) and another line in the control region (C) in the first 10 to 15 minutes. Do not interpret after 15 minutes.

IgM and IgG Reagent Test: Formation of two red lines in the test region (M = 1 and G = 2) and a red line in the control (C) in the first 10 to 15 minutes. Do not interpret after 15 minutes.

Non-Reagent Test: Formation of a red line in the region control (C) and complete absence of a red line in the regions tests (M = 1 and G = 2). Do not interpret after 15 minutes.

Invalid Result: Complete absence of a line in the control region (C) with or without a red line in the test regions (M = 1 and G = 2). Test the sample again.

Notes:
1- The line in the control region can appear before 10 minutes incubation of the test, and this does not mean that results do not reagents can be interpreted ahead of time.
2- The results should not be interpreted after 15 minutes to avoid incorrect results.
3- The results provided by this kit must be interpreted responsible medical professional, and it is not the only criterion for determining the diagnosis and/or treatment of the patient.

CROSS REACTIVITY
1349 samples negative for COVID-19 and positive for influenza, rhinovirus and/or respiratory syncytial virus, were subjected to the cross-reactivity test. 2 false positive results were observed for influenza, 2 for rhinovirus and 2 for respiratory syncytial virus. Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL
The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

LIMITATIONS OF THE PROCESS
1- The COVID-19 IgG/IgM BIO Bioclin kit is a qualitative test and should only be used for in vitro diagnostic use.
2- The COVID-19 IgG/IgM BIO Bioclin Kit determines only the presence of antibodies against the COVID-19 virus and should not be used alone as a diagnostic criterion. Reagent results should be evaluated with other available clinical information. Isolated test results cannot be used for definitive diagnosis.
3- A non-reactive result does not exclude the possibility of infection by the COVID-19 virus. Scientific studies reveal that people infected with the COVID-19 virus may have an immunological window for the detection of IgM and IgG.
4- Results of immunosuppressed patients should be evaluated with caution.
5- Like any diagnostic test, the results must be interpreted together with other clinical data of the patient.
6- If the non-reactive result remains and the clinical symptoms persist, an additional test should be performed using another method. A non-reactive result does not exclude, at any time, the possibility of infection by the COVID-19 virus.

PRODUCT PERFORMANCE
QUALITY CONTROL
Accuracy
METHODS COMPARISON, METHODOLOGICAL SPECIFICITY AND SENSITIVITY
To establish the sensitivity and specificity of the kit COVID-19 IgG/IgM BIO compared to another test for COVID-19, 547 known clinical samples were studied. Another qualitative kit commercially available was used for comparison with the kit and determination of relative sensitivity and specificity. The clinical sensitivity found was 96,3% and the specificity was 99,2%.

Precision
REPEATABILITY
The repeatability was calculated from 10 successive determinations, using 5 different samples, obtaining the following results:

Samples	Repetition Nº	Expected Result	Found Result
Human Blood 01	10	100% Reagent IgG	100% Reagent IgG
Human Blood 02	10	100% Non-Reagent	100% Non-Reagent
Human Serum 01	10	100% Non-Reagent	100% Non-Reagent
Human Serum 02	10	100% Reagent IgM	100% Reagent IgM
Human Plasma 01	10	100% Reagent IgM/IgG	100% Reagent IgM/IgG

REPRODUCIBILITY
The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Samples	Nº of Days	Repetition Nº	Expected Result	Found Result
Human Blood 01	03	10	100% Non-Reagent	100% Non-Reagent
Human Blood 02	03	10	100% Reagent IgM	100% Reagent IgM
Human Serum 01	03	10	100% Reagent IgM/IgG	100% Reagent IgM/IgG
Human Serum 02	03	10	100% Non-Reagent	100% Non-Reagent
Human Plasma 01	03	10	100% Reagent IgM/IgG	100% Reagent IgM/IgG

DIAGNOSTIC SIGNIFICANCE
Coronaviruses (CoV) belong to the Coronaviridae family, and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are fever, cough, dyspnoea, myalgia or fatigue. It can cause diseases ranging from the common cold to more serious illnesses, such as the Middle East Respiratory Syndrome (MERS-CoV) and the Severe Acute Respiratory Syndrome (SARS-CoV). In December 2019, a series of pneumonia cases of unknown cause appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples it was shown to be Coronavirus infection, which was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents, cases of infection have already been reported.

BIBLIOGRAPHIC REFERENCES
1- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
4- ANVISA. RDC Nº107, DE 5 DE SETEMBRO DE 2016. DOU nº172, 06 de setembro de 2016. Ministério da Saúde.
5- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE
Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE
Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for COVID-19 IgG/IgM BIO kit: 10269360322

Review: June/2020 (1)

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED