CREATININA

REF K016

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Creatinine. Colorimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Jaffe modified

End point colorimetric: the Creatinine reacts with the Picric Acid, forming a complex of reddish-yellow. This pH occurs in a maximum formation of the Picric-Creatinine color complex, and, also with other plasmatic elements. With the adding of the Acid Reagent, the pH is reduced and the correct color creatinine is undone, remaining the cromogens color. For difference between the readings obtained in the acid pH, obtains the real value of Creatinine.

Fixed time Kinetic: the variation on the alkaline Picric formation speed, without acidification of the products formed, its gained through the spectrophotometers reading on the first minutes. The readings obtained are free from reaction with cromogens pricrate, allowing, the creatinine determination.

REAGENTS

Number 1 - Alkaline Reagent - Store between 15 and 30°C. In low temperatures the reagent can present turbidity, to discard it heat up to 37°C. Mix before use. Contains: Sodium Hydroxide 110 mmol/L, Sodium Carbonate 75 mmol/L and surfactant.

Number 2 - Picric Acid - Store between 15 and 30°C. Contains: Picric Acid 60 mmol/L.

Number 3 - Acid Reagent - Store between 15 and 30°C. Contains: Acetic Acid 12,25 mol/L.

Number 4 - Standard - Store between 15 and 30°C. After handling keep in refrigerator between 2 and 8°C to avoid evaporation. Contains: Creatinine 3 mg/dL

PRESENTATION

Reagent	Presentation 1	Presentation 2
Nº 1	200 mL	100 0 mL
Nº 2	50 mL	250 mL
Nº 3	10 mL	50 mL
Nº 4	10 mL	30 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, water-bath at 37°C, watches and stopwatches, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized in Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- The water level in water bath must by higher than the reagent level in test tubes.

6- The reagent Nº 1 is caustic, therefore you should avoid akin contact. The reagent Nº 3 must be manipulated with care, because it is irritating to skin and mucous membrane, can cause minor burns.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum or plasma collected with Heparin, Oxalate, Fluoride, Citrate or EDTA Bioclin. The analyte is stable for 7 days between 2 and 8°C.

Isolated urine: The urine should be collected and preserved in refrigerator until the time of dosage.

24 hours Urine: The urine should be collected within a period of 24 hours and preserved in refrigerator until the time of dosage.

PROCESS DESCRIPTION

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

TECHNIQUE I - END POINT

Mark 3 test tubes with: B (Blank), P (Standard) and A (Sample), and proceed as follows:

	Blank	Standard	Sample
Reagent Nº 1	2,0 mL	2,0 mL	2,0 mL
Distilled Water	250 µL	---	---
Sample	---	---	250 µL
Reagent Nº4	---	250 µL	---
Reagent Nº 2	500 µL	500 µL	500 µL

Homogenize and incubate in water-bath at 37°C for 10 minutes. Read the absorbance of the Sample and the Standard at 510 nm (500 - 540 nm), hitting the zero with the Blank. The absorbance of the Samples is A1 and the Standard is P.

And then add:

	Blank	Standard	Sample
Reagent Nº 3	100 µL	---	100 µL

Homogenize and wait for 5 minutes between 15 e 30°C. Read the absorbance A2 of the Sample at 510 nm (500 - 540 nm), hitting the zero with the Blank. The color of the reaction is stable for 30 minutes.

URINE DOSAGE

Take note of the collected volume in mL, centrifuge a aliquot of the urine sample to be dosed in proceed the dilution of 1:25 (0,1 mL of urine and 2,4 mL of distilled or deionized water). Then, making the determination as the technique described above, omitting the second part (acidification). Multiple the obtained results for 25.

CALCULATIONS

Creatinine (mg/dL) = $\frac{A1 - A2}{\text{Standard Absorbance}} \times 3$

According to the Lambert-Beer Law, the Calibration Factor can be used.

Calibration Factor = $\frac{\text{Standard Concentration (3 mg/dL)}}{\text{Standard Absorbance}}$

Creatinine (mg/dL) = (A1 - A2) x Calibraton Factor

Creatinine in Urine (mg/24h) = $\frac{\text{mg/dL} \times \text{Volume in 24h (mL)}}{100}$

TECHNIQUE II – FIXED POINT KINETIC METHOD

WORKING REAGENT PREPARATION

Mixture a share of reagent Nº 2 (Picric Acid) with four shares of reagent Nº 1 (Alkaline Reagent). Store in a amber flask. Stable for 24 hours between 15 and 30°C. Homogenized the Working Reagent before the technique begins.

REACTION

Add 0,1 mL of Standard or Sample (serum, plasma or urine) in 1 mL of Working Reagent (homogenized previously), homogenize and transfer immediately for thermostated cuvette at 37°C. Measure absorbance of Standard and Sample at 510 nm (500 - 540 nm) in 30 and 90 seconds.

CALCULATIONS

Absorbance from Standard = Abs. 90 seconds – Abs. 30 seconds
or Sample

Calibration Factor = $\frac{\text{Standard Concentration (3 mg/dL)}}{\text{Standard Absorbance}}$

Creatinine (mg/dL) = Sample Absorbance x Calibraton Factor

DEPURATION

Hydrate the patient with minimum 600 mL of water, and then, orientation to empty the bladder. The depuration of Creatinine must be carried out using the total volume of collected urine 4, 12 or 24 hours, according to standard procedure in the laboratory. In any moment of the chosen interval, collect a blood sample. Measure the total volume of the urine and calculate the volume per minute (VM).

VM = $\frac{\text{Determined time urinary volume (mL)}}{\text{Determined time (in minutes)}}$

Dose the serum Creatinine and urine according to the proposed methodologies.

Calculate the Creatinine Depuration:

U = Urine Creatinine (mg/dL)

S = Serum Creatinine (mg/dL)

VM = Volume/minute

A = Body surface value of the patient in m²

Creatinine depuration = $\frac{U \times VM}{S}$
(mL/min)

Attention: This value can be corrected to the patients body surface. Use a chronogram of weight-height correlation. Multiply the depuration value by 1,73 and divide by the patient body surface.

Creatinine depuration = $\left(\frac{U \times VM \times 1.73}{S \times A} \right)$
(mL/min./1,73m²)

PROCEDURE LIMITATIONS

The kinetic process can only be carried out in a spectrophotometer with thermostated cuvette at 37°C, due to the small time interval of the reaction, The methodologies described shall not be used in extremely lipemic samples.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 914, making it equivalent to the IDMS method (Isotope Dilution Mass Spectrometry), as recommended by the NKDEP (National Kidney Disease Education Program).

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of Creatinine in healthy populations of male and female.

Serum / Plasma	Adult	0,4 - 1,4 mg/dL
	Newborn	0,31 - 0,92 mg/dL
	2 weeks - 1 year	0,16 - 0,39 mg/dL
	1 - 3 years	0,17 - 0,35 mg/dL
	3 - 5 years	0,26 - 0,42 mg/dL
	5 - 7 years	0,29 - 0,48 mg/dL
	7 - 9 years	0,34 - 0,55 mg/dL
	9 - 11 years	0,32 - 0,64 mg/dL
	11 - 13 years	0,42 - 0,71 mg/dL
Urine	Adult (Male)	1500 - 2500 mg/24h or 21 - 26 mg/kg/24h
	Adult (Female)	800 - 1500 mg/24h or 16 - 22 mg/kg/24h
	Children 2 - 3 anos	6 - 22 mg/kg/24h
	Children > 3 anos	12 - 30 mg/kg/24h
Creatinine Clearance	Adult (Male)	97 - 137 mL/min/1,73m²
	Adult (Female)	88 - 128 mL/min/1,73m²
	Children	70 - 140 mL/min/1,73m²

mg/Kg weight = mg/24h divided by body weight
To convert the values of mg/dL in mmol/L (SI), multiplyby 0,0884.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and / or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY
Creatinine kit was compared with other commercially available methods for measurement of Creatinine. 07 analyzes were performed and the results were evaluated. The linear equation obtained was Y = 1,072X - 0,056, with correlation coefficient 0,999. With these results, we can conclude that the kit shows good methodological specificity.

Precision
REPEATABILITY
The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0,679	1,128	1,441
Standard Deviation (mg/dL)	0,008	0,008	0,008
Coefficient of Variation (%)	1,198	0,697	0,573

REPRODUCIBILITY
The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0,679	1,129	1,440
Standard Deviation (mg/dL)	0,001	0,001	0,001
Coefficient of Variation (%)	0,085	0,089	0,080

Sensitivity
The sensitivity was calculated from 20 determinations of the sample free of Creatinine. The average found was 0,0377 mg/dL, with standard deviation of 0,0006 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,0395 mg/dL.

Linearity
The color reaction is linear up to 10 mg/dL. For higher values, dilute the sample with Sodium Chloride 0,85% and repeat the determination. Multiply the obtained value by dilution factor.

DIAGNOSTIC SIGNIFICANCE
The Creatinine, synthesized mainly in the liver and in kidney, is produced and excrement in a constant rhythm, independently of diet, hydration level and protein metabolism. With the kidney function failure, the laboratorial data proportionate secure indexes as the kidney's ability to excrement, re-absorption and secrete.
The Creatinine is raised in the acute and chronic kidney insufficiency, in the urine tract obstruction, congestive heart insufficiency, dehydration, chock and mellitus diabetes.
Lower levels can be observed on: muscle dystrophy, reduction by innutrition of the muscle mass, severe kidney disease.

NUMBER OF TESTS
K016-1 ... 100 Tests / 250 µL of sample / 2 mL of reagent
200 Tests / 125 µL of sample / 1 mL of reagent

K016-2 ... 500 Tests / 250 µL os sample / 2 mL of reagent
1000 Tests / 125 µL of sample / 1 mL of reagent

BIBLIOGRAPHIC REFERENCES
1 - JAFFE, M.: J. Phisiol. Chem. 10:381, 1886.
2 - SLOT, C.: Sand. J. Clin. Lab. Invest. 17:381, 1965.
3 - HENRY, R. J. -: Clinical Chemistry - Principles and Technics, Ed.,Harper and Row. New York, 1964.
4 - TODD, SANFORD, DAVIDSOHN: Diagnósticos Clínicos, 16 Ed.,Editora Manole, 1983.
5 - CARL, A. B. and EDWARD, R.A.: Tietz Textbook of Clinical Chem 2 nd ed., 1994, 1531 - 1539.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration number for Creatinine kit: 10269360089

Review: January/2018

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 BATCH CODE

 DATE OF MANUFACTURE

 USED BY
(last day of month)

 TEMPERATURE LIMITATION
(store at)

 CONTAINS SUFFICIENT
FOR <N> TESTS

 CONSULT INSTRUCTIONS
FOR USE

 IN VITRO DIAGNOSTIC DEVICE

 EUROPEAN AUTHORIZED
REPRESENTATIVE

 KEEP AWAY
FROM SUNLIGHT

 MANUFACTURED BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 INFLAMMABLE

 CORROSIVE

 POISON

 CE MARK

 DO NOT USE IF
PACKAGE IS
DAMAGED