CREATININA ENZIMÁTICA

REF K161

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FINALIDADE

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REF K161
USAGE INSTRUCTIONS


FUNCTION

Method for determination of Creatinine in human samples of serum, plasma and urine. Enzymatic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric Enzymatic
The creatinine present in the sample is converted into Creatine through the action of the enzyme Amidohydrolase Creatinine. The creatine is hydrolyzed by the action of the enzyme Creatine Amidohydrolase into Sarcosine and Urea. The enzyme Sarcosine Oxidase produces Glycine, Formaldehyde and Hydrogen Peroxide, through the oxidative demethylation of the Sarcosine. The Hydrogen Peroxide is then reacted with TOPS and 4-aminoantipyrine by Peroxidase to produce Quinoneimine. The maximum absorbance of this formation is 546nm and the intensity of color of the product is proportional to the concentration of creatinine in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Buffer < 100 mmol/L, Ascorbate Oxidase < 10 KU/L, Sarcosine Oxidase < 20 KU/L, Creatine Amidinohydrolase < 60 KU/L, TOPS < 30 g/L, Catalase < 400 KU/L, chelator, surfactant, stabilizer and preservative.
Number 2 - Reagent color - Store between 2 and 8°C. Contains: Creatinine Amidohydrolase < 600 KU/L Peroxidase < 100 KU/L 4-Aminoantipyrine < 2 g/L, buffer < 100 mmol/L, surfactant and preservative.
Number 3 - Standard - Store between 2 and 8°C. Contains: Creatinine 3,0 mg/dL, stabilizer and preservative.

PRESENTATION

Presentation	Reagent N°1	Reagent N°2	Reagent N°3
1	2 x 54 mL	2 x 18 mL	1 x 3 mL
2	2 x 30 mL	2 x 10 mL	1 x 3 mL
3	3 x 30 mL	3 x 10 mL	1 x 3 mL
4	4 x 30 mL	4 x 10 mL	1 x 3 mL
5	5 x 30 mL	5 x 10 mL	1 x 3 mL
6	1 x 60 mL	1 x 20 mL	1 x 3 mL
7	2 x 60 mL	2 x 20 mL	1 x 3 mL
8	3 x 60 mL	3 x 20 mL	1 x 3 mL
9	4 x 60 mL	4 x 20 mL	1 x 3 mL
10	5 x 60 mL	5 x 20 mL	1 x 3 mL
11	1 x 90 mL	1 x 30 mL	1 x 3 mL
12	2 x 90 mL	2 x 30 mL	1 x 3 mL
13	3 x 90 mL	3 x 30 mL	1 x 3 mL
14	4 x 90 mL	4 x 30 mL	1 x 3 mL
15	5 x 90 mL	5 x 30 mL	1 x 3 mL
16	1 x 54 mL	1 x 18 mL	1 x 3 mL
17	1 x 30 mL	1 x 10 mL	1 x 3 mL
18	3 x 54 mL	3 x 18 mL	1 x 3 mL
19	4 x 54 mL	4 x 18 mL	1 x 3 mL
20	5 x 54 mL	5 x 18 mL	1 x 3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, stopwatches, test tubes, Biocal Bioclin, Biocontrol N and Biocontrol P Bioclin. These itens are found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature shall be between 2 and 8°C and transport temperatures shall be up to 30°C, not exceeding 5 days. **Do not freeze.** Protect from light and avoid moisture.

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
2- Strictly follow the methodology proposed to obtain exact results.
3- Water used in material cleaning must to be recent and free of contaminants.
4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

- 6- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
7- Do not use the product in case of damaged packaging.
8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma collected with Heparin, Oxalate, Fluoride, Citrate or EDTA Bioclin. The compound is stable for 7 days between 2 and 8°C, and 3 months at -20°C.⁵
The urine is stable for 6 days between 2 and 8°C, and 6 months at -20°C.⁵
Isolated urine: The urine should be collected and preserved in refrigerator until the time of dosage.
24 hours Urine: The urine should be collected within a period of 24 hours and preserved in refrigerator until the time of dosage.

PROCESS DESCRIPTION

The calibration stability of the Creatinine Enzymatic kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends usign the Kit Biocal Bioclin as calibrator and as control serum, Biocontrol N and P Bioclin Kits.

This kit is for use in automated biochemical analyzers only.
Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.

CALCULATIONS

24 HOURS URINE

Creatinine in Urine (mg/24h) = $\frac{\text{mg/dL} \times \text{Volume in 24h (mL)}}{100}$

CLEARANCE

Hydrate the patient, with at least 600 mL of water and then guide you to empty the bladder. Creatinine clearance should be performed on a sample taken at 4, 12 or 24 hours. In any time of the chosen time interval, a sample of blood. Measure the total volume of urine and calculate the volume per minute (VM).

VM = $\frac{\text{Urinary volume in determined period of time (mL)}}{\text{Determined time (in minutes)}}$

Dose the serum and urine creatinine following the methodologies proposed.

To calculate creatinine clearance:
U = Urinary Creatinine (mg/dL)
S = Serum Creatinine (mg/dL)
VM = Volume/minute
A = Value of the patient's body surface area in m²

Creatine Clearance = $\frac{U}{S} \times VM$
(mL/min)

Note: This value should be corrected for patient's body surface. Using the nomogram of height-weight correlation. Multiply the value of clearance by 1,73 and divide by the surface of patient's body.

Creatine Clearance = $\frac{U}{S} \times VM \times \frac{1,73}{A}$
(mL/min/1,73 m²)

Example: Clearance (24-hour urine)
U = Urinary Creatinine = 110 mg/dL
S = Serum Creatinine = 0,9 mg/dL
24-hour urine volume = 1,480 mL in 24 hours (1,440 minutes)

Volume/minute = VM = $\frac{1,480}{1,440} = 1,0277$

Patient's weight = 60 kg
Patient height = 165 cm
Body surface area = 1,66 m²

Clearance = $\frac{110}{0,9} \times 1,0277 = 125,60$
(mL/min)

Clearance = $\frac{110}{0,9} \times 1,0277 \times \frac{1,73}{1,66} = 130,889$
(mL/min./1,73 m²)

INTERFERING

Triglycerides concentrations up to 1000 mg/dL, bilirubin (free and conjugated) up to 16 mg/dL, Hemoglobin up to 400 mg/dL, Ascorbic Acid up to 40 mg/dL and Creatine up to 20 mg/dL does not interfere in the reaction.
Samples with Bilirubin, Haemoglobin and triglycerides at levels greater than the above should be diluted in NaCl 150 mmol/L (0.85%) before testing. The result obtained should be multiplied by the dilution factor.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard and the BIOCAL calibrator are traceable to the reference material NIST (National Institute of Standards and Technology) SRM 914, making it equivalent to the IDMS method (Isotope Dilution Mass Spectrometry), as recommended by the NKDEP (National Kidney Disease Education Program).

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of creatinine in healthy populations of male and female.

Serum/Plasma	Adult	0,4 - 1,4 mg/dL
	Newborn	0,31 - 0,92 mg/dL
	2 weeks - 1 year	0,16 - 0,39 mg/dL
	1 - 3 years	0,17 - 0,35 mg/dL
	3 - 5 years	0,26 - 0,42 mg/dL
	5 - 7 years	0,29 - 0,48 mg/dL
	7 - 9 years	0,34 - 0,55 mg/dL
	9 - 11 years	0,32 - 0,64 mg/dL
	11 - 13 years	0,42 - 0,71 mg/dL
Urine	Adult (Male)	1500 - 2500 mg/24h or 21 - 26 mg/kg/24h
	Adult (Female)	800 - 1500 mg/24h or 16 - 22 mg/kg/24h
	Children 2 - 3 anos	6 - 22 mg/kg/24h
	Children > 3 anos	12 - 30 mg/kg/24h
Creatinine clearance	Adult (Male)	97 - 137 mL/min/1,73m ²
	Adult (Female)	88 - 128 mL/min/1,73m ²
	Children	70 - 140 mL/min/1,73m ²

mg/kg weight = mg/24h divided by body weight

To convert the values to mg/dL to mmol/L (SI), multiply by 0,0884.

These values should be used as guidelines, each laboratory should create their own range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS
The Creatinine Enzymatic kit was compared with another method commercially available for measurement of Creatinine. 42 analyses were performed and the results were evaluated. The linear equation was obtained $Y = 0,944X + 0,225$ and correlation coefficient 0,999. With these results, we can conclude that the kit shows good methodology specificity.

Precision
REPEATABILITY
The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Samples		
	1	2	3
Average Concentration (mg/dL)	0,80	1,45	5,66
Standard Deviation (mg/dL)	0,01	0,01	0,03
Variation Coefficient (%)	1,39	0,66	0,45

REPRODUCIBILITY
The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Samples		
	1	2	3
Average Concentration	0,78	1,46	5,67
Standard Deviation	0,04	0,01	0,01
Variation Coefficient (%)	4,59	0,86	0,12

Sensitivity
Sensitivity was calculated from 40 determinations of a sampel free of Creatinine. The average found was 0,035 mg/dL, with a standard deviation 0,005 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 0,050 mg/dL.

Linearity
The reaction color is linear to the concentration of 150 mg/dL. For higher values dilute the sample with sodium chloride 0,85% and repeat the determination. Multiply the value obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE
Creatinine is synthesized in the liver and kidneys, is produced and excreted at a constant rate independent of diet, degree of hydration and protein metabolism. With renal failure, the laboratory data provide insurance rates and the ability kidneys to excrete, secrete and reabsorb creatinine. Creatinine is elevated in acute and chronic renal failure in urinary tract obstruction, congestive heart failure in dehydration, shock and in diabetes *mellitus*. Decreased levels may be observed in muscular dystrophies, malnutrition, reduced muscle mass and severe renal disease.

BIBLIOGRAPHIC REFERENCES
1- SLOT, C.; Sand. J. Clin. Lab. Invest. 17:381, 1965.
2- HENRY, R. J.: Clinical Chemistry - Principles and Technics, 4 Ed., Harper and Row. New York, 1964.
3- TODD,SANFORD,DAVIDSOHN;Diagnósticos Clínicos, 16 Ed., Editora Manole, 1983.
4- CARL, A. B. and EDWARD, R. A.: Tietz Textbook of Clinical Chem 2nd ed., 1994, 1531-1539.
5 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC	REP
----	-----

OBELIS S.A.
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454.
E-mail: sac@bioclin.com.br

ANVISA registration for Creatinine Enzymatic kit: 10269360260

Review: March/2020

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED