DENGUE NS1

REF K130

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FINALIDADE

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1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.
9. Dussart P, et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
10. Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
11. Bioclin - Dados de Arquivo.

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REFERENCIAS BIBLIOGRÁFICAS
1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48 .
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.
9. Dussart P., et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
10. Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
11. Bioclin - Dados de Arquivo.

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USAGE INSTRUCTIONS

FUNCTION

Screening immunochromatographic assay for the qualitative determination of NS1 antigen for Dengue virus in whole blood, serum or plasma.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
The Dengue NS1 kit is an immunochromatographic test for the detection of NS1 antigen of Dengue virus. The kit has a “T” test line with immobilized anti-NS1 antibodies and a “C” control line with anti-goat antibodies. When a sample is added, it interacts with anti-NS1 antibodies conjugated to colloidal gold particles. Next, the sample migrates by capillarity through the membrane, and interacts with anti-NS1 antibodies present in the “T” line. If the sample has NS1 antigen, a colored band appears in the “T” test region, indicating a reagent result. The sample continues migrating through the membrane until it reaches the “C” control region, where it interacts with anti-goat antibodies, generating a colored line that indicates the correct functioning of the product.

REAGENTS

Reagent Number 1 (R1) Cassette - Store between 15 and 30°C. **Do not freeze.** Each sachet contains:

- A)** 01 plastic cassette with 01 test strip containing:
- Control region: rabbit anti-goat antibody IgG;
 - Test region T: Goat anti-NS1 antibody;
 - Conjugate: Anti-NS1 goat antibody conjugated to colloidal gold particles;

- B)** 01 sachet of silica

Comes with disposable plastic pipette.

PRESENTATION

Presentation	Cassette per packing
1	1
2	10
3	20
4	25
5	40
6	50
7	80
8	100

EQUIPMENTS AND OPERATIONAL INPUTS

Watch or stopwatch. Can be found in specialized markets for articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The temperature of storage and transportation must be 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain accurate results.

3- The cassettes do not show contaminants, but how the samples are biological, the safety routines should be followed accurately and the use of disposable gloves is essential.

- 4- Use clean and dry containers for storage of samples.
- 5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 6- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 7- Do not use the product in case of damaged packaging.
- 8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

We recommend the sample collection between the 1st and 4th day of the first symptoms to carry out the test.

Whole Blood: Whole blood should be collected with EDTA. Do not use fingerstick. The samples should be stored between 2 and 8°C and tested within 24 hours after collection.

Serum / Plasma: Must be separated from the blood interface immediately after the centrifugation, that must occur after the sample collection. The serum or plasma can be stored for up to 3 days when kept between 2 and 8°C and for 06 months at -20°C.

DESCRIPTION OF PROCESS


- 1- The sample should be in room temperature before starting the test.
- 2- Remove the cassette from its protective packaging and appropriately identify it so.
- 3- Transfer 3 drops, using a disposable plastic pipette that came with the kit, or 100 µL of Whole Blood, Serum or Plasma, within the sample window.

- 4- Wait for the formation of lines after resting for 15 to 20 minutes. **Do not interpret it after 20 minutes.**

INTERPRETATION OF RESULT

Reagent Test: Formation of two red lines, one in the test region (T) and one in the control region (C) in the first 15 minutes. Any intensity of the red line in the test region (T), in the presence of the control line, represents a reagent result. That is, the presence of a weak test line, already identifies a reagent result. Do not interpret after 20 minutes.

Test Non-reagent: Formation of one red line in the control region (C) with no red line in the test region (T). Do not interpret after 20 minutes.

Inadequate Test: Complete absence of line in the control region (C) with or without a red line in the test region (T). Test the sample again.

Notes:

- 1- The line on the control region can appear before 15 minutes incubation test, it does not indicate that non-reagent results can be interpreted before the time.
- 2- The results must not be interpreted after 20 minutes to avoid incorrect results.
- 3- The presence of red line in the test region (T) in any intensity, weak or strong, when the control line is also present, identifies a reagent result.

INTERFERENCES

Cross reaction with other virus from the flavivirus could be common, according to bibliographic references, however were not observed during the development tests of this product.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PROCEDURE LIMITATIONS

- 1- The Dengue NS1 kit is a qualitative test and must be used only for *in vitro* diagnostic.
- 2- The Dengue NS1 kit is a screening test to identify the presence of NS1 antigen from the Dengue virus.
- 3- As every diagnostic test, the results must be interpreted along with other clinical data of the patient.
- 4- The isolated result of the test can not be used for a definitive diagnostic.
- 5- If the non-reagent result continues and the clinical symptoms persist, an additional test should be carried out using a different methodology. A non-reagent result does not exclude, in any circumstance, the possibility of infection by the Dengue virus.
- 6- Reagent results can be confirmed with another test and available clinical data.

7- The fingerstick is not recommended because the amount of blood may be insufficient or clot, hindering the correct performance of the test.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS, METHODOLOGY AND SPECIFICITY
To establish the sensitivity and specificity of the Dengue NS1 kit accordingly with another Dengue test, 309 known clinical samples were studied. The clinical sensitivity found was of 94,39% were only 6 from the 107 reagent samples were not identified by both methods. The methodological sensitivity was of 100%. The specificity found was of 99,00%.

Precision
REPEATABILITY
The repeatability was calculated from 20 successive determinations, using 5 different samples, obtaining the following results:

Samples	Nº of repetitions	Expected Result	Found Result
Human Blood 01	20	100% Non-reagent	100% Non-reagent
Human Blood 02	20	100% Reagent	100% Reagent
Human Serum 01	20	100% Non-reagent	100% Non-reagent
Human Serum 02	20	100% Reagent	100% Reagent
Human Plasma 01	20	100% Non-reagent	100% Non-reagent

REPRODUCIBILITY
The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Samples	Nº of Days	Nº of Repetitions	Expected Result	Found Result
Human Blood 01	03	20	100% Reagent	100% Reagent
Human Blood 02	03	20	100% Non-reagent	100% Non-reagent
Human Serum 01	03	20	100% Reagent	100% Reagent
Human Serum 02	03	20	100% Non-reagent	100% Non-reagent
Human Plasma 01	03	20	100% Reagent	100% Reagent

DIAGNOSTIC SIGNIFICANCE
The Dengue virus is a Flavivirus transmitted by the mosquito *Aedes aegypti*. It is distributed by the tropical and subtropical areas of the world, causing more than 100 millions of infections annually. The classic Dengue infection is characterized by fever (normally between 38 and 40°C) out of the suddenly, unwellness,

anorexia (small appetite), headache, muscular and eye sore. Hemorrhagic Dengue can cause gum hemorrhage and epistaxis, internal bleeding and disseminated intravenous coagulation, with many organs damage, leading to death. The NS1 is a glycoprotein highly conserved present in high concentrations in the serum of patients infected with Dengue. The Ns1 antigen is found from the 1st until the 9th days after the appearance of the fever in patients samples with primary and secondary infections. However, its concentration decreases over days of infection, significantly reducing the chances of detection of antigen. Articles report 50% less sensitivity for detection after the 4th day of infection by several diagnostic methods. The primary dengue infection causes an increase in the IgM antibodies after 3 to 5 days from the beginning of the fever. The IgM antibodies normally remain in the system from 30 to 90 days. Patients in endemic areas can present secondary infections, which result in elevated levels of the IgG antibodies, isolated or simultaneously with the IgM response. In case of secondary infection, the NS1 antigen may not be detected due to the possibility of lower viremia.

BIBLIOGRAPHIC REFERENCES
1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.
9. Dussart P, et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
10. Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
11. Bioclin - Dados de Arquivo.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC

REF

OBELIS S.A.

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1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for Dengue NS1 kit: 10269360203

Review: June/2019

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

EC

REP

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED