FERRO SÉRICO

REF K017

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Iron Serum. Colorimetric test, only for *in vitro* diagnostic.

PRINCIPLE OF ACTION

Methodology: Modified Goodwin
Iron is released from Transferrin in acidic and low to its ferrous state by the action of Hydroxylamine. Subsequently, reacts with Ferrozine leading to the formation of a violet color complex.

REAGENTS

- Number 1 - Reduction Buffer** - Store between 15 and 30°C. Contains: Succinic Acid < 500 mmol/L, Hydroxylamine Hydrochloride < 300 mmol/L, surfactant and stabilizer.
- Number 2 - Color Reagent** - Store between 15 and 30°C. Contains: Ferrozine < 14 mmol/L and Glacial Acetic Acid < 2 mol/L.
- Number 3 - Standard** - Store between 15 and 30°C. Contains: Iron 100 µg/dL, Hydroxylamine Hydrochloride < 350 mmol/L and preservative.

PRESENTATION

Reagent	Volumen
Nº 1	40 mL
Nº 2	2,5 mL
Nº 3	5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch and stopwatch, pipettes test tubes, Biocontrol N and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature shall be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- 2-** Strictly follow the methodology proposed to obtain exact results.
- 3-** Water used in material cleaning must to be recent and free of contaminants.
- 4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5-** Water level in water bath should be above the level of reagents in test tubes.
- 6-** Every material used for analyses must be previously emerged in a recent HCL (7%) solution for 30 minutes. Rinse it with deionized water.

- 7-** Hemolysis, even if discrete, can interfere with dosage.
- 8-** Do not use plasma.
- 9-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 10-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.
- 11-** Do not use the product in case of damaged packaging.
- 12-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Sample obtained free of hemolysis, stable for 07 days between 2 and 8°C and 1 year at -20°C ⁶. Therapeutic control samples should always be collected at the same time.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes with letters B (Blank), P (Standard), A (Sample) and proceed as follows:

	Blank	Standard	Sample
Reagent Nº 1	2,0 mL	2,0 mL	2,0 mL
Distilled Water	500 µL	---	---
Sample	---	---	500 µL
Reagent Nº 3	---	500 µL	---

Smoothly homogenize and determine the absorbance of the sample at 560 nm (540 - 580 nm), hitting the zero with the Blank. This will be the absorbance **A₁**.

And then add:

	Blank	Standard	Sample
Reagent Nº 2	1 drop	1 drop	1 drop

Smoothly homogenize and incubate in water-bath at 37°C, per 10 minutes. Determine the absorbancies of Sample and the Standard in 560 nm (540 - 580 nm), matching the zero with the Blank. This will be the absorbance of the sample **A₂**.

CALCULATIONS

Iron Serum (µg/dL) =
$$\frac{(A_2 - A_1)}{\text{Standard Abs.}} \times 100$$

The color reaction follows the Lambert Beer Law, it can use the Calibration Factor:

Calibration Factor =
$$\frac{100}{\text{Standard Absorbance}}$$

Iron Serum (µg/dL) = (A₂ - A₁) x Calibration Factor

The results are expressed in µg/dL.

PROCEDURE LIMITATIONS

The use of water contaminated by Iron increase Blank color intensity, changing the results.

INTERFERENCES

Some drugs can raise levels of Serum Iron such as: Chloramphenicol, Estrogens, Ethanol, Methotrexate and Anovulatory. Other drugs may decrease the levels of Serum Iron such as: Acetylsalicylic Acid (high doses), Allopurinol, Cortisone and Corticotrophin.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 937.

REFERENCE VALUES

The reference values in µg/dL for this method were obtained through the determination of Iron in healthy populations of male and female.

Iron Serum (µg/dL)		
Newborn		100 - 250
Infant		40 - 100
Children		50 - 120
Adult	Male	65 - 170
	Female	50 - 170

To convert the values of $\mu\text{g/dL}$ to $\mu\text{mol/L}$, multiply by 0,179.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
RECOVERY

The recovery analysis was made with 5 sample determinations. The accurate was calculated and was found in a good parity with the reference values, obtaining a recovery between 94% and 104%.

COMPARISON OF METHODSAND METHODOLOGICAL SPECIFICITY

Iron Serum kit was compared with other commercially available methods for measurement of Iron. 07 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1,011X - 1,534$, and the correlation coefficient 0,999. With theses results, we can conclude that the kit shows good methodological specificity.

Precision
REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration ($\mu\text{g/dL}$)	57,95	85,00	192,95
Standard Deviation ($\mu\text{g/dL}$)	0,83	0,79	0,89
Coefficient of Variation (%)	1,42	0,93	0,46

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration ($\mu\text{g/dL}$)	57,92	85,05	193,03
Standard Deviation ($\mu\text{g/dL}$)	0,25	0,05	0,08
Coefficient of Variation (%)	0,43	0,06	0,04

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Iron. The average found was 1,1506 $\mu\text{g/dL}$ and the standard deviation of 0,0005 $\mu\text{g/dL}$. The sensitivity, that indicates the methods detection limit, corresponds the average plus 3 times the standard deviation and is equal to 1,1521 $\mu\text{g/dL}$.

Linearity

The reaction is linear up to 600 $\mu\text{g/dL}$. For larger values than 600 $\mu\text{g/dL}$, dilute the serum with distilled water, repeat dosage and multiply the result by a factor of dilution.

DIAGNOSTIC SIGNIFICANCE

Iron is the most important transition metal and abundant in the human body. It is essential to the synthesis of hemoglobin and myoglobin, besides being present in several biochemical reactions of metabolism.

The serum iron is elevated in sideroblastic anemia, hemolytic anemias, hereditary anemias (Thalassemia Major), acute hepatitis, chronic liver diseases, hemochromatosis and massive blood transfusions.

Low serum iron levels are found in chronic infection diseases, neoplasms, parasites, menstrual period, large bleeding, pregnancy and dietary deficiencies.

NUMBER OF TESTS

20 Tests / 500 μL of Sample / 2mL of Reagent
40 Tests / 250 μL of Sample / 1mL of Reagent

BIBLIOGRAPHIC REFERENCES

1 - PETERS T., et coll.: J. Lab and Clin Med. 48:280,1956.
2 - HENRY, R. J., Cannon. D. C. and Winkleman, J. W., eds.: Clinical Chemistry: Principles and Technics, Haper and Row Publishers, 1974.
3 - ZAK, B.: Clin Chim. Acta 3: 328, 1958.
4 - BRITTENHAM, G.: Clin Chim. Acta. 91: 203, 1979.
5 - TONKS, D. B.: Chim. Acta 2: 393, 1983.
6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
7 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
email bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

English 2/2

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
Email: sac@bioclin.com.br

ANVISA registration for the Iron Serum kit: 10269360110

Review: February/2020

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED