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REF K070

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Iron Serum. Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Modified Goodwin
Iron is released from Transferrin in an acid enviroment and lowered to its ferrous state by the action of Hydroxylamine. Subsequently, reacts with Ferrozine leading to the formation of a violet colored complex.

REAGENTS

- Number 1 - Reduction Buffer** - Store between 15 and 30°C. Contains: Succinic Acid 250 mmol/L Hydroxylamine 300 mmol/L, surfactant and stabilizer.
Number 2 - Color Reagent - Store between 15 and 30°C. Contains: Ferrozine 12 mmol/L.
Number 3 - Standard - Store between 15 and 30°C. Contains: Iron 100 µg/dL, stabilizer.

PRESENTATION

Reagent	Volume
Nº 1	80 mL
Nº 2	2 mL
Nº 3	3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The transport and storage temperatures must not exceed 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
2- Strictly follow the methodology proposed to obtain exact results.
3- Water used in material cleaning must to be recent and free of contaminants.
4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
5- Hemolysis, even if discrete, can interfere with dosage.
6- Do not use plasma.
7- Every material used for analyses must be previously emerged in a recent HCL (7%) solution for 30 minutes. Rinse it with deionized water.
8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

- 10-** Do not use the product in case of damaged packaging.
11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, stable for 07 days if kept in temperature between 2 and 8°C and 1 year at -20°C ⁷. The samples for therapeutic control must be collected always at the same time.

PROCESS DESCRIPTION

WORKING COLOR REAGENT PREPARATION

The Reagent N° 2 (Color Reagent) is not ready for use. The Working Color Reagent should be prepared by mixing 2 mL of Reagent N° 1 with 0,5 mL of Reagent N° 2, or proportional amounts.
The Working Color Reagent, after preparation, is stable for 30 days when stored at temperatures between 15 and 30°C and protected from light.
To use the kit in automated equipment, the R1 is Reduction Buffer and the R2 is Working Color Reagent prepared as described above.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

PROCEDURE LIMITATIONS

The use of water contaminated by Iron increase white color intensity, changing the results. Some drugs can raise levels of Serum Iron such as Chloramphenicol, Estrogens, Ethanol, Methotrexate and Anovulatory. Other drugs may decrease the levels of Serum Iron such as: Acetylsalicylic Acid (high doses), Allopurinol, Cortisone and Corticotrophin.
Samples with high Fibrin content may show falsely elevated results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 937.

REFERENCE VALUES

The reference values in µg/dL, for this method, were obtained through the determination of Iron in healthy populations of male and female.

Iron Serum (µg/dL)		
Newborn		100 - 250
Infant		40 - 100
Children		50 - 120
Adult	Male	65 - 170
	Female	50 - 170

To convert the values from µg/dL to µmol/L, multiply by 0,179.
These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.
The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
RECOVERY

The recovery analysis was performed with 05 determinations samples. Accuracies were calculated and found in good agreement with the reference values, obtaining a recovery between 95% and 104%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Iron Serum Automated kit was compared with other method for determination of Iron dosage commercially available. 07 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 0,997X + 0,095$ and correlation coefficient 0,999. With these results we can conclude that the kit shows good methodological specificity.

Accuracy

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration Promedio ($\mu\text{g/dL}$)	66,99	113,30	196,00
Standard Deviation ($\mu\text{g/dL}$)	1,45	1,66	1,30
Coefficiente de Variación (%)	2,16	1,46	0,66

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration Promedio ($\mu\text{g/dL}$)	66,70	113,62	196,25
Standard Deviation ($\mu\text{g/dL}$)	0,23	0,39	0,35
Coefficiente de Variación (%)	0,34	0,34	0,18

Sensitivity

The sensitivity was calculated from 20 determinations from a sample free of iron. The average found was $3,3 \mu\text{g/dL}$ with standard deviation of $0,1 \mu\text{g/dL}$. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to $3,6 \mu\text{g/dL}$.

Linearity

The reaction is linear up to $600 \mu\text{g/dL}$. For larger values than $600 \mu\text{g/dL}$, dilute the serum with distilled water, repeat dosage and multiply the result by a factor of dilution.

DIAGNOSTIC SIGNIFICANCE

Iron is the most important transition metal and abundant in the human body. It is essential to the synthesis of hemoglobin and myoglobin, besides being present in several biochemical reactions of metabolism.

The Serum Iron is elevated in sideroblastic anemia, hemolytic anemias, hereditary anemias (Thalassemia Major), acute hepatitis, chronic liver diseases, hemochromatosis and massive blood transfusions. Low Serum Iron levels are found in chronic infection diseases, neoplasms, parasites, menstrual period, large bleeding, pregnancy and dietary deficiencies.

BIBLIOGRAPHIC REFERENCES

- 1 - PETERS T., et coll.: J. Lab and Clin Med. 48:280,1956.
- 2 - HENRY, R. J., Cannon. D. C. and Winkelman, J. W., eds.: Clinical Chemistry: Principles and Technics, Haper and Row Publishers, 1974.
- 3 - ZAK, B.: Clin Chim. Acta 3: 328, 1958.
- 4 - BRITTENHAM, G.: Clin Chim. Acta. 91: 203, 1979.
- 5 - TONKS, D. B.: Chim. Acta 2: 393, 1983.
- 6 - BURTIS, C.A., ASHWOOD, E.R., Tietz Textbook of Clinical Chemistry, 2ª Edição, Philadelphia: W.B. Saunders, 1986:2175-2211.
- 7 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Iron Serum Automated kit:
10269360121

Review: August/2018

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED