FOSFATASE ALCALINA

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FINALIDADE

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of Alkaline Phosphatase. Colorimetric test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

**Methodology:** Modified Roy.

The Alkaline Phosphatase present in the sample hydrolyzes Substrate Thymolphthalein Monophosphate, releasing Thymolphthalein and Inorganic Phosphate in a predetermined time of 10 minutes. The enzyme activity is proportional to the amount Thymolphthalein formed, which is measured colorimetrically. With the addition of Color Reagent, the enzymatic reaction is inactivated and acquires color Thymolphthalein blue. As the Substrate and the reaction product present different colors, in an alkaline environment, the final color is established by a mixture of colors.

REAGENTS

**Number 1 - Substrate** - Store between 15 and 30°C. After it is used store between 2 and 8°C in well sealed flasks to avoid evaporation. Reagent can present turbidity or precipitation, fact that does not interfere in its quality. Agitate before using it. Contains: Thymolphthalein Monophosphate 40 mmol/L, solubilizer and stabilizer.

**Number 2 - Buffer** - Store between 15 and 30°C. After its use seal the flask to avoid contamination with CO<sub>2</sub> from the air. Diethanolamine 0,3 mol/L pH 9,9, Sodium Citrate 10 mmol/L surfactant and activator.

**Number 3 - Color Reagent** - Store between 15 and 30°C. Contains: Sodium Carbonate 150 mmol/L and Sodium Hydroxide 100 mmol/L.

**Number 4 - Standard 40 U/L** - Store between 15 and 30°C. After it is used store between 2 and 8°C in well sealed flasks to avoid evaporation. Contains: Thymolphthalein 0,4 mmol/L (40U/L) and solubilizer.

PRESENTATION

| Reagents     | Volume |
|--------------|--------|
| Reagent Nº 1 | 5 mL   |
| Reagent Nº 2 | 50 mL  |
| Reagent Nº 3 | 200 mL |
| Reagent Nº 4 | 3 mL   |

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter with reading at 340 nm with a bath at 37°C, watch or stopwatch, pipettes,

test tubes, Calibrator, Bicontrol N and Biocontrol P Bicolin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The transport at temperatures between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must to be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Water level at water bath must be superior to the level of the reagents in the test tubes.
- 6- Reagent Nº 3 must be handled carefully, for it presents corrosive action.
- 7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Advisory Service Customer) of Quibasa.
- 9- Do not use the product in case of damaged packaging.

PROCESS DESCRIPTION

SAMPLES

Serum or plasma (heparin) obtained free of hemolysis. The analyte is stable for 07 days if kept between 2 and 8°C and 30 days if kept at -10°C.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark three test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

|                                                         | Blank  | Standard | Sample |
|---------------------------------------------------------|--------|----------|--------|
| Reagent Nº1                                             | 50 µL  | 50 µL    | 50 µL  |
| Reagent Nº2                                             | 500 µL | 500 µL   | 500 µL |
| Reagent Nº4                                             | --     | 50 µL    | --     |
| Incubate at 37°C for 2 minutes.                         |        |          |        |
| Sample                                                  | --     | --       | 50 µL  |
| Homogenize and Incubate at 37°C for exactly 10 minutes. |        |          |        |
| Reagent Nº3                                             | 2,0 mL | 2,0 mL   | 2,0 mL |

Homogenize well and determine the absorbance of sample and Standard at 578 nm (570-590), correcting the zero with the blank. The color is stable for 60 minutes.

CALCULATIONS

Alkaline Phosphatase =  $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 40$  (U/L)

As the reaction follows the Beer-Lambert Law, the calibration factor can be used.

Calibration =  $\frac{\text{Standard Concentration (40 U/L)}}{\text{Factor}} \times \text{Standard Absorbance}$

U/L = Sample Absorbance x Calibration Factor

Results are expressed in U/L.

PROCEDURE LIMITATIONS

This methodology can not be adapted for use in automation equipment.

INTERFERENCES

Do not use hemolyzed or lipemic samples, because such samples can falsely elevate the results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

The reference values in U/L, for this method, were obtained through the determination of Alkaline Phosphatase in healthy populations of male and female.

Adults: 12 to 43 U/L

Children up to 12 years of age: 56 to 156 U/L

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

## PRODUCT PERFORMANCE

### QUALITY CONTROL

#### Accuracy

#### COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Alkaline Phosphatase was compared with other commercially available methods of measurement of phosphatase. For this test, We analyzed seven (07) clinical samples. The obtained results allow the construction of a linear equation  $Y = 1,009X - 0,117$ , with correlation coefficient equal 0,999. With these results, we can conclude that the kit presents good methodological specificity.

#### Precision

#### REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

|                                     | Sample 1 | Sample 2 | Sample 3 |
|-------------------------------------|----------|----------|----------|
| <b>Average Concentration (U/L)</b>  | 17,35    | 41,55    | 101,60   |
| <b>Standard Deviation (U/L)</b>     | 0,75     | 0,76     | 0,75     |
| <b>Coefficient of Variation (%)</b> | 4,29     | 1,83     | 0,74     |

#### REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

|                                     | Sample 1 | Sample 2 | Sample 3 |
|-------------------------------------|----------|----------|----------|
| <b>Average Concentration (U/L)</b>  | 17,60    | 41,38    | 101,32   |
| <b>Standard Deviation (U/L)</b>     | 0,22     | 0,14     | 0,25     |
| <b>Coefficient of Variation (%)</b> | 1,24     | 0,35     | 0,24     |

#### Sensitivity

The sensitivity indicates the detection limit of the method. It was calculated from twenty (20) determinations of phosphatase, a sample absent of analytical. An average value of 0,419 U/L, with standard deviation of 0,004 U/L. The sensitivity corresponds to the sum of the average found with three times the standard deviations: 0,429 U/L.

#### Linearity

The reaction is linear up to 500 U/L. For samples with higher values dilute the serum using Sodium Chloride 0,85%, repeat the dosage and multiply the result by the dilution factor.

#### DIAGNOSTIC SIGNIFICANCE

The Alkaline Phosphatase is present in high concentrations in bones, liver, intestine and placenta. Constitutes a useful indicator of liver and bone diseases associated osteoblastic hyperactivity.

Conditions in which Alkaline Phosphatase is increased: in liver diseases and biliary tract in metastasis liver and bone metastasis, in acromegaly, in hyperthyroidism in rickets in infectious mononucleosis, Paget's disease (osteoblastic activity) and physiological bone growth.

Decreased values of Alkaline Phosphatase can be observed in hypothyroidism, growth retardation in children, hypophosphatasia (inborn metabolic error), serious malnutrition.

#### NUMBER OF TESTS

100 tests / 50  $\mu$ L of sample / 2,5 mL of Reagent

200 tests / 25  $\mu$ L of sample / 1,25 mL of Reagent

#### BIBLIOGRAPHIC REFERENCES

1 - COLEMAN, C. M.; STROJE, R. C.: Clin. Chem. Acta 13: 401, 1966.

2 - ROY, A. V.: Clin. Chem. 16:431, 1970.

3 - ZAK, B.: Clin Chim. Acta 3: 328, 1958.

4 - BRITTENHAM, G.: Clin Chim. Acta. 91: 203, 1979.

5 - TONKS, D. B.: Chim. Acta 2: 393, 1983.

6 - CARL, A. B. and EDWARD, R. A.: Tietz Textbook of Clinical Chem 2nd ed., 1994, 1942-1943.

#### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

#### **OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium

#### CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454.

E-mail: sac@bioclin.com.br

ANVISA registration for Alkaline Phosphatase Kit:  
10269360104.

**Review:** July/2017

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY  
(last day of month)

TEMPERATURE LIMITATION  
(store at)

CONTAINS SUFFICIENT  
FOR <N> TESTS

CONSULT INSTRUCTIONS  
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED  
REPRESENTATIVE

KEEP AWAY  
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF  
PACKAGE IS  
DAMAGED