FOSFATASE ALCALINA CINÉTICA

REF K021

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FINALIDADE

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ALKALINE PHOSPHATASE KINETIC
REF K021
USAGE INSTRUCTIONS

FUNCTION

Method for determination of Alkaline Phosphatase. Kinetic test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Kinetic IFCC

The Alkaline Phosphatase catalyzes the transfer of Phosphate group from the substrate p-Nitrophenylphosphate (pNPP) to 2-Amino-2-Methyl-1-Propanol (AMP), releasing p-Nitrophenol according to the equation below:

The rate of p-Nitrophenol releasing, which has high absorbance at 405 nm, is proportional of Alkaline Phosphatase enzymatic activity of the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: 2-Amino 2-Methyl Propanol < 1,0 mol/L, activators, stabilizer and preservative.

Number 2 - Substrate - Store between 2 and 8°C. Contains: p-NFF Substrate (p-Nitrophenylphosphate) < 100 mmol/L, stabilizer and preservative.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2
1	108 mL	12 mL
2	1 x 40 mL	1 x 10 mL
3	2 x 40 mL	2 x 10 mL
4	4 x 40 mL	4 x 10 mL
5	2 x 40 mL	1 x 20 mL
6	4 x 40 mL	2 x 20 mL
7	4 x 27 mL	1 x 12 mL
8	6 x 36 mL	6 x 4 mL
9	6 x 36 mL	2 x 12 mL
10	3 x 36 mL	3 x 4 mL
11	3 x 36 mL	1 x 12 mL
12	6 x 36 mL	1 x 24 mL
13	1 x 20 mL	1 x 5 mL
14	2 x 20 mL	2 x 5 mL
15	3 x 20 mL	3 x 5 mL
16	4 x 20 mL	4 x 5 mL
17	5 x 20 mL	5 x 5 mL
18	6 x 20 mL	6 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatch, test tubes, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must to be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers, that may alter the results significantly.
- 5- It is important, for the good development of the test, a rigorous control of time and temperature.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma (Heparin) free o hemolysis. The analyte is stable for 07 days if kept in temperature between 2 and 8°C and 30 days if kept at 10°C negative.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Mix nine (09) parts of Reagent Nº 1 with one (01) part of Reagent Nº 2. Working Reagent is stable during 14 days, if kept between 2 and 8°C. Store away from light.

REACTION CONDITIONS

It is indispensable conditions to use a thermostated cuvette at 37°C, 1 cm optical path and reading at 405 nm.

TECHNIQUE

Bioclin recommends using the kit Biocal Bioclin as calibrator and as control serum, Biocontrol N and P Bioclin Kits.

Add 20 µL of Sample to 1,0 mL of Working Reagent, mix and transfer to a thermostated cuvette at 37°C and wait for 1 minute. Make the initial reading simultaneously starting the stopwatch. Repeat readings after 1, 2 and 3 minutes. Calculate the mean differences in absorbance per minute (ΔA/min.) and use it to calculate the result.

CALCULATIONS

ALP (U/L) = ΔA/min. x 2757
Results are expressed as U/L.

LIMITATIONS OF THE PROCEDURE

The kinetic method is based on the molar absorptivity. By this reason, the readings must be conducted in a spectrophotometer that satisfies the following conditions: Wavelength 405 nm
Semi trajectory of the pass band 10 nm
Stray light less than 0,5%
1cm thermostated cuvette

INTERFERENCES

Do not use hemolyzed or lipemic samples, for the results obtained might be falsely elevated.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the p-Nitrophenol, or through the BIOCAL calibrator. Bioclin recommends the usage of the BIOCAL calibrator which is also traceable to the molar absorptivity of the p-Nitrophenol.

REFERENCE VALUES

The reference values in U/L, for this method were obtained through the determination of Alkaline Phosphatase in healthy populations.
Adults: 27 to 100 U/L
Children: 75 to 390 U/L
(Including teenager up to 15 years old)

These values should be use as reference, being that every laboratory must create their own range of values for reference, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS AND SPECIFICITY
METHODOLOGY
The Alkaline Phosphatase Kinetic was compared with other commercially available methods for measurement of Phosphatase. 42 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 0,9979X - 0,486$ and the correlation coefficient 0,9996. With these results, we can conclude that the kit shows good methodological specificity.

Precision
REPEATABILITY
The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	105,80	278,40	196,93
Standard Deviation (U/L)	0,76	1,88	2,63
Coefficient of Variation (%)	0,72	0,67	1,34

REPRODUCIBILITY
The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	105,93	277,02	195,05
Standard Deviation (U/L)	2,70	1,71	1,76
Coefficient of Variation (%)	2,55	0,62	0,90

Sensitivity
The sensitivity was calculated from 40 determinations of a sample free the Phosphatase. The average found was 1,58 U/L and the standard deviation of 0,50 U/L. The sensitivity, that indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 3,08 U/L.

Linearity
The reaction is linear up to the concentration of 1500 U/L. For sample with value above 1500 U/L, we recomend diluting the sample with NaCl 0,85%. Multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE
Alkaline Phosphatase is present in high concentrations bone, liver, intestine and placenta. Constitutes a useful indicator of liver and bone diseases associated with osteoblastic hyperactivity. The conditions under which Alkaline Phosphatase is increased are: the liver and biliary tract, metastasis liver and bone metastasis, in acromegaly, in hyperthyroidism, in rickets, infectious mononudeosis, Paget's disease (osteoblastic activity) and in physiologic bone growth. Decreased values of Alkaline Phosphatase may be observed in hypothyroidism, growth retardation in children, hypophosphatasia (inborn metabolic error) and severe malnutrition.

BIBLIOGRAPHIC REFERENCES
1 - International Federation of Clinical Chemistry (1978). Section Nº 2. Provisional recommendations on the theory of reference values. Part. I. The concept of reference values. Expert panel on the Theory of Reference Values – Clin. Chim. Acta 87/3: 459F (1978).
2 - Spandrio, L. et al – Considerazioni sulla determinazione della fosfatesi alcalina del siero com método cinético – LAB2/2:123 (1975).
3 - Gardner, M. D. et al – Frecuency distribution and 'reference values' of plasma alkaline phosphatase (EC3.1.3.1) activity in the adult population of Scotish new town – J. Clin. Pathol. 31/12:1202 (1978).
4 - Demaría, I; Setta, F.; Lorenzo, L. – VIII Congreso Argentino de Bioquímica, Octubre 1990 – Revista de la Asociación Bioquímica Argentina 54/3 (1990).

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC

REP

OBELIS S.A.

English 2/2

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Alkaline Phosphatase Kinetic kit: 10269360072

Review: February/2019

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

EC

REP

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED