FÓSFORO

REF K020

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Phosphorus. Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Colorimetric

The Inorganic Phosphorus reacts with the Ammonium Molybdate (Reagent N° 2) forming Ammonium Phosphomolybdate, that later on is reduced to Blue of Molybdate, which intensity of color developed is proportional to the concentration of Phosphorus present in the sample.

REAGENTS

Number 1 - Reductor - Store between 15 and 30°C. Contains: Ascorbic Acid < 100 mmol/L and stabilizer.

Number 2 - Molybdate - Store between 15 and 30°C. Contains: Ammonium Molybdate < 50 mmol/L and stabilizer.

Number 3 - Alkaline Reagent - Store between 15 and 30°C. Contains: Sodium Hydroxide < 15 mol/L.

Number 4 - Standard - Store between 15 and 30°C. Contains: Phosphorus 5,0 mg/dL and preservatives.

PRESENTATION

Reagent	Volume
N° 1	7 mL
N° 2	7 mL
N° 3	20 mL
N° 4	4 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, whatch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTAION AND STORAGE CONDITIONS

The storage and transport temperatures shall be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Do not interchange bottle caps from reagents.

6- Use dropper in a vertical position, in order to drop over the reagent.

7- In places at high temperatures, it is advisable to make the determinations with the tubes in water bath.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, plasma harvested with Heparin and urine of 24 hours. The sample is stable for 4 days at 2 to 8°C and 1 year at -20°C.⁷ Serum or plasma should be separated within one hour after collection, due to release of hematic phosphorus. Acidification with HCl preserves the urine for 15 days at a temperature between 2 and 8°C.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Mark 3 tubes: B (Blank), A (sample), P (Standard) and proceed as following:

	Blank	Standard	Sample
Distilled Water	2,5 mL	2,5 mL	2,5 mL
Sample	---	---	100 µL
Reagent N° 4	---	100 µL	---
Reagent N° 1	1 drop	1 drop	1 drop
Reagent N° 2	1 drop	1 drop	1 drop
Homogenize well and wait 2 minutes. Then, add:			
Reagent N° 3	2 drops	2 drops	2 drops

Homogenize and let it rest for 5 minutes. Make the readings in 650 nm (640 - 670 nm), adjusting the zero with tube “B” (Blank). Color reaction is stable for 5 minutes.

For urine dosage: homogenize, separate 10 mL and adjust the pH between 1 and 3 with concentrated HCl. When performing the technique, dilute the urine 1:10 with distilled or deionized water, multiply the results by 10.

PHOSPHOLÍPIDS

The Phospholipids are extract selectively with Isopropyl Alcohol and later dosed according to the methodology for Phosphorus doses in serum. In a tube marked with the letter “A” (Sample), proceed as follow:

Isopropyl Alcohol	2,7 mL
Serum	300 µL

Add the serum drop by drop. Roughly shake and let it rest for 5 minutes. Centrifuge at 3000 rpm during 5 minutes. Transfer 1,0 mL of translucent supernatant for a test tube marked “A” (Sample). Place it in boiling water-bath until the complete evaporation of Isopropyl Alcohol. Heat the gas tip until the mixture is totally clear (5 minutes). Cool it between 15 and 30°C. This tube will be used for the color reaction. Carry out dosage according to the technique for serum, hiding the step of coloring of sample.

CALCULATIONS

Serum

Phosphorus (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}}$ x 5

Since the reaction follows the Lambert-Beer Law, the Calibration Factor can be used.

Calibration Factor = $\frac{\text{Standard Concentration (5 mg/dL)}}{\text{Standard Absorbance}}$

Phosphorus (mg/dL) = Sample Absorbance x Calibration Factor

Urine

Phosphorus (mg/24h) = $\frac{(\text{mg/dL}) \times \text{Urine Volume (mL)}}{100}$

Phospholipids = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}}$ x 125

The Phospholipids values are expressed in mg/dL of Lecithin.

INTERFERENTS

Hemolysis, lipemia or jaundice are factors that lead to values falsely elevated.

Citrated, fluoridated, oxalated or EDTA plasma produce falsely diminished results.

PROCEDURE LIMITATIONS

The development of a bluish color in Blank tube indicates that the water used is of poor quality, must be analyzed with criteria.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 186G.

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of the phosphorus in healthy populations of male and female.

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of Phosphorus. An average value of 0,1458 mg/dL, with standard deviation of 0,0005 mg/dL. The sensitivity that indicates the method detection limit, corresponds the average plus three times the standard deviation, and is equal to 0,1474 mg/dL.

Linearity

The color of the reaction is linear up to the concentration of 13 mg/dL. For higher values, dilute the sample with distilled or deionized water, perform a new determination and multiply the value obtained by the dilution factor.

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Phosphorus kit: 10269360117

EC

REP

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

English 2/2

To convert the values of mg/dL to mmol/L (SI) multiply by 0,322. These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS
The Phosphorus kit was compared with other commercially available methods for measurement of Phosphorus. 7 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 0,9708X + 0,1384$ and the correlation coefficient 0,9889. With these results, we can conclude that the kit shows good methodological specificity.

Precision
REPEATABILITY
The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	2,85	3,53	5,36
Standard Deviation (mg/dL)	0,01	0,01	0,01
Coefficient of Variation (%)	0,30	0,30	0,24

REPRODUCIBILITY
The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	2,86	3,55	5,36
Standard Deviation (mg/dL)	0,01	0,02	0,01
Coefficient of Variation (%)	0,35	0,43	0,19

DIAGNOSTIC SIGNIFICANCE
Phosphorus homeostasis is maintained primarily by renal function, intestinal absorption, bone metabolism and function parathyroid gland. Phosphorus has a role in bone composition, the acid-base balance, metabolism of fats and carbohydrates, besides being present the formation of compounds such as nucleic acids. We can observe the increased values of Phosphorus in the following conditions: renal failure (renal tuberculosis, chronic glomerulonephritis, hydronephrosis, pyelonephritis) hypoparathyroidism, excessive vitamin D. The decreased level of Phosphorus comes from many causes as vitamin D deficiency (rickets, osteomalacia) intake of antacids, liver disease, syndrome Fanconi, hyperparathyroidism, use of thiazide diuretics.

NUMBER OF TESTS
140 Tests / 100 µL of sample

BIBLIOGRAPHIC REFERENCES
1 - BAGINSKI, E. S.: Amer J. Med. Tec., 35: 475, 1969.
2 - GODWIN, J. F., Clin. Chem., 16: 776, 1970.
3 - TONKS, D. B., Clin. Chem., 9 – 217, 1963.
4 - CHRISTIAN, D. G.: Amer. J. Clin. Pathol. 54: 118, 1979.
5 - HENRY, R. J.: Clinical Chemistry – Principles and Technics, 2 Ed., Haper-Row, York, 1974.
6 - TODD; SANFORD; DAVIDSOHN, Diagnósticos Clínicos,1983.
7 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
8 - Bioclin – Dados de arquivo.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

EC

REP

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Phosphorus kit: 10269360117

Review: March/2019

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

EC

REP

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED