# Bioclin

## HbA1c BIRREAGENTE

REF K162

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#### FINALIDADE

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- 1- The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- 2- LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985-92.
- 3- JEPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.
- 4- HOEZEL W, WEYKAMP C et al, Clin. Chem., (2004), 50:1:166-74.
- 5- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 6- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for quantitative determination of A1c Hemoglobin in total blood, in photometric systems. For *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Immune turbidimetric

Hb1Ac is directly determined without the measurement of total hemoglobin. The total hemoglobin and the Hb1Ac in hemolyzed blood combine with latex particles with similar affinity. The number of connections is proportional to the relative concentration of both substances in blood.

Mouse monoclonal antibodies anti-human HbA1c bind to particles linked to Hb1Ac. Goat polyclonal antibodies anti mouse IgG react with the monoclonal anti-human Hb1Ac and agglutination occurs. The agglutination is proportional to the amount of HbA1c absorbed on the surface of latex particles. The agglutination is measured in absorbance. The HbA1c value is obtained by calibration curve.

#### REAGENTS

**Number 1 - Latex Buffer** - Store between 2 and 8°C. Contains: Latex < 5%, Tris buffer < 100 mmol/L, stabilizers and preservative.

**Number 2 - Anti-Hb1Ac Antibody** - Store between 2 and 8°C. Contains: Tris buffer < 100 mmol/L, mouse Monoclonal Antibody anti-Human HbA1c < 10 mg/dL, goat Polyclonal Antibody anti-mouse IgG < 10 mg/dL, stabilizers and preservative.

**Number 3 - Hemolyzing Solution** - Store between 2 and 8°C. Contains: Hemolyzing solution and preservative.

#### PRESENTATION

Presentation	Reagent		
	1	2	3
Nº 1	15 mL	5 mL	100 mL
Nº 2	2 x 15 mL	2 x 5 mL	2 x 100 mL
Nº 3	3 x 15 mL	3 x 5 mL	3 x 100 mL
Nº 4	4 x 15 mL	4 x 5 mL	4 x 100 mL
Nº 5	5 x 15 mL	5 x 5 mL	5 x 100 mL
Nº 6	30 mL	10 mL	2 x 100 mL
Nº 7	2 x 30 mL	2 x 10 mL	2 x 200 mL
Nº 8	2 x 30 mL	2 x 10 mL	4 x 100 mL
Nº 9	3 x 30 mL	3 x 10 mL	3 x 200 mL
Nº 10	3 x 30 mL	3 x 10 mL	6 x 100 mL
Nº 11	60 mL	20 mL	4 x 100 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, HbA1c Calibrator and HbA1c Control Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

**2- Strictly follow the methodology proposed to obtain exact results.**  
**3- Water used in material cleaning must be recent and free of contaminants.**

**4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.**

**5- Make a new calibration curve when necessary and every batch of product.**

**6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.**

**7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.**

**8- Do not use the product in case of damaged packaging.**

**9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.**

#### SAMPLES

Total Blood. Blood must be collected through venipuncture, avoiding prolonged tourniquet. Blood must be collected with EDTA.

#### Stability of sample

Total Blood: 1 week between 2 and 8°C.  
 Hemolyzed: 10 hours between 15 and 30°C and 10 days between 2 and 8°C.

#### PROCESS DESCRIPTION

The calibration stability of the HbA1c Bireagent kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using controls.

#### SAMPLE PREPARATION

Technique	
Sample	10 µL
Hemolyzing Solution	500 µL
Homogenize for 5 minutes until complete Hemolysis	

#### TECHNIQUE

The reagents are ready to use. For the use of the kit, use as calibrator the HbA1c Calibrator Bioclin kit and as control the HbA1c Control Bioclin kit.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC.**

#### Calibration

To perform the calibration curve, the **Calibrator HbA1c - K092** kit should be used.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### REFERENCE VALUES

Diagnostic criteria for diabetes *mellitus* recommended by the Brazilian Society of Diabetes.

Normal	< 5,7%
Pre-diabetes	5,7 – 6,4%
Diabetes	> 6,5%

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### PRODUCT PERFORMANCE QUALITY CONTROL

#### Accuracy

#### COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The HbA1c Bireagent kit was compared with another method commercially available to measure Glycated Hemoglobin dosage. 42 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 0,9742X - 0,0351$  and correlation coefficient 0,9645. With these results we can conclude that the kit shows good methodological specificity.

#### Precision

#### REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (%)	9,81	4,58	7,27
Standard Deviation (%)	0,06	0,04	0,05
Coefficient of Variation (%)	0,56	0,96	0,62

#### REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (%)	9,96	4,71	7,39
Standard Deviation (%)	0,15	0,12	0,10
Coefficient of Variation (%)	1,46	2,49	1,31

**Sensitivity**

Minimum limit of detection is of 1.5% Hb1Ac.

**Linearity**

This test determines the concentration of Hb1Ac in the range of 1.5% up to the highest point of the calibration curve.

**DIAGNOSTIC SIGNIFICANCE**

HbA1c is formed by the non enzymatic reaction of glucose with native hemoglobin. This process occurs continuously throughout red cells life (average life 100 - 120 days). Reason from glucosylation is directly proportional to the concentration of blood glucose. The level of HbA1c in blood represents the average level of glucose 6 to 8 weeks earlier. Therefore, it is convenient to monitor the concentration of glucose in the blood in a great period in individuals with mellitus diabetes. Clinical studies show that low levels of HbA1c may help prevent or slow the rate of complications in late diabetes.

**BIBLIOGRAPHIC REFERENCES**

- 1- The Diabetes Control and Complications Trial Research Group, N. Engl. J. Med., (1993), 329:977-86.
- 2- LITTLE PR, ROHLFING CL, WIEDMEYER HM, MYERS GL et al. Clin. Chem., (2001), 47:1985-92.
- 3- JEPPSSON JO, KOBOLD U, BARR J, FINKE A et al. Clin. Chem. Lab. Med., (2002), 40:78-79.
- 4- HOEZEL W, WEYKAMP C et al, Cin. Chem., (2004), 50:1:166-74.
- 5- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 6- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil


**CUSTOMER SERVICE**

Customer Advisory Service  
Phone: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for HbA1c Bireagent kit: 10269360270

**Review:** November/2019

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED