HBsAg
REF K145

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Immunochromatographic test for rapid determination of qualitative Hepatitis B virus surface antigen (HBsAg / ad and ay subtypes) in serum, plasma or whole blood.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
The HBsAg kit is an immunochromatographic test that contains IgG conjugate (anti HBs) - colloidal gold that binds to the HBsAg presented in the sample and forms a complex antigen / antibody / colloidal gold. The complex will migrate through the membrane by capillarity action in the direction of the anti-HBsAg antibody in the (T) test region. When the antibody presente in the test region has been bound, a line will appear to determinate the reagent reaction. In the absense of the HBsAg surface antigen, the line will not appear. The complex will continue to migrate until it reaches the (C) control region. The conjugate not linked to the antigen binds to the anti-IgG antibody of the C region forming a colored line indicating that the test is working and that the result is valid.

REAGENTS:

Reagent Number 1 (R1): Cassette - Store between 2 and 30°C. **Do not freeze.** Contains:

- A)** 01 plastic cassette with 01 test strip containing:
- Control Region (C): Anti-IgG antibody;
- Test Region (T): Goat anti-HBsAg antibody;
- Conjugate: Conjugate IgG (anti-HBs) - colloidal gold.

- B)** 01 sachet of silica

Reagent Number 2 (R2): Diluent - Store between 2 and 30°C. **Do not freeze.** Contains: Buffer Solution < 200 mmol/L and preservative.

Collection Accessories:

- Disposable plastic pipette
- Disposable lancet (Sterile)

PRESENTATION

Reagent	Number Cassettes per Packaging	Diluent
Cassette + Diluent (Includes pipette) *	5	1,0 mL
	10	1,5 mL
	15	2,0 mL
	20	3,0 mL
	25	3,5 mL
	30	4,0 mL
	40	5,0 mL
	50	6,0 mL
	100	2 x 6,0 mL
Cassette + Diluent (Includes lancet and pipettes) *	5	1,0 mL
	10	1,5 mL
	15	2,0 mL
	20	3,0 mL
	25	3,5 mL
	30	4,0 mL
	40	5,0 mL
	50	6,0 mL
	100	2 x 6,0 mL

* The number of accessories is proportional to the numbe of cassettes.

EQUIPMENTS AND OPERATIONAL INPUTS

Pipettes, tips, watches and stopwatches, alcohol 70% (p/p), cotton and waste for disposal. They can be found at markets specialized in laboratories of clinical analysis.

TRANSPORTATION AND STORAGE CONDITIONS

Transport can be carried out at temperatures up to 45°C by maximum, 45 days. The storage temperature must be 2 at 30°C. Avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
2- Strictly follow the methodology proposed to obtain accurate results.
3- The cassettes (test devices) do not contain contaminating substances, but since the samples are biological, the safety routines must be followed strictly and the use of disposable gloves is essential.
4- Use clean dry recipients to settle the samples.
5- Open the aluminum envelop, to remove the cassette (test device), only at the time of use.
6- The kit should not be frozen or exposed to high temperatures (above 45 ° C) as this will cause deterioration of the kit.
7- The reagents, as well as the sample, should reach room temperature before starting the test.
8- Do not mix kit reagents with the same lot or with different lot numbers.
9- Do not use kit beyond its expiration date. The expiry date is stated on the product label.
10- The cassette (test device), the pipette and the lancet are disposable and can not be reused.

- 11-** Kit materials that have had contact with the patient's biological material should have a special destination. Disposal of lancets and pipettes should be disposed of in hazardous waste. The cassettes (test devices) and diluent vials should be discarded in biological waste.
12- Cartons, Instructions for Use, aluminum envelope and sachet of silica may be disposed of in regular waste.
13- We recommend the application of the local, state and federal standards of environmental protection so that the disposal of the reagents and the biological material is done in accordance with the current legislation.
14- In order to obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or at the request of the SAC Customer Assistance) of Quibasa.
15- Do not use the product in case of damage to the packaging.
16- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES


Serum or Plasma: Specimens may be refrigerat at 2 a 8°C for a maximum period of 3 days. If samples can not be assayed within 3 days, they can be stored for up to 30 days at a temperature of -20°C (freezer). To obtain the plasma, collect EDTA, Heparin or Citrate as anticoagulant.

Whole Blood: Collect whole blood by venipuncture using the appropriate anticoagulant (EDTA, Heparin or Citrate) or by digital puncture with the aid of a lancet. Whole blood collected by lancet should be used for immediate testing. Whole blood collected by venipuncture with anticoagulant may be stored under refrigeration at 2 to 8°C for a maximum period of 3 days. **Do not freeze.**


DESCRIPTION OF PROCESS

TECHNIQUE - VENOPUNTION

- 1-** The sample should be at room temperature before starting the test.
2- Remove the cassette (test drive) from the protective packaging, place it on a clean and level surface and appropriately identify it so.
3- For Whole Blood: Transfer 75 µL of whole blood into the sample well. Hold the bottle vertically and apply 2 drops of Diluent (Reagente Nº 2) in the sample well.


- 4- For Serum or Plasma:** Transfer 80 - 100 µL of sérum or plasma into the sample well.

- 5-** Wait for the formation of lines. Interpret results between 15 and 30 minutes. Do not interpret after 30 minutes.

TECHNIQUE - PUNCTURE DIGITAL / LANCET

- 1-** Remove the cassette (test device) from the protective packaging, put it on a clean and level surface and identify it in a way proper.
2- Wash and dry the patient's hands. Passing alcohol 70% (p/p) at the tip of the finger for asepsis of the area used.
3- Press the tip of the finger that will be pierced by the lancet for accumulation of blood in this region.
4- Remove the protective cap from the lancet.
5- Position and press the lancet firmly in the area to be lancet. Then the blood will leave the perforated area.
6- Collect the blood with the disposable plastic pipette that comes with the kit. Put the pipette on the drop of blood and let the blood flow by capillarity, **without** pressing the bulb, to the marked trace in the disposable plastic pipette (equivalent to 75 µL of blood).
7- In the sample well, dispense the blood by pressing the bulb of the pipette.
8- Hold the Diluent bottle (Reagent Nº 2) vertically and apply 2 drops (50 µL) of Diluent in the sample well.
9- Wait for the lines to form. Interpret the results between 15 and 30 minutes. Do not interpret after 30 minutes.

INTERPRETATION OF RESULTS


Reagent Test: Formation of two red lines, one in the control region (C) and one in the test region (T) after 15 minutes. Do not interpret after 30 minutes.

Non-Reagent Test: Formation of a red line in the control region (C) and complete absence of line in the test region (T) after 15 minutes. Do not interpret after 30 minutes.

Invalid Testing: The absence of line formation in the control region (C) indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

PROCEDURE LIMITATIONS

- 1- The line in the control region will appear before 15 minutes test run, and this does not mean that the results can be interpreted ahead of schedule.
- 2- It is important to use the correct sample volume because lower or higher volumes can determine results mistakes.
- 3- The results provided by this kit should be interpreted responsible professional, and is not the sole criterion for the determination of the diagnosis and/or treatment of the patient.
- 4- A non-reactive result does not exclude the possibility of infection with the Hepatitis B virus (HBV). The concentration of HBsAg in the sample may vary according to the time of the presence of anti-HBs antibodies, and be less than the limit of analytical detection of the HBsAg kit.
- 5- If a suspicion of HBV infection persists, a new sample should be collected 30 days after the date of collection of the first sample to perform a new test.
- 6- This kit was validated only for serum samples, plasma and whole blood. Do not use other biological fluids or *pool* of the sample.
- 7- Do not use *post-mortem* samples.
- 8- The kit allows the detection of HBsAg in all phases of the infection, acute or chronic, provided that at detectable levels.
- 9- The isolated test results can not be used for definitive diagnosis. For confirmation of results reagents, follow the "Manual Técnico para o Diagnóstico das Hepatites Virais" approved by Ordinance nº 25 of 01 December 2015.
- 10- A non-reactive result does not exclude the possibility of infection with the hepatitis B virus (HBV). The concentration of HBsAg in the sample may vary according to disease stages at time of testing with the presence of anti-HBs antibodies, and in patients on antiviral treatment that cause viral load reduction. These conditions may decrease the patient's HBsAg levels lower than the kit's analytical detection limit.

INTERFERENCES

Hemolyses with hemoglobin values up to 10 g/L do not interfere in the HBsAg dosage. No significant interference were found in up to 30 g/L of lipids, 0,2 g/L of Bilirubin and 60 g/L of proteins. Interferentes above values can lead to false reagent results or false non-reagent. There was no evidence of interference with Rheumatoid Factor up to 1083 IU/mL, C Reactive Protein up to 41,5 mg/dL and Anti Streptolysin O up to 1046 IU/mL.

INTERNAL QUALITY CONTROL

The Clinical laboratory must have an internal quality control program where procedures, regulations, limits, and tolerance to variations are clearly established. It is important that every measurement system present a analytical variability, that must be monitored by the laboratories. Therefore, it is recommended the use of controls that allow the dosages accuracy and precision evaluation.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The HBsAg kit was compared with other commercially available HBsAg dosing kit. 86 analyzes were performed and the results were evaluated. The agreement found was > 99,9%. With these results, we can conclude that the kit shows good methodological specificity.

	Reference Kit	HBsAg
Sample Reagent	32	32
Sample No-Reagent	54	54
Total samples	86	

Agreement: > 99,9% (86/86)

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 3 different samples, obtaining the following results. The results presented a 100% correlation between the expected and found results.

Samples	Nº of Repetitions	Expected Results	Found Results
01 (Serum)	10	Reagent	Reagent
02 (Plasma)	10	Reagent	Reagent
03 (Blood Total)	10	Reagent	Reagent

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 different samples, obtaining the following results. The results presented a 100% correlation between the expected and found results.

Samples	Nº of Days	Nº of Repetitions	Expected Results	Found Results
01 (Serum)	03	10	Reagent	Reagent
02 (Plasma)	03	10	Reagent	Reagent
03 (Blood Total)	03	10	Reagent	Reagent

Clinical Sensitivity and Spesificity

The HBsAg kit was tested using clinically known samples. For this test, 1365 samples were analyzed, of which 349 are reagent for HBsAg and 1016 no-reagent for HBsAg. The following results were found:

	Expected Results	HBsAg
Sample Reagent	349	349
Sample No-Reagent	1016	1014
Total samples	1365	

The found results suggest the following values:

Clinical Sensitivity: > 99,9% (349/349)

Clinical Specificity: 99,8% (1014/1016)

Analytical Sensitivity

The analytical sensitivity is 2 ng/mL.

High Dose Hook Effect

Hook effect wasn't observed with a high concentration of HBsAg sample up to 800 µg/mL.

DIAGNOSTIC SIGNIFICANCE

Hepatitis B is an infectious disease caused by the virus HBV that affects the liver. The HBV virus is transmitted through sexual contact, blood exposure (including sharing of sharps) and transmission from mother to child via the placenta or during childbirth.

HBV infection results in an apparent number of serological markers and one of the first of these markers is the Hepatitis B Surface Antigen (HBsAg). The main subtypes of HBsAg includes ad and ay, all sharing the common determinant 'a'.

HBsAg appears between 1 and 10 weeks after infection and before of biochemical evidence of liver disease or jaundice. However, for most people, due to its high production, the HBsAg can be detected after 30 days of infection. This one period is known as the diagnostic window. In the state of chronic carrier, HBsAg persists for more than 6 months. The Antibodies against HBV antigens appear between 30 and 60 days after infection, varying between individuals. This is the period known as the immunological window.

HBV infection can lead to acute auto - limiting hepatitis, fulminant hepatitis, chronic hepatitis with progression for cirrhosis and hepatic insufficiency, and carrier status chronic asymptomatic. The HBV virus in infected people persists for the rest of their lives and can be passed on to other people. Thus, the hepatitis B has become a public health problem. Therefore, screening for HBsAg is recommended for all donors, pregnant women and people in high-risk groups.

NUMBER OF TESTS

The number of tests corresponds to the number of cassettes (test device) in each presentation.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for HBsAg kit: 10269360307

Review: August/2019

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED