HCV

REF K144

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FINALIDADE

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1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Four-antigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
5. Bioclin – Dados de arquivos
5. Bioclin – Dados de arquivo.

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1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lellie. Confirmation of Hepatitis C Virus Infection by New Four-antigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
5. Bioclin – Datos de archivos
5. Bioclin – Datos de archivo.

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USAGE INSTRUCTIONS

FUNCTION

Immunochromatographic method for rapid determination of qualitative anti HCV in serum or whole blood.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The HCV kit is an immunochromatographic test that contains colloidal gold conjugate that binds to the anti HCV antibody presented in the sample and forms a complex that will migrate through the membrane by capillarity action in the direction of the HCV antigen (Core, NS3, NS4 and NS5) immobilized in the (T) Test region, when it is linked occurs the appearance of the strip that determinates the positive reaction. In the absence of the anti HCV, the strip will not appear. The complex will continue to migrate until it reaches the (C) Control region. The conjugate not linked to the antigen binds to the C region forming a colored strip indicating that the test is working as planned and that the result is valid.

REAGENTS:

Reagent N° 1: Cassette: Store between 15 and 30°C. **Do not freeze.** Contains:

- Control Region (C): antibody anti-IgG.
- Region Test (T):HCV antigens.
- Conjugate: Colloidal gold conjugate.

Reagent N° 2 Diluent: Store between 15 and 30°C. **Do not freeze.** Contains: Tris Buffer 100 mmol/L pH 7,4, preservative.

PRESENTATION

Reagent	Presentation	Number of Cassettes for Package	Diluent
Casete + Diluent	1	1	1 mL
	2	5	1 mL
	3	10	1 mL
	4	15	2,5 mL
	5	20	2,5 mL
	6	25	2,5 mL
	7	30	2 x 2,5 mL
	8	40	2 x 2,5 mL
	9	50	2 x 2,5 mL
	10	100	4 x 2,5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Watch or stopwatch. Can be found in specialized markets for articles for Clinical Analysis Laboratories.

TRANSPORTAION AND STORAGE CONDITIONS

The temperature of storage and transportation must be 15 to 30°C. Avoid moisture and protect from light.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain accurate results.

3- The cassettes do not show contaminants, but how the samples are biological safety routines should be followed accurately and the use of disposable gloves is essential.

4- Use clean and dry containers for storage of samples;

5- Treat all samples as potentially infecting.

6- Reagent N°2 contains Sodium Azide, irritating to mucous and skin. Handle with caution.

7- The kit should not be frozen or exposed to high temperatures, as this will cause deterioration.

8- The reagents and the sample must be equilibrated environment before starting the test.

9- Open the envelope of aluminum, to remove the cassette, only in time of use.

10- It is interesting to conduct the test immediately after sample collection.

11- The kit must not be used after expiration date.

12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

14- Do not use the product in case of damaged packaging.

15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES


Serum or whole blood (without anticoagulants).

The samples can be stored under refrigeration, between 2 and 8°C for a maximum period of 3 days. If the samples can not be analyzed within those three days, they can be stored at -20°C for up to 30 days.

DESCRIPTION OF PROCESS

TECHNIQUE

- 1- The sample should be at room ambient temperature before starting the test.
- 2- Remove the cassette from its protective packaging, place it on a clean and leveled surface and appropriately identify it so.
- 3- Transfer 10 µL of sample in the sample well.
- 4- Hold the vial vertically and apply 2 drops (100 µL) of Diluent in the sample well.

5- Wait for the formation of lines after resting from 15 minutes. Interpret results between 15 and 20 minutes.

6- Do not take a reading 20 minutes past the reaction.

INTERPRETATION OF RESULTS


Testing Positive: Formation of two lines in the (C) control region and (T) test region after 15 minutes. Do no interpret after 20 minutes.

Test Negative: Formation of one line in the (C) control region and complete absence of a line in the (T) test region after 15 minutes. Do not interpret after 20 minutes.

Inadequate Testing: Complete absence of line on the (C) Control region with or without a line in the (T) Test region. Repeat the sample test.

Note:

- 1- The line in the Control region will appear before the 15 minutes of the test incubation, and that does not mean that the negative results can be interpreted before the time.
- 2- The positive result will not change after 20 minutes, but the result must not be interpreted after this time to avoid incorrect results.
- 3- The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

INTERFERENCES

Hemolyses with hemoglobin values up to 10 g/L do not interfere in the HCV dosage. No significant interference were found in up to 30 g/L of lipids, 0,2 g/L of bilirubin and 60 g/L of proteins.

INTERNAL QUALITY CONTROL

The Clinical laboratory must have an internal quality control program where procedures, regulations, limits, and tolerance to variations are clearly established. It is important that every measurement system present a analytical variability, that must be monitored by the laboratories. Therefore, it is recommended the use of controls that allow the dosages accuracy and precision evaluation.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS, METHODOLOGY AND SPECIFICITY SENSITIVITY

The Bioclin HCV Kit was compared to other commercially available kits for the anti HCV testing. For this test two hundred and forty (240) clinical samples were tested. The following results were found:

HCV Bioclin		Reference Kit		Subtotal
		+	-	
Test Kit	+	119	1	120
	-	0	120	120
Subtotal		119	121	240

The found results suggest the following values:

Sensitivity: > 99,9% (119/119)

Specificity: 99,17% (120/121)

Accuracy

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 5 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Found Results
01	20	100% Positive	100% Positive
02	20	100% Negative	100% Negative
03	20	100% Positive	100% Positive
04	20	100% Positive	100% Positive
05	20	100% Negative	100% Negative

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Samples	Nº of Days	Nº of Repetitions	Expected Results	Found Results
01	03	20	100% Positive	100% Positive
02	03	20	100% Negative	100% Negative
03	03	20	100% Positive	100% Positive
04	03	20	100% Positive	100% Positive
05	03	20	100% Negative	100% Negative

DIAGNOSTIC SIGNIFICANCE

The Hepatitis C Virus has a small envelope and it is the main cause of the parenteral transmission of the Hepatitis not-A and not-B. The infection by HCV cause a vast variety of chronic hepatic diseases, cirrhoses and hepatic cancer. The virus main transmissions are through blood transfusion and hemoderivatives, organs transplant, and contaminated needles and syringes shared. Antibody against HCV are found in more than 80% of patients with Hepatitis non-A or non-B. Cloning the viral genome had become possible the development of the tests by serum that use recombined antigens.

NUMBER OF TESTS

Reagents	Presentation	Number of Tests
Cassettes	1	1
	2	5
	3	10
	4	15
	5	20
	6	25
	7	30
	8	40
	9	50
	10	100

BIBLIOGRAPHIC REFERENCES

1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. *Science*. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. *Science*. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Four-antigen Recombinant Immunoblot Assay. *Lancet*. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. *J. Clinical Immunoassay*. 1993;16:204.
5. Bioclin – Dados de arquivos

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for HCV: 10269360231

Review: March/2017

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED