# Bioclin

## HDL DIRETO

REF K071

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#### FINALIDADE

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- 1 - WILLIAM, P., ROBINSON, D., BAILY, A.. High Density Lipoprotein and Coronary Risk Factor. Lancet, 1:72 (1979).
- 2 - BADIMON, J.J., BADIMON, L., FUESTER, V.. Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma. Fraction in the Cholesterol-Fed Rabbit. Journal of Clinical Investigation, 1990;85:1234-41.
- 3 - BARR, D.P., RUSS, E.M., EDER, H.A.. Protein-Lipid Relationships in Human Plasma. Am. J. Med., 11:480, 1951.
- 4 - GOTTO, A.M.. Lipoprotein Metabolism and Etiology of Hyperlipidemia. Hospital Practice, 23: Suppl. 1.4 1988.
- 5 - Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 7- Executive Summary of the Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA 2001: 285:2486-97.
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of HDL Cholesterol in samples of serum or plasma (Heparin or EDTA). Colorimetric Enzymatic test, for *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Enzymatic - Colorimetric

The HDL Direct kit allows a direct determination without need for prior precipitation or sample treatment. Reagents are used to carry out selective dosage of cholesterol bound to HDL. The surfaces of low density lipoproteins (LDL), very low density (VLDL) and chylomicrons are stabilized by adsorption of polyanion and do not undergo the action of modified enzymes present in reagent. HDL, however, solubilize under the action of a detergent, allowing the enzymatic action on cholesterol attached to it. Thus, the intensity of color formed is proportional to the concentration of HDL cholesterol in the sample.

#### REAGENTS

**Number 1 - Polyanion** - Store between 2 and 8°C. Contains: Polyanion 0.2% and 4-Aminophenazone 2.0 mmol/L.

**Number 2 - Enzymatic Reagent** - Store between 2 and 8°C. Contains: Buffer 100 mmol/L (pH 7.0), DSBT 24 mmol/L, Cholesterol Esterase < 300 U/L, Cholesterol Oxidase < 300 U/L, Peroxidase < 1000 U/L and detergent.

**Number 3 - Calibrator** - Store between 2 and 8°C. Contains: Lyophilized Serum with HDL value set. **Potentially infectious.**

**Warning: The concentration varies with the lot - see bottle label.**

#### PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2	Reagent N° 3
1	60 mL	20 mL	2 mL
2	240 mL	80 mL	2 mL
3	30 mL	10 mL	2 mL
4	2 x 30 mL	2 x 10 mL	2 mL
5	4 x 30 mL	4 x 10 mL	2 mL
6	6 x 30 mL	6 x 10 mL	2 mL
7	2 x 60 mL	2 x 20 mL	2 mL
8	3 x 60 mL	3 x 20 mL	2 mL
9	4 x 60 mL	4 x 20 mL	2 mL
23	60 mL	20 mL	5 mL
24	240 mL	80 mL	5 mL
25	30 mL	10 mL	5 mL
26	2 x 30 mL	2 x 10 mL	5 mL
27	4 x 30 mL	4 x 10 mL	5 mL
28	6 x 30 mL	6 x 10 mL	5 mL
29	2 x 60 mL	2 x 20 mL	5 mL
30	3 x 60 mL	3 x 20 mL	5 mL
31	4 x 60 mL	4 x 20 mL	5 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, watch or stopwatches, pipettes, test tubes, water bath at 37°C, Biocontrol N Bioclin and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum obtained free from hemolysis, plasma collected with Heparin or EDTA. After collection, the sample should be centrifuged in up to 3 hours.

HDL cholesterol is stable in serum for up to 7 days at 2 to 8°C<sup>6</sup> and for 30 days at -20°C.

Samples stored at temperatures between 15 and 30°C have stabilities of at most 14 hours.

#### PROCESS DESCRIPTION

The calibration stability of the HDL Direct kit installed on refrigerated equipment is at least 21 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

#### REAGENTS PREPARATION

**Before reconstitute the Calibrator should be noted the volume indicated on the bottle and proceed as described below.**

**Presentation 1 to 9:** Open carefully the bottle and add exactly **2 mL of distilled water**, as volume described in the label of the bottle.

**Presentation 23 to 31:** Open carefully the bottle and add exactly **5 mL of distilled water**, as volume described in the label of the bottle.

Mix. Leave the bottle resting for 30 minutes. Every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid formation of foam.

After reconstitution, the calibrator is stable for 7 days at 2 to 8°C or 30 days at -20°C.

#### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), C (Calibrator), S (Sample) and proceed as follows:

	Blank	Calibrator	Sample
<b>Reagent N° 1</b>	300 µL	300 µL	300 µL
<b>Calibrator</b>	---	3 µL	---
<b>Sample</b>	---	---	3 µL

Homogenize and place it in a water bath at 37°C for 5 minutes. Read the absorbance from Sample (S1) and Calibrator (C1) in 550 nm (540 - 560 nm), adjusting the zero with the Blank.

Reagent N° 2	100 µL	100 µL	100 µL

Homogenize and place it in a water bath at 37°C for 5 minutes. Read the absorbance from Sample (S2) and Calibrator (C2) in 550 nm (540 - 560 nm), adjusting the zero with the Blank. Color is stable for 60 minutes.

#### CALCULATIONS

HDL Cholesterol =  $\frac{(A_2 - A_1)}{(C_2 - C_1)} \times \text{Calibrator Concentration}$

Results are expressed in mg/dL.

#### PROCEDURE LIMITATIONS

Do not use anticoagulants containing citrate. Do not expose the reactive to light.

#### INTERFERENCES

The Ascorbic Acid up to 1000 mg/dL, Bilirubin up to 60 mg/dL, Hemoglobin up to 1000 mg/dL, Gammaglobulin up to 5000 mg/dL and Triglycerides up to 1800 mg/dL did not alter the results. N-acetylcysteine samples should not be use.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's calibrator is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 1951.

**REFERENCE VALUES <sup>7</sup>**

The reference values in mg/dL, for this method, were obtained through the determination of HDL cholesterol in healthy populations of male and female.

Reference values below are for individuals with 12 hours fasting and follow Adult Treatment Panel III (ATP III).

<b>Adults (≥ 20 years)</b>	<b>Desirable</b>	> 60 mg/dL
	<b>Limit</b>	40 - 60 mg/dL
	<b>Low</b>	< 40 mg/dL

<b>Childrens and Adolescents (2 to 19 years)</b>	<b>Desirable</b>	> 45 mg/dL
--	------------------	------------

To convert the data in International System (IS), multiply the results: mg/dL x 0.02586 = mmol/L

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE****QUALITY CONTROL****Accuracy  
RECOVERY**

The recovery analysis was performed with 05 determinations samples. Accuracies were calculated and found themselves in good agreement with the reference values, obtaining a recovery between 95% and 103%.

**COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY**

The HDL Direct kit was compared with other method for HDL Cholesterol dosage commercially available. 42 analyses were performed and the results were evaluated. The linear equation obtained was  $Y = 0.989X + 0.3611$  and the correlation coefficient 0.9987. With these results we can conclude that the kit shows good methodological specificity.

**Precision  
REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average concentration (mg/dL)</b>	16.98	26.18	21.05
<b>Standard Deviation (mg/dL)</b>	0.16	0.38	0.22
<b>Coefficient of Variation (%)</b>	0.93	1.47	1.05

**REPRODUCIBILITY**

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average concentration (mg/dL)</b>	16.70	26.16	21.01
<b>Standard Deviation (mg/dL)</b>	0.52	0.03	0.49
<b>Coefficient of Variation (%)</b>	3.12	0.11	2.33

**Sensitivity**

The sensitivity was calculated from 40 determinations of a sample free of HDL Cholesterol. The average found was 0.66 mg/dL with a standard deviation of 0.08 mg/dL. The sensitivity, which indicates the method detection limit, correspond the average plus 3 times the standard deviation, and is equal to 0.91 mg/dL.

**Linearity**

Reaction is linear from 3 to 200 mg/dL. For values higher than 200 mg/dL, repeat the determination diluting the sample with NaCl 0.85%. Multiply the results obtained by the dilution factor.

**DIAGNOSTIC SIGNIFICANCE**

The plasmatic lipoproteins are spherical particles consisting of Cholesterol, Triglycerides, and Phospholipids and Proteins. These particles are responsible for transportation of the cholesterol from peripheral cells and liver. The relative proportion between proteins and lipids determine the density of lipoproteins - chylomicrons, VLDL, LDL and HDL. Monitoring HDL Cholesterol in serum is of great importance due to inverse correlation between plasmatic concentrations and the risk coronary disease.

**BIBLIOGRAPHIC REFERENCES**

- 1 - WILLIAM, P., ROBINSON, D., BAILY, A.. High Density Lipoprotein and Coronary Risk Factor. Lancet, 1:72 (1979).
- 2 - BADIMON, J.J., BADIMON, L., FUESTER, V. Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma. Fraction in the Cholesterol-Fed Rabbit. Journal of Clinical Investigation, 1990; 85: 1234-41.
- 3 - BARR, D.P., RUSS, E.M., EDER, H.A.. Protein-Lipid Relationships in Human Plasma. Am. J. Med., 11:480, 1951.
- 4 - GOTTO, A.M.. Lipoprotein Metabolism and Etiology of Hyperlipidemia. Hospital Practice, 23: Suppl. 1,4 1988.
- 5 - Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 7- Executive Summary of the Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA 2001; 285:2486-97.
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**  
Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium


**CUSTOMER SERVICE**

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for HDL Direct kit: 10269360122

**Review:** June/2020

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED