HEMOGLOBINA

REF K023

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Hemoglobin. Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Cyanomethaemoglobin

The method is based in oxidizing the Iron atom (Iron II) of the Hemoglobin molecule by Potassium Ironicyanate in pH low alkaline, forming the Methaemoglobin which is converted in Cyanomethaemoglobin after the reaction with the Potassium Cyanate. The reddish color is proportional to the Hemoglobin concentration presented in the sample.

REAGENTS

Number 1 - Stock Color Reagent - Store between 15 and 30°C. Contains: Phosphate buffer 200 mmol/L, Potassium Ferricyanide 120 mmol/L, Potassium Cyanide 150 mmol/L and surfactant.

Handle this reagent very carefully, toxic reagent.

PRESENTATION

Reagent	Volume
Nº1	2 x 10 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, watch or stopwatch, pipettes, test tubes, volumetric balloon and **Bioclin** Standard Hemoglobin Kit. They can be found at markets specialized in laboratories of clinical analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperatures should be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in cleaning must be to recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- When preparing the Working Reagent, bubbles can emerge, and it is totally normal. However, should be avoided, because it difficult the meniscus accuracy.
- It is of great importance that the flask storing the Working Reagent do not contain acid, base or organic solvents such as detergents.
- The Working Reagent and the residual reactions can never be left in presence of acid reagents, because they lead to the formation of toxic gases.
- The total blood, when collected in oxalate or its mixtures demands a larger time interval, so the reading can be carried out.
- The blood sample can be homogenized properly before pipette.
- The tips of automatic pipettes must be clean, dry and in physical use conditions.
- The color reagent can not be used if show change in the color.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Total blood obtained free of hemolysis, collected with EDTA, citrate or oxalate.

PROCESS DESCRIPTION

WORKING REAGENT PREPARATION

In a 1 liter volumetric balloon place approximately 500 mL of distilled or deionized water. Transfer the whole content of the Stock Color Reagent flasks, washing it 3 times with distilled or deionized water, transferring the washing waters to the same balloon. Complete the volume of 1.000 mL with distilled or deionized water. Store it in a amber flask to the environment temperature, protect from direct light and moisture. Stable for 6 months between 15 and 30°C.

TECHNIQUE

Mark test tube "S" (Sample) and proceed as follows:

	Sample
Working Reagent	2,5 mL
Whole Blood	10 µL

Homogenized and wait 5 minutes at environment temperature, carry out the reading at 540 nm (520 - 550 nm) matching the zero with the working reagent. The color reaction is stable for 60 minutes.

CALCULATIONS

To calculate the value in g/dL see the instructions contained in **Bioclin** Standard Hemoglobin kit.

Calibration Factor = $\frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$

Hemoglobin = Sample Absorbance x Calibration Factor (g/dL)

PROCEDURE LIMITATIONS

The Color Reagent contains Cyanate. Handle carefully. **Toxic.**

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

The reference values in g/dL, for this method, were obtained through the determination of the hemoglobin in healthy populations of male and female.
Men12,5 to 17,5 g/dL
Women11,5 to 15,5 g/dL

To convert g/dL to mmol/L values, multiply it by 1,61. These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The Hemoglobin kit was compared with other commercially available methods for measurement of hemoglobin. 07 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1,008X - 0,107$, with correlation coefficient 0,999. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (g/dL)	9,92	13,12	17,05
Standard Deviation (g/dL)	0,01	0,01	0,01
Coefficient of Variation (%)	0,08	0,05	0,08

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (g/dL)	9,92	13,12	17,05
Standard Deviation (g/dL)	0,01	0,01	0,01
Coefficient of Variation (%)	0,06	0,04	0,03

Sensitivity

Sensitivity was calculated from 20 determinations of a sample free of hemoglobin. The average found was 0,329 g/dL and the standard deviation of 0,001 g/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,330 g/dL.

Linearity

The reaction is linear up to 20 g/dL.

DIAGNOSTIC SIGNIFICANCE

Hemoglobin is a protein present inside the red blood cells and has as its main function the transport of oxygen from the lungs to the cells of the body. Its determination is important for the diagnosis and control of anemia, polycythemia, dehydration, hemorrhages, nutritional deficiencies, erythematous lupus, among others.

NUMBER OF TESTS

800 Tests / 10 μ L of Sample / 2,5 mL of Reagent

BIBLIOGRAPHIC REFERENCES

- 1 - HENRY, R. J., Clinical Chemistry – Principles and Technics – Hoeber Division, Haper and Row, New York.
- 2 - TONKS, D. B., Clin. Chem., 9:217, 1983.
- 3 - DRABKIN, D. L., Austin, J. H., J. Biol. Chem. , 112:51, 1935.
- 4 - BARON, D. N., Bioquímica de la Patología Médica Ediciones Toray S/A – Barcelona.
- 5 - Internacional Committe for Standartization in Haematology: Recommendations for Haemoglobinometry in Human Blood: Brit., J. Haemat., 13 (suppl.) 71, 1967.
- 6 - Word Health Organization, International Biological Reference Preparations 1968. W. H. O. Techn. Rep. Ser. 384, 85, 1968.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP

OBELIS S.A.

English 2/2

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454.
E-mail: sac@bioclin.com.br

ANVISA registration for Hemoglobin kit: 10269360097

Review: May/2017

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED