HIV TRI LINE

REF K087

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Immunochromatographic method for rapid and qualitative determination antibodies (IgM, IgG, IgA and IgE) against HIV 1, including subtype O and anti-HIV 2. For in vitro diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The 3rd generation Tri Line HIV kit, manufactured with antigens both the solid and conjugate phases, is capable of to detect total antibodies (IgM, IgG, IgA and IgE) to HIV 1, including subtype O, and HIV 2. Recombinant antigens HIV 1 and HIV 2 react with antibodies present in serum, plasma and whole blood samples. At samples move through a chromatographic membrane by capillary action. Reagent samples for HIV 1 will form a red line in the region where recombinant HIV antigen 1 is immobilized. Reagent samples for HIV 2 will form a second red line in the region corresponding to the antigen HIV 2 immobilized. Samples continue to be absorbed by the membrane to the control region, with the formation of another line red, confirming the correct processing of the test.

REAGENTS

Reagent Number 1 (R1) Cassette - Store between 2 and 30°C. **Do not freeze.** Contains:

- A)** 01 plastic cassette with 01 test strip containing:
- Control region (C): Polyclonal goat anti-IgG antibody;
 - Test region (T): HIV-1 and HIV-2 recombinant antigens;
 - Conjugate: Recombinant HIV-1 and HIV-2 antigens conjugated with colloidal gold;

B) 01 sachet of silica;

Reagent Number 2 (R2) Diluent - Store between 2 and 30°C. **Do not freeze.** Contains: Buffer Tris (pH 8,0).

Collection Accessories:

- Disposable plastic pipette.
- Disposable lancet

PRESENTATION

The kit possess the following presentations:

Reagent	Presentation	Number Cassette per Packaging	Diluent
Cassette + Diluent (Includes pipettes)*	2	5	5,0 mL
	3	10	5,0 mL
	4	20	5,0 mL
	5	25	5,0 mL
	6	30	5,0 mL
	7	40	2 X 5,0 mL
	8	50	2 X 5,0 mL
	9	100	3 X 5,0 mL
Cassette + Diluent (Includes lancets and pipettes)*	10	5	5,0 mL
	11	10	5,0 mL
	12	20	5,0 mL
	13	25	5,0 mL
	14	30	5,0 mL
	15	40	2 X 5,0 mL
	16	50	2 X 5,0 mL
	17	100	3 X 5,0 mL

* These presentations accompany lancet and/or plastic pipette disposables. The number of accessories is proportional to the number of cassettes.

EQUIPMENTS AND OPERATIONAL INPUTS

Pipettes, tips, watches and stopwatches, alcohol 70% (p/p), cotton and waste for disposal. They can be found at markets specialized in laboratories of clinical analysis.

CONDITIONS OF STORAGE AND TRANSPORTATION

Transport can be carried out at temperatures up to 45°C by maximum, 3 weeks. The storage temperature must be 2 at 30°C. Avoid moisture. Do not freeze.

WARNINGS

- 1- For *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- The used device has potential infectious risk. The safety routine must be followed strictly and the use of disposable gloves is indispensable.
- Use clean and dry recipients to settle the samples.
- Open the aluminum envelop, to take the cassette off, only for its immediate use.
- The kit must not be frozen or exposed to temperatures above 45°C this exopside may cause its deterioration.
- The reagents, as well as the sample, must reach environment temperature before the test begins.
- The disposal of diagnostic devices must follow the local infectious disposal lawor Laboratory regulation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE


Serum/Plasma: the serum and plasma samples are stable for 3 days when kept between 2 and 8°C and for 3 months at -20°C, it can only be defrosted for use once. Hemolyzed and lipemic samples shall be avoided, because they can cause a wrongly positive result. Serum or plasma containing precipitate can supply inconsistent results. To obtain the plasma, collect EDTA, Heparin or Citrate as anticoagulant.

Whole Blood: whole blood may be collected by the aid of a lancet (digital puncture) or by venipuncture, using the appropriate anticoagulant (EDTA, Heparin or Citrate) or Aid of a lancet (digital puncture). The whole blood obtained by digital puncture should be used for immediate testing. Whole blood obtained by venipuncture with anticoagulant may be stored between 2 and 8°C for up to 2 days. **Do not freeze the whole blood.**

PROCESS DESCRIPTION

TECHNIQUE – VENOPUNTION






- 1- Allow the samples equilibrate to room temperature before start the test.
- 2- Remove the cassette from the protective packaging, place it on a clean and level surface and identify it in a proper.
- 3- For Serum, Plasma or Whole Blood: Transfer 10 µL of serum, plasma or whole blood in well A.
- 4- Hold the Diluent bottle (Reagent N° 2) vertically and apply 2 drops (70 µL) of Diluent to well B.
- 5- Wait for the lines to form. Interpret the results between 15 and 30 minutes. Do not interpret after 30 minutes.

TECHNIQUE - PUNCTURE DIGITAL / LANCET

- 1- Remove the cassette from the protective packaging, put it on a clean and level surface and identify it in a way proper.
- 2- Wash and dry the patient's hands. Passing alcohol 70% (p/p) at the tip of the finger for asepsis of the area used.
- 3- Press the finger tip to be pierced by the lancet for accumulation of blood in this region.
- 4- Remove the protective cap from the lancet.
- 5- Position and press the lancet firmly in the area to be lancet. Then the blood will leave the perforated area.
- 6- Collect blood with disposable plastic pipette that accompanies the kit. Aspirate the blood until the trace marked on the disposable plastic pipette.
- 7- Dispense blood in to well A by pressing the pipette.

- 8- Hold the Diluent bottle (Reagent N° 2) vertically and apply 2 drops (70 µL) of Diluent to well B.
- 9- Wait for the lines to form. Interpret the results between 15 and 30 minutes. Do not interpret after 30 minutes.

INTERPRETATION OF RESULTS

Reagent Test: Formation of control line and one or more two red lines after 15 minutes in the test region. No interpret after 30 minutes.

HIV 1 Reagent: In addition to the red line in the Control region (C), another line is formed in the test region 1 (T1).


HIV 2 Reagent: In addition to the red line in the Control region (C), another line is formed in the Test region 2 (T2).

HIV 1 and 2 Reagents: In addition to the red line in the control (C), two other lines are formed in the Test regions 1 and 2 (T1 and T2). If the color intensity of line T1 is greater than line T2, it may be consider the result as HIV-1 reagent.

If the color intensity of the line T2 is greater than the line T1 can be consider the result as HIV-2 reagent.

Non-Reagent Test: Formation of only one red line, after 15 minutes, in the Control region. Do not interpret after 30 minutes. Note: See Notes.

Invalid Testing: The absence of line formation in the control (C) indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

Notes:

- 1-** The line in the control region will appear before 15 minutes test incubation, and this does not mean that the results can be interpreted ahead of schedule.
- 2-** It is important to use the correct sample volume because lower or higher volumes can determine results mistakes.
- 3-** HIV-2 has 60% homology to HIV-1. Considering that HIV-1 and HIV-2 co-infection is rare, if the test shows HIV-1 and HIV-2 reactive lines, it is likely that this result was caused by cross-reactivity between them. Therefore, if the HIV Tri Line kit is HIV-2 reactive, follow the guidelines in the Technical Manual for the diagnosis of HIV infection.
- 4-** Isolated test results can not be used for definitive diagnosis. For confirmation of results reagents, follow the Ordinance nº 29 of December 2013, approving the Technical Manual for the Diagnosis of HIV infection.
- 5-** The literature reports that, on average, 95% of those infected present seroconversion within 30 days. However, this period may vary according to the immune response of each individual. There may be very rare cases of individuals who have a prolonged course of seroconversion, presenting an immune window up to 6 months.
- 6-** Considering that the HIV Tri Line kit is 3rd generation, that is, manufactured with recombinant antigens both in the solid phase and in the conjugate, it is able to detect total antibodies (IgM, IgG, IgA and IgE) for HIV 1, including subtype O, and HIV 2. A literature describes that third generation kits are capable of detect antibodies on average 22 to 25 days after infection, for being able to detect IgM antibodies.

7- If a suspected HIV infection continues, a new sample should be collected 30 days after the date of collection of the first sample to perform a new test.

8- False non-reagents results may occur if the test be used in samples of patients known who are receiving or who have received antiretroviral therapy.

9- The results provided by this kit should be interpreted responsible professional, and is not the sole criterion for the determination of the diagnosis and / or treatment of the patient.

INTERFERENT

Hemolysis with hemoglobin values up to 10 g/L does not interfere with the dosage of HIV. No significant interferences were found up to 30 g/L of lipids, 0,2 g/L of bilirubin and 60 g/L of proteins. No cross-reactivity was observed with HBsAg, HCV, and Syphilis.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

METHODS COMPARISON, METHODOLOGICAL SPECIFI-CITY AND SENSITIVITY

To establish the sensitivity and specificity of the HIV Tri Line kit, other qualitative kit was compared to other methods available in the market.

A concordance of 100% was found. Comparative studies using 192 truly reagent samples and 1023 truly non-reagent samples were tested, and demonstrate that the HIV Tri Line kit has > 99,9% sensitivity (192/192) and 99,9% of specificity (1022/1023).

Cross reactive test were carried out using positive anti-HCV serum, anti-TP (syphilis) and anti-HBc, and no reactivity was observed.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 5 different samples, obtaining the following results:

Samples	Repetition Nº	Expected Result	Found result
Human Blood 01	20	100% Reagent	100% Reagent
Human Blood 02	20	100% Non-Reagent	100% Non-Reagent
Human Serum 01	20	100% Reagent	100% Reagent
Human Serum 02	20	100% Reagent	100% Reagent
Human Plasma 01	20	100% Non-Reagent	100% Non-Reagent

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Samples	Nº of Days	Repetition Nº	Expected Result	Found result
Human Blood 01	03	20	100% Reagent	100% Reagent
Human Blood 02	03	20	100% Non-Reagent	100% Non-Reagent
Human Serum 01	03	20	100% Reagent	100% Reagent
Human Serum 02	03	20	100% Reagent	100% Reagent
Human Plasma 01	03	20	100% Non-Reagent	100% Non-Reagent

DIAGNOSTIC SIGNIFICANCE

The Human Immunodeficiency Virus (HIV) is a retrovirus, identified with etiologic agent of the Acquired Immunodeficiency Syndrome (AIDS). AIDS is characterized by changes in T-lymphocytes population, that play an important rule in the immunological system. In the infected individual the virus causes a reduction of the T-cells subpopulation called T "helper", that leave these patients susceptible to opportunist infections and some malignances. The main transmission ways are: unprotected sexual act, blood contamination or blood derivatives and the transmission from mother to child during birth. The HIV is constitute of a RNA molecule, protected by a capsid and an envelop. There are two kinds of HIV virus: HIV 1 and HIV 2. The envelop of HIV is the major target of the immune response. The presence of the virus makes the immune system of the patients produce anti-HIV antibodies. The detection of this antibodies can be used as a diagnosis tool.

NUMBER OF TESTS

The number of tests corresponds to the number of cassettes in each presentation.

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QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE
Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for HIV Tri Line kit: 10269360148

Review: April/2018

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

MANUFACTURED BY

CONTROL

CONTROL +

CONTROL -

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED