IMUNOGLOBULINA M

REF K063

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method of determination of IgM. Immunoturbidimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry

The reagent allows the quantitative determination of IgM in human serum by immunoprecipitation, in the presence of an activator polymer that increases the sensitivity and speed of the immunoturbidimetric assay. The Immunoglobulins form with the specific antiserum an insoluble complex, producing turbidity, whose intensity increases the absorbance intensity is proportional to the concentration of IgM in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: PEG 6.000 < 15%, buffer and preservative.

Number 2 - Antiserum - Store between 2 and 8°C. Contains: Goat antiserum for Immunoglobulin M < 8 g/L, buffer and preservative.

PRESENTATION

Presentations	Reagent	
	Nº 1	Nº 2
1	25 mL	2 mL
2	2 x 25 mL	2 x 2 mL
3	3 x 25 mL	3 x 2 mL
4	4 x 25 mL	4 x 2 mL
7	20 mL	5 mL
8	2 x 20 mL	2 x 5 mL
9	3 x 20 mL	3 x 5 mL
10	4 x 20 mL	4 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the Biochemistry Automated Equipment, Multical kit and Multicontrol Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport temperature up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Reagents nº 1 and 2 contain Sodium Azide, irritant to skin and mucous. Handle carefully.
- Do not freeze the reagents, the freezing of antibody or buffer can affect their performance.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Fresh serum or plasma obtained with heparin or EDTA. The analyte is stable for 2 days at 2 to 8°C. If you need to store for a longer period, the sample should be frozen.

PROCESS DESCRIPTION

PREPARATION OF WORK REAGENT

Reagents are ready to use.

TECHNIQUE

For calibration and control of the reaction, utilize the Multical kits and Multicontrol Bioclin.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

PROCESS LIMITATIONS

Hemolysis, icterus and lipemia (low or moderate) do not interfere in test performance, using the recommended dilution for the sample. Gross lipemia and turbid samples should be processed with a greater dilution.

Turbidimetric assays of Immunoglobulin cannot discriminate monoclonal, polyclonal and oligoclonal elevations. However, monotypic Immunoglobulins reach excess Antigen at concentrations much lower than in the case of Immunoglobulin polyclonal normal. It is recommended that testing of these proteins be held in serum dilutions which allow readings in lower limits of the calibration curve. It is suggested that such cases results are verified by immunoelectrophoresis. Usually the phenomenon of excess antigen is characterized by reduced reaction coefficient of Immunoprecipitation.

It is recommended that tests are performed in duplicate when is not possible to control by immunoelectrophoresis. Perform proposed dilution of 1:10 and also 1:20. Calculate the concentration of 1:20 dilution found by multiplying the results by 1.961. The result cannot make a significant difference (± 15%) of dilution 1:10.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

Age	Unit mg/dL
Newborns	9 to 35
1 to 3 months	10 to 87
3 to 5 months	15 to 120
6 to 12 months	16 to 105
1 to 2 years	20 to 190
2 to 4 years	20 to 200
4 to 6 years	20 to 210
6 to 9 years	20 to 175
9 to 11 years	30 to 155
11 to 15 years	30 to 150
Adults	38 to 280

These values should be used as a guideline, and each laboratory should establish its range of reference values, according with the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

RECOVERY

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 94 and 103%.

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

Immunoglobulin M kit was compared with other commercially available methods for determination of IgM dosage. 40 tests were performed and the results were evaluated. The linear equation obtained was $Y = 0,948X - 1,169$ and correlation coefficient 0,991. With these results we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (mg/dL)	80,00	47,78
Standard Deviation (mg/dL)	0,42	0,12
Coefficient of Variation (%)	0,53	0,26

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (mg/dL)	80,06	47,81
Standard Deviation (mg/dL)	0,18	0,30
Coefficient of Variation (%)	0,22	0,63

Sensitivity

Sensitivity was calculated from 10 determinations of a sample free of IgM. The average found was 0,23 mg/dL with a standard deviation of 0,05 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,37 mg/dL.

Linearity

The reaction is linear up to the level of the highest calibrator concentration.

DIAGNOSTIC SIGNIFICANCE

IgM constitutes approximately 10% of all immunoglobulins. It is found mainly intravascular, being in a class of "precocious" antibodies present in the initial phases of diseases that trigger the humoral response. For this reason, it is useful in differential diagnosis of viral or acute parasitic infections. It is a protein that does not cross the placenta and is the only summarizes the newborn. Increased serum levels in newborns indicate stimulation of the immune system in utero by virus of rubella, cytomegalovirus, toxoplasmosis, or syphilis. It is also found on the surface of B lymphocytes, performing the function of receptor for antigens. As IgM is pentameric, is more efficient to connect antigens and the complement system. The levels of IgM can be increased in Macroglobulinemy of Waldenströrm, Acute Trypanosomiasis, Acute Actinomycosis, Carrion's Disease (Bartonella), Malaria, Infectious Mononucleosis, Lupus Erythomatosus, Rheumatoid Arthritis, and Desgamaglobulinemy. Decreased levels are found in Agammaglobulinemy, Lymphoproliferative Disorders, Lymphoid Aplasia, Myeloma IgA or IgG, Dysgammaglobulinemia and Chronic Lymphoblastic Leukemia.

NUMBER OF TESTS

Variable according to the automated equipment used. Check the number of tests on schedule.

BIBLIOGRAPHIC REFERENCES

1 - KILLINGSWORTH, L. M; SAVORY, J.: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971, 17, 936.
2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry; C. V. Mosby Company, 1987.
3 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.
4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.
5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26, 1459.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Immunoglobulin M kit: 10269360154

Review: February/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED