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FINALIDADE

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- Wu, A.H.B., ed. Tietz clinical guide to laboratory tests, 4th edition, p. 880. W.B. Saunders Company, St. Louis (2006).
- Bergmeyer, H.U., Gawehn, K., and Grassl, M. (1974) in *Methods of Enzymatic Analysis*. Second Edition, Volume I, 509-510, Academic Press, Inc., New York.
- M.N. Berry,R. D. Mazzachi, M. Pejakovlc,and M. J. Peake Enzymatic Determination of Potassium in Serum. *CLIN. CHEM.* 35/5, 817-820 (1989).

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REF K131
USAGE INSTRUCTIONS

FUNCTION

Method for the quantitative determination of Potassium in human serum. Enzymatic test, for *In vitro* diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic.
The Potassium is spectrophotometrically determinate through a kinetic assay using the Pyruvate Kinase dependent of potassium. The Pyruvate generated is converted into Lactate accompanying a conversion of NADH to NAD⁺. The correspondent reduction of the optical density of 380 nm is proportional to the concentration of the Potassium in the serum.

REAGENTS

Number 1: Buffer - Store between 2 and 8°C. Contains: LDH < 50 KU/L, NADH < 10 mmol/L, substrate, stabilizers and preservative.
Number 2: Enzymatic Reagent - Store between 2 and 8°C. Contains: Pyruvate Kinase < 50 KU/L, stabilizers and preservative.
Number 3: Standard Level 1 - Store between 2 and 8°C. Contains: Potassium Chloride (concentration varies to each lot), buffers, surfactant, stabilizers and preservatives.
Number 4: Standard Level 2 - Store between 2 and 8°C. Contains: Potassium Chloride (concentration varies to each lot), buffers, stabilizers and preservatives.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2	Reagent Nº 3	Reagent Nº 4
5	40 mL	10 mL	3 mL	3 mL
6	2 x 40 mL	2 x 10 mL	3 mL	3 mL
7	3 x 40 mL	3 x 10 mL	3 mL	3 mL
8	4 x 40 mL	4 x 10 mL	3 mL	3 mL
9	30 mL	10 mL	3 mL	3 mL
10	2 x 30 mL	2 x 10 mL	3 mL	3 mL
11	3 x 30 mL	3 x 10 mL	3 mL	3 mL
12	4 x 30 mL	4 x 10 mL	3 mL	3 mL
13	60 mL	20 mL	3 mL	3 mL
14	2 x 60 mL	2 x 20 mL	3 mL	3 mL
15	3 x 60 mL	3 x 20 mL	3 mL	3 mL
16	4 x 60 mL	4 x 20 mL	3 mL	3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
2- Strictly follow the methodology proposed to obtain exact results.
3- Water used in material cleaning must to be recent and free of contaminants.
4- Do not ingest. Avoid contact with eyes and skin.
5- The samples must be treated as potentially infectious.
6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
8- Do not use the product in case of damaged packaging.
9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum.
Do not use hemolyzed serum. The analyte is stable for 5 days between 2 and 8°C and 14 days at -20°C.

PROCESS DESCRIPTION

The calibration stability of the Potassium Ions Automated kit installed on refrigerated equipment is at least 03 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

REAGENT PREPARATION

The reagents are ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

The kit is indicated for biochemistry automated analyzers only.
Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.

RESULTS

The results are expressed in mmol/L.

PROCEDURE LIMITATIONS

Periodically test the calibration to detect possible changes in the response of the instrument.

INTERFERENCE

The test does not suffer interference from the substances with the top concentrations of:
Na⁺ 150 mmol/L
NH₄⁺ 0,5 mmol/L
Ca²⁺ 7,5 mmol/L
P_i 2,0 mmol/L

Ascorbic Acid 10,0 mmol/L
Zn²⁺ 0,5 mmol/L
Fe³⁺ 0,5mmol/L
Cu²⁺ 0,5 mmol/L
Triglycerides 1000 mg/dL
Hemoglobin 500 mg/dL
Bilirubin conjugated 20 mg/dL
Bilirubin non conjugated 15 mg/dL
It should only be dosed in serum, free of anticoagulants. The sample must be processed as soon as possible after collection, with a maximum dosage of 5 days for serum frozen after centrifugation. In the case of serum freezing is no possibility that unknown substances outside the listed substances may cause interferences in dosages.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standards are traceable to the reference material NIST (National Institute of Standards and Technology) SRM 918B.

REFERENCE VALUES

The reference values in mmol/L, for this method, were obtained through the determination of potassium in healthy populations of male and female.

Serum..... 3,5 a 5,1 mmol/L
13,7 - 19,9 mg/dL

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.
The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY
The Potassium Ions Automated was compared with another method commercially available to measure Potassium. 42 analysis were conducted and the results were evaluated. The linear equation was obtained Y = 0,9802X + 0,0506, with a correlation coefficient equal to 0,9825. With these results we can conclude that the kit has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mmol/L)	3,09	5,03	7,54
Standard Deviation (mmol/L)	0,08	0,08	0,10
Coefficient of Variation (%)	2,56	1,58	1,33

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

	Sample 1	Sample 2	Sample 3
Average Concentration (mmol/L)	3,09	5,03	7,54
Standard Deviation (mmol/L)	0,08	0,08	0,10
Coefficient of Variation (%)	2,56	1,58	1,33

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mmol/L)	3,13	5,05	7,63
Standard Deviation (mmol/L)	0,05	0,04	0,13
Coefficient of Variation (%)	1,74	0,70	1,65

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Potassium. The average found was of -0,055 mmol/L with standard deviation of 0,028 mmol/L. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,030 mmol/L.

Linearity

The Potassium Ions Automated kit is linear between 2 and 8 mmol/L.

DIAGNOSTIC SIGNIFICANCE

In health subjects, the level of the liquid extra cellular potassium is regulated to be maintained at 3,5 - 5,1 mmol/L. Small deviation in the normal levels may cause severe consequences in one's health. The monitoring of the potassium serum concentration is important to both routine and emergency routines.

BIBLIOGRAPHIC REFERENCES

1. Wu, A.H.B., ed. Tietz clinical guide to laboratory tests, 4th edition, p. 880. W.B. Saunders Company, St. Louis (2006).
2. Bergmeyer, H.U., Gawehn, K., and Grassl, M. (1974) in *Methods of Enzymatic Analysis*. Second Edition, Volume I, 509-510, Academic Press, Inc., New York.
3. M.N. Berry,R. D. Mazzachi, M. Pejakovic,and M. J. Peake
Enzymatic Determination of Potassium in Serum. *CLIN. CHEM.* 35/5, 817-820 (1989).

 QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC

REP

 OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Potassium Ions Automated kit: 10269360209

Review: January/2020

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 BATCH CODE

 DATE OF MANUFACTURE

 USED BY
(last day of month)

 TEMPERATURE LIMITATION
(store at)

 CONTAINS SUFFICIENT
FOR <N> TESTS

 CONSULT INSTRUCTIONS
FOR USE

 IN VITRO DIAGNOSTIC DEVICE

 EUROPEAN AUTHORIZED
REPRESENTATIVE

 MANUFACTURED BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 INFLAMMABLE

 CORROSIVE

 POISON

 CE MARK