FINALIDADE

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of Lactate in plasma or human cephalorachidian liquid. UV Enzymatic test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: UV Enzymatic using the Dehydrogenase Lactate.

Lactic Dehydrogenase (LDH) catalyzes the oxidation of L-Lactate to Pyruvate, with consequent reduction of NAD^+ to NADH. The L-Lactate concentration is measured by the increase in absorbance at 340 nm, from the NADH formed.

REAGENTS

Number 1 - Enzyme Reagent - Store between 2 and 8°C. Contains: LDH < 100 KU/L, buffer, stabilizers and preservative.

Number 2 - Coenzyme - Store between 2 and 8°C. Contains: NAD < 20 mmol/L, buffers and preservative.

Number 3 - Standard - Store between 2 and 8°C. Contains: Lithium L-Lactate 30 mg/dL and preservative.

PRESENTATION

Presentation	Reagent n°1	Reagent n°2	Reagent n°3
1	1 x 20 mL	1 x 5 mL	1 x 3 mL
2	2 x 20 mL	2 x 5 mL	1 x 3 mL
3	3 x 20 mL	3 x 5 mL	1 x 3 mL
4	4 x 20 mL	4 x 5 mL	1 x 3 mL
5	1 x 40 mL	1 x 10 mL	1 x 3 mL
6	2 x 40 mL	2 x 10 mL	1 x 3 mL
7	4 x 40 mL	4 x 10 mL	1 x 3 mL
8	2 x 40 mL	1 x 20 mL	1 x 3 mL
9	4 x 40 mL	2 x 20 mL	1 x 3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N Bioclin and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature is between 2 and 8°C and the transport temperature shall be up to 30°C, not exceeding 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- The water level in water-bath must be superior to the reagents level in the test tubes.

6- Do not use jaundice or hemolyzed samples.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Costumer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used ar properly calibrated and subjected to periodic maintenance.

SAMPLE

Plasma or cephalorachidian liquid. Do not use serum. As an anticoagulant use glycolytic inhibitors, example, fluoride/oxalate or fluoride/heparin. Do not use jaundice or hemolyzed samples. Plasma stability: 8 hours up to 20 and 25°C or 14 days up to 2 and 8°C.

PROCESS DESCRIPTION

The calibration stability of the Lactate kit installed on refrigerated equipment is at least 25 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

REAGENT PREPARATION

Mixture 4 parts of Reagent N°1 with 1 part of Reagent N°2 (Ex.: 20 mL R1 + 5 mL R2) = Monoreagent
Stability: 2 weeks up to 2 and 8°C.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes with letters B (Blank), P (Standard), A (Sample) and proceed as follows:

	Blank	Standard	Sample
Standard	---	10 μL	---
Water	10 μL	---	---
Sample	---	---	10 μL
Monoreagent	1.0 mL	1.0 mL	1.0 mL

Smoothly homogenize, incubate for 5 minutes at a 37°C $\pm$ 2°C and determine the Sample absorbance at 340 nm, matching the zero with the Blank. Read the absorbance in, no longer, than 30 minutes.

CALCULATIONS

$$\text{Lactate (mg/dL)} = \frac{\text{Sample Abs.} \times \text{Standard Concentration}}{\text{Standard Absorbance}}$$

Following the Lambert-Beer law, we can use the Calibration Factor:

$$\text{Calibration Factor} = \frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$$

$$\text{Lactate (mg/dL)} = \text{Sample absorbance} \times \text{Calibration factor}$$

The values are expressed in mg/dL. To convert the values to mmol/L, multiply the obtained value by 0.1109:

$$\text{Lactate (mmol/L)} = \text{Lactate (mg/dL)} \times 0.1109$$

INTERFERENTS

Ascorbic Acid above 30 mg/dL, direct and indirect Bilirubin above 60 mg/dL, Lipemia above 2000 mg/dL of Triglycerides, Hemoglobin above 1000 mg/dL, Dopamine above 10 mg/dL, L-Dopamine above 20 mg/L, Methyl dopamine above 10 mg/L and Glycolic Acid above 1200 mg/L can interfere in Lactate dosage.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of Lactate in healthy populations of male and female.

PLASMA		
	mg/dL	mmol/L
Vein	4.5 - 19.8	0.5 - 2.2
Arterial	4.5 - 14.4	0.5 - 1.6
CEPHALORACHIDIAN LIQUID		
	mg/dL	mmol/L
Adults	10 - 22	1.1 - 2.4
Newborns	10 - 60	1.1 - 6.7
3 - 10 days	10 - 40	1.1 - 4.4
> 10 days	10 - 25	1.1 - 2.8

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy RECOVERY

The recovery analysis was made with 05 sample determinations. The accuracy was calculated and was found in a good concordance with the reference values, obtaining a recovery between 94 and 104%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Lactate kit was compared with other method commercially available for dosage of Lactate. 42 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1.008X - 0.6958$ and the correlation coefficient 0.9991. With these results, we can conclude that the kit shows good methodological specificity.

Precision REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	30.34	20.14	40.11
Standard Deviation (mg/dL)	0.16	0.23	0.65
Coefficient of Variation (%)	0.52	1.16	1.62

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	30.11	19.30	40.26
Standard Deviation (mg/dL)	0.21	0.73	0.59
Coefficient of Variation (%)	0.68	3.77	1.48

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Lactate. The average found was -0.438 mg/dL with standard deviation of 0.726 mg/dL. The sensitivity, that indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 1.740 mg/dL.

Linearity

The reaction is linear up to 120 mg/dL. For values above 120 mg/dL, dilute the sample in Sodium Chloride solution 9 g/L, repeat the dosage and multiply the obtained result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Lactate is the final product of the anaerobic glycolysis and it is used as an indicative of the oxygen state in cellular tissues.

The elevated Lactate levels in blood occur in anoxia due to shock, congestive cardiac stroke, intoxication and thiamine's deficiency. That is why the Lactate is dosed with extreme caution.

The determination of muscular Lactate is used to evaluate the metabolism variation in athletes training.

BIBLIOGRAPHIC REFERENCES

1. Thomas L. Clinical Laboratory Diagnostics 1st ed. Frankfurt: THBooks Verlagsgesellschaft; 1998. p. 160-166.
2. David B Sacks, M. B., Ch.B. F.A.C.P. Carbohydrates In: Burtis CA, Ashwood ER, editors. Tietz Book of Clinical Chemistry. 3 ed. Philadelphia: W.B Saunders Company; 1999. p. 787-790.
3. Westgard J.O., Lahmeyer B.L., Bimbaum M.L. Use of the Du Pont "Automatic Clinical Analyser" in Direct Determination of Lactic Acid in Plasma Stabilized with sodium Fluoride. Clin chem 1972; 18:1334-8.
4. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Lactate kit: 10269360143

Review: June/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED