FINALIDADE

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USAGE INSTRUCTIONS



FUNCTION

Method for quantitative determination of Lipase in samples of human serum or plasma heparinized. Colorimetric test, for *in vitro* diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric
The Lipase chromogenic substrate 1.2 - O - Dilauryl - Rac - Glicero - 3 - Acid - Glutaric - (6 - Methylresorufin) - Ester is cleaved by the catalytic action of alkaline lipase forming 1.2 - O - Dilauryl - Rac - Glycerol and a intermediate unstable, Glutaric Acid - (6 - Methylresorufin) Ester. This spontaneously decompose to form Glutaric Acid in alkaline solution and Methylresorufin. The intensity of red color is formed directly proportional to the Lipase activity and can be determined photometrically.

REAGENTS

Number 1 – Buffer Reagent – Store between 2 and 8°C. Contains: Buffer BICIN (N,N-Bis(2-Hydroxyethyl)-Glycine) 50 mmol/L (pH 8.0), Colipase > 0.9 mg/L, Na-Deoxycolate 1.6 mmol/L, Calcium Chloride 10 mmol/L, detergent and preservative.
Number 2 – Substrato – Store between 2 and 8°C. Contains: Buffer Tartrate 10 mmol/L (pH 4.0), 1.2-O-Dilauryl-Rac-Glycero-3-Ácido Glutaric-(6-Methylresorufin) Ester 0.27 mmol/L, Taurodeoxicolato 8.8 mmol/L, detergent and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 20 mL	1 x 5 mL
2	2 x 20 mL	2 x 5 mL
3	3 x 20 mL	3 x 5 mL
4	4 x 20 mL	4 x 5 mL
5	5 x 20 mL	5 x 5 mL
6	1 x 40 mL	1 x 10 mL
7	2 x 40 mL	2 x 10 mL
8	4 x 40 mL	4 x 10 mL
9	2 x 40 mL	1 x 20 mL
10	4 x 40 mL	2 x 20 mL
11	1 x 15 mL	1 x 10 mL
12	2 x 15 mL	2 x 10 mL
13	3 x 15 mL	3 x 10 mL
14	4 x 15 mL	4 x 10 mL
15	1 x 30 mL	1 x 20 mL
16	2 x 30 mL	2 x 20 mL
17	3 x 30 mL	3 x 20 mL
18	4 x 30 mL	4 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes and tips, test tubes, watch or stopwatch, Biocal Bioclin, Biocontrol N and Biocontrol P Bioclin. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transportation temperature should be between 2 and 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid humidity.
Do not freeze.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in the cleansing of the material should be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Reagent N° 2 must be protected from light.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or heparinized plasma can be uses.
Serum samples are stable for 7 days between 2 and 30°C and 12 months when frozen at -20°C.
Plasma samples are stable for 7 days between 2 and 30°C and 12 months when frozen at -20°C.

PROCESS DESCRIPTION

The calibration stability of the Lipase Automated kit installed on refrigerated equipment is at least 28 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

For the use of Lipase Automated, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol N and Biocontrol P Bioclin kits.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

INTERFERENT

Studies have been performed to determine the level o interference of different compounds: Conjugate Bilirubin and non conjugate do not produce interference up to 60 mg/dL; Hemoglobin does not produce significant interference up to 500 mg/dL; Lipids do not produce significant interference up to 1000 mg/dL; Rare cases of gammopathy, particularly the IgM type (Waldenström macroglobulinemia) may cause unreliable results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The calibration kit must be carried on using the BIOCAL calibrator, which is traceable to the reference material BCR-694 (Community Bureau of Reference).

REFERENCE VALUES

Adults: < 60 U/L

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS, METHODOLOGICAL SPECIFICITY AND SENSITIVITY

Lipase Automated kit was compared with other method commercially available for dosage of Lipase. 42 analyzer were performed and the results were evaluated. The linear equation obtained was $Y = 0.8219X + 3.2388$ and correlation coefficient 0.9978. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	50.30	139.20	11.45
Standard Deviation (U/L)	0.85	1.16	0.81
Coefficient of Variation (%)	1.70	0.83	7.12

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	51.28	144.25	12.78
Standard Deviation (U/L)	0.91	6.93	1.17
Coefficient of Variation (%)	1.77	4.80	9.17

Sensibility

Sensitivity of Lipase Automated kit was calculated from 40 determinations of a sample free of Lipase. The average found was of 2.063 U/L with a standard deviation of 0.046 U/L. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 2.201 U/L.

Linearity

Test is linear up to 300 U/L.

DIAGNOSTIC SIGNIFICANCE

Lipase is an enzyme produced almost exclusively by pancreas. In the presence of pancreatitis, the lipase is released at high concentrations in the blood, persisting for much longer than the Amylase, allowing confirmation of late diagnosis disease. The values of Lipase in pancreatitis can reach from 5 to 10 times the upper limit of reference, while the other diseases Lipase levels are always below 3 times the values of reference. Lipase increases in acute pancreatitis, pancreatitis disease (cancer of pancreas), obstruction of the pancreatic ducts by calculation or malignancies. Since the excretion of Lipase is via the kidneys, kidney failure in general, Lipase levels are high. In intestinal infarction, in acute cholangitis and in small bowel obstruction increased values of Lipase can also be observed. The lipase may be reduced physiologically during pregnancy, in the course of tuberculosis and other infectious diseases.

BIBLIOGRAPHIC REFERENCES

- 1 - Steinberg WM, Goldstein SS, Davies ND et al. Diagnostics assays in acute pancreatitis [Review]. Ann Intern Med 1985; 102:576-580.
- 2 - Panteghini M et al. Diagnostic value of four assays for lipase determination in serum: A comparative reevaluation. Clin Biochem 1991; 24:497-503.
- 3 - TIETZ, N. W., Chlinical Guide to Laboratory Tests., 3ª edição. Philadelphia, PA: WB Saunders, 1995: 865.
- 4 - Neumann U et al. New substrates for the optical determination of lipase. EP 207252 (1987).
- 5 - Gargouri Y, Julien R, Bois A et al. Studies on the detergent inhibition of pancreatic lipase activity. J of Lipase Research 1983; 24:1336-1342.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA resgistration for Lipase Automated kit: 10269360158

Review: July/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED